TRUE™ DILATATION

Balloon Valvuloplasty Catheter

ENGLISH

INSTRUCTIONS FOR USE

CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

Device Description

The TRUE™ DILATATION Balloon Valvuloplasty Catheter is an over-the-wire coaxial catheter with a balloon fixed at the tip. The catheter is available in 110cm and 55cm lengths, and has two lumens: one lumen is used to inflate and deflate the balloon and the other permits the use of a guidewire to position the catheter. The balloon inflation luer-lock hub (angled) connects to a syringe inflation device to deliver radiopaque contrast media for inflation. The guidewire luer-lock hub (straight) connects to the guidewire lumen. The balloon is non-compliant and is designed to reach a known diameter and length when inflated within the specified pressure range. Two radiopaque marker bands are provided for fluoroscopic positioning of the device across the aortic valve. These bands are positioned at the proximal and distal balloon shoulders. Balloon catheter dimensions, balloon nominal pressure, maximum inflation pressure, recommended introducer size, and recommended guidewire size are indicated on the package label.

Indications for Use

The TRUE™ DILATATION Balloon Valvuloplasty Catheter is indicated for balloon aortic valvuloplasty.

Contraindications

The TRUE™ DILATATION Balloon Valvuloplasty Catheter is contraindicated for use in patients with annular dimensions < 18 mm.

Warnings

1. Contents supplied **STERILE** using ethylene oxide (EO). Non-pyrogenic. Do not use if sterile barrier is opened or damaged. Single patient use only. Do not reuse, reprocess or re-sterilize.
2. This device has been designed for single use only. Reusing this medical device bears the risk of cross-patient contamination as medical devices – particularly those with long and small lumina, joints, and/or crevices between components – are difficult or impossible to clean once body fluids or tissues with potential pyrogenic or microbial contamination have had contact with the medical device for an indeterminate amount of time. The residue of biological material can promote the contamination of the device with pyrogens or microorganisms which may lead to infectious complications.
3. Do not resterilize. After reesterilization, the sterility of the product is not guaranteed because of an indeterminate degree of potential pyrogenic or microbial contamination which may lead to infectious complications. Cleaning, reprocessing and/or reesterilization of the present medical device increases the probability that the device will malfunction due to potential adverse effects on components that are influenced by thermal and/or mechanical changes.
4. Catheter balloon inflation diameter must be carefully considered in selecting a particular size for any patient. It is critical to perform a clinical diagnostic determination of valve anatomical dimensions prior to use; imaging modalities such as transthoracic echocardiogram (TTE), computerized tomography (CT), angiography, and/or transesophageal echocardiogram (TEE) should be considered. The inflated balloon diameter should not be significantly greater than valvular diameter.
5. When the catheter is exposed to the vascular system, it should be manipulated while under high-quality fluoroscopic observation. Do not advance or retract the catheter unless the balloon is fully deflated. If resistance is met during manipulation, determine the cause of the resistance before proceeding. Applying excessive force to the catheter can result in tip breakage or balloon separation, or cause injury to the patient (such as vessel perforation).
6. If flow through catheter becomes restricted, do not attempt to clear catheter lumen by infusion. Doing so may cause catheter to rupture, resulting in vessel trauma. Remove and replace catheter.
7. Do not exceed the RBP recommended for this device. Balloon rupture may occur if the RBP rating is exceeded. To prevent over-pressurization, use of a pressure monitoring device is recommended.
8. After use, this product may be a potential biohazard. Handle and dispose of in accordance with acceptable medical practices and applicable local, state, and federal laws and regulations.
9. If using device to support Transcatheter Aortic Valve Implantation (TAVI), consult TAVI system's Instructions for Use for any additional procedural instructions related to selection and use of valvuloplasty balloon.

Precautions

1. Carefully inspect the catheter prior to use to verify that catheter has not been damaged during shipment and that its size, shape and condition are suitable for the procedure for which it is to be used. Do not use if product damage is evident.

2. The catheter should only be used by physicians trained in the performance of percutaneous transluminal valvuloplasty.
3. The minimal acceptable French size is printed on the package label. Do not attempt to pass the catheter through a smaller size sheath introducer than indicated on the label.
4. Use the recommended balloon inflation medium of 1/3 to 2/3 contrast to saline ratio. Never use air or other gaseous medium to inflate the balloon.
5. If resistance is felt during post procedure withdrawal of the catheter through the introducer sheath, determine if contrast is trapped in the balloon with fluoroscopy. If contrast is present, push the balloon out of the sheath and then completely evacuate the contrast before proceeding to withdraw the balloon.
6. If resistance is still felt during post procedure withdrawal of the catheter, it is recommended to remove the balloon catheter and guidewire/introducer sheath as a single unit.
7. In the very unlikely event of balloon burst or rupture, balloon could be more difficult to remove through the sheath and could require introducer sheath removal.
8. Do not torque, excessively bend catheter or continue to use if the shaft has been bent or kinked.
9. Prior to re-insertion through the introducer sheath, the balloon should be wiped clean with gauze and rinsed with sterile normal saline.
10. Do not remove guidewire from catheter during procedure.
11. Dilation procedures should be conducted under high-quality fluoroscopic guidance.
12. Careful attention must be paid to the maintenance of tight catheter connections. Aspirate before proceeding to avoid air introduction into the system.
13. If inflating balloon in patient to facilitate re-folding, ensure balloon is positioned so that it can be inflated safely.

Potential Adverse Reactions

The complications which may result from a percutaneous transluminal valvuloplasty procedure include:

- Additional intervention
- Allergic reaction to drugs or contrast medium
- Aneurysm or pseudoaneurysm
- Arrhythmias
- Cardiovascular injury
- Conduction system injury
- Embolization
- Hematoma
- Hemorrhage, including bleeding at the puncture site
- Hypotension/hypertension
- Inflammation
- Occlusion
- Pain or tenderness
- Pneumothorax or hemothorax
- Sepsis/infection
- Shock
- Short term hemodynamic deterioration
- Stroke
- Thrombosis
- Valvular tearing or trauma
- Vessel dissection, perforation, rupture, or spasm

Directions for Use

Handling & Storage

Store in a cool, dry, dark place. Do not store near radiation or ultra-violet light sources.

Rotate inventory so that the catheters and other dated products are used prior to the "Use By" date.

Do not use if packaging is damaged or opened.

Equipment for Use

- Contrast medium
- Sterile saline solution
- Luer lock syringe/inflation device with manometer (50 ml or larger)
- Appropriate introducer sheath and dilator set
- .035" guidewire
- 3-way high pressure stopcock

Dilatation Catheter Preparation

1. Remove catheter from package. Verify the balloon size is suitable for the procedure and the selected accessories accommodate the catheter as labeled.
2. Keep balloon guard on balloon until removing it immediately prior to inserting catheter into introducer.

3. Mix contrast / saline with a ratio of 1/3 contrast, 2/3 saline by volume.
4. Prior to use, the air in the balloon catheter should be removed. To facilitate purging, select a syringe or inflation device with a 50 ml or larger capacity and fill approximately half of it with the appropriate balloon inflation medium. Do not use air or any gaseous medium to inflate the balloon.
5. Connect a stopcock to the balloon inflation female luer hub on the dilatation catheter.
6. Connect the syringe to the stopcock.
7. Hold the syringe with the nozzle pointing downward, open the stopcock and aspirate for approximately 15 seconds. Release the plunger.
8. Repeat step #7 two more times or until bubbles no longer appear during aspiration (negative pressure). Once completed, evacuate all air from the barrel of the syringe/inflation device.
9. Prepare the wire lumen of the catheter by attaching a syringe to the wire lumen hub and flushing with sterile saline solution.

Use of the TRUE™ DILATATION Balloon Valvuloplasty Catheter

1. Prepare patient and place introducer and guidewire into access site using standard percutaneous catheterization techniques.
2. Remove balloon guard distally, ensuring guard is no longer on catheter.
3. Backload the distal tip of the TRUE™ DILATATION Balloon Valvuloplasty Catheter over the prepositioned guidewire and advance the tip to the introduction site.
4. Using fluoroscopic guidance, advance the catheter through the introducer sheath and over the wire to the site of inflation. Discontinue product use if marker bands are not visible under fluoroscopy.
5. Position the balloon within the relevant area of the aortic valve to be dilated, ensure the guidewire is in place, and while ensuring the balloon is held in a static position, inflate the balloon to a pressure not greater than RBP.
6. Apply negative pressure to fully evacuate fluid from the balloon. Confirm that the balloon is fully deflated under fluoroscopy.
7. While maintaining negative pressure and the position of the guidewire, withdraw deflated catheter over the guidewire and through the introducer sheath. Use of a gentle clockwise twisting motion may be used to help facilitate catheter removal through the introducer sheath.
8. If unusual resistance is met when attempting to withdraw balloon, position balloon in anatomical position in which it can be inflated safely. Inflate balloon and then deflate it. Balloon re-folding can be observed under fluoroscopy. Use of recommended contrast concentration will facilitate fluoroscopic visualization of balloon re-folding.

Balloon Reinsertion

Precaution: Do not continue to use the catheter if the shaft is bent or kinked.

Precaution: Prior to re-insertion through the introducer sheath, the balloon should be wiped clean with gauze and rinsed with sterile normal saline.

1. Load the catheter onto a guidewire.
2. Advance the catheter over the pre-positioned guidewire to the introduction site and through the introducer sheath. If resistance is encountered, replace the previously used balloon catheter with a new balloon.
3. Continue the procedure according to the "Use of the TRUE™ DILATATION Balloon Valvuloplasty Catheter" section herein.

Catalogue Number	0184511/ T184511	0204511/ T204511	0214512/ T214512	0224512/ T224512	0234512/ T234512	0244512/ T244512	0254513/ T254513	0264513/ T264513	0284514/ T284514
Balloon Description	18 mm x 4.5 cm	20 mm x 4.5 cm	21 mm x 4.5 cm	22 mm x 4.5 cm	23 mm x 4.5 cm	24 mm x 4.5 cm	25 mm x 4.5 cm	26 mm x 4.5 cm	28 mm x 4.5 cm
Actual Diameter (mm)	17.9 mm - 18.2 mm	20.1 mm - 20.4 mm	21.0 mm - 21.2 mm	22.1 mm - 22.4 mm	22.9 mm - 23.2 mm	24.1 mm - 24.4 mm	25.0 mm - 25.3 mm	26.1 mm - 26.4 mm	28.0 mm - 28.3 mm
Over Indicated Pressure Range (3atm-6atm)									

Clinical Reference

The following is a summary of clinical data collected between January and August 2012 for cases performed in Europe:

- For 69 BAV procedures performed using 69 TRUE™ DILATATION balloons during this time period, the involved physicians completed data collection forms reporting on procedure information, device performance, and device-related adverse events. Sixty-eight (98.5%) of these BAV procedures were performed as part of a transcatheter aortic valve implantation (TAVI) procedure. There were no (0) device-related adverse events reported.

Warning: After use, this product may be a potential biohazard. Handle and dispose of in accordance with acceptable medical practices and applicable local, state and federal laws and regulations.

An issue or revision date and revision number for these instructions are included for the user's information on the last page of this booklet.

In the event 36 months have elapsed between this date and product use, the user should contact Bard Peripheral Vascular to see if additional product information is available.

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EC REP

Authorised Representative in the European Community

Représentant autorisé au sein de la Communauté européenne
Bevollmächtigter in der Europäischen Gemeinschaft
Rappresentante autorizzato nella Comunità Europea
Representante autorizado en la Comunidad Europea
Gemachtigde vertegenwoordiger binnen de Europese Gemeenschap
Representante autorizado na Comunidade Europeia
Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα
Autoriseret repræsentant i Det Europæiske Fællesskab
Auktoriserad representant i EG
Valtuutettu edustaja Euroopan yhteisössä
Autorisert representant i EU
Autoryzowany przedstawiciel na terenie Unii Europejskiej
Bejegyzett képviselő az Európai Közösségben
Autorizovaný zástupce pro Evropskou unii
Avrupa Topluluğu'ndaki Yetkili Temsilci
歐洲共同體的授權代表
유럽 공동체 내의 인증받은 대리업체
Уполномоченный представитель в Европейском сообществе

Protect From Heat

Protéger de la chaleur
Vor Hitzeinwirkung schützen
Proteggere dal calore
Proteger del calor
Tegen hitte beschermen
Proteger do calor
Προφύλασσετε το προϊόν από τη θερμότητα
Beskyttes mod varme
Skyddas från värme
Suojattava kuumuudelta
Beskyttes mot varme
Chronic przed wysokimi temperaturami
Hőtől védendő
Chraňte před teplem
Isidan Koryununuz
不要受熱
열로부터 보호하십시오
Защищать от нагревания

Do Not Use if the Product Sterile Barrier System or its Packaging is Compromised

Ne pas utiliser si le système de barrière stérile du produit ou son emballage sont ou paraissent endommagés.
Bei beeinträchtigtem Sterilschutzsystem oder beschädigter Verpackung das Produkt nicht verwenden
Non utilizzare se il sistema della barriera sterile del prodotto o la confezione sono compromessi
No utilizar si el sistema de barrera estéril del producto o su envase están dañados
Niet gebruiken wanneer het steriele-barrièresysteem van het product of de verpakking is aangetast
Não utilizar se o sistema de barreira estéril do produto ou respectiva embalagem estiverem comprometidos
Μη χρησιμοποιείτε το προϊόν εάν έχει παραβιαστεί το σύστημα φραγμού αποστείρωσης ή η συσκευασία του
Må ikke anvendes, hvis produktets sterilbarrieresystem eller emballagen er kompromitteret
Använd inte om produktens sterila barriär eller förpackning är skadad
Ei saa käyttää, jos tuotteen steriilisuojaus tai pakkaus on vaurioitunut
Skal ikke brukes hvis produktets sterile barrieresystem eller pakningen er forringet
Nie używać, jeśli naruszono sterylność produktu lub jego opakowanie
Ne használnia, ha a termék védő steril lezárási rendszer vagy a csomagolás sérült!
Pokud je narušena sterilní bariéra nebo poškozen obal, výrobek nepoužívejte
Ürünün Steril Bariyer Sistemi ya da Ambalajı Bozulmuşsa Ürünü Kullanmayın
如果產品消毒屏障系統或包裝受損，請勿使用
제품의 멸균 장벽 시스템 또는 포장에 손상된 경우 사용하지 마십시오.
Запрещается применять изделие, если его стерильная упаковка или внешняя упаковка повреждены

SL

Shaft Length

Longueur de la gaine
Schaftlänge
Lunghezza del corpo
Longitud del cuerpo
Hulslänge
Comprimento do eixo
Μήκος στελέχους
Skaftlængde
Skaftlängd
Varren pituus
Skaftlengde
Długość trzonu
Tengely hossza
Délka dráku
Mil Uzunluğu
軸長度
샤프트 (shaft) 길이
Длина стержня

Not Made with Natural Rubber Latex

Fabriqué sans latex naturel
Nicht aus Naturgummlatex hergestellt.
Non prodotto con lattice di gomma naturale
Este producto no se fabrica con látex de caucho natural
Niet gemaakt met latex van natuurlijk rubber
Não contém látex de borracha natural
Δεν κατασκευάζεται από φυσικό ελαστικό λάτεξ
Ikke fremstillet med naturlig gummlatex
Ej tillverkad med naturgummlatex
Valmistuksessa ei ole käytetty luonnonkumilateksia
Dette produktet er ikke produsert med naturlig gummlateks
Produkt nie jest wytwarzany z lateksu kauczuku naturalnego
Nem tartalmaz természetes gumilatexet.
Není vyroben z přírodního pryžového latexu
Doğal Kauçuk Lateksten Üretimimştir
製造未採用天然膠乳
천연 고무 라텍스를 사용해 제조되지 않음
Не содержит натуральный каучуковый латекс

Bard and True are trademarks and/or registered trademarks of C. R. Bard, Inc. or an affiliate

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Copyright ©2014 C. R. Bard, Inc. All rights reserved. Printed in the USA.

Copyright ©2014 C. R. Bard, Inc. Tous droits réservés. Imprimé aux États-Unis.
Copyright ©2014 C. R. Bard, Inc. Alle Rechte vorbehalten. Gedruckt in den USA.
Copyright ©2014 C. R. Bard, Inc. Tutti i diritti riservati. Stampato negli Stati Uniti.
Copyright ©2014, C. R. Bard, Inc. Reservados todos los derechos. Impreso en EE.UU.
Copyright ©2014 C. R. Bard, Inc. Alle rechten voorbehouden. Gedrukt in de VS.
Copyright ©2014, C. R. Bard, Inc. Todos os direitos reservados. Impreso nos E.U.A.
Πνευματικά Δικαιώματα ©2014 C. R. Bard, Inc. Με την επιφύλαξη παντός δικαιώματος. Εκτυπώθηκε στις Η.Π.Α.
Copyright ©2014 C. R. Bard, Inc. Alle rettigheder forbeholdes. Trykt i USA.
Copyright ©2014 C. R. Bard, Inc. Alla rättigheter förbehållna. Tryckt i USA.
Copyright ©2014 C. R. Bard, Inc. Kaikki oikeudet pidätetään. Painettu Yhdysvalloissa.
Copyright ©2014 C. R. Bard, Inc. Alle rettigheter forbeholdt. Trykt i USA.
Copyright ©2014 C. R. Bard, Inc. Wszelkie prawa zastrzeżone. Wydrukowano w USA.
Copyright ©2014, C. R. Bard, Inc. Minden jog fenntartva. Nyomtatva az USA-ban.
Copyright ©2014 C. R. Bard, Inc. Všechna práva vyhrazena. Vytlačeno v USA.
Telif Hakkı ©2014 C. R. Bard, Inc. Tüm hakları saklıdır. ABD'de basılmıştır.
Copyright © 2014 C. R. Bard, Inc. 保留所有權利。在美國印製。
저작권 © 2014 C. R. Bard, Inc. 版權 소유. 미국에서 인쇄.
Авторские права © 2014 Компания C. R. Bard, Inc. Все права защищены. Напечатано в США.

Manufacturer:
Bard Peripheral Vascular, Inc.
1625 West 3rd Street
Tempe, AZ 85281
USA

Tel: 1-480-894-9515
1-800-321-4254
Fax: 1-480-966-7062
1-800-440-5376
www.bardpv.com

**Authorised Representative
in the European Community**
Bard Limited
Forest House
Tilgate Forest Business Park
Brighton Road, Crawley
West Sussex, RH11 9BP, UK

BARB | PERIPHERAL
VASCULAR