ASAHI® FUBUKI®

Neurovascular Guide Catheter Neurovascular Guide Catheter Dilator Kit

INSTRUCTIONS FOR USE

SYMBOLS

Legal manufacturer

Do not reuse

Do not resterilize

Do not use if package is damaged

Caution, consult accompanying documents

Consult instructions for use

Keep dry

Keep away from sunlight

Sterilized using ethylene oxide


Guide catheter's outer diameter

Catheter's inner diameter

Maximum injection pressure

Maximum guide wire O.D.

Caution: Federal (USA) Law restricts this device to sale by or on the order of a physician


Use by

Catalogue number

LOT number

Guide catheter's effective length

AMK-DT501 Ver.1.00 / 18TS061

Caution: Federal (USA) Law restricts this device to sale by or on the order of a physician.

Read these instructions carefully before use and observe the **Indications for Use, Contraindications, Warnings, Precautions, Malfunctions and Adverse Effects and How to Use** sections in the **Instructions for Use**. Failure to do so may result in complications, including serious injury to the patient or death.


These Instructions for Use apply to the **ASAHI FUBUKI Neurovascular Guide Catheter and ASAHI FUBUKI Neurovascular Guide Catheter Dilator Kit** (hereafter called "this product"). For specifications of this product, refer to this product label.

Descriptions

This product is composed of a guide catheter (hereafter called "this guide catheter") and a dilator (hereafter called "the supplied dilator") packed in a sterile pack. This sterile pack is packed in an outer box with these Instructions for Use (some types of this product are not supplied with the dilator.). The outer surface of this guide catheter is coated with a hydrophilic coating to provide high lubricity when the surface is wet. The inner lumen of the shaft (excluding the connector portion) is lined with a fluoropolymer layer to facilitate the manipulation of the guide wire and other devices. The shaft is entirely radiopaque. This guide catheter can withstand a pulling force of 3.43N (350gf) towards the direction of the distal end.

Shape and Structure

●Guide catheter

●Dilator

Indications for Use

This product is intended to be used to guide interventional devices for Neurovascular therapy to a lesion or a procedural site for a percutaneous intravascular procedure in the Neurovasculature. This product is also intended to be used for injection of contrast media.

Do not use this product other than for use in the Neurovasculature.

Contraindications

None

Warnings

- This guide catheter and package are presterilized with ethylene oxide gas (EOG). Therefore, be sure not to open the package until just prior to use, and use the aseptic technique in handling this guide catheter. Do not resterilize this guide catheter. If resterilized, the performance or quality of this guide catheter may be compromised and there is a risk of complications including infection.
- This guide catheter is intended for a single use only. Do not reuse it. Reusing this guide catheter may damage it, and/or cause blood vessel damage.
- Check that this guide catheter and/or package is not damaged during transportation before use. Do not use if this guide catheter and/or package is suspected to be damaged.
- Do not use this guide catheter after the expiration date indicated on the label. Discard any guide catheter that exceeds the expiration date.
- This guide catheter must be used only by a physician who is fully trained in neuroradiology treatment.
- This guide catheter must be used in an institution where emergency surgical operation can be performed immediately.
- Do not modify this guide catheter for any reason. Never make new/additional side holes on this guide catheter shaft. (Using a modified guide catheter may damage this guide catheter and cause blood vessel damage. Making a side hole by any tool available in the hospital may cause blood clot formation in this guide catheter shaft and/or deteriorate this guide catheter performance of, or damage to this guide catheter shaft.)
- Do not use this guide catheter or the supplied dilator alone when inserting into a blood vessel. Doing so may damage this guide catheter, and cause to damage to the wall of the vessel by the tip of the supplied dilator or of this guide catheter (size or stiffness).
- Operate this guide catheter carefully, and if any resistance is felt, stop the manipulation and identify the cause of the resistance under X-ray fluoroscopy. Note that when using this guide catheter in combination with other devices, especially when the clearance between the lumen of this guide catheter and the combined device is small, the combined device may be stuck in this guide catheter. (Otherwise, this guide catheter may be damaged, causing the blood vessel to be damaged.)
- Do not insert the guide wire by force or advance it rapidly when this guide catheter is bent or twisted. (The guide wire may cause perforation or damage of this guide catheter, or result in blood vessel damage.)
- When applying torque manipulation to this guide catheter that is placed along the tortuous course of the vessel, take extra care and handle with caution because the manipulation may lead to kink or the like of this guide catheter, and may result in damage of the blood vessel or damage to this guide catheter.
- During procedures, use appropriate anticoagulant therapy by taking the patient's condition into consideration.
- Always refer to the insert and/or Instructions for Use attached to the medical devices to be used together with this guide catheter.
- Do not use ethiodol and/or lipiodol agents, or any agent containing these drugs. (These agents may cause the erosion of the connector resin and result in damage.)

- Do not use this guide catheter for a patient with blood coagulation disorder or severe change in coagulation ability due to any cause. Doing so may worsen the symptoms.
- Do not use this guide catheter for a patient who cannot adapt to surgical operations.
- Do not use this guide catheter for a patient who has had any severe and distinct reaction against the agents or the like necessary for the intended procedure. Doing so may cause side effects such as allergic symptoms or shock disease.
- Do not use this guide catheter for a patient who is or could be pregnant. (The fetus may be affected by X-rays under fluoroscopy.)
- Do not use this guide catheter for a patient who has had undue peripheral vascular disease which inhibits the insertion of a sheath introducer in the appropriate size.
- Do not use this guide catheter for a patient who is judged unsuitable for the procedure by the physician.

Precautions

- Prior to use, confirm that all devices including this guide catheter function normally.
- Check that the specifications of this guide catheter are suitable for the purpose of use and procedures.
- In case this guide catheter is inserted over an angiographic guide wire with hydrophilic polymer coating, advance this guide catheter gently, paying careful attention not to abrade the hydrophilic coating of the guide wire by the tip of the supplied dilator.
- Before use, make sure that this guide catheter is appropriate in size and shape, and compatible with other medical devices to be used in combination. Because the supplied dilator is dedicated to this guide catheter, do not use it with other Asahi Intecc guide catheters of different Fr sizes or other manufacturers' guide catheters.
- Since the maximum guide wire diameter is as follows, use this guide catheter in combination with a guide wire having a diameter at or below the maximum guide wire diameter.

Maximum guide wire diameter	
When the supplied dilator is used.	When the supplied dilator is not used.
0.89 mm (0.035")	0.97 mm (0.038")

- Since the tip of this guide catheter is not tapered, it may cause occlusion in small vessels. Use care to avoid total blockage of blood flow.
- If the supplied dilator is kept inserted into this guide catheter, it may adversely affect the shape-memory property of this guide catheter tip. Therefore, do not insert the supplied dilator into this guide catheter until immediately before inserting this guide catheter into the patient.
- In case the clearance between the device to be inserted in and the lumen of this guide catheter is small, fully open the hemostatic valve and insert the device slowly in order not to cause air embolism.
- When taking this guide catheter out of the package, keep this guide catheter attached to the protective thick paper and take out the protective paper together with this guide catheter so as to avoid damage to this guide catheter or tip of this guide catheter.
- When withdrawing this guide catheter from the blood vessel, it is recommended that the angiographic guide wire be inserted into the guide catheter first, and that the dilator be reinserted into this guide catheter completely and locked to the guide catheter prior to withdrawal. If this guide catheter is the type not supplied with a dilator, be sure to insert the angiographic guide wire into this guide catheter's lumen before withdrawing this guide catheter. Manipulating this guide catheter alone may cause damage of the blood vessel or damage to the guide catheter.
- Before injection of contrast media, check the guide catheter for twist, crush or kink, and reconfirm this guide catheter is not occluded. Injection of contrast media into an occluded guide catheter may further damage this guide catheter.
- This guide catheter is coated with a hydrophilic coating and has a high level of lubricity. Consequently, always confirm the location of the tip under this guide catheter X-ray fluoroscopy and manipulate this guide catheter carefully.
- Operate slowly when inserting the supplied dilator into this guide catheter. (Advancing it rapidly may apply excessive load onto the tip of this guide catheter and result in damage.)
- Do not wipe the surface of this guide catheter with gauze, absorbent cotton, etc. soaked with alcohol, chlorhexidine gluconate solution, etc. (It may damage this guide catheter, or reduce the lubricity of the surface.)
- The guide catheter has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of this guide catheter in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.
- Since this guide catheter is not designed for use in combination with a power injector, do not use it in combination with a power injector.
- Take preventive measures against infection after use. Discard this guide catheter as medical waste.

Malfunctions and Adverse Effects

When using this guide catheter, the following malfunctions and adverse effects may occur. A serious adverse effect may lead to severe complication or death. However, such malfunctions and adverse effects are not limited to these.

- Malfunctions**
 - Damage
 - Kink
 - Separation
 - Withdrawal difficulty
 - Insertion difficulty
 - Obstruction
 - Bend
 - Damage to the hydrophilic coating
- Adverse Effects**
 - Death
 - Infection
 - Vessel dissection
 - Damage to a vessel, including possible vessel perforation
 - Bleeding complications
 - Distal embolism
 - Thrombus
 - Infarction
 - Residue
 - Ischemia
 - Arrhythmia
 - Vasospasm
 - Blood pressure reduction
 - Allergic reaction

How to Use

Insert this guide catheter into vascular system percutaneously.

○ Insertion via a sheath

■ Preparation of guide catheter

- Take out this guide catheter with the protective thick paper from the package.
- Remove this guide catheter from the protective thick paper.
- After taking it out, check that this guide catheter is not damaged.
- Flush this guide catheter with heparinized and sterilized saline using a syringe.

■ How to insert

- Prepare the patient for neuroradiology with the usual process.
- Insert a puncture needle with a plastic outer tube into the intended puncture site, and withdraw the puncture needle while leaving the plastic tube remaining in the vessel.

- Insert the guide wire for sheath into the plastic tube, and withdraw the plastic tube.
- Insert the dilator of the sheath into the sheath completely and lock the connectors, then insert the sheath and the dilator into the vessel along the guide wire for sheath. The combinations of this guide catheter and sheath are as follows:
 - 6Fr guide catheter ●●● 6Fr sheath
 - 7Fr guide catheter ●●● 7Fr sheath
 - 8Fr guide catheter ●●● 8Fr sheath
- Withdraw the dilator of the sheath and the guide wire for sheath, leaving the sheath inserted in the vessel.
- Insert an angiographic guide wire into the sheath.
- Insert this guide catheter gently and carefully along the angiographic guide wire into the sheath. For smooth insertion of this guide catheter, add drops of heparinized and sterilized saline with a syringe to the surface of this guide catheter to enhance the lubricity of the hydrophilic coating.
- Pull out the angiographic guide wire, and flush the inside of this guide catheter by aspirating the air inside the guide catheter with a syringe.
- Advance this guide catheter to the target site and secure it.

■ Withdrawal/replacement of the guide catheter

- Insert the angiographic guide wire until it protrudes by approximately 5 cm from the tip of this guide catheter.
- Pull out this guide catheter gently and carefully along the angiographic guide wire.
- To replace this guide catheter, insert a new guide catheter by following step 7) in the "How to insert" section above.

○ Insertion not via a sheath

■ Preparation of guide catheter

- Take out this guide catheter with the protective thick paper from the package.
- Remove this guide catheter and the supplied dilator from the protective thick paper.
- After taking them out, check that the guide catheter and the supplied dilator are not damaged.
- Flush this guide catheter and the supplied dilator with heparinized and sterilized saline using a syringe.

■ How to insert

- Prepare the patient for neuroradiology with the usual process.
- Insert a puncture needle with plastic outer tube into the intended puncture site, and withdraw the puncture needle while leaving the plastic tube remaining in the vessel.
- Insert an angiographic guide wire into the outer plastic tube and advance the guide wire to near the intended site, and then withdraw the plastic tube.
- Insert the supplied dilator into this guide catheter completely, and lock this guide catheter with the connector of the supplied dilator.
- Insert the guide catheter combined with the supplied dilator gently and carefully along the angiographic guide wire. For smooth insertion of this guide catheter, add drops of heparinized and sterilized saline with a syringe to the surface of this guide catheter to enhance the lubricity of the hydrophilic coating. If any resistance is felt, make a small incision in the skin using a surgical knife.
- When this guide catheter has almost reached the intended site, separate the supplied dilator from the guide catheter, and then pull out the dilator.
- Pull out the angiographic guide wire, and flush the inside of this guide catheter by aspirating the air inside the guide catheter with a syringe.
- Advance this guide catheter to the target site and secure it.

■ Withdrawal/replacement of the guide catheter

- Insert the angiographic guide wire until it protrudes by approximately 5 cm from the tip of this guide catheter.
- When the tip of this guide catheter is withdrawn back along the angiographic guide wire proximal to the targeted site, insert the supplied dilator completely along the angiographic guide wire, and lock the connector of the supplied dilator with this guide catheter.
- While keeping this guide catheter connected with the supplied dilator, pull out this guide catheter along the angiographic guide wire.
- To replace this guide catheter, insert a new guide catheter by following step 4) in the "How to insert" section above.

Storage method

Do not keep the product in a bent or heavily-loaded condition. This product must be kept out of water. Store in a cool, dark and dry place.

Expiration date

The expiration date is indicated on the label of the product package.

Contents

One set/bag

Liability Disclaimer

By no means shall "ASAHI INTECC CO., LTD. and its affiliated companies" (hereinafter referred to as the "Company") be liable for accidents, personal injuries, adverse effects due to any improper use of the product(s) or any other use inconsistent with these instructions. In no event shall the Company be liable for any damages either (i) arising out of storage of the product(s) after the shipment from the Company or (ii) due to selection of patients, surgery techniques, or any other medical activities by the medical institution that uses the product(s).

"ASAHI" and "FUBUKI" are trademarks or registered trademarks of ASAHI INTECC CO., LTD. in Japan and other countries.

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© 2018 ASAHI INTECC CO., LTD.
Micro Therapeutics, Inc.