

VENOUS WALLSTENT™ SELF-EXPANDING STENT

Clinically Proven Case After Case

Backed by a decade of clinical data, Venous Wallstent Self-Expanding Stent delivers proven technology for the most challenging venous applications.

Unmatched Clinical Experience

The **safety and durability** of the Wallstent Self-Expanding Stent have been repeatedly validated.

This proven stent has demonstrated consistently low complication rates across all indicated uses.

It is indicated for improving luminal diameter in the iliofemoral veins for the treatment of symptomatic venous outflow obstruction.



37

Publications reporting clinical results including Wallstent*



650K

Total stents sold



4,500

Patients with reported outcomes


*From 2000-2018

Dynamic, Innovative Design



HIGH DURABILITY

Venous Wallstent is specifically engineered to provide **fracture resistance**. It is designed to withstand the forces common in venous anatomy.



BRAIDED ARCHITECTURE

Braided design with Elgiloy® material provides **compression resistance** in venous anatomy. In addition, its reconstrainability eases placement.



UNIQUE CONFORMABILITY


The stent is designed for optimal **conformability** so it can be deployed effectively in curved and tapered vessels.

Comprehensive Range of Sizes


Diameters from 10 mm to 20 mm*

VENOUS WALLSTENT SELF-EXPANDING STENT | SKUs

		SIZE / UPN	
Stent Diameter (cm)	10	10X20X75CM	H74912044102070
		10X42X75CM	H74912044104270
		10X68X75CM	H74912044106870
		10X94X75CM	H74912044109470
	12	12X20X75CM	H74912044122070
		12X40X135CM	H74912044124010
		12X40X75CM	H74912044124070
		12X60X135CM	H74912044126010
		12X60X75CM	H74912044126070
		12X90X135CM	H74912044129010
	14	12X90X75CM	H74912044129070
		14X20X75CM	H74912044142070
		14X40X75CM	H74912044144070
		14X60X75CM	H74912044146070
	16	14X90X75CM	H74912044149070
		16X20X75CM	H74912044162070
		16X40X75CM	H74912044164070
		16X60X75CM	H74912044166070
	18	16X90X75CM	H74912044169070
		18X40X75CM	H74912044184070
20	18X90X75CM	H74912044189070	
	20X40X75CM	H74912044204070	
	20X55X75CM	H74912044205570	
		20X80X75CM	H74912044208070



VENOUS WALLSTENT SELF-EXPANDING STENT



VENOUS WALLSTENT™

CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician. Rx only. Prior to use, please see the complete "Directions for Use" for more information on Indications, Contraindications, Warnings, Precautions, Adverse Events, and Operator's Instructions. **INDICATIONS FOR USE/INTENDED USE:** The VENOUS WALLSTENT is indicated for improving central venous luminal diameter following unsuccessful angioplasty in patients on chronic hemodialysis with stenosis of the venous outflow tract. Unsuccessful angioplasty is defined as residual stenosis $\geq 30\%$ for a vein ≤ 10 mm in diameter or $\geq 50\%$ for a vein > 10 mm in diameter, a tear which interrupts the integrity of the intima or lumen, abrupt lesion site occlusion, or refractory spasm. The vessels that can be treated with the VENOUS WALLSTENT are the innominate and subclavian veins, ranging from 8 mm to 15 mm in diameter. The VENOUS WALLSTENT is also indicated for improving luminal diameter in the iliofemoral veins for the treatment of symptomatic venous outflow obstruction. **CONTRAINDICATIONS:** Patients with uncorrected bleeding disorders. **WARNINGS:** Subsequent restenosis may require repeat dilation of the vessel segment containing the stent. A stent cannot be repositioned or removed after the deployment threshold has been exceeded. **PRECAUTIONS:** Do not advance a partially ($\leq 50\%$) deployed stent. **MAGNETIC RESONANCE IMAGING (MRI) SAFETY INFORMATION:** Non-clinical testing has demonstrated that the VENOUS WALLSTENT™ system is MR Conditional for single and overlapping lengths up to 120 mm. A patient with this stent can be scanned safely, immediately after placement, under the following conditions: Static magnetic field of 1.5 Tesla or 3 Tesla • Highest spatial gradient magnetic field of 19 Tesla/m (1900 Gauss/cm) or less • Maximum MR system reported, whole body averaged specific absorption rate (SAR) of ≤ 1 W/kg for patient landmarks above the umbilicus (patient navel) and ≤ 2 W/kg (Normal Operating Mode) for patient landmarks below the umbilicus. **RF Heating:** Under the scan conditions defined above, VENOUS WALLSTENT is expected to produce a maximum in-vivo temperature rise of 3.1°C after 15 minutes of continuous scanning. **Image Artifact:** In non-clinical testing, the image artifact caused by the device extends approximately 13 mm from the stent when imaged with a gradient echo pulse sequence and a 3 Tesla MRI system. The artifact does obscure the lumen. **Recommendations:** It is recommended that patients register the conditions under which the implant can safely be scanned with the MedicalAlert Foundation (www.medicalert.org) or an equivalent organization. **ADVERSE EVENTS:** Adverse Events associated with use of stents for Vascular applications may include, but are not limited to, the following: Allergic reactions (drug, contrast, device, or other) • Angina • Arteriovenous fistula • Cerebrovascular accident/stroke/Transient Ischemic Attack • Death • Embolism (air, plaque, thrombus, device, or other) • Fever • Hematoma • Hemorrhage/Bleeding • Ischemia • Hypotension/hypertension • Myocardial infarction/Ischemia • Need for urgent intervention or surgery • Pain • Pulmonary embolism • Renal insufficiency or failure • Restenosis of stented vessel • Sepsis/infection • Stent fracture • Stent migration • Stent/Vessel occlusion • Thrombus/thrombosis • Vasospasm • Venous congestion • Vessel injury (perforation, trauma, rupture, dissection, pseudoaneurysm or other) 90960659 AC

WALLSTENT is a registered or unregistered trademark of Boston Scientific Corporation or its affiliates. All other trademarks are property of their respective owners.

Boston Scientific
Advancing science for life™

Peripheral Interventions

300 Boston Scientific Way
Marlborough, MA 01752-1234
www.bostonscientific.com

To order product or for more information
contact customer service at 1.888.272.1001.

© 2020 Boston Scientific Corporation
or its affiliates. All rights reserved.

PI-781701-AA