ARGON MEDICAL DEVICES, INC.
1445 Flat Creek Road Athens, Texas 75751 USA
Tel: 800-927-4669;
Tel: +1 903-675-9321
www.argonmedical.com

Argon Medical Devices UK Ltd
Eastgate Business Centre
Eastern Avenue
Burton-on-Trent
Staffordshire
DE13 0AT

X9585730/0321A

EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands

Symbols • Symboles • Símbolos • Simboli • Symbole • Símbolos • Símbolos • Symboler • Symbolit • Symboler • Symboler • Symboler • Símbolen • Semboller • Σύμβολα • Symbole • Symboly • Symboly • Sýmboľid • Szimbólumok • Simboliai • Simboluri • Simboli الرموز • Білу тэраг •

∅	Snare Diameter • Диаметър на примката • Diámetro do laço • Průměr smyčky • Snarediameter • Snare-Durchmesser • Διάμετρος βρόχου • Diámetro del lazo • Silmuse läbimõõt • Silmukan halkaisija • Diamètre du collet • Snare átmérő • Diámetro del filtro • Kilpų skersmuo • Diameter snaar • Snarediameter • Średnica chwytaka • Diámetro da alça • Diametruł anseı • Priemer slučky • Premer zanke • Snarans diameter • Kement Çapı • Đường Kinh Thòng Lọng • "Snare" قطر المنطار
⊙	Catheter Diameter • Диаметър на катетъра • Diámetro do cateter • Průměr katétru • Kateterdiameter • Katheterdurchmesser • Διάμετρος καθετήρα • Diámetro del catéter • Kateetri läbimõõt • Katetrin halkaisija • Diamètre du cathéter • Katéterátmérő • Diámetro del catetere • Kateterio skersmuo • Diameter catheter • Kateterdiameter • Średnica cewnika • Diámetro do cateter • Diametruł cateterului • Priemer katétra • Premer katetra • Kateterns diameter • Kateter Çapı • Đường Kinh Ống Thông • قطر الصغرة

<https://www.argonmedical.com/resources/product-information>

English	The symbols glossary is located electronically at www.argonmedical.com/symbols
Български	Речникът на символите е наличен в електронен вид на адрес www.argonmedical.com/symbols
Potuguês do Brasil	O glossário dos símbolos está localizado eletronicamente em www.argonmedical.com/symbols
Česky	Slovník symbolů v elektronické formě je umístěn na www.argonmedical.com/symbols
Dansk	Symbolforklaringen findes i elektronisk form på www.argonmedical.com/symbols
Deutsch	Das Glossar der Symbole kann unter www.argonmedical.com/symbols abgerufen werden
Ελληνικά	Το γλωσσάρι των συμβόλων παρέχεται ηλεκτρονικά στη διεύθυνση www.argonmedical.com/symbols
Español	El glosario de símbolos puede consultarse en línea en la dirección www.argonmedical.com/symbols
Eesti	Sümbolite seletuste kogu on elektrooniliselt saadaval aadressil www.argonmedical.com/symbols
Suomi	Elektroninen symbolisanasto on osoitteessa www.argonmedical.com/symbols .
Français	Le glossaire des symboles est disponible au format électronique à l'adresse www.argonmedical.com/symbols
Magyar	A szimbólumok összefoglalása elektronikusan itt található: www.argonmedical.com/symbols
Italiano	Il glossario dei simboli è disponibile in formato elettronico all'indirizzo www.argonmedical.com/symbols
Lietuvių k.	Simbolių žodyno elektroninė forma pateikta www.argonmedical.com/symbols
Nederlands	Uitleg over de symbolen is elektronisch beschikbaar op www.argonmedical.com/symbols
Norsk	Symbolforklaringen befinner seg elektronisk på www.argonmedical.com/symbols
Polski	Glosariusz symboli w postaci elektronicznej dostępny jest pod adresem www.argonmedical.com/symbols
Português	O glossário de símbolos está localizado eletronicamente no endereço www.argonmedical.com/symbols
Română	Glosarul simbolurilor este stocat electronic la www.argonmedical.com/symbols
Slovenčina	Glosár symbolov je v elektronickej podobe k dispozícii na stránkach www.argonmedical.com/symbols
Slovensko	Slovar simbolov se nahaja v elektronski obliki na www.argonmedical.com/symbols
Svenska	Symbolforklaringen hittas elektroniskt på www.argonmedical.com/symbols
Türkçe	Semboller sözlüğü elektronik ortamda www.argonmedical.com/symbols adresinde yer almaktadır
Tiếng Việt	Bảng chú giải biểu tượng được hiển thị dưới dạng điện tử tại www.argonmedical.com/symbols
Arabic	يوفر مسرد الرموز إلكترونياً على www.argonmedical.com/symbols

Atrieve™ Vascular Snare Kit

Intended Use/Purpose:

The Atrieve™ Vascular Snare Kit is intended for use in the cardiovascular system or hollow viscous to retrieve and manipulate foreign objects.

Indication for Use:

The Atrieve™ Vascular Snare Kit is intended for use in the cardiovascular system or hollow viscous to retrieve and manipulate foreign objects. Manipulation procedures include indwelling venous catheter repositioning, indwelling venous catheter fibrin sheath stripping, and central venous access veni-puncture procedure assistance.

Description:

The Atrieve™ Vascular Snare Kit comprises a Snare device with three preformed loops constructed of super-elastic nitinol and platinum. The loops are non-interwoven and can slide relative to one another. The loops of the Snare can be introduced into a delivery catheter without the risk of permanent device deformation. The Atrieve™ Vascular Snare Kit also includes a Delivery Catheter, a torque handle, and an Introducer. The Delivery Catheter has a radiopaque marker band at the distal tip. The product comes packaged in a protective hoop holder.

The Atrieve™ Vascular Snare Kit is available in the following sizes:

Snare Diameter Range	Snare Length	Delivery Catheter Size	Delivery Catheter Length
2-4 mm	175 cm	3.2 F	150 cm
4-8 mm	175 cm	3.2 F	150 cm
6-10 mm	120 cm	6 F	100 cm
9-15 mm	120 cm	6 F	100 cm
12-20 mm	120 cm	6 F	100 cm
18-30 mm	120 cm	7 F	100 cm
27-45 mm	120 cm	7 F	100 cm

Warnings:

- This device is not intended for the removal of foreign objects entrapped by tissue growth.
- This device should not be used for fibrin sheath stripping in the presence of septal defects of Persistent Foramen Ovale.
- This device is not intended for removal of implanted pacing leads.
- Pull forces applied to catheters during fibrin sheath stripping may damage, stretch, or break indwelling catheters 6 French or smaller in diameter. Do not use excessive pull force when attempting fibrin sheath stripping of catheters 6 French or smaller in diameter.
- Do not use excessive force when manipulating the catheter through an introducer. Excessive force may damage the catheter.
- This device was designed, tested and manufactured for single patient use only. Reuse or reprocessing of this device may lead to its failure and subsequent patient injury. Reprocessing and/or resterilization of this device may create the risk of contamination and patient infection. Do not reuse, reprocess or sterilize this device.
- Inspect the package integrity before use.
- Do not use if package is open or damaged or if the expiry date has been exceeded.
- Do not continue to use if any of the component are damaged during the procedure.
- Nitinol is a nickel titanium alloy. Possible reaction may occur for those patients who exhibit sensitivity to nickel.
- Possible reaction may occur for those patients who exhibit sensitivity to platinum.

Precautions:

- This device is intended for use by a qualified physician.
- Care should be observed when using this device for removal of a large fibrin sheath in order to minimize risk of pulmonary embolism.
- Bleeding may occur through the delivery catheter hub during some procedures. A hemostasis valve (not included) can be attached to the delivery catheter to reduce the bleeding as needed or per physician's judgment.

Potential Complications:

- Potential complications associated with foreign body retrieval devices in arterial vasculature include, but are not limited to:
 - Embolization
 - Stroke
 - Myocardial infarction (depending upon placement)
- Potential complications associated with foreign body retrieval devices in venous vasculature include, but are not limited to:
 - Pulmonary embolism
- Other potential complications associated with foreign body retrieval devices include, but are not limited to:
 - Vessel perforation
 - Device entrapment

Catheter damage can occur when attempting fibrin sheath stripping on small French size diameter catheters (See WARNINGS). Incidence of pulmonary embolism after fibrin sheath stripping may occur (See PRECAUTION).

Prepare the Snare:

Select the appropriate Snare diameter range for the site in which the foreign body is located. The Snare diameter range should approximate the size of the vessel in which it will be used.

- Remove the Snare and Delivery Catheter components of the kit from their hoop holders and inspect for any damage.

- Remove the Introducer and Torque Handle device from the proximal end of the Snare shaft.
- Load the Snare into the Delivery Catheter by inserting the proximal end of the Snare into the distal (non-hubbed) end of the Delivery Catheter, until the proximal end of the Snare shaft exits the hub and the loops can be retracted into the distal end of the Delivery Catheter.
- Test and inspect the Snare device by extending and retracting the loops through the distal end of the Catheter 2-3 times, while carefully examining the Delivery Catheter, radiopaque band and the device for any damage or defects.
- When appropriate, the Snare and Catheter can be advanced to the desired site as a single unit assembled as described above.

Alternative Preparation of the Snare

If the Delivery Catheter is already positioned within the vasculature, the provided Introducer (located on the proximal end of the Snare and just distal to the Torque Handle) may be used to position the Snare in the indwelling Delivery Catheter.

- Remove the Snare from the hoop holder and inspect for any damage.
- Move the provided Introducer (located on the proximal end of the Snare, just distal to the Torque Handle) distally until the loops of the Snare are enclosed within the tubing portion of the Introducer.
- Insert the distal end of the Introducer into the hub of the indwelling Delivery Catheter until resistance is felt. This will indicate the tip of the Introducer is properly aligned with the inner lumen.
- Hold the Introducer as straight as possible, grasp the shaft of the Snare just proximal to the hub of the Introducer and advance the Snare until it is secure within the lumen of the Delivery Catheter. The Introducer can be removed by first removing the Torque Handle and pulling the Introducer off the proximal end of the Snare's shaft. After removing the Introducer, re-apply the Torque Handle onto the shaft.

Snare Assisted Retrieval and Manipulation Suggestions:

- If present, remove the indwelling delivery catheter.
- If a guidewire is in a patient at the location of a foreign body, advance the Delivery Catheter over the guidewire to the desired location. Then remove the guidewire and advance the Snare through the Delivery Catheter. Alternatively, cinch one loop of the Snare over the proximal end of the guidewire and advance the Snare and Catheter assembly into a guide catheter or introducer sheath until the distal end of the Delivery Catheter is positioned proximal to the foreign body.
- If a guidewire is not present, pull the Snare into the distal end of the Delivery Catheter and advance through a guide catheter or introducer sheath (not included) until it is positioned proximal to the foreign body. Alternatively, collapse the Snare loops by pulling the device into the distal end of the Introducer. Place the tapered end of the Introducer into the proximal (hub) end of the Delivery Catheter, guide catheter or sheath and advance the Snare forward maintaining constant contact between the Introducer and Delivery Catheter hub.
NOTE: When attempting to utilize guide catheters or sheaths not specifically manufactured for use with the Atrieve™ Vascular Snare, it is important to test product compatibility prior to use.
- Gently push the Snare shaft forward to completely open the loops outside of the Delivery Catheter. The loops are then slowly advanced forward, and may be rotated if desired, around the proximal end of the foreign body. Alternatively, the Snare may be advanced beyond the target location and the loops brought back around the distal end of the foreign body.
- By advancing the Delivery Catheter, the loops of the device are closed to capture the foreign body. (Note that attempting to close the loops by pulling the Snare into the Delivery Catheter will move the loops from their position around the foreign body.)
- To manipulate a foreign body, maintain tension on the Delivery Catheter to retain the hold on the foreign body, and move the Snare and Delivery Catheter together to manipulate a foreign body to the desired position.
- To retrieve a foreign body, maintain tension on the Delivery Catheter and move the Snare and Delivery Catheter assembly together proximally to, or into a guide catheter or sheath. The foreign body is then withdrawn through or together with the guiding catheter or introducer sheath. Withdrawal of large foreign bodies may require the insertion of larger sheaths, guiding catheters, or a cut-down at the peripheral site.

Snare Assisted Removal of Fibrin Sheaths from Indwelling Catheters:

- Using standard technique, prepare a femoral vein approach, advance the selected Snare to the inferior vena cava or right atrium.
- Advance a 0.035" guidewire through the end port (distal or venous port if more than one lumen) of the indwelling catheter and into the inferior vena cava or right atrium.
- Position one of the Snare loops around the guidewire.
- Advance the Snare over the distal end of the indwelling catheter to a position proximal to the fibrin sheath.
- Close the Snare around the indwelling catheter and continue applying light traction while gently pulling the Snare down toward the distal end of the indwelling catheter over the end ports.
- Repeat steps 4 & 5 until the indwelling catheter is free of fibrin sheath.

Snare Assisted Venous Canalization:

- Introduce the Snare at a patent venous access site and position in the vasculature at the desired site.
- Open the Snare loops to provide a target to guide an entry needle into the desired venous access site.
- Introduce a guidewire through the needle and through the Snare loops. Remove the needle.
- Close the Snare over the guidewire by advancing the Delivery Catheter.
- Pull the guidewire into the desired location.

How Supplied

The Atrieve™ Vascular Snare Kit is sterilized by ethylene oxide gas. Content is sterile and non-pyrogenic in unopened and undamaged package. Do not use device if package has been damaged or has been opened.

Storage

Store in controlled room temperature.

Disposal

After use the Atrieve™ Vascular Snare Kit may be potential biohazard. Handle and dispose of in accordance with accepted medical practice and all applicable laws and regulations.

NOTE: In the event a serious incident related to this device occurs, the event should be reported to Argon Medical at quality.regulatory@argonmedical.com as well as to the competent health authority where the user/patient resides.

Atrieve™ Vascular Snare is a trademark of Argon Medical Devices, Inc.
©2013 Argon Medical Devices, Inc.. All rights reserved.

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©2013 Argon Medical Devices, Inc. All rights reserved.