HALO ONE™

Thin-Walled Guiding Sheath

ENGLISH

INSTRUCTIONS FOR USE

Caution: Federal (U.S.A.) Law restricts this device to sale by or on the order of a physician.

DEVICE DESCRIPTION

The HALO ONE™ Thin-Walled Guiding Sheath is designed to perform as both a guiding sheath and an introducer sheath.

The HALO ONE™ Thin-Walled Guiding Sheath consists of a thin-walled (up to 1F reduction in outer diameter compared to standard sheaths of equivalent French size) sheath made from braided single-lumen tubing, fitted with a female luer hub at the proximal end and a formed atraumatic distal tip (A) (Figure 1). The thin-wall design reduces the thickness of the sheath wall to help facilitate intravascular access from access sites including but not limited to radial, femoral, popliteal, tibial, and pedal. (Figure 2).

A detachable hemostasis valve (C), employing a crosscut silicone membrane and incorporating a side arm terminating in a 3-way stopcock (D), is connected to the sheath luer hub (E). The sheath is supplied with a compatible vessel dilator (B) that snaps securely into the hemostasis valve hub (C). The sheath has a strain relief feature located at the luer hub and a radiopaque platinum-iridium marker (F) located close to the distal tip (A). Refer to product label for product specifications.

Sheaths are provided both with and without hydrophilic coating as labeled. Hydrophilic coating where labeled is applied over the distal portion of the sheath to provide a lubricious surface to ease insertion.


Figure 1


Figure 2: Thin-wall design (A) of 5F HALO ONE™ sheath in relation to a standard 4F (B) and 5F sheath (C)

INDICATION FOR USE

The HALO ONE™ Thin-Walled Guiding Sheath is indicated for use in peripheral arterial and venous procedures requiring percutaneous introduction of intravascular devices. The HALO ONE™ Thin-Walled Guiding Sheath is NOT indicated for use in the neurovasculature or the coronary vasculature.

CONTRAINDICATIONS

There are no known contraindications for the HALO ONE™ Thin-Walled Guiding Sheath.

WARNINGS

1. Contents supplied **STERILE** using ethylene oxide (EtO). **Non-pyrogenic. Do not re-use, reprocess or re-sterilize.** This device is intended for single use only.
2. **Do not resterilize.** After resterilization, the sterility of the product is not guaranteed because of an indeterminable degree of potential pyrogenic or microbial contamination which may lead to infectious complications. Cleaning, reprocessing

and/or resterilization of the present medical device increases the probability that the device will malfunction due to potential adverse effects on components that are influenced by thermal and/or mechanical changes.

3. This device has been designed for single use only. Reusing this medical device bears the risk of cross-patient contamination as medical devices – particularly those with long and small lumina, joints, and/or crevices between components – are difficult or impossible to clean once body fluids or tissues with potential pyrogenic or microbial contamination have had contact with the medical device for an indeterminable period of time. The residue of biological material can promote the contamination of the device with pyrogens or microorganisms which may lead to infectious complications.
4. Visually inspect the packaging to verify that the sterile barrier is intact. Do not use if the sterile barrier is opened or damaged.
5. Use the sheath prior to the “Use By” date specified on the package.
6. Do not advance the guidewire, sheath/dilator, procedural device, or any component if resistance is met, without first determining the cause and taking remedial action.
7. Do not use a power injector through the sideport or the three-way stopcock
8. The HALO ONE™ Thin-Walled Guiding Sheath has not been evaluated for use in the neurovasculature or the coronary vasculature.
9. After use, this product may be a potential biohazard. Handle and dispose of in accordance with acceptable medical practices and applicable local, state and federal laws and regulations
10. Only advance or retract the sheath with the dilator inserted and only advance or retract the sheath and dilator while placed over a properly sized guidewire.
11. Failure to deactivate the procedural device prior to removal through the sheath may cause damage to the sheath and may result in patient injury.

PRECAUTIONS

1. The HALO ONE™ Thin-Walled Guiding Sheath shall only be used by trained physicians. Access procedures should be conducted under fluoroscopic guidance with appropriate x-ray equipment and / or ultrasound.
2. Prior to beginning radial artery access, an assessment such as the Allen or Barbeau test should be performed to assess the presence/adequacy of dual arterial circulation to the hand. Radial artery access is not recommended for patients with abnormal Allen or Barbeau test results or radial pulse, or insufficient dual arterial supply.
3. Prior to beginning pedal access, physicians should assess the vascular anatomy to assure there is adequate antegrade flow at the level of the ankle.
4. The minimum acceptable sheath French size is printed on the package label. Do not attempt to pass devices through a smaller size sheath introducer than indicated on the device label.
5. The pouch should be inspected prior to opening to ensure the sterile barrier is not compromised. The device should be carefully removed and placed in the sterile field. The entire procedure from skin puncture or incision to sheath withdrawal must be carried out aseptically.
6. Carefully inspect the sheath prior to use to verify that the sheath has not been damaged during shipment and that its size, shape and condition are suitable for the procedure for which it is to be used. Do not use if product damage is evident.
7. Careful attention must be paid to the maintenance of tight valve connections for duration of procedure to avoid blood leakage or the introduction of air into the system. Take remedial action if any excessive blood leakage is observed.
8. Insert dilator into the center of the sheath valve. Forced insertion of the dilator which misses the center of the valve may cause damage and result in blood leakage.
9. Advance or withdraw the sheath slowly. If resistance is met do not advance or withdraw until the cause of resistance is determined.
10. When inserting, manipulating or withdrawing a device through the introducer always maintain the introducer position.
11. Remove the dilator from the sheath slowly to avoid incomplete closing of the valve resulting in blood leakage.
12. When using procedural devices close to the tip of the sheath care must be taken to ensure the active mechanism portion of the procedural device (e.g., balloon, stent zone, material removal section of atherectomy device) is not within the tip of the sheath. The radiopaque marker is located within 5 mm of the end of the tip but does not mark the true distal tip of the sheath.
13. Before removing or inserting the interventional/diagnostic device through the sheath, aspirate blood from the 3-way stopcock to remove any fibrin deposition which may have accumulated in or on the tip of the sheath.
14. Ensure to deactivate the procedural device prior to removal through the sheath.
15. Take care when back loading the tip of the dilator over the guidewire to avoid damage to the dilator.

16. Ensure the dilator is securely connected with the sheath prior to advancing otherwise only the sheath may advance into the vessel and the sheath tip may cause damage to the vessel.
17. Ensure the dilator is in place within the sheath before advancing the sheath as otherwise damage may be caused to the vessel.
18. Do not place sutures on the sheath tubing since this may restrict access/flow through the sheath. When puncturing, suturing or incising near the sheath be careful not to damage the sheath. Proper functioning of the sheath depends on its integrity. Care should be used when handling the sheath.
19. If using fluid injection through the 3 way stopcock, ensure the dilator or procedural device is not in place at the same time.
20. If resistance is felt during post-procedure withdrawal of the procedural device, it is recommended to remove the procedural device, guidewire, and sheath as a single unit.
21. In order to activate the hydrophilic coating, it is recommended to wet the HALO ONE™ Thin-Walled Guiding Sheath with heparinized saline solution immediately prior to its insertion in the body. Failure to activate the coating may lead to sub-optimal trackability of the sheath. To maintain lubricity this surface must be kept completely wet.
22. Proper functioning of the HALO ONE™ Thin-Walled Sheath depends on its integrity. Care should be used when handling the sheath. Damage may result from kinking, stretching, or forceful wiping of the HALO ONE™ Thin-Walled Guiding Sheath. Do not continue to use the sheath if the shaft has been bent or kinked.

POTENTIAL ADVERSE EFFECTS

Potential adverse effects that may result from a percutaneous vascular procedure (directly or indirectly associated with the device) may include, but are not limited to:

- Air embolism
- Aneurysm or pseudoaneurysm
- Arteriovenous fistula
- Compartment Syndrome
- Death
- Embolism
- Endocarditis
- Hematoma
- Hemorrhage, including bleeding at the puncture site
- Intimal tear
- Radial artery occlusion/spasm
- Sepsis/infection/inflammation
- Tissue necrosis
- Thrombus formation
- Vessel spasm, perforation or dissection
- Potential systemic indirect/inherent adverse effects related to general endovascular procedures may include, but are not limited to:
 - Arrhythmias
 - Drug reactions, allergic reaction to contrast media
 - Hypotension/hypertension
 - Pain and tenderness

STORAGE

Store in a cool, dry, dark place. Rotate inventory so that the catheter and other dated products are used prior to the "Use By" date. Do not use if packaging is damaged or opened.

DIRECTIONS FOR USE

Equipment Required

- Sterile saline solution (heparinized)
- Sterile saline solution (heparinized) / Contrast medium (2:1 ratio recommended)
- Luer lock syringe/inflation device with manometer (10 ml or larger)
- 0.035" (0.89 mm) guidewire or 0.018" (0.46 mm) guidewire as labeled per device

HALO ONE™ Thin-Walled Guiding Sheath Preparation

1. Verify the French size is suitable for the procedure and can accommodate the required procedural devices as labeled (Figure 3). Remove sheath and dilator from package.
2. Prior to use, the air in the sheath and dilator should be removed. To facilitate purging, the device is packaged with the dilator inside the sheath in a reverse direction allowing both to be flushed simultaneously. Select a syringe or inflation device with a 10 ml or larger capacity and fill approximately half of it with sterile saline solution. Prepare the lumen by connecting the syringe or inflation device containing the solution into the stopcock on the side-port of the

sheath and flush with the sterile heparinized saline solution (Figure 4). Close the stopcock to maintain the air tightness following flushing.

3. Prior to use the reverse loaded dilator must be removed from the distal end of the sheath. If not already completed at step 2, the air in the dilator lumen should be removed. To facilitate purging, select a syringe or inflation device with a 10 ml or larger capacity and fill approximately half of it with sterile saline solution. Prepare the lumen by connecting the syringe or inflation device containing the solution into the luer connector of the dilator hub and flushing with the sterile heparinized saline solution (Figure 7).
4. Insert the provided vessel dilator through the hemostatic valve and click the dilator hub into place in the valve housing (Figures 5 and 6).
5. In order to activate the hydrophilic coating, where labeled, it is recommended to wet the HALO ONE™ Thin-Walled Guiding Sheath with sterile saline solution immediately prior to its insertion in the body.

Use of the HALO ONE™ Thin-Walled Guiding Sheath

6. Identify the insertion site, including but not limited to radial, femoral, popliteal, tibial, or pedal access sites and prepare the site using proper aseptic technique and local anesthesia as required.
7. At the operator's discretion, make a small skin incision at the puncture site with a surgical scalpel. Recommended for scar tissue at the access site.
8. Backload the distal tip of the HALO ONE™ Thin-Walled Guiding Sheath dilator over the pre-positioned guidewire and advance the tip to the introduction site.
9. Advance the dilator and the sheath through the skin and into the vessel. Grasp the sheath and dilator assembly close to the skin while advancing to avoid buckling, employ a rotating motion while advancing as required. (Note: If using a hydrophilic guidewire, ensure that it is kept hydrated with sterile heparinized saline, at all times.)
10. Carefully advance the dilator and the sheath over the wire to the site required within vessel. The radiopaque marker identifies the sheath tip location under fluoroscopy.
11. Disconnect the dilator hub from the valve by bending to the side to unsnap from the valve cap (Figure 8). Remove the dilator slowly while holding the sheath in place and ensuring the guidewire remains in place (Figure 9) as required.
12. Load the procedural device over the pre-positioned guidewire.
13. Advance the procedural device carefully into the centre of the valve diaphragm and through the sheath to the treatment site (Figure 10) while maintaining the sheath position.
14. Position the procedural device relative to the lesion to be treated (Figure 11) ensuring the active mechanism portion of the procedural device is not within the sheath.
15. Following use ensure the procedural device is fully deactivated before carefully withdrawing the procedural device through the sheath while maintaining the sheath position.
16. After the procedure is completed insert the dilator over the wire into the sheath
17. Slowly withdraw the sheath and dilator as a unit and then remove the guide wire.
18. Apply hemostasis at the access site according to standard practice.
19. Dispose of the single-use device.

Warning: After use, this product may be a potential biohazard. Handle and dispose of in accordance with acceptable medical practices and applicable local, state and federal laws and regulations.

Limited Product Warranty

ClearStream Technologies Ltd ("ClearStream") warrants to the original purchaser that this product will be free from defects in material and workmanship for a period of one (1) year from the date of purchase. If this product proves to be so defective, purchaser may return same to ClearStream for repair, replacement, refund or credit at ClearStream's option. All returns must be authorized in advance in accordance with ClearStream's Returned Goods Policy found in its then current Price List. The liability of ClearStream under this limited product warranty does not extend to any abuse, misuse, improper storage, alteration, further manufacture, packaging or processing of this product or its repair by anyone other than an authorized ClearStream representative.

THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING, WITHOUT LIMITATION ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.

THE LIABILITY AND REMEDY STATED IN THIS LIMITED PRODUCT WARRANTY WILL BE THE SOLE LIABILITY OF ClearStream AND REMEDY AVAILABLE TO PURCHASER FOR THIS PRODUCT, WHETHER IN CONTRACT, TORT (INCLUDING NEGLIGENCE) OR OTHERWISE, AND ClearStream WILL NOT BE LIABLE TO PURCHASERS FOR ANY INDIRECT, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES ARISING OUT OF ITS HANDLING OR USE.

Some states/countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your state/country.

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Copyright © 2019 BD. All rights reserved.
Copyright © 2019 BD. Tous droits réservés.
Copyright © 2019 BD. Alle Rechte vorbehalten.
Copyright © 2019 BD. Tutti i diritti riservati.
Copyright © 2019 BD. Todos los derechos reservados.
Copyright © 2019 BD. Alle rechten voorbehouden.