ENGLISH

CROSSER™ Recanalization System CROSSER™ Generator - Manual of Operations

Instructions for Use

Warning! Carefully read all instructions prior to use. Observe all warnings and precautions noted throughout these instructions. Failure to do so may result in complications.

Caution! Federal (USA) law restricts this device to sale by or on the order of a physician

Description

The CROSSER™ Recanalization System consists of a re-usable electronic CROSSER™ Generator, Foot Switch, high frequency Transducer, the FLOWMATE® Injector (optional), and single-use CROSSER™ Catheter. A CROSSER™ Catheter is connected to the electronic CROSSER™ Generator through the Transducer. The Foot Switch is used to activate the system. The CROSSER™ Generator and Transducer convert AC power into high frequency mechanical vibrations which are propagated through a core wire to the metal tip of the CROSSER™ Catheter. The main body of the CROSSER™ Catheter is constructed from Pebax and a hydrophilic coating covers the distal portion of the CROSSER™ Catheter.

The CROSSER™ Generator is packaged with the Transducer and Foot Switch. A timer on the front of the CROSSER™ Generator displays the total activation time.

The FLOWMATE® Injector is a sterile saline injector for use as an accessory to the CROSSER™ Recanalization System. The FLOWMATE® Injector is sold separately by Bard.

Indications

The CROSSER™ Recanalization System is indicated to facilitate the intra-luminal placement of conventional guidewires beyond chronic total occlusions.

The CROSSER™ Generator is intended for use only with the CROSSER™ Catheter. Refer to the CROSSER™ Catheter Instructions for Use.


Figure 1 – The CROSSER™ Recanalization System

Figure 2- Front of CROSSER™ Generator

Figure 3- Rear of CROSSER™ Generator

Specifications

Model	GEN200
Weight (all components)	18 lbs (7 kg)
Dimensions (CROSSER™ Generator)	15 x 12 x 4 in. (38 x 30 x 10 cm)
Electrical	100-120/220-240 V ~ at 50/60 Hz 2/1 A
Power Cord	3 conductor, 6A, 250 V, 1.0 mm ² , HAR, VDE or OVE cord
Safety	Meets IEC60601-1, UL 60601-1, and CAN/CSA C22.2 No. 601.1
EMC	IEC60601-1-2: 2007
Ingress of Fluids	CROSSER™ Generator/Transducer – IPX0, ordinary, Foot Switch – IPX8 
Protection Against Electrical Shock	Type CF  Applied Part ; Degree: Class 1
Fuses	Two fuses rated at 2 X T2AL, 250V
Operating Frequency Range	20.5 KHz ± 500 Hz
Activation Time Limits	Single depress of foot switch 30 seconds maximum Max per CROSSER™ Catheter 5 minutes. Time accuracy is ± 1 second.
Transport and Storage Environment	Temperature -10 °C (14 °F) to +55 °C (131 °F) Relative Humidity 10% to 85% Non-condensing Atmospheric Pressure: 50-106 kPa
Operating Environment	Temperature +10 °C (50 °F) to +40 °C (104 °F) Relative Humidity 30% to 75% Non-condensing Atmospheric Pressure: 70-106 kPa
Mode of Operation	Continuous

This equipment has been tested and found to comply with the EMC limits for the Medical Device Directive 93/42/EEC (IEC 60601-1-2: 2007). The limits are designed to provide reasonable protection against harmful interference in a typical medical installation. The equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with these instructions, may cause harmful interference to other devices in the vicinity. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference with other devices, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving device
- Increase the separation between the equipment
- Connect the equipment into an outlet on a circuit different from that to which the other device(s) is connected
- Consult the manufacturer or field service technician for help

Electromagnetic Environmental Guidance

Electromagnetic Emissions		
The CROSSER™ Generator is intended for use in the electromagnetic environment specified below. The customer or the user of the CROSSER™ Generator should assure that it is used in such an environment.		
Emissions test	Compliance	Environmental guidance
RF emissions CISPR 11	Group 1	The equipment uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class A	The equipment is suitable for use in all establishments other than domestic, and may be used in domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes, provided the following warning is heeded:
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	Warning: The CROSSER™ Generator is intended for use by health care professionals only. The CROSSER™ Generator may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as re-orienting or relocating the equipment or shielding the location.
The CROSSER™ Generator is not suitable for interconnection with other equipment with the exception of the FLOWMATE® Injector.		
Portable and mobile RF communications equipment can affect operation of the CROSSER™ Generator.		
The use of accessories and cables other than those specified, with the exception of cables qualified and sold by Bard for use with the CROSSER™ Generator may result in increased emissions or decreased immunity of the equipment and may cause the equipment to be non-compliant with requirements of IEC60601-1-2:2007.		
The CROSSER™ Generator should not be used adjacent to or stacked with equipment other than the FLOWMATE® Injector. If adjacent or stacked use is necessary, the equipment should be observed to verify normal operation in the configuration in which it will be used.		

Electromagnetic Immunity			
The CROSSER™ Generator is intended for use in the electromagnetic environment specified below. The customer or the user of the CROSSER™ Generator should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Environmental guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV line(s) to line(s) ±2 kV line(s) to earth	±1 kV line(s) to line(s) ±2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11 U _T = 230 Vac	<5 % U _T (>95 % dip in U _T) for 0,5 cycle 40 % U _T (60 % dip in U _T) for 5 cycles 70 % U _T (30 % dip in U _T) for 25 cycles <5 % U _T (>95 % dip in U _T) for 5 sec	<5 % U _T (>95 % dip in U _T) for 0,5 cycle 40 % U _T (60 % dip in U _T) for 5 cycles 70 % U _T (30 % dip in U _T) for 25 cycles No anomalies 95% dip meets requirements.	Mains power quality should be that of a typical commercial or hospital environment. If the user of the equipment requires continued operation during power mains interruptions, it is recommended that the equipment be powered from an uninterruptible power supply or a battery.
Power frequency (50 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Electromagnetic Immunity			
The CROSSER™ Generator is intended for use in the electromagnetic environment specified below. The customer or the user of the CROSSER™ Generator should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Environmental guidance
Conducted RF IEC 61000-4-6	3 V _{rms} 150 kHz to 80 MHz	3 V _{rms}	Portable and mobile RF communications equipment should be used no closer to any part of the equipment including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = \left[\frac{3.5}{V_1} \right] \sqrt{P}$ $d = \left[\frac{3.5}{E_1} \right] \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = \left[\frac{7}{E_1} \right] \sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of known RF transmitting devices and equipment marked with the following symbol: 
NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.			
NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the equipment is used exceeds the applicable RF compliance level above, the equipment should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the equipment.			
b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

Recommended separation distances between portable and mobile RF communications equipment and the CROSSER™ Generator			
The CROSSER™ Generator is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the CROSSER™ Generator can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the CROSSER™ Generator as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter		
	150 kHz to 80 MHz $d = \left[\frac{3.5}{V_1} \right] \sqrt{P}$	80 MHz to 800 MHz $d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{E_1} \right] \sqrt{P}$
0.01	0.12m	0.12m	0.23m
0.1	0.37m	0.37m	0.95m
1	1.17m	1.17m	2.33m
10	3.69m	3.69m	9.49m
100	11.67m	11.67m	23.33m
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

Warnings and Precautions

- **The CROSSER™ Recanalization System should only be used by individuals trained in percutaneous transluminal angioplasty (PTA or PTCA).**
- **Warning: To avoid risk of electric shock, this equipment must only be connected to a supply mains with protective earth.**
- **Warning: No modification of this equipment is allowed.**
- **Warning: Risk of injury or death to operator or patient can result from unauthorized modification to CROSSER™ Generator.**
- **Note: Position CROSSER™ Generator to allow easy access to separable power cord during use.**
- **Prior to use, the CROSSER™ Generator, Foot Switch and Transducer should be inspected for signs of damage. Never use damaged components. Contact manufacturer for a replacement.**
- **The CROSSER™ Generator is not protected against the ingress of fluids (IPX0). Do not place fluids above or on the CROSSER™ Generator as damage to the CROSSER™ Generator or shock may result.**
- **Explosion Hazard: Do not use in the presence of flammable anesthetics.**
- **Electrical Shock Hazard: Do not attempt to service CROSSER™ Generator, Foot Switch or Transducer. Contact manufacturer for servicing requirements.**
- **In the event of a power outage, power interruption, or system reset, replace CROSSER™ Catheter before restarting recanalization procedure.**
- **This device generates, uses, and can radiate radio frequency energy and should only be used and installed according to these instructions for use.**
- **DO NOT activate the CROSSER™ Recanalization System without proper irrigation. Make sure to establish proper irrigation prior to introduction into guide catheter.**
- **Do not exceed 5 minutes of activation time as CROSSER™ Catheter malfunction may occur. If the 5 minutes of activation time is achieved, exchange the device for a new CROSSER™ Catheter before resetting the CROSSER™ Generator.**
- **Position Foot Switch and cable to minimize potential tripping hazard.**
- **Ensure CROSSER™ Generator is securely mounted to a 5 caster IV pole with minimum caster diameter of 7.62cm (3 inches) and a minimum base diameter of 55.9cm (22 inches) no higher than 1.5 meters off the ground.**
- **Should high frequency vibration fail to stop when Foot Switch is released, power CROSSER™ Generator 'OFF' or unplug from power receptacle.**
- **Keep out of reach of patients. The CROSSER™ Generator must be at least 1.5 meters away from the patient.**
- **Caution! The unit is equipped with an exhaust vent. Do not block the openings or internal heat build-up may occur, causing heat build-up or damage to the equipment.**

Required for Use

A CROSSER™ Catheter is required for use with the CROSSER™ Generator.

Set Up

Note: This symbol (IEC 60417-5021) denotes the Potential Equalization Conductor that is defined in IEC 60601-1 3rd edition as,

Conductor other than a Protective Earth Conductor or a neutral conductor, providing a direct connection between electrical equipment and the potential equalization busbar of the electrical installation.

Medically used rooms in most countries have no facilities for the use of detachable Potential Equalization Conductors. If however the CROSSER™ Generator is installed in locations where Potential Equalization Conductors are used, the appropriate requirements have been satisfied.

1. Set up CROSSER™ Generator on a 5 caster IV pole not more than 1.5 meters above ground. Position Foot Switch and cable to minimize potential tripping hazard.

Note: Refer to other sections of this operations manual for proper use and set up, cleaning and maintenance of the CROSSER™ Generator, Transducer and Foot Switch.

2. Connect power cord to back of CROSSER™ Generator and then to the wall receptacle. Make sure the power cord is positioned to minimize tripping hazard.

Note: To ensure grounding efficiency is maintained only connect the CROSSER™ Generator to a 'Hospital Grade' or 'Hospital Use' receptacle.

3. Connect Foot Switch to CROSSER™ Generator.

Warning! Do not switch on CROSSER™ Generator without first connecting the Foot Switch.

4. Connect CROSSER™ Catheter to transducer.

Note: See CROSSER™ Catheter IFU for details on CROSSER™ Catheter attachment.

Note: This connection to the CROSSER™ Catheter is considered a connection to an applied part.

5. Switch on CROSSER™ Generator using power switch on back of CROSSER™ Generator.

Operational Use

Activation

Activation occurs when the Foot Switch is depressed. An audible tone may be heard whenever the CROSSER™ Recanalization System is activated.

Timer Display

The timer will record CROSSER™ Catheter activation time in minutes and seconds. The CROSSER™ Generator will stop functioning once 5 minutes of activation has been achieved.

Warning! Do not exceed 5 minutes of activation time. If 5 minutes of activation time is achieved, exchange for a second CROSSER™ Catheter and press the RESET button on the face of the CROSSER™ Generator to reset the timer. Using a CROSSER™ Catheter for more than 5 minutes of activation time may cause CROSSER™ Catheter malfunction and is NOT permitted.

Warning! Only reset the CROSSER™ generator when a new CROSSER™ catheter is connected.

Note: CROSSER™ Generator will only deliver energy for a maximum of 30 seconds with one continuous depression of the Foot Switch and will not deliver more than 5 minutes total time without resetting the CROSSER™ Generator.

Shutdown Procedure

1. Refer to CROSSER™ Catheter IFU for catheter removal from patient and Transducer.
2. Switch the CROSSER™ Generator 'off'. The power switch is located at the mains connection on the back panel.

Cleaning

Warning! Only use a mild detergent solution to disinfect the unit.

Caution! Always disconnect the power cord before cleaning.

1. Clean the Transducer and Foot Switch with a cloth that has been dampened in the solution.

Caution! Do not immerse the Transducer in the solution.

2. Clean the CROSSER™ Generator by wiping with a cloth that has been dampened in the solution. Do not use any abrasive cleaners.

Caution! Do not immerse the CROSSER™ Generator. The CROSSER™ Generator unit is not sealed against the ingress of fluids. Damage to the CROSSER™ Generator and electrical shock could occur.

Maintenance

Performance of preventive maintenance is not essential. Routine inspection can, however, identify potential problems before they become serious, thus enhancing the instrument's reliability and extending its useful operating life. All National requirements for periodically inspecting the general condition of the device and other checks shall be followed by the user.

Warning! To reduce the risk of electrical shock, do not open the equipment's inner housing. Refer servicing to qualified personnel only. Removal of panels by unauthorized personnel will void the unit's warranty.

Note: Follow local governing ordinances and recycling plans regarding disposal or recycling of device components.

How Supplied

The CROSSER™ Generator, Foot Switch and Transducer are supplied ready for use.

Troubleshooting

Symptom	Cause	Solution
Power switch is ON but front panel display is OFF	Power cord not properly connected to wall receptacle or CROSSER™ Generator	Check connections
	Fuses blown	Replace fuses
	Wall receptacle defective or not supplying correct voltage	Check wall receptacle Return to manufacturer for servicing or replacement
Power switch is ON, front panel display is ON, Transducer not activating when Foot Switch is depressed	Foot Switch not connected or not connected properly	Check connection at CROSSER™ Generator
	Foot Switch is defective	Return to manufacturer for servicing or replacement
	Transducer is defective	Return to manufacturer for servicing or replacement
	Power switch turned off/on too quickly or unit unplugged during last use without first switching off	Turn off, wait 5 seconds, turn on
Excessive vibration noise coming from Transducer	Transducer/CROSSER™ Catheter connection loose	Tighten connection
Poor performance	CROSSER™ Catheter is not tightened properly to Transducer	Refer to CROSSER™ Catheter Instructions for Use for proper connection procedure
	Locking collar is not engaged	Slide locking collar on Transducer distally over proximal hub of the CROSSER™ Catheter
	CROSSER™ Catheter kinked or broken	Replace CROSSER™ Catheter
	Transducer defective	Return to manufacturer for servicing or replacement

Warranty and Servicing

Bard Peripheral Vascular warrants to the first purchaser of this product, that this product will be free from defects in materials and workmanship for a period of one year from the date of first purchase and liability under this limited product warranty will be limited, to repair or replacement of the defective product, in Bard Peripheral Vascular's sole discretion, or refunding your net price paid. Wear and tear from normal use or defects resulting from misuse of this product are not covered by this limited warranty.

TO THE EXTENT ALLOWABLE BY APPLICABLE LAW, THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL BARD PERIPHERAL VASCULAR BE LIABLE TO YOU FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM YOUR HANDLING OR USE OF THIS PRODUCT.

Some countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your country.

An issue or revision date and revision number for these instructions are included for the user's information on the last page of this booklet. In the event 36 months have elapsed between this date and product use, the user should contact Bard Peripheral Vascular to see if additional product information is available.

Contact Bard for servicing information.

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