ENGLISH
CROSSER™ Recanalization System
FLOWMATE® Injector - Manual of Operations

Instructions for Use

Warning! Carefully read all instructions prior to use. Observe all warnings and precautions noted throughout these instructions. Failure to do so may result in a hazardous situation.

Caution! Federal law restricts this device to sale by or on the order of a physician.

Description

The CROSSER™ Recanalization System consists of a re-usable electronic CROSSER™ Generator, Foot Switch, high frequency Transducer, the FLOWMATE® Injector (optional) and single-use CROSSER™ Catheter. A CROSSER™ Catheter is connected to the electronic CROSSER™ Generator through the Transducer. The Foot Switch is used to activate the system. The CROSSER™ Generator and Transducer convert AC power into high frequency mechanical vibrations which are propagated through a core wire to the metal tip of the CROSSER™ Catheter. During use (or activation of the CROSSER™) saline is required to flush the CROSSER™ Catheter. A FLOWMATE® Injector or other injection system is required.

The FLOWMATE® Injector is a sterile saline injector for use as an accessory to the CROSSER™ Recanalization System. The FLOWMATE® Injector is designed to infuse saline at a rate of 0.3ml/sec at a maximum of 200 PSI. The Injector is operated by a Foot Switch. A jumper cable is supplied with the FLOWMATE® Injector to allow connection to the CROSSER™ Generator. The FLOWMATE® Injector is programmed to first start the FLOWMATE® Injector, and then the CROSSER™ Generator each time the Foot Switch is depressed.

Figure 1 – The CROSSER™ Recanalization System


Figure 2 - FLOWMATE® Injector Drawings


Figure 3- Control Panels / Buttons

Step 1

SET-UP MODE

FLASHING LIGHT, Set-up Required
STEADY LIGHT, Ready to Fill Syringe

Step 2

RUN MODE

STEADY LIGHT, Ready to Activate

ALARM

Press to Reset Alarm

Withdraw Button

Advance Button

Indications

The FLOWMATE® Injector is indicated for use as a sterile saline injector for use as an accessory to the CROSSER™ Recanalization System.

Refer to the CROSSER™ Generator and the CROSSER™ Catheter Instructions for Use for complete system operation and interventional use procedure.

Specifications	
Model No.	FG1050-02
Catalog No.	INJ100
Weight (all components)	17.5 lbs (8.0kg)
Dimensions (Injector alone)	20 x 12 x 9 in. (51 x 30 x 23cm)
Electrical	100-120/220-240 V~at 50/60 Hz 2/1 A
Power Cord	European Cord: 3 conductor, 6A, 250 V, 1.0 mm ² , <HAR> cord North American Cord: 3 conductor, 6A, 250 V, 18AWG, CSA/UL cord
Mode of Operation	Continuous operation
Injector Speed	Setup: forward or reverse 5 ml/sec +/- 11% Run: 0.3 ml/sec +/- 5%
Pressure Limit	200 psi (1379 kPa)
Alarms	Overpressure (at 200psi) – yellow light, three pulse beep Low pressure (at approximately 6 psi) – yellow light, three pulse beep Empty syringe - blue light turns off, single beep Missing syringe – yellow light, three pulse beep
Safety	Meets IEC60601-1, 2 nd edition; UL 60601-1, 2 nd edition; CAN/CSA C22.2 No. 601.1 CSA Certified to above Canadian, US and International Medical Standards
EMC	See tables below
Ingress of Fluids	Injector – IPX0, ordinary, Foot Switch – IPX8
Protection Against Electrical Shock	Type CF Applied Part ; Degree: Class I
Fuses	Two fuses rated at 2 X T2AL, 250V
Activation Time Limits	Single depress of Foot Switch 3 second delay to start infusion plus 30 seconds maximum of activation of the CROSSER™ Generator. Time accuracy is ± 1 second.
Transport and Storage Environment	Temperature: -10 °C (14 °F) to +55 °C (131 °F) Relative Humidity: 10% to 85% Non-condensing Atmospheric Pressure: 50-106 kPa
Operating Environment	Temperature: +10 °C (50 °F) to +40 °C (104 °F) Relative Humidity: 30% to 75% Non-condensing Atmospheric Pressure: 70-106 kPa

All operating parameters are specified at typical room temperatures (70 degrees F or 21 degrees C)

This equipment has been tested and found to comply with the EMC limits for IEC 60601-1-2:2007. These limits are designed to provide reasonable protection against harmful interference in a typical medical installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with these instructions, may cause harmful interference to other devices in the vicinity. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference with other devices, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving device
- Increase the separation between the equipment
- Connect the equipment into an outlet on a circuit different from that to which the other device(s) is connected
- Consult the manufacturer or field service technician for help.

Electromagnetic Emissions		
The FLOWMATE® Injector is intended for use in the electromagnetic environment specified below. The customer or the user of the FLOWMATE® Injector should assure that it is used in such an environment.		
Emissions Test	Compliance	Environmental Guidance
RF emissions CISPR 11	Class 1	The FLOWMATE® Injector uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class A	The FLOWMATE® Injector is suitable for use in all establishments other than domestic, and may be used in domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes, provided the following warning is heeded: Warning: This equipment/system is intended for use by healthcare professionals only. This equipment/ system may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as re-orienting or relocating the FLOWMATE® Injector or shielding the location.
Harmonic emissions IEC 61000 3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	
The FLOWMATE® Injector is not suitable for interconnection with other equipment with the exception of the CROSSER™ Generator		
Portable and mobile RF communications equipment can affect operation of the FLOWMATE® Injector.		
The use of accessories and cables other than those specified, with the exception of cables qualified and sold by Bard for use with the FLOWMATE® Injector may result in increased emissions or decreased immunity of the equipment and may cause the equipment to be non-compliant with requirements of IEC60601-1-2:2007.		
The FLOWMATE® Injector should not be used adjacent to or stacked with equipment other than the CROSSER™ Generator. If adjacent or stacked use is necessary, the equipment should be observed to verify normal operation in the configuration in which it will be used.		

Electromagnetic Immunity			
The FLOWMATE® Injector is intended for use in the electromagnetic environment specified below. The customer or the user of the FLOWMATE® Injector should assure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Environmental Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines, ±1 kV for input/output lines	±2 kV for power supply lines, ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV for line(s) to line(s), ±2 kV for line (s) to earth	±1 kV for line(s) to line(s), ±2 kV for line (s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11 U _T =230 VAC	<5% U _T (>95% dip in U _T) for 0.5 cycle 40% U _T (60% dip in U _T) for 5 cycles 70% U _T (30% dip in U _T) for 25 cycles <5% U _T (>95% dip in U _T) for 5 sec	<5% U _T (>95% dip in U _T) for 0.5 cycle 40% U _T (60% dip in U _T) for 5 cycles 70% U _T (30% dip in U _T) for 25 cycles <5% U _T (>95% dip in U _T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the FLOWMATE® Injector requires continued operation during power mains interruptions, it is recommended that the FLOWMATE® Injector be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Electromagnetic Immunity			
The FLOWMATE® Injector is intended for use in the electromagnetic environment specified below. The customer or the user of the FLOWMATE® Injector should assure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Environmental Guidance
Conducted RF IEC 61000-4-6	3 V _{rms} 150 kHz to 80 MHz	3V	Portable and mobile RF communications equipment should be used no closer to any part of the FLOWMATE® Injector, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Separation distance according to frequency of transmitter $d = \left[\frac{3.5}{V_1} \right] \sqrt{P}$ $d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$ $d = \left[\frac{7}{E_1} \right] \sqrt{P}$ Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (mm). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3V/m	
<p>NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.</p> <p>NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			
<p>^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the FLOWMATE® Injector is used exceeds the applicable RF compliance level above, the FLOWMATE® Injector should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the FLOWMATE® Injector.</p> <p>^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.</p>			

Recommended separation distances between portable and mobile RF communications equipment and the FLOWMATE® Injector

The FLOWMATE® Injector is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the FLOWMATE® Injector can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the FLOWMATE® Injector as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter		
	150 kHz to 80 MHz $d = \left[\frac{3.5}{V_1} \right] \sqrt{P}$	80 MHz to 800 MHz $d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{E_1} \right] \sqrt{P}$
0.01	0.12m	0.12m	0.23m
0.1	0.37m	0.37m	0.95m
1	1.17m	1.17m	2.33m
10	3.69m	3.69m	9.49m
100	11.67m	11.67m	23.33m

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Warnings and Precautions

- The FLOWMATE® Injector should only be used by individuals trained in percutaneous transluminal angioplasty (PTA or PTCA).
- **WARNING: To avoid the risk of electric shock, this equipment must only be connected to a supply mains with protective earth.**
- **WARNING: No modification of this equipment is allowed.**
- **WARNING: Risk of injury or death to operator or patient can result from unauthorized modification to FLOWMATE® Injector.**
- **Note: Position FLOWMATE® Injector to allow easy access to separable power cord during use.**
- Prior to use, the FLOWMATE® Injector, Foot Switch and all connection to power and the CROSSER™ Generator should be inspected for signs of damage. Never use damaged components. Contact the manufacturer for a replacement or servicing.
- The FLOWMATE® Injector is intended to be used with sterile normal saline only. Do not use the FLOWMATE® Injector with contrast media.
- **Explosion Hazard: Do not use in the presence of flammable anesthetics.**
- **Electrical Shock Hazard: Do not attempt to service FLOWMATE® Injector or Foot Switch. Contact manufacturer for servicing requirements.**
- The Grounding reliability can only be achieved when the equipment is connected to an equivalent receptacle marked Hospital Only or Hospital Grade. Grounding reliability cannot be achieved by connecting the system to a Multiple Portable Socket-Outlets (MPSO – power strip) or extension cord.
- In the event of a power outage or power interruption, run through the set-up procedures. The FLOWMATE® Injector is designed to prevent operation following such events until the set-up procedures are repeated.
- This FLOWMATE® Injector generates, uses, and can radiate radio frequency energy and should only be used and installed according to these instructions for use.
- Position Foot Switch and cable to minimize potential tripping hazard.
- Ensure FLOWMATE® Injector is securely mounted to a 5 caster IV pole with minimum caster diameter of 7.62cm (3 inches) and a minimum base diameter of 55.9 cm (22 inches), no higher than 1.3 meters (51 inches) off the ground. The FLOWMATE® Injector should not be mounted above the CROSSER™ Generator.
- Should high frequency vibration or saline infusion fail to stop when Foot Switch is released, power FLOWMATE® Injector/ CROSSER™ Generator 'OFF' or unplug from power receptacle.
- Do not reuse disposable items. Dispose of per normal hospital procedures.
- **Air Embolism Risk: to minimize air embolism risk, make certain the operator is trained on proper syringe loading techniques. Before injecting, be sure to clear air trapped in the syringe, infusion line and the CROSSER™ Catheter. Do not introduce the CROSSER™ Catheter into the patient until all air has been purged from the infusion line and the CROSSER™ Catheter. A steady stream of saline exiting from the tip of the CROSSER™ Catheter will be observed once all air has been purged.**

Required for Use

- A MedRad (part#150-FT-Q), Coeur (part #C853-0155), or DeRoyal (part #77-400277) 150ml power injection syringe (these part numbers are provided as a reference only, please verify the part numbers with the manufacturer when placing any orders)
- A minimum of 48" (122cm) of sterile contrast infusion line
- Normal Sterile Saline (recommended to be refrigerated to 45°F [7°C] prior to loading. See CROSSER™ Catheter Instructions for Use).

Out of Box Set Up / Installation

Note: This symbol (IEC 60417-5021) denotes the Potential Equalization Conductor that is defined in IEC 60601-1, 3rd Edition, as

Conductor other than a Protective Earth Conductor or a neutral conductor, providing a direct connection between electrical equipment and the potential equalization busbar of the electrical installation.

Medically used rooms in most countries have no facilities for the use of detachable Potential Equalization Conductors. If, however, the FLOWMATE® Injector is installed in locations where Potential Equalization Conductors are used, the appropriate requirements have been satisfied.

1. Set up FLOWMATE® Injector on a 5 caster IV pole not more than 1.3 meters above ground using the IV pole clamps on the back of the FLOWMATE® Injector. Position Foot Switch and cable to minimize potential tripping hazard and tighten the clamps firmly to prevent the FLOWMATE® Injector from sliding down the pole.
2. Connect power cord to back of FLOWMATE® Injector and then to the wall receptacle. Make sure the power cord is positioned to minimize tripping hazard.

Note: To ensure grounding efficiency is maintained only connect the FLOWMATE® Injector to a 'Hospital Grade' or 'Hospital Use' receptacle.

3. Connect the Foot Switch Jumper Cable from the FLOWMATE® Injector (see figure 2B) to the CROSSER™ Generator (the Foot Switch connection; Z).
4. Connect Foot Switch (supplied with the CROSSER™ Generator) to the FLOWMATE® Injector, the connection is labeled Z (see figure 2C).
5. The FLOWMATE® Injector is now ready for operational use.

Note: It is recommended to test the FLOWMATE® Injector by following the Operational Use instructions below during training of the staff and prior to use in an actual procedure.

Operational Use

The FLOWMATE® Injector has two modes of operation, Set Up and Run. In Set Up Mode the Turret is placed in the upward position. This is the position where the syringe may be loaded, saline may be withdrawn into the syringe and air is purged prior to placing the FLOWMATE® Injector into Run Mode. When the Turret is rotated into the downward position the FLOWMATE® Injector is placed into Run Mode. Placing the Turret in between the two mode locations results in an idle configuration and a flashing White setup light. No motion of the pump will be allowed in this state.

CAUTION! Do not rotate the Turret into Run Mode until proper Set Up has been completed.


Set Up Mode

1. The FLOWMATE® Injector is in Set Up mode whenever the turret is in the 'up' position. Make sure the FLOWMATE® Injector is in the upward position before initiating the start up procedure.
2. Switch the unit 'on'. The power switch is located by the ratings label of the stationary portion of the FLOWMATE® Injector. The White light next to the Set Up icon should be flashing.
3. Open Syringe Locking Mechanism and insert a compatible 150ml syringe. See Required For Use section above for a list of compatible syringes.

Make sure the syringe is seated fully into the turret, with the flat edge on the syringe flange aligned with the corresponding edge on the Turret.

Note: The FLOWMATE® Injector is equipped with a syringe sensor. A compatible syringe must be seated correctly for the FLOWMATE® Injector to run properly.

Warning! Do not attempt to fill the syringe by removing the plunger. The FLOWMATE® Injector is designed to be set up and used only with a new syringe.

- Close the Syringe Locking Mechanism. The **White** light next to the Set Up icon should be steady on. If **FLASHING**, ensure the Turret is rotated fully into the Set Up position and the syringe Locking mechanism is fully closed.
- Press and hold the Advance button to move the plunger all the way forward until the plunger head stops.
- Connect the saline supply and withdraw saline into the FLOWMATE® Injector by pressing and holding the Withdraw button. Keep pressing until the plunger is withdrawn all the way (the plunger will stop once the plunger is completely withdrawn).


Note: During the saline loading process a vacuum may be created within the syringe. It is acceptable to momentarily release the Withdraw button to allow the incoming saline to equalize.

- Gently tap on the side of the syringe neck / syringe retainer to ensure that no air bubbles remain in the syringe.
- Once the saline has been loaded disconnect the saline supply and connect the infusion line leading to the CROSSER™ Catheter. It is recommended that an assistant off the sterile field press the Advance button to infuse saline through the infusion line to provide a fluid to fluid connection to the CROSSER™ Catheter.

Note: This saline connection to the CROSSER™ Catheter is considered a connection to an applied part.

Warning! Do not introduce the CROSSER™ Catheter into the patient until all air has been purged from the infusion line and the CROSSER™ Catheter. A steady stream of saline exiting from the tip of the CROSSER™ Catheter will be observed once all air has been purged.

- Purge the air from the system by momentarily pressing the Advance button multiple times each for 3-5 second duration. Advance the plunger at least until its sealing edge aligns with the mark around the clear syringe retainer (see drawing below).

Note: The FLOWMATE® Injector is equipped with an overpressure alarm that is intended to sense a closed system. In the event there are kinks in the infusion line, the plunger is advanced too rapidly in SET UP mode with a CROSSER™ Catheter attached, or the tip of the CROSSER™ Catheter is plugged the overpressure alarm will sound. If the alarm sounds check the fluid path and then reset the alarm condition by depressing the Alarm Reset button. While the FLOWMATE® Injector is in alarm mode the Foot Switch and Injector motor are both disabled.


- Once all air has been purged, the FLOWMATE® Injector is ready to be placed into **Run Mode**.

Run Mode

- To place the FLOWMATE® Injector into Run Mode, rotate the turret completely so that the syringe barrel is pointing downward (approximately 20 degrees from horizontal). The **Green** light next to the Run Mode icon should be steady on and the **Blue** light on the top of the unit will be flashing.

Note: The FLOWMATE® Injector will not be functional unless the Turret is fully rotated into either the Set Up or Run Mode.

Warning! Once the FLOWMATE® Injector is placed in Run Mode the Foot Switch is active and when depressed will activate the FLOWMATE® Injector and then the CROSSER™ Generator.

- The FLOWMATE® Injector is now armed and may be operated by depressing the Foot Switch or the Advance button (the Withdraw button is disabled during Run Mode).
- To activate the CROSSER™ Recanalization System depress the Foot Switch. The FLOWMATE® Injector will infuse saline for a period of 3 seconds while emitting 3 audible tone bursts at 1 second intervals, and then the CROSSER™ Generator will be activated. The FLOWMATE® Injector will continue to infuse saline as long as the Foot Switch is depressed or until the saline reservoir in the syringe is depleted. Whenever the Injector motor is actively infusing the **Blue** light will be steady on.

Note: The CROSSER™ Generator will continue to deliver energy while the Foot Switch is depressed for a maximum of 30 seconds following which the audible tone from the CROSSER™ Generator will cease. Refer to the CROSSER™ Generator Manual of Operations for more detail.

Note: The Advance button will operate the FLOWMATE® Injector in Run Mode (and emit 3 audible tone bursts at 1 second intervals), but will not activate the CROSSER™ Generator. The CROSSER™ Generator can only be activated by the Foot Switch.

Note: The FLOWMATE® Injector is equipped with an overpressure alarm that is intended to sense a closed system. In the event there are kinks or an obstruction in the infusion line the overpressure alarm will sound. If the alarm sounds check the fluid path for possible obstructions. The alarm may be reset by depressing the Alarm Reset button, should the alarm fail to reset see the Troubleshooting section of this manual.

Note: While the FLOWMATE® Injector is in alarm mode, the Foot Switch and Injector motor are both disabled.

- The FLOWMATE® Injector will continue to operate until all saline has been depleted.

Shutdown Procedure

- Refer to CROSSER™ Catheter IFU for catheter removal from patient.
- Disconnect saline connection to CROSSER™ Catheter.
- Place FLOWMATE® Injector into Set Up mode.
- Withdraw the plunger by pressing and holding the Withdraw button. Keep pressing until the plunger is withdrawn all the way (the plunger will stop once the plunger is completely withdrawn).
- Disconnect the saline supply from the end of the syringe. Open the Syringe Locking Mechanism and remove the compatible 150 ml syringe.
- Switch the unit off. The power switch is located at the mains connection on the back panel.
- Remove the 150ml power injection syringe and infusion line. Dispose of per normal hospital procedures.

Cleaning

Warning! Only use a mild detergent solution to disinfect the unit.

Caution! Always disconnect the power cord before cleaning.

- Clean the Foot Switch with a cloth that has been dampened in the solution.
- Clean the FLOWMATE® Injector by wiping with a cloth that has been dampened in the solution. Do not use any abrasive cleaners.

Caution! Do not immerse the FLOWMATE® Injector. The FLOWMATE® Injector is not sealed against the ingress of fluids. Damage to the FLOWMATE® Injector and electrical shock could occur.

Maintenance

Performance of preventive maintenance is not essential. Regular maintenance may, however, contribute to identifying potential problems before they become serious, thus enhancing the instrument's reliability and extending its useful operating life. All local requirements for periodically inspecting the general condition of the FLOWMATE® Injector and other checks shall be followed by the user.

Warning! To reduce the risk of electrical shock, do not open the equipment's inner housing. Refer servicing to manufacturer only. Removal of panels by unauthorized personnel will void the unit's warranty.

Note: Follow local governing ordinances and recycling plans regarding disposal or recycling of FLOWMATE® Injector (INJ100) and its components at the end of its useful life.

How Supplied

The FLOWMATE® Injector is supplied ready for use following the simple Out of Box Set Up Instructions above.

Troubleshooting		
Symptom	Cause	Solution
Power switch is ON but <u>White</u> light next to the Set Up icon (or Run Mode Icon) is OFF	Power cord not properly connected to wall receptacle or FLOWMATE® Injector	Check connections
	Fuses blown	Replace fuses
	Wall receptacle defective or not supplying correct voltage	Check wall receptacle Return to manufacturer for servicing or replacement
Power switch is ON and <u>White</u> light next to the Set Up icon is FLASHING	FLOWMATE® Injector needs to be placed into Set Up mode with syringe installed	Rotate Turret fully into Set Up mode and verify syringe is installed correctly.
Power switch is ON, <u>Green</u> light next to Run Mode Icon is ON, <u>Blue</u> Light is FLASHING, but FLOWMATE® Injector, CROSSER™ Generator, or both are not activating when Foot Switch is depressed	Foot Switch not connected or not connected properly	Check connection at FLOWMATE® Injector
	Foot Switch jumper cable not connected or not connected properly	Check connection at CROSSER™ Generator
	Foot Switch is defective	Return to manufacturer for servicing or replacement
	CROSSER™ Generator/Transducer is defective	Return to manufacturer for servicing or replacement
	Power switch turned off/on too quickly or unit unplugged during last use without first switching off	Turn FLOWMATE® Injector and CROSSER™ Generator off, wait 5 seconds, turn on; Run FLOWMATE® Injector through set up procedure
Malfunctioning FLOWMATE® Injector	Return to manufacturer for servicing or replacement	
Post-procedure, the syringe door will not open, and there is difficulty removing the syringe.	The syringe plunger prevents the syringe door from opening.	Tap the Advance button to move the plunger forward just slightly, and then open the syringe door. Total advancement should be ½ mm maximum.
	Incompatible syringe	Only use a compatible 150ml syringe* Return to manufacturer for servicing or replacement
Plunger will not withdraw when Withdraw Button is pressed	Syringe not properly loaded	Check seating of syringe
	FLOWMATE® Injector is not in Set Up mode	Rotate turret to Set Up position. Withdraw button is deactivated in Run Mode
	Incompatible syringe	Only use a compatible 150ml syringe* Return to manufacturer for servicing or replacement
Plunger will not advance or withdraw when Advance or Withdraw Button is pressed	Syringe not properly loaded	Check seating of syringe
	Malfunctioning FLOWMATE® Injector	Return to manufacturer for servicing or replacement
	Incompatible syringe	Only use a compatible 150ml syringe*
Alarm sound (Three pulses, repeating) sounds when Foot Switch is depressed	Over pressure	Check system for blockages, reset alarm
Alarm sound (Three pulses, repeating) when FLOWMATE® Injector is first powered on	System failure	Return to manufacturer for servicing or replacement
CROSSER™ Catheter feels hot to the touch during use	Blocked irrigation lumen	Dispose of CROSSER™ catheter, get a new CROSSER™ catheter, notify manufacturer
	Irrigation equipment not turned 'ON'	Check irrigation system
	Irrigation equipment not set properly or malfunctioning irrigation equipment	Check irrigation system
	Leak in irrigation system	Check all luer connections
	Improper hand/finger placement on the CROSSER™ Catheter	Avoid creating a focal bend on the CROSSER™ catheter shaft. Use only two fingers to advance the CROSSER™ Catheter.

*See required for use section for list of compatible syringes.

Warranty and Servicing

Bard Peripheral Vascular warrants to the first purchaser of this product, that this product will be free from defects in materials and workmanship for a period of one year from the date of first purchase and liability under this limited product warranty will be limited, to repair or replacement of the defective product, in Bard Peripheral Vascular's sole discretion, or refunding your net price paid. Wear and tear from normal use or defects resulting from misuse of this product are not covered by this limited warranty.

TO THE EXTENT ALLOWABLE BY APPLICABLE LAW, THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL BARD PERIPHERAL VASCULAR BE LIABLE TO YOU FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM YOUR HANDLING OR USE OF THIS PRODUCT.

Some countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your country.

An issue or revision date and revision number for these instructions are included for the user's information on the last page of this booklet. In the event 36 months have elapsed between this date and product use, the user should contact Bard Peripheral Vascular to see if additional product information is available.

Contact Bard for servicing information.

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