# IMPRA®

## ePTFE Vascular Grafts

# Instructions for Use

### Instructions for use

IMPRA® ePTFE Vascular Grafts

IMPRA CARBOFLO® ePTFE Vascular Grafts

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# INSTRUCTIONS FOR USE

Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

## Device Description, Indications, Contraindications, Warnings, Precautions and Adverse Reactions

### Device Description

IMPRA® ePTFE vascular grafts are constructed of expanded polytetrafluoroethylene (ePTFE). IMPRA CARBOFLO® ePTFE vascular grafts also contain carbon impregnated into the inner portions of the graft walls.

### Indications for Use

IMPRA® ePTFE grafts are indicated for use as vascular prostheses.

**Straight, Tapered, Short Tapered, Stepped, CENTERFLEX™ Graft and ENDFLEX™ Graft Configurations** are intended for use as subcutaneous arteriovenous conduits for blood access, bypass, or reconstruction of peripheral arterial blood vessels. Tapered, short tapered, and stepped configurations may help minimize the risk of steal syndrome and high cardiac output. **CENTERFLEX™** graft and **ENDFLEX™** graft configurations have a non-removable external spiral support (beading) and can be used where resistance to compression or kinking is desired.

**IMPRA FLEX™ Graft Configurations** are intended for bypass or reconstruction of peripheral arterial blood vessels and have removable spiral support (beading) over the entire graft. These grafts can be used where resistance to compression or kink is desired.

Insufficient clinical data are available on which to base any conclusions regarding the use of Thinwall grafts in blood access or to support the use of IMPRA® vascular grafts for applications involving; pulmonary arteries, cerebral arteries, coronary arteries, brachiocephalic trunk, cardiac vein, pulmonary veins, or the inferior or superior vena cava.

### Contraindications

None known.

### Warnings

1. All IMPRA® ePTFE vascular grafts are supplied sterile and non-pyrogenic unless the package is open or damaged. IMPRA® ePTFE vascular grafts are sterilized by ethylene oxide. Each graft is intended for single patient use only. **DO NOT RESTERILIZE.**
2. This device has been designed for single use only. Reusing this medical device bears the risk of cross-patient contamination as medical devices – particularly those with long and small lumina, joints, and/or crevices between components – are difficult or impossible to clean once body fluids or tissues with potential pyrogenic or microbial contamination have had contact with the medical device for an indeterminable period of time. The residue of biological material can promote the contamination of the device with pyrogens or microorganisms which may lead to infectious complications.
3. Do not resterilize. After resterilization, the sterility of the product is not guaranteed because of an indeterminable degree of potential pyrogenic or microbial contamination which may lead to infectious complications. Cleaning, reprocessing and/or resterilization of the present medical device increases the probability that the device will malfunction due to potential adverse effects on components that are influenced by thermal and/or mechanical changes.
4. Do not use after expiration date printed on the label.
4. Anastomotic or graft disruption has been associated with Axillofemoral, Femoral Femoral, or Axillobifemoral bypass procedures if implanted improperly. Refer to Specific Operative Procedures (Extra-Anatomic Bypass Procedures) for further instructions. Thinwall and IMPRA FLEX™ Thinwall grafts are NOT recommended for these types of bypass procedures.
6. For Extra-Anatomic procedures, (e.g., Axillofemoral, Femoral Femoral, or Axillobifemoral Bypass) the patient should be cautioned that sudden, extreme or strenuous movements should be totally avoided for a period of at least six to eight weeks to allow for proper stabilization of the graft. Routine activities such as raising the arms above the shoulders, reaching out in front, extended reaching, throwing, pulling, striding or twisting should be avoided.
7. IMPRA® ePTFE grafts do not stretch (are non-elastic) in the longitudinal direction. The correct graft length for each procedure must be determined by considering the patient's body weight, posture, and the range of motions across the anatomical area of graft implantation. Failure to cut the grafts to an appropriate length may result in anastomotic or graft disruption, leading to excessive bleeding, and loss of limb or limb function, and/or death.
8. Aggressive and/or excessive graft manipulation when tunneling, or placement within a too tight or too small tunnel, may lead to separation of the spiral beading and/or graft breakage.
9. When embolectomy or balloon angioplasty catheters are used within the lumen of the graft, the inflated balloon size must match the inner diameter of the graft. Over-inflation of the balloon or use of an inappropriately sized balloon may dilate or damage the graft.
10. Do not remove the external spiral support (beading) from any CENTERFLEX™ graft or ENDFLEX™ graft. Attempts to remove the beading may damage the graft wall. If damage occurs, discard the graft.
11. Do not cannulate IMPRA FLEX™ grafts with external support over the full length of the graft or the externally supported portion of CENTERFLEX™ or ENDFLEX™ grafts. Cannulation at these sites may lead to beading embolization and/or pseudoaneurysm formation.
12. Avoid repeated or excessive clamping at the same location on the graft. If clamping is necessary, use only atraumatic or appropriate vascular smooth jawed clamps to avoid damage to the graft wall.

13. Exposure to solutions (e.g., alcohol, oil, aqueous solutions, etc.) may result in loss of the graft's hydrophobic properties. Loss of the hydrophobic barrier may result in graft wall leakage. Preclotting of this graft is unnecessary.
14. Avoid excessive graft manipulation after exposure to blood or body fluids. Do not forcibly inject any solution through the lumen of the graft, or fill the graft with fluid prior to pulling it through the tunnel as loss of the graft's hydrophobic properties may occur. Loss of the hydrophobic barrier may result in graft wall leakage.
15. Do NOT expose IMPRA® ePTFE grafts to temperatures greater than 500°F (260°C). PTFE decomposes at elevated temperatures producing highly toxic decomposition products <sup>1</sup>.
16. After use, the product may be a potential biohazard. Handle and dispose of in accordance with acceptable medical practice and applicable local, state and federal laws and regulations.
17. During tunneling create a tunnel that closely approximates the outer diameter of the graft. A tunnel that is too loose may result in delayed healing and may also lead to perigraft seroma formation.

### Precautions

1. Only physicians qualified in vascular surgery techniques should use this prosthesis. The healthcare provider is responsible for all appropriate postoperative care instructions to the patient.
2. The healthcare provider must observe aseptic technique during implantation and postoperatively.
3. When removing the external spiral support (beading) of IMPRA FLEX™ grafts, the beading must be removed slowly and at a 90° angle to the graft. Rapid unwinding and/or removal at less than a 90° angle may result in graft damage. Do not use surgical blades or sharp, pointed instruments to remove the beading as this may damage the graft wall. If damage occurs, that segment of the graft should not be used. (Reference Figure 1.)

**Note: Do not remove spiral support beading from CENTERFLEX™ or ENDFLEX™ grafts.**

Figure 1

4. When suturing, avoid excessive tension on the suture line, inappropriate suture spacing and bites, and gaps between the graft and host vessel. Failure to follow correct suturing techniques may result in suture hole elongation, suture pull-out, anastomotic bleeding and/or disruption. Refer to "Suturing" for further instructions.
5. To minimize lymph collection around the graft in Extra-Anatomic bypass procedures or in peripheral reconstructive procedures, the lymphatics should be carefully ligated and sealed, especially in the groin area.
6. Consider intraoperative and postoperative patient anticoagulation therapy for each patient as appropriate.

### Adverse Reactions

Potential complications which may occur with any surgical procedure involving a vascular prosthesis include, but are not limited to: disruption or tearing of the suture line, graft, and/or host vessel; suture hole bleeding; graft redundancy; thrombosis, embolic events, occlusion or stenosis; ultrafiltration; seroma formation; swelling of the implanted limb; formation of hematoma or pseudoaneurysm; infection; aneurysm/dilation; blood leakage; hemorrhage; steal syndrome; and/or skin erosion.

# DIRECTIONS FOR USE

### Equipment Required

Tunneler, suture, atraumatic clamp, and/or scissors.

### Opening the Package

Hold the outer tray in one hand. Peel back the lid. Remove the inner tray. Peel back the inner tray lid slowly and carefully remove the graft using sterile atraumatic instruments or sterile gloves. Protect the graft against damage from sharp or heavy instruments.

## General Operative Techniques

### Tunneling Techniques

Prior to utilizing a sheath tunneler, verify that the graft outer diameter dimensions fit the sheath internal diameter dimensions.

Use of a sheath tunneler is recommended as it will minimize graft handling, protect the external reinforcement (**Flex** Beading) and help maintain graft integrity.

Always follow the instructions for use for the specific tunneler utilized to place the graft.

Create a tunnel that closely approximates the outer diameter of the graft. A tunnel that is too loose may result in delayed healing and may lead to perigraft seroma formation. Reference **WARNING #3 and #17.**

**Anastomotic Preparation of Grafts with Flex Beading**

Use sterile technique to remove the **Flex** Beading. Hold the graft firmly in one gloved hand. With the other gloved hand, gently grasp and lift the end of the beading with gloved fingers. Slowly unwind the beading by peeling at a 90° angle to the graft. (Reference **Figure 1** and **PRECAUTION #3**). Remove and trim off enough beading to cut the graft to length.

**Suturing**

Size the graft appropriately to minimize excessive tension at the suture line. Use a tapered, noncutting needle with a nonabsorbable monofilament suture approximately the same size as the needle. Take 2mm suture bites in the graft following the curve of the needle and gently pull the suture at a 90° angle. Proper sizing of the graft length prior to implant will minimize suture hole elongation caused by excessive tension. Reference **PRECAUTION #3** and **WARNING #7**.

**Thrombectomy**

Techniques for declotting IMPRA® grafts include but are not limited to the use of balloon catheters. Reference **WARNING #9**.

**Longitudinal Incision:** Place stay sutures before introducing the embolectomy catheter. For IMPRA FLEX™ products, cut through the spiral beading and base graft. The spiral beading will realign itself after closure. Place a longitudinal incision in the graft that is long enough to accommodate the extraction of a fully dilated thrombectomy catheter balloon. A patch may be considered as an aid to graft closure.

**Transverse Incision:** No stay sutures are necessary. A horizontal mattress suture is recommended for graft closure.

During the early postoperative period, the natural progression of healing renders the graft translucent in appearance. In this state, a longitudinal incision with stay sutures is recommended. If a transverse incision is performed, a horizontal mattress suture technique and PTFE pledgets may aid in closure.

**Angiography**

Should angiography be performed at the time of procedure, the artery proximal to the graft should be used for injection if possible.

## Specific Operative Procedures

**Extra-Anatomic Bypass Procedures (e.g., Axillofemoral, Femoral Femoral, and Axillobifemoral)**

For Extra-Anatomic bypass procedures, careful attention must be given to the following techniques. Failure to follow these technical considerations may result in suture hole elongation, mechanical disruption, or tearing of the graft, suture line or host vessel, thrombosis, extreme blood loss, loss of limb function, loss of limb, or death. Reference **WARNING #7** and **PRECAUTION #4**.

- To avoid extreme stress on the anastomosis and the graft, include the patient's weight and range of limb motion when determining graft length, tunnel length and location.
- To determine the correct graft length, drape the patient to allow full movement of the arm, shoulder girdle or legs.
- Avoid protracted hyperabduction of the arm, during the surgical procedure. Prolonged hyperabduction may lead to brachial plexus injury.
- Allow sufficient graft length to avoid stressing of axillary or femoral anastomosis throughout the full range of movement of the arm, shoulder girdle, or legs. The graft should be placed under both the pectoralis major and pectoralis minor. (Reference **Figure 2**.)

**Figure 2**

- Cutting the graft slightly longer than necessary has been reported by some surgeons to reduce further the risk of stressing the graft or the anastomosis.
- Correctly bevel the axillary anastomosis. Stress on the graft is minimized when the graft is placed perpendicular to the axillary artery. Therefore, the anastomotic angle should be as small as possible and should not exceed 25° relative to the cut edge of the graft.
- Place the graft anastomosis close to the rib cage on the first portion of the axillary artery<sup>2</sup>. Do not place the anastomosis on the third portion of the axillary artery.
- Notify the patient that sudden, extreme or strenuous movements of the arm, shoulder or leg should be totally avoided for a period of at least six to eight weeks to allow for proper stabilization of the graft. Routine activities such as raising the arms above the shoulders, reaching out in front, extended reaching, throwing, pulling, striding, or twisting should be avoided. Reference **WARNING #5, #6, and #7**.

**Blood Access Procedures**

Leave the graft in place for approximately two weeks prior to use. There is an increased risk of hematoma formation if the graft is punctured prior to complete healing.

1. Insert the blood access needle at a 20° to 45° angle, with the bevel up. When the graft is penetrated, advance the needle parallel to the graft.
2. ROTATE (CHANGE) THE CANNULATION SITES.

3. Do not repeat cannulation in the same area. Repeat cannulation may lead to formation of a hematoma or a pseudoaneurysm. Do **NOT** cannulate within the dialysis needle's length of the proximal and distal anastomosis.
4. Strictly adhere to aseptic technique to minimize infection.
5. Apply moderate digital pressure to the cannulation site after needle withdrawal. This compression assists in hemostasis.

**Note: There should always be a pulse or thrill of near equal intensity proximal and distal to the area of compression.**

6. Instruct the patient as to proper postoperative care. Reference **WARNING #10** and **#11**.

**References**

1. *Guide to the Safe Handling of Fluoropolymer Resins*, 4th Edition, The Fluoropolymers Division of the Society of The Plastics Industry, Inc.
2. Victor M. Bernhard, M.D. and Jonathan B. Towne, M.D., Editors, *Complications in Vascular Surgery*, Second Edition, Grune and Stratton, Inc. (Harcourt Brace Jovanovich; Publishers), Orlando, 1985, 56.
3. "Perigraft Seromas, Complicating Arterial Grafts", Robert M. Blumenberg, M.D., et al, *Surgery*, Vol. 97, No. 2, February 1985.

**Warranty**

Bard Peripheral Vascular warrants to the first purchaser of this product that this product will be free from defects in materials and workmanship for a period of one year from the date of first purchase and liability under this limited product warranty will be limited to repair or replacement of the defective product, in Bard Peripheral Vascular's sole discretion or refunding your net price paid. Wear and tear from normal use or defects resulting from misuse of this product are not covered by this limited warranty.

**TO THE EXTENT ALLOWABLE BY APPLICABLE LAW, THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL BARD PERIPHERAL VASCULAR BE LIABLE TO YOU FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM YOUR HANDLING OR USE OF THIS PRODUCT.**

Some countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your country.

An issue or revision date and a revision number for these instructions are included for the user's information on the last page of this booklet. In the event 36 months have elapsed between this date and product use, the user should contact Bard to see if additional product information is available.

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**Literaturstellen:**

1. [Guide to the Safe Handling of Fluoropolymer Resins](#), 4th Edition, The Fluoropolymers Division of the Society of The Plastics Industry, Inc.
2. Victor M. Bernhard, M.D. and Jonathan B. Towne, M.D., Editors, [Complications in Vascular Surgery](#), Second Edition, Grune and Stratton, Inc. (Harcourt Brace Jovanovich; Publishers), Orlando, 1985, 56.
3. "Perigraft Seromas, Complicating Arterial Grafts", Robert M. Blumenberg, M.D., et al, *Surgery*, Vol. 97, No. 2, February 1985.

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Referencias

1. Guide to the Safe Handling of Fluoropolymer Resins, 4th Edition, The Fluoropolymers Division of the Society of The Plastics Industry, Inc.
2. Victor M. Bernhard, M.D. and Jonathan B. Towne, M.D., Editors, Complications in Vascular Surgery, Second Edition, Grune and Stratton, Inc. (Harcourt Brace Jovanovich; Publishers), Orlando, 1985, 56.
3. "Perigraft Seromas, Complicating Arterial Grafts", Robert M. Blumenberg, M.D., et al, Surgery, Vol. 97, No. 2, February 1985.

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- Guide to the Safe Handling of Fluoropolymer Resins, 4th Edition, The Fluoropolymers Division of the Society of The Plastics Industry, Inc.
- Victor M. Bernhard, M.D. and Jonathan B. Towne, M.D., Editors, **Complications in Vascular Surgery**, Second Edition, Grune and Stratton, Inc. (Harcourt Brace Jovanovich; Publishers), Orlando, 1985, 56.
- "Perigraft Seromas, Complicating Arterial Grafts", Robert M. Blumenberg, M.D., et al, *Surgery*, Vol. 97, No. 2, February 1985.

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### Referenser

- Guide to the Safe Handling of Fluoropolymer Resins*, 4th Edition, The Fluoropolymers Division of the Society of The Plastics Industry, Inc.
- Victor M. Bernhard, M.D. and Jonathan B. Towne, M.D., Editors, *Complications in Vascular Surgery*, Second Edition, Grune and Stratton, Inc. (Harcourt Brace Jovanovich; Publishers), Orlando, 1985, 56.
- "Perigraft Seromas, Complicating Arterial Grafts", Robert M. Blumenberg, M.D., et al, *Surgery*, Vol. 97, No. 2, February 1985.

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- Guide to the Safe Handling of Fluoropolymer Resins, 4th Edition, The Fluoropolymers Division of the Society of The Plastics Industry, Inc.
- Victor M. Bernhard, M.D. and Jonathan B. Towne, M.D., Editors, *Complications in Vascular Surgery*, Second Edition, Grune and Stratton, Inc. (Harcourt Brace Jovanovich; Publishers), Orlando, 1985, 56.
- "Perigraft Seromas, Complicating Arterial Grafts", Robert M. Blumenberg, M.D., et al, *Surgery*, Vol. 97, No. 2, February 1985.

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**Irodalom**

- Guide to the Safe Handling of Fluoropolymer Resins, 4th Edition, The Fluoropolymers Division of the Society of The Plastics Industry, Inc.
- Victor M. Bernhard, M.D. and Jonathan B. Towne, M.D., Editors, [Complications in Vascular Surgery](#), Second Edition, Grune and Stratton, Inc. (Harcourt Brace Jovanovich; Publishers), Orlando, 1985, 56.
- "Perigraft Seromas, Complicating Arterial Grafts", Robert M. Blumenberg, M.D., et al, Surgery, Vol. 97, No. 2, February 1985.

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**Literatura**

1. **Guide to the Safe Handling of Fluoropolymer Resins**, 4th Edition, The Fluoropolymers Division of the Society of The Plastics Industry, Inc.
2. Victor M. Bernhard, M.D. and Jonathan B. Towne, M.D., Editors, **Complications in Vascular Surgery**, Second Edition, Grune and Stratton, Inc. (Harcourt Brace Jovanovich; Publishers), Orlando, 1985, 56.
3. "Perigraft Seromas, Complicating Arterial Grafts", Robert M. Blumenberg, M.D., et al, *Surgery*, Vol. 97, No. 2, February 1985.

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**Referanslar**

- [Guide to the Safe Handling of Fluoropolymer Resins](#), 4th Edition, The Fluoropolymers Division of the Society of The Plastics Industry, Inc.
- Victor M. Bernhard, M.D. and Jonathan B. Towne, M.D., Editors, [Complications in Vascular Surgery](#), Second Edition, Grune and Stratton, Inc. (Harcourt Brace Jovanovich; Publishers), Orlando, 1985, 56.
- "Perigraft Seromas, Complicating Arterial Grafts", Robert M. Blumenberg, M.D., et al, [Surgery](#), Vol. 97, No. 2, February 1985.

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**U.S. Patent: 6,436,135. Other U.S. and Foreign Patents Pending.****Copyright   2015 C. R. Bard, Inc. All rights reserved. Printed in the U.S.A.**
**Bard, Impra, Impra Carboflo, CenterFlex, EndFlex, Impra Flex and the blue lines design are trademarks and/or registered trademarks of C. R. Bard, Inc. or an affiliate.**
**Bard Peripheral Vascular, Inc.**