# VENAFLO® II

## ePTFE Vascular Graft

# Instructions for Use

**VENAFLO® II ePTFE Vascular Graft**  
**Instructions For Use**

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# INSTRUCTIONS FOR USE

Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

## Device Description, Indications, Contraindications, Warnings, Precautions and Adverse Reactions

### Device Description

VENAFLO® II ePTFE Vascular Grafts are constructed of expanded polytetrafluoroethylene (ePTFE) with carbon impregnated into the inner portions of the graft wall. This graft is also modified with a proprietary cuff with trim lines at the venous end.

Stepped configurations may help minimize the risk of arterial steal syndrome and high cardiac output. CENTERFLEX™ graft configurations have an external spiral bead support on the center of the graft and can be used where resistance to compression or kinking is desired.

### Indications for Use

VENAFLO® II Vascular Grafts in **Straight, Stepped, CENTERFLEX™, and Stepped CENTERFLEX™** configurations are indicated for use as subcutaneous arteriovenous conduits for blood access only.

### Contraindications

None known.

### Warnings

- All VENAFLO® II ePTFE vascular grafts are supplied sterile and non-pyrogenic unless the package is opened or damaged. VENAFLO® II grafts are sterilized by ethylene oxide. Each graft is intended for single patient use only. **DO NOT RESTERILIZE.**
- This device has been designed for single use only. Reusing this medical device bears the risk of cross-patient contamination as medical devices – particularly those with long and small lumina, joints, and/or crevices between components – are difficult or impossible to clean once body fluids or tissues with potential pyrogenic or microbial contamination have had contact with the medical device for an indeterminable period of time. The residue of biological material can promote the contamination of the device with pyrogens or microorganisms which may lead to infectious complications.
- Do not resterilize. After resterilization, the sterility of the product is not guaranteed because of an indeterminable degree of potential pyrogenic or microbial contamination which may lead to infectious complications. Cleaning, reprocessing and/or resterilization of the present medical device increases the probability that the device will malfunction due to potential adverse effects on components that are influenced by thermal and/or mechanical changes.
- Do not use after expiration date printed on the label.
- Do NOT cannulate the externally supported portion of any VENAFLO® II ePTFE vascular graft due to the risk of beading embolization and/or pseudoaneurysm.
- DO NOT REMOVE THE EXTERNAL SPIRAL BEAD SUPPORT FROM ANY CENTERFLEX™ CONFIGURED GRAFT.** Attempts to remove the beading may damage the graft. If damage occurs, discard the graft.
- VENAFLO® II grafts do not stretch (are non-elastic) in the longitudinal direction. Failure to cut the graft to an appropriate length may result in anastomotic or graft disruption, leading to excessive bleeding, loss of limb or limb function, and/or death.
- Aggressive and/or excessive graft manipulation when tunneling, or placement within a too tight or too small tunnel, may lead to separation of the spiral beading and/or graft breakage. **DO NOT** pass the cuff portion (venous end) of the VENAFLO® II graft through a tunneler sheath or the tissue tunnel, as this could lead to separation of the spiral beading and/or graft breakage.
- When embolectomy or balloon angioplasty catheters are used within the lumen of the graft, the inflated balloon size must match the inner diameter of the graft. Over-inflation of the balloon or use of an inappropriately sized balloon may dilate or damage the graft.
- Avoid repeated or excessive clamping at the same location on the graft. If clamping is necessary, use only atraumatic or appropriate vascular smooth jawed clamps to avoid damage to the graft wall. Do not clamp the cuffed portion of the graft.
- Exposure to solutions (e.g., alcohol, oil, aqueous solutions, etc.) may result in loss of the graft's hydrophobic properties. Loss of the hydrophobic barrier may result in graft wall leakage. Preclotting of this graft is unnecessary.<sup>1</sup>
- Avoid excessive graft manipulation after exposure to blood or body fluids. Do not forcibly inject any solution through the lumen of the graft, or fill the graft with fluid prior to pulling it through the tunnel as loss of the graft's hydrophobic properties may occur. Loss of the hydrophobic barrier may result in graft wall leakage.
- Do NOT expose VENAFLO® II grafts to temperatures greater than 500°F (260°C). PTFE decomposes at elevated temperatures, producing highly toxic decomposition products.<sup>2</sup>
- After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable laws and regulations.
- During tunneling, be sure to create a tunnel that closely approximates the outer diameter of the graft. A tunnel that is too loose may result in delayed healing and may also lead to perigraft seroma formation.
- Care should be taken when cannulating the graft for dialysis access; avoid the cuff area and the CENTERFLEX™ external support, properly rotate cannulation sites, and use proper compression to achieve hemostasis during post-cannulation care.

### Precautions

- Only physicians qualified in vascular surgery techniques should implant this prosthesis. The healthcare provider is responsible for all appropriate postoperative care instructions to the patient.
- The healthcare provider must observe aseptic technique during implantation, postoperatively, and during cannulation.
- When suturing, avoid excessive tension on the suture line, inappropriate suture spacing and bites, and gaps between the graft and host vessel. Failure to follow correct suturing techniques may result in suture-hole elongation, suture pull-out, anastomotic bleeding and/or disruption. Refer to "Suturing" for further instructions.
- Do not cannulate the cuffed portion of the graft. For further instructions, refer to the "Blood Access" section.
- Consider intraoperative and postoperative anticoagulation therapy for each patient as appropriate.

### Adverse Reactions

Potential complications that may occur with any surgical procedure involving a vascular prosthesis include, but are not limited to: Disruption or tearing of the suture line, graft, and/or host vessel; suture-hole bleeding; graft redundancy; thrombosis; embolic events; occlusion or stenosis; ultrafiltration; seroma formation; swelling of the implanted limb; formation of hematomas or pseudoaneurysm; infection; arterial steal syndrome; aneurysm/dilation; blood leakage; hemorrhage; and/or skin erosion.

# DIRECTIONS FOR USE

### Equipment Required

Tunneler, double armed suture, atraumatic clamp, scalpel and/or scissors.

### Opening the Package

Hold the outer tray in one hand. Peel back the lid. Remove the inner tray. Peel back the inner tray lid slowly and carefully remove the graft using sterile atraumatic instruments or gloves. Protect the graft against damage from sharp or heavy instruments.

## General Operative Techniques

### Sizing


Use the following guidelines when sizing the cuffed end of all VENAFLO® II ePTFE vascular grafts for the venous anastomosis.

Size of Vein (External Diameter)	Cuff Sizing (See Figure 1)
>8mm	Do not trim
6-8mm	Trim on the line labeled 6-8
3-5mm	Trim on the line labeled 3-5

### Trimming Recommendations

Hold the graft with the cuff in an open configuration, as shown in Figure 1. Using curved scissors, start trimming the cuff on the dashed line located nearest the heel. (See Figure 2.) Continue trimming on the dashed lines following the curvature suggested.

**Note:** Do not fold the graft cuff in half when trimming. Do not trim material from the heel.

### Tunneling Recommendations

Create a tunnel that closely approximates the outer diameter of the graft. Always follow the Instructions for Use for the specific tunneler utilized to place the VENAFLO® II ePTFE vascular graft.

**WARNING:** A tunnel that is too loose may result in delayed healing and may also lead to perigraft seroma formation. See Warning #15.

**Note:** Prior to utilizing a sheath tunneler, verify that the VENAFLO® II graft O.D. dimensions fit the sheath I.D. dimensions.

VENAFLO® II ePTFE Vascular Grafts in **Straight Anatomical Configurations:** After selecting the appropriate tunneler, introduce the tunneler through the arterial site. After creating the tunnel, feed the graft from the venous end to the arterial end. (See Figure 3.) Ensure the cuffed end does not pass through the tunnel. Only the arterial end of the VENAFLO® II ePTFE vascular graft should be sized for length.

Figure 3

**VENAFLO® II ePTFE Vascular Grafts in Looped Anatomical Configurations:** Perform an incision to expose the artery and vein to be anastomosed and a counter-incision at the site where the graft will curve. For CENTERFLEX™ grafts, this is where the flexed portion of the graft should be placed. Position the graft overlying the forearm to optimize graft placement prior to tunneling. Introduce the tunneler at the counter-incision site. After creating the first tunnel, replace the bullet tip (if necessary), attach the arterial end of the graft to the tunnel, and pull the graft through the tunnel leaving the cuffed end exposed near the vein to be anastomosed. For the second tunnel, repeat the tunneling procedure, initiating the tunnel at the anastomotic site and advancing it down to the counter incision. (See Figure 4.) Once again, ensure that the venous (cuffed) graft end does not pass into the tunnel.

Figure 4

Only the arterial end of the VENAFLO® II ePTFE vascular graft should be sized for length.

#### **Suturing**

Size the graft appropriately to minimize excessive tension at the suture line. Use a tapered, non-cutting needle with a non-absorbable monofilament suture approximately the same size as the needle. Take 2 mm suture bites in the graft following the curve of the needle and gently pull the suture at a 90° angle. Proper sizing of the graft length prior to implant will minimize suture-hole elongation caused by excessive tension. See **Precaution #3**.

#### **Suggested Suturing Technique**

After sizing the cuff to the vein diameter, perform a venotomy five to eight millimeters shorter than the length of the cuff. Using double-armed suture, tack the heel of the cuff, aligning the toe in the direction of venous blood flow. Suture continuously from the heel to just beyond the midpoint of the cuff on one side of the venotomy. Using the other arm of the heel suture, continue suturing to just beyond the midpoint of the graft on the other side of the venotomy. Extend the venotomy to the exact cuff length. (See Figure 5.)

Figure 5

Complete the "cuffed" anastomosis by suturing with one of the suture arms, around the toe, to the other side of the cuff. Always tie off on the side of the cuff (See Figure 6).

Figure 6

Perform the arterial anastomosis using standard surgical techniques. In order to retain a cuffed configuration, only the arterial end of the VENAFLO® II ePTFE vascular graft should be sized for length.

#### **Thrombectomy**

Techniques for declogging VENAFLO® II ePTFE vascular grafts include but are not limited to the use of balloon catheters. See **Warning #9**.

#### **Longitudinal Incision:**

Place stay sutures before introducing the embolectomy catheter. (If the graft is a CENTERFLEX™ product, cut through the spiral beading and base tube. After closure, the spiral beading will realign itself.) Place a longitudinal incision in the graft that is long enough to accommodate the extraction of a fully dilated thrombectomy catheter balloon. A patch may be considered as an aid to graft closure.

#### **Transverse Incision:**

No stay sutures are necessary. A horizontal mattress suture is recommended for graft closure. During the early postoperative period, the natural progression of wound healing renders the graft translucent in appearance. In this state, a longitudinal incision with stay sutures is recommended. If a transverse incision is performed, a horizontal mattress suture technique and PTFE pledgets may aid in closure.

#### **Angiography**

Should angiography be performed at the time of the procedure, the artery proximal to the graft should be used for injection, if possible.

#### **Blood Access**

Best results are achieved by:

- 1 Leaving the graft in place for approximately two weeks prior to use. Immediate use may increase the risk of hematoma formation.

2. Inserting blood access needle at a 20° to 45° angle, with the bevel up until the graft is penetrated, and then advancing the needle parallel to the graft. Routine flipping of needles is not recommended.
3. ROTATING THE CANNULATION SITES. Repeated cannulation in the same area may lead to formation of hematoma or pseudoaneurysm.
4. Cannulating greater than a dialysis needle's length of the proximal anastomosis. See **Precaution #4**.
5. Strictly adhering to aseptic technique to minimize infection.
6. Applying moderate digital pressure to the cannulation site after needle withdrawal. This compression aids in hemostasis.

**Note: There should always be a pulse or thrill near equal intensity proximal and distal to the area of compression. Ensure that hemostasis is maintained with the least amount of pressure required.**

#### **References**

1. "Perigraft Seromas, Complicating Arterial Grafts", Robert M. Blumenberg, M.D., et al, Surgery, Vol. 97, No. 2, February 1985.
2. Guide to the Safe Handling of Fluoropolymer Resins, 4th Edition, The Fluoropolymers Division of the Society of the Plastics Industry, Inc.

#### **Warranty**

Bard Peripheral Vascular warrants to the first purchaser of this product that this product will be free from defects in materials and workmanship for a period of one year from the date of first purchase and liability under this limited product warranty will be limited to repair or replacement of the defective product, in Bard Peripheral Vascular's sole discretion or refunding your net price paid. Wear and tear from normal use or defects resulting from misuse of this product are not covered by this limited warranty.

**TO THE EXTENT ALLOWABLE BY APPLICABLE LAW, THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL BARD PERIPHERAL VASCULAR BE LIABLE TO YOU FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM YOUR HANDLING OR USE OF THIS PRODUCT.**

Some countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your country.

An issue or revision date and a revision number for these instructions are included for the user's information on the last page of this booklet. In the event 36 months have elapsed between this date and product use, the user should contact Bard Peripheral Vascular to see if additional product information is available.

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1. "Perigraft Seromas, Complicating Arterial Grafts", Robert M. Blumenberg, M.D., et al, Surgery, Vol. 97, No. 2, February 1985.
2. **Guide to the Safe Handling of Fluoropolymer Resins**, 4th Edition, The Fluoropolymers Division of the Society of the Plastics Industry, Inc.

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**EC REP****Authorised Representative in the European Community**

Représentant autorisé au sein de la Communauté européenne  
 Bevollmächtigter in der Europäischen Gemeinschaft  
 Rappresentante autorizzato nella Comunità Europea  
 Representante autorizado en la Comunidad Europea  
 Gemachtigde binnen de Europese Gemeenschap  
 Representante autorizado na Comunidade Europeia  
 Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα  
 Autoriseret repræsentant i Det Europæiske Fællesskab  
 Auktoniserad representant inom EG  
 Valtuutettu edustaja Euroopan yhteisössä  
 Autoriseret representant i EU  
 Autoryzowany przedstawiciel na terenie Unii Europejskiej  
 Meghatalmazott képviselő az Európai Közösségben  
 Autorizovaný zástupce v Evropské unii  
 Avrupa Topluluğu Yetkili Temsilcisi  
 Уполномоченный представитель в Европейском сообществе  
 歐共同体授權代表  
 유럽 공동체 내의 인증받은 대리업체

**Do Not Use If Package Is Damaged**

Ne pas utiliser si l'emballage est endommagé  
 Nicht verwenden, wenn die Verpackung beschädigt ist  
 Non utilizzare se la confezione è danneggiata  
 Non usar si el envase está dañado  
 Niet gebruiken wanneer de verpakking beschadigd is  
 Não utilizar se a embalagem estiver danificada  
 Μη χρησιμοποιείτε το προϊόν εάν η συσκευασία του έχει υποστεί ζημιά  
 Må ikke bruges, hvis emballagen er beskadiget  
 Använd inte produkten om förpackningen skadats  
 Ei saa käyttää, jos pakkaus on vaurioitunut  
 Skal ikke brukes hvis pakningen er skadet  
 Nie stosować, jeśli opakowanie jest uszkodzone  
 Ne használja fel, ha a csomagolás sérült  
 Nepoužívejte, pokud je obal poškozen  
 Paket Hasarlıysa Kullanmayınız  
 Не использовать, если упаковка повреждена  
 如果包裝已受損，請勿使用  
 포장에 손상되어 있다면 사용하지 마십시오

**Non-Pyrogenic**

Apyrogène  
Pyrogenfrei  
Apirogeno  
Apirógeno  
Niet-pyrogeen  
Apirogénico  
Μη πυρογόνο  
Pyrogenfri  
Pyrogenfri  
Pyrogeenitön  
Pyrogenfri  
Apirogenny  
Pirogèmentes  
Apyrogenni  
Pirojenik deǵildir  
Аπιρογενно  
無熱原  
비발열성

**Not Made with Natural Rubber Latex**

Fabriqué sans latex naturel  
Nicht aus Naturgummlatex hergestellt.  
Non prodotto con lattice di gomma naturale  
Este producto no se fabrica con látex de caucho natural  
Niet gemaakt met latex van natuurlijk rubber  
Não contém látex de borracha natural  
Δεν κατασκευάζεται από φυσικό ελαστικό λάτεξ  
Ikke fremstillet med naturlig gummlatex  
Ej tillverkad med naturgummlatex  
valmistuksessa ei ole käytetty luonnonkumilateksia  
Dette produktet er ikke produsert med naturlig gummlatex  
Produkt nie jest wytwarzany z lateksu kauczuku naturalnego  
Nem tartalmaz természetes gummlatexet.  
Není vyroben z přírodního prýžového latexu  
Doğal Kauçuk Lateksten Üretilmemiştir  
Не содержит натуральный каучуковый латекс  
製造未採用天然膠乳  
천연 고무 라텍스를 사용해 제조되지 않음

**U.S. Patent: 6,436,135. Other U.S. and Foreign Patents Pending. Copyright © 2015 C. R. Bard, Inc. All rights reserved. Printed in the U.S.A.**

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**Bard, Venaflō, CenterFlex and the blue lines design are trademarks and/or registered trademarks of C. R. Bard, Inc. or an affiliate.**

Bard, Venaflō, CenterFlex et la conception des lignes bleues sont des marques déposées de C. R. Bard, Inc., ou d'une filiale.  
Bard, Venaflō, CenterFlex und das blaue Linien-Design sind eingetragene Warenzeichen von C. R. Bard, Inc., oder einer Tochtergesellschaft.  
Bard, Venaflō, CenterFlex e il design blue lines sono marchi registrati della C. R. Bard, Inc. o di un'azienda affiliata.  
Bard, Venaflō, CenterFlex y el diseño con líneas azules son marcas comerciales registradas de C. R. Bard, Inc. o una filial.  
Bard, Venaflō, CenterFlex er het ontwerp met blauwe lijnen zijn gedeponeerde handelsmerken van C. R. Bard, Inc. of een geaffilieerde onderneming.  
Bard, Venaflō, CenterFlex e o desenho das linhas azuis são marcas registradas da C. R. Bard, Inc., ou de uma empresa afiliada.  
Οι ονομασίες Bard, Venaflō, CenterFlex και το σχέδιο με τις μπλε γραμμές είναι εμπορικά σήματα ή και σήματα κατατεθέντα της C. R. Bard, Inc. ή θυγατρικής της.  
Bard, Venaflō, CenterFlex og den blå linie designen er registrerede varemærker tilhørende C. R. Bard eller en filial.  
Bard, Venaflō, CenterFlex och utformningen med blå linjer är inregistrerade varumärken som ägs av C. R. Bard, Inc. eller ett med dem associerat företag.  
Bard, Venaflō, CenterFlex ja sinisten viivojen kuvio ovat C. R. Bard, Inc:n tai sen tytäryhtiön rekisteröityjä tavaramerkkejä, tytäryhtiön rekisteröityjä tavaramerkkejä.  
Bard, Venaflō, CenterFlex og det blå strekmønsteret er varemærker og/eller registrerte varemærker som tilhører C.R. Bard, Inc., eller et tilknyttet selskap.  
Bard, Venaflō, CenterFlex i logo z niebieskimi liniami są zarejestrowanymi znakami towarowymi firmy C. R. Bard, Inc. lub jej spółek stowarzyszonych.  
A Bard, a Venaflō, a CenterFlex és a blue lines design a C.R. Bard, Inc., vagy valamelyik fiókállalatának bejegyzett védjegyei.  
Bard, Venaflō, CenterFlex a the blue lines design jsou registrovanými ochrannými známkami společnosti C.R. Bard, Inc. nebo její pobočky.  
Bard, Venaflō, CenterFlex ve blue lines design (mavi çizgi tasarımı), C.R. Bard, Inc. veya bağlı kuruluşunun tescilli ticari markalarıdır.  
Bard, Venaflō, CenterFlex и дизайн в виде голубых линий являются товарными знаками и/или зарегистрированными товарными знаками компании C. R. Bard, Inc. или ее дочерних компаний.  
Bard · Venaflō · CenterFlex 和藍線設計是 C. R. Bard, Inc. 或其附屬機構的商標和/或註冊商標。  
Bard, Venaflō, CenterFlex 및 파란선 디자인은 C. R. Bard, Inc. 또는 계열사의 상표 및/또는 등록상표입니다.

**Manufacturer:**

**Bard Peripheral Vascular, Inc.**

1625 West 3rd Street  
Tempe, AZ 85281  
USA

TEL: 1-480-894-9515  
1-800-321-4254  
FAX: 1-480-966-7062  
1-800-440-5376  
[www.bardpv.com](http://www.bardpv.com)

**Authorized Representative  
in the European Community  
BD Switzerland Sàrl**  
Terre Bonne Park – A4  
Route de Crassier 17  
1262 Eysins – Switzerland

**BARD** | PERIPHERAL  
VASCULAR

PK1139600

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Copyright © 2015 C. R. Bard, Inc. Tous droits réservés. Imprimé aux États-Unis.
Copyright © 2015 C. R. Bard, Inc. Alle Rechte vorbehalten. Gedruckt in den USA.
Copyright © 2015 C. R. Bard, Inc. Tutti i diritti riservati. Stampato negli Stati Uniti.
Copyright © 2015, C. R. Bard, Inc. Reservados todos los derechos. Impreso en EE.UU.
Copyright © 2015 C. R. Bard, Inc. Alle rechten voorbehouden. Gedrukt in de VS.
Copyright © 2015, C. R. Bard, Inc. Todos os direitos reservados. Impresso nos E.U.A.