Instructions for Use
MYNX CONTROL™ Vascular Closure Device

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ENGLISH

Instructions for Use

MX5060 MYNX CONTROL™ Vascular Closure Device 5F
MX6760 MYNX CONTROL™ Vascular Closure Device 6F/7F

To ensure proper deployment and use of this device and to prevent injury to patients, read all information contained in these instructions for use.

DEVICE DESCRIPTION

The MYNX CONTROL Vascular Closure Device (VCD) is designed to achieve femoral artery hemostasis via delivery of the GRIP TECHNOLOGY™ sealant, an extravascular, water-soluble synthetic hydrogel, using a balloon catheter in conjunction with a standard procedural sheath. The GRIP TECHNOLOGY™ sealant is made of a polyethylene glycol (PEG) material which expands upon contact with subcutaneous fluids to seal the arteriotomy. The sealant is resorbed by the body within 30 days.

The MYNX CONTROL VCD is supplied with a 10 ml locking syringe used for balloon inflation and deflation.

The device contains no components manufactured from latex rubber. The 5F MYNX CONTROL VCD device has a gray sheath catch. The 6F/7F MYNX CONTROL VCD device has a green sheath catch. Refer to Figure 1 in the Procedure and Device Preparation Section for the MYNX CONTROL VCD components.

INDICATIONS FOR USE

The MYNX CONTROL VCD is indicated for use to seal femoral arterial access sites while reducing times to hemostasis and ambulation in patients who have undergone diagnostic or interventional endovascular procedures utilizing a 5F, 6F or 7F procedural sheath.

CONTRAINDICATIONS

There are no known contraindications for the MYNX CONTROL VCD.

WARNINGS

Do not use if components or packaging appear to be damaged or defective or if any portion of the packaging has been previously opened.

DO NOT REUSE OR RESTERILIZE. The MYNX CONTROL VCD is for single use only. The catheter is loaded with a single Hydrogel sealant. Reuse of the device would result in no delivery of Hydrogel sealant.

Do not use the MYNX CONTROL VCD if the puncture site is located above the most inferior border of the inferior epigastric artery (IEA) and/or above the inguinal ligament based upon bony landmarks, since such a puncture site may result in a retroperitoneal hematoma/bleed. Perform a femoral angiogram to verify the location of the puncture site.

Do not use the MYNX CONTROL VCD if the puncture is through the posterior wall or if there are multiple punctures, as such punctures may result in a retroperitoneal hematoma/bleed.

PRECAUTIONS

The MYNX CONTROL VCD should only be used by a trained licensed physician or healthcare professional.

The MYNX CONTROL VCD should not be used in patients with a known allergy to PEG.

Cardinal Health
5452 Betsy Ross Drive
Santa Clara, CA 95054 USA
Phone: +1 650 864 5473
Fax: +1 650 903 1018

© 2018 Cardinal Health. All Rights Reserved. CARDINAL HEALTH, the Cardinal Health LOGO, ESSENTIAL TO CARE, MYNX, MYNX CONTROL, and GRIP TECHNOLOGY are trademarks of Cardinal Health and may be registered in the US and/or in other countries. All other marks are the property of their respective owners. Patent www.cardinalhealth.com/patents.

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The MYNX CONTROL VCD should not be used with sheaths longer than 12 cm effective length or incompatible sheaths listed in Table 9.

SPECIAL PATIENT POPULATIONS

The safety and effectiveness of the MYNX CONTROL VCD have not been established in the following patient populations:

- ◆ Pediatric patients or others with small common femoral arteries (< 5 mm in diameter)
- ◆ Patients with clinically significant peripheral vascular disease in the vicinity of the puncture
- ◆ Patients with prior surgical procedure, PTA, stent placement, or vascular graft in the common femoral artery
- ◆ Patients with bleeding disorders such as thrombocytopenia (platelet count < 100,000/mm³), hemophilia, von Willebrand's disease or anemia (Hgb < 10 g/dL, Hct < 30%)
- ◆ Patients with uncontrolled hypertension (systolic BP > 180 mm Hg)
- ◆ Patients with morbid obesity (BMI > 40 kg/m²)
- ◆ Patients who are pregnant or lactating
- ◆ Patients with documented INR > 1.5 or patients currently receiving glycoprotein IIb/IIIa platelet inhibitors

ADVERSE EVENTS

The MYNX® VCD was evaluated in a prospective multi-center, non-randomized clinical trial involving 190 patients to achieve femoral arterial access site hemostasis following diagnostic angiography (n=95) or interventional procedures (n=95). The control arm (standard compression) from the MATRIX VSG™ System Clinical Trial (MATRIX) was used as a historical control group. Table 1 is a report of the total number and the rate of major and minor complications in the MYNX VCD Trial and the standard compression arm of the MATRIX Trial. Clinical data collected on the MYNX VCD is directly applicable to the MYNX CONTROL VCD.

THIS PAGE INTENTIONALLY LEFT BLANK

**Table 1: Reported Major and Minor Complications (Intent-to-Treat)
 All Patients**

Major Complications by Event All Patients	MYNX* (n=190)	Standard Compression (Historical Control) (n=164)	p-value [†]
Vascular Repair	0.0% (0/190)	0.0% (0/164)	N/A
Permanent Access Site-Related Nerve Injury	0.0% (0/190)	0.0% (0/164)	N/A
Surgery For Access Site-Related Nerve Injury	0.0% (0/190)	0.0% (0/164)	N/A
Access Site-Related Bleeding Requiring Transfusion	0.5% (1/190)	0.0% (0/164)	0.0021
New Ipsilateral Lower Extremity Ischemia Requiring Invasive/Non Invasive Intervention	0.0% (0/190)	0.0% (0/164)	N/A
Access Site-Related Infection – Major	0.0% (0/190)	0.0% (0/164)	N/A
Local Access Site Inflammatory Reaction – Major	0.0% (0/190)	0.0% (0/164)	N/A
Generalized Infection	0.0% (0/190)	0.0% (0/164)	N/A
Any Major Complication	0.5% (1/190)	0.0% (0/164)	0.0021
Minor Complications by Event All Patients	MYNX* (n=190)	Standard Compression (Historical Control) (n=164)	p-value [†]
Pseudoaneurysm – Treated With Thrombin Injection	0.5% (1/190)	0.0% (0/164)	0.0021
Pseudoaneurysm – Not Requiring Treatment [‡]	2.6% (5/190)	0.0% (0/164)	0.0831
AV Fistula	0.0% (0/190)*	0.0% (0/164)	N/A
Hematoma ≥ 6 cm	3.2% (6/190)*§	0.6% (1/164)	0.0853
Access Site-Related Bleeding Requiring > 30 min to Achieve Hemostasis	0.0% (0/190)	0.6% (1/164)	0.0002
Late Access Site-Related Bleeding (Following Hospital Discharge)	0.0% (0/190)	0.0% (0/164)	N/A
Ipsilateral Lower Extremity Arterial Emboli	0.0% (0/190)	0.0% (0/164)	N/A
Transient Loss of Ipsilateral Lower Extremity Pulse	0.0% (0/190)	0.0% (0/164)	N/A
Ipsilateral Deep Vein Thrombosis	0.0% (0/190)*	0.0% (0/164)	N/A
Transient Access Site-Related Nerve Injury	0.0% (0/190)	0.0% (0/164)	N/A
Access Site-Related Vessel Laceration	0.0% (0/190)	0.0% (0/164)	N/A
Access Site Wound Dehiscence	0.0% (0/190)	0.0% (0/164)	N/A
Local Access Site Infection – Minor	0.0% (0/190)	0.0% (0/164)	N/A
Local Access Site Inflammatory Reaction – Minor	0.0% (0/190)	0.0% (0/164)	N/A
Any Minor Complication	3.7% (7/190)	1.2% (2/164)	0.0921

* Does not include four (4) non-device related events consisting of: AV Fistula (n=2), Ipsilateral Deep Vein Thrombosis (n=1), and Hematoma (n=1).

§ Does not include 1 patient with a pre-existing hematoma prior to deployment of MYNX.

‡ In a subset of MYNX patients, ultrasound assessment was performed prior to patient discharge which revealed some patients had developed small pseudoaneurysms which resolved spontaneously. The events are reported in the table above, but are not included in the overall rate of minor complications.

† p-value based on test for non-inferiority.

THIS PAGE INTENTIONALLY LEFT BLANK

Table 2: Reported Major and Minor Complications (Intent-to-Treat) Diagnostic Patients

Major Complications by Event Diagnostic Patients	MYNX* Diagnostic (n=95)	Standard Compression (Historical Control) Diagnostic (n=83)	p-value †
Vascular Repair	0.0% (0/95)	0.0% (0/83)	N/A
Permanent Access Site-Related Nerve Injury	0.0% (0/95)	0.0% (0/83)	N/A
Surgery For Access Site-Related Nerve Injury	0.0% (0/95)	0.0% (0/83)	N/A
Access Site-Related Bleeding Requiring Transfusion	0.0% (0/95)	0.0% (0/83)	N/A
New Ipsilateral Lower Extremity Ischemia Requiring Invasive/Non Invasive Intervention	0.0% (0/95)	0.0% (0/83)	N/A
Access Site-Related Infection – Major	0.0% (0/95)	0.0% (0/83)	N/A
Local Access Site Inflammatory Reaction – Major	0.0% (0/95)	0.0% (0/83)	N/A
Generalized Infection	0.0% (0/95)	0.0% (0/83)	N/A
Any Major Complication	0.0% (0/95)	0.0% (0/83)	N/A
Minor Complications by Event Diagnostic Patients	MYNX* Diagnostic (n=95)	Standard Compression (Historical Control) Diagnostic (n=83)	p-value †
Pseudoaneurysm – Treated With Thrombin Injection	1.1% (1/95)	0.0% (0/83)	0.0459
Pseudoaneurysm – Not Requiring Treatment‡	1.1% (1/95)	0.0% (0/83)	0.0459
AV Fistula	0.0% (0/95)*	0.0% (0/83)	N/A
Hematoma ≥ 6 cm	2.1% (2/95)*§	1.2% (1/83)	0.0599
Access Site-Related Bleeding Requiring > 30 min to Achieve Hemostasis	0.0% (0/95)	0.0% (0/83)	N/A
Late Access Site-Related Bleeding (Following Hospital Discharge)	0.0% (0/95)	0.0% (0/83)	N/A
Ipsilateral Lower Extremity Arterial Emboli	0.0% (0/95)	0.0% (0/83)	N/A
Transient Loss of Ipsilateral Lower Extremity Pulse	0.0% (0/95)	0.0% (0/83)	N/A
Ipsilateral Deep Vein Thrombosis	0.0% (0/95)	0.0% (0/83)	N/A
Transient Access Site-Related Nerve Injury	0.0% (0/95)	0.0% (0/83)	N/A
Access Site-Related Vessel Laceration	0.0% (0/95)	0.0% (0/83)	N/A
Access Site Wound Dehiscence	0.0% (0/95)	0.0% (0/83)	N/A
Local Access Site Infection – Minor	0.0% (0/95)	0.0% (0/83)	N/A
Local Access Site Inflammatory Reaction – Minor	0.0% (0/95)	0.0% (0/83)	N/A
Any Minor Complication	3.2% (3/95)	1.2% (1/83)	0.1458

* Excludes non-device related events consisting of: AV Fistula (n=2) and Hematoma (n=1).

§ Excludes 1 patient with pre-existing hematoma prior to deployment of MYNX.

‡ In a subset of MYNX patients, ultrasound assessment was performed prior to patient discharge which revealed some patients had developed small pseudoaneurysms which resolved spontaneously. The events are reported in the table above, but are not included in the overall rate of minor complications.

† p-value based on test for non-inferiority.

THIS PAGE INTENTIONALLY LEFT BLANK

Table 3: Reported Major and Minor Complications (Intent-to-Treat) Interventional Patients

Major Complications by Event Interventional Patients	MYNX* Interventional (n=95)	Standard Compression (Historical Control) Interventional (n=81)	p-value†
Vascular Repair	0.0% (0/95)	0.0% (0/81)	N/A
Permanent Access Site-Related Nerve Injury	0.0% (0/95)	0.0% (0/81)	N/A
Surgery For Access Site-Related Nerve Injury	0.0% (0/95)	0.0% (0/81)	N/A
Access Site-Related Bleeding Requiring Transfusion	1.1% (1/95)	0.0% (0/81)	0.0459
New Ipsilateral Lower Extremity Ischemia Requiring Invasive/Non Invasive Intervention	0.0% (0/95)	0.0% (0/81)	N/A
Access Site-Related Infection – Major	0.0% (0/95)	0.0% (0/81)	N/A
Local Access Site Inflammatory Reaction – Major	0.0% (0/95)	0.0% (0/81)	N/A
Generalized Infection	0.0% (0/95)	0.0% (0/81)	N/A
Any Major Complication	1.1% (1/95)	0.0% (0/81)	0.0459
Minor Complications by Event Interventional Patients	MYNX* Interventional (n=95)	Standard Compression (Historical Control) Interventional (n=81)	p-value†
Pseudoaneurysm – Treated With Thrombin Injection	0.0% (0/95)	0.0% (0/81)	N/A
Pseudoaneurysm – Not Requiring Treatment‡	4.2% (4/95)	0.0% (0/81)	0.4819
AV Fistula	0.0% (0/95)	0.0% (0/81)	N/A
Hematoma ≥ 6 cm	4.2% (4/95)	0.0% (0/81)	0.4819
Access Site-Related Bleeding Requiring > 30 min to Achieve Hemostasis	0.0% (0/95)	1.2% (1/81)	0.0060
Late Access Site-Related Bleeding (Following Hospital Discharge)	0.0% (0/95)	0.0% (0/81)	N/A
Ipsilateral Lower Extremity Arterial Emboli	0.0% (0/95)	0.0% (0/81)	N/A
Transient Loss of Ipsilateral Lower Extremity Pulse	0.0% (0/95)	0.0% (0/81)	N/A
Ipsilateral Deep Vein Thrombosis	0.0% (0/95)*	0.0% (0/81)	N/A
Transient Access Site-Related Nerve Injury	0.0% (0/95)	0.0% (0/81)	N/A
Access Site-Related Vessel Laceration	0.0% (0/95)	0.0% (0/81)	N/A
Access Site Wound Dehiscence	0.0% (0/95)	0.0% (0/81)	N/A
Local Access Site Infection – Minor	0.0% (0/95)	0.0% (0/81)	N/A
Local Access Site Inflammatory Reaction – Minor	0.0% (0/95)	0.0% (0/81)	N/A
Any Minor Complication	4.2% (4/95)	1.2% (1/81)	0.2950

* Excludes non-device related event consisting of: Ipsilateral Deep Vein Thrombosis (n=1).

‡ In a subset of MYNX patients, ultrasound assessment was performed prior to patient discharge which revealed some patients had developed small pseudoaneurysms which resolved spontaneously. The events are reported in the table above, but are not included in the overall rate of minor complications.

† p-value based on test for non-inferiority.

THIS PAGE INTENTIONALLY LEFT BLANK

THIS PAGE INTENTIONALLY LEFT BLANK

POTENTIAL ADVERSE EVENTS

In addition to the complications noted in the MYNX VCD clinical trial, the following potential complications, which may be related to the endovascular procedure or the vascular closure, may occur: allergic reaction, ecchymosis, foreign body/local reaction, retroperitoneal bleed, vessel occlusion, pulmonary embolism or death.

CLINICAL TRIAL

The MYNX VCD was evaluated in a prospective, multi-center, non-randomized clinical trial designed to evaluate its safety and effectiveness in sealing femoral arterial access sites in patients following diagnostic or interventional catheterization procedures. Clinical data collected on the MYNX VCD is directly applicable to the MYNX CONTROL VCD. The study was conducted in Germany at five institutions involving 190 patients. Patients eligible for participation included candidates for early ambulation and patients who were clinically indicated for a diagnostic or interventional procedure involving access through the femoral artery using a 5F, 6F or 7F sheath with an overall sheath length ≤ 15.7 cm.

The objective for the safety endpoints of the study was to demonstrate non-inferiority to a historical control group, and the objective for the effectiveness endpoints of the study was to demonstrate superiority to the historical group. The historical control group consisted of patients treated with standard compression who were enrolled in the control arm of the MATRIX Trial conducted under IDE# G030182.

Patients were required to be at least 18 years of age, to have signed an Informed Consent Form, and to have undergone a diagnostic or interventional procedure through the common femoral artery. Patients were excluded if they presented with clinically significant peripheral vascular disease, prior procedure in the ipsilateral common femoral artery ≤ 30 days before the MYNX study catheterization procedure, ipsilateral closure device, known allergy to contrast medium or device materials, a myocardial infarction with elevated ST segment ≤ 24 hours prior to procedure, uncontrolled hypertension, existing bleeding disorder, evidence of infection or local inflammation, chronic corticosteroid therapy ≥ 1 month duration, common femoral artery diameter < 5 mm, pre-existing bleeding around the procedural sheath, preexisting hematoma, intraluminal thrombus, pseudoaneurysm, AV fistula, any type of dissection, fibrotic, calcified, or $> 50\%$ stenotic femoral artery, arterial puncture outside the common femoral artery, ipsilateral venous sheath, multiple arterial sticks, suspected posterior femoral arterial wall puncture, antegrade puncture, intra-aortic balloon pump, or planned extended hospitalization.

The majority of the patients were male ($n=133$) vs. female ($n=57$) with all patients' ages ranging from 40 to 85 years. Of the 190 patients enrolled, 50% were diagnostic patients and 50% were interventional patients. The average activated clotting time (ACT) for interventional patients was 223 ± 84 seconds with a range of 108-634 seconds, and for diagnostic patients, 187 ± 81 seconds with a range of 133-280 seconds.

EFFECTIVENESS RESULTS

The study objective for the primary effectiveness endpoints was to demonstrate that the MYNX VCD patient results were superior to the historical control arm (standard compression patients from the MATRIX Trial). The results of the effectiveness measures are summarized in Table 4 for all patients and in Table 5 for diagnostic and interventional patients.

**Table 4: Effectiveness and Secondary Endpoints
 All Patients**

Effectiveness All Patients	MYNX	Standard Compression (Historical Control)	p-value [†]
Time to Hemostasis (minutes)*			
mean ± standard deviation (n)	1.3 ± 2.3 (183)	25.4 ± 16.2 (161)	< 0.0001
median (Q1, Q3) (n)	0.5 (0.0, 2.0) (183)	20.0 (15.0, 30.0) (161)	< 0.0001
range (min, max)	(0.0, 22.5)	(6.0, 120.0)	N/A
Time to Ambulation (hours)*			
mean ± standard deviation (n)	2.6 ± 2.6 (181)	7.4 ± 4.8 (160)	< 0.0001
median (Q1, Q3) (n)	2.0 (1.8, 2.2) (181)	6.0 (4.5, 7.4) (160)	< 0.0001
range (min, max)	(1.3, 20.0)	(1.6, 26.9)	N/A
Secondary Endpoints	MYNX	Standard Compression (Historical Control)	p-value[†]
Device Success	93.2% (177/190)	N/A	N/A
Procedure Success	99.5% (189/190)	100.0% (164/164)	1.0000

* The number of patients used to calculate effectiveness measures differ from overall study sample size due to missing values.

† p-value based on test for superiority.

Note: Times to ambulation varied across investigational sites due to individual hospital practices and protocol.

THIS PAGE INTENTIONALLY LEFT BLANK

THIS PAGE INTENTIONALLY LEFT BLANK

Table 5: Effectiveness and Secondary Endpoints Diagnostic and Interventional Patients

Effectiveness	MYNX Diagnostic	Standard Compression (Historical Control) Diagnostic	p-value†	MYNX Interventional	Standard Compression (Historical Control) Interventional	p-value†
Time to Hemostasis (minutes)*						
mean ± standard deviation (n)	1.0 ± 1.3 (92)	23.6 ± 17.1 (83)	< 0.0001	1.5 ± 2.9 (91)	27.3 ± 15.2 (78)	< 0.0001
median (Q1, Q3) (n)	0.5 (0.0, 1.6) (92)	19.0 (14.0, 25.0) (83)	< 0.0001	0.6 (0.0, 2.0) (91)	25.0 (19.0, 30.0) (78)	< 0.0001
range (min, max)	(0.0, 6.0)	(6.0, 120.0)	N/A	(0.0, 22.5)	(10.0, 120.0)	N/A
Time to Ambulation (hours)*						
mean ± standard deviation (n)	2.5 ± 2.1 (91)	5.4 ± 2.7 (82)	< 0.0001	2.8 ± 3.0 (90)	9.4 ± 5.6 (78)	< 0.0001
median (Q1, Q3) (n)	2.0 (1.9, 2.3) (91)	5.2 (4.3, 6.1) (82)	< 0.0001	1.9 (1.8, 2.2) (90)	7.1 (5.5, 11.8) (78)	< 0.0001
range (min, max)	(1.4, 19.6)	(1.6, 26.8)	N/A	(1.3, 20.0)	(2.5, 22.3)	N/A
Secondary Endpoints	MYNX Diagnostic	Standard Compression (Historical Control) Diagnostic	p-value†	MYNX Interventional	Standard Compression (Historical Control) Interventional	p-value†
Device Success	93.7% (89/95)	N/A	N/A	92.6% (88/95)	N/A	N/A
Procedure Success	100% (95/95)	100% (83/83)	N/A	98.9% (94/95)	100% (81/81)	1.0000

* The number of patients used to calculate effectiveness measures differ from overall study sample size due to missing values.

† p-value based on test for superiority.

Note: Times to ambulation varied across investigational sites due to individual hospital practices and protocol.

Procedure success was defined as successfully achieving hemostasis using any method with freedom from major complications. Device success was defined as the ability to deploy the delivery system, deliver the sealant, and achieve hemostasis with MYNX VCD at the femoral artery puncture site. The procedure success rate was 99.5% demonstrating successful hemostasis in all patients with only one major complication reported. The device success rate for MYNX VCD was 93.2%.

Time to discharge was a secondary endpoint for the diagnostic patients only. The time to discharge was defined as the time from advancer tube removal to the time a patient was discharged from the hospital. The mean time to discharge was 32.3 ± 55.6 hours for the MYNX VCD group compared to 20.1 ± 36.1 hours for the historical control standard compression group (p=0.02). In the subset of diagnostic patients, the mean time to discharge for MYNX VCD patients was 35.3 ± 71.4 hours compared to 17.3 ± 47.6 hours for the historical control standard compression group (p=0.06). The mean time to discharge for MYNX VCD interventional patients was 29.0 ± 31.1 hours compared to 23.0 ± 16.9 hours for the interventional historical control standard compression group (p=0.15). The differences in time to discharge may be due to the differences in standard of care for patients treated in Europe (MYNX VCD patients) compared to patients treated in the United States (standard compression historical control group).

**Table 6:
 Cumulative Time to Hemostasis
 All Patients***

Time to Hemostasis (minutes)	MYNX (n=183)	Standard Compression (Historical Control) (n=161)
1	59.0% (108)	0.0% (0)
2	83.6% (153)	0.0% (0)
3	91.8% (168)	0.0% (0)
4	94.5% (173)	0.0% (0)
5	96.2% (176)	0.0% (0)
10	98.9% (181)	5.0% (8)
15	99.4% (182)	26.7% (43)
20	99.4% (182)	53.4% (86)
25	100% (183)	64.6% (104)
> 30	100% (183)	100% (161)

**Table 7:
 Cumulative Time to Ambulation
 All Patients***

Time to Ambulation (hours)	MYNX (n=181)	Standard Compression (Historical Control) (n=160)
2	50.8% (92)	1.3% (2)
3	90.0% (163)	3.8% (6)
4	92.3% (167)	12.5% (20)
5	93.9% (170)	35.6% (57)
10	97.2% (176)	83.8% (134)
15	98.3% (178)	89.4% (143)
20	100% (181)	95.6% (153)
25	100% (181)	99.4% (159)
> 30	100% (181)	100% (160)

THIS PAGE INTENTIONALLY LEFT BLANK

*The number of patients differs from overall study sample size due to missing values

**Table 8: Cumulative Time to Discharge
 Diagnostic Patients***

Time to Discharge (hours)	MYNX Diagnostic (n=79)	Standard Compression Diagnostic (n=82)
2	0.0% (0)	0.0% (0)
3	2.5% (2)	1.2% (1)
4	16.5% (13)	4.9% (4)
5	19.0% (15)	18.3% (15)
10	25.3% (20)	81.7% (67)
15	29.1% (23)	82.9% (68)
20	38.0% (30)	84.1% (69)
25	67.1% (53)	89.0% (73)
30	82.3% (65)	91.5% (75)
> 30	100% (79)	100% (82)

*The number of patients differs from overall study sample size due to missing values

THIS PAGE INTENTIONALLY LEFT BLANK

Overall, the primary effectiveness endpoints were successfully met in the MYNX study. Patients treated with the MYNX VCD had reduced times to hemostasis and ambulation compared to patients treated with standard compression. In addition, the MYNX VCD procedures demonstrated a high level of both procedural and device success.

CONCLUSIONS

The results from this clinical trial demonstrate that patients who have undergone diagnostic or interventional procedures utilizing a 5F, 6F or 7F procedural sheath and treated with the MYNX VCD have superior times to hemostasis and ambulation compared to patients in the standard compression control arm of the Matrix Trial.

SHEATH INTRODUCERS INCOMPATIBLE WITH MYNX CONTROL VASCULAR CLOSURE DEVICE (VCD)

Appropriate training is required prior to use. Read the instructions packaged with the vascular sheath introducer carefully. USE ONLY WITH A STANDARD SHEATH INTRODUCER with up to 12 cm effective length. Cardinal Health test results have shown that the sheath introducers listed in Table 9 are incompatible with the MYNX CONTROL VCD Closure device due to their sheath hub and hemostasis valve design.

Table 9

MANUFACTURER	DESCRIPTION
Medtronic	Input® Introducer Sheaths 11 cm
Cook	Check-Flo® Performer® Introducer

Table 9: List of incompatible sheath introducers when used with corresponding French size specific to VCD. The third-party trademarks used herein are trademarks of their respective owners. **Note:** Testing conducted on sheath introducers available as of January 2017.

The safety and effectiveness of the MYNX CONTROL VCD used with a reinforced sheath introducer that may have a larger outer diameter than that of a standard sheath introducer has not been assessed.

PROCEDURE AND DEVICE PREPARATION

The techniques and procedures described in these Instructions for Use do not represent all medically acceptable protocols, nor are they intended as a substitute for the physician's experience and judgment in treating any specific patient.

HOW SUPPLIED

The MYNX CONTROL VCD is supplied sterile. Do not use if the MYNX CONTROL VCD components or packaging appear to be damaged or defective or if any portion of the packaging has been previously opened.

DO NOT REUSE OR RESTERILIZE. The MYNX CONTROL VCD is for single use only.


The MYNX CONTROL VCD includes (**Figure 1**):

- (1) MYNX CONTROL VCD including balloon catheter and integrated sealant
- (1) 10 ml locking syringe

THIS PAGE INTENTIONALLY LEFT BLANK

THIS PAGE INTENTIONALLY LEFT BLANK

Figure 1: MYNX CONTROL Vascular Closure Device

The 5F MYNX CONTROL VCD has a gray sheath catch. The 6F/7F MYNX CONTROL VCD has a green sheath catch.

PROCEDURE PREPARATION

When using the 5F MYNX CONTROL VCD confirm that the procedural sheath is 5F with effective working length not exceeding 12 cm. When using the 6F/7F MYNX CONTROL VCD confirm that the procedural sheath is 6F or 7F with effective working length not exceeding 12 cm.

Do not attempt to use the MYNX CONTROL VCD to close access sites where a larger procedural sheath than indicated was used.

NOTE:

- ◆ If a patient has had a procedural sheath left in place for an extended period of time, consideration should be given to the use of prophylactic antibiotics before inserting the MYNX CONTROL VCD.
- ◆ If the puncture is at or below the femoral bifurcation, or is an antegrade puncture, the balloon may be prepped with a diluted contrast solution (50% contrast / 50% saline), in place of 100% saline in order to visualize the balloon while pulling back to the arteriotomy and to ensure that the balloon properly abuts the arteriotomy. Do not use 100% contrast solution as this will impact balloon inflation / deflation.

Confirm via femoral arteriogram prior to using the MYNX CONTROL VCD:

- ◆ Common femoral artery single wall puncture.
- ◆ Evidence of adequate flow.
- ◆ No evidence of significant PVD in the vicinity of the puncture.

DEVICE PREPARATION AND POSITIONING

- ◆ Flush the procedural sheath with sterile heparinized saline.
- ◆ Remove device components from packaging
- ◆ Prepare balloon
 - Fill locking syringe with 2 to 3 ml of sterile saline, attach to stopcock and draw vacuum.
 - Check Luer connector and tighten if necessary.
 - Inflate the balloon until the black marker on the inflation indicator is fully visible.
 - Check for leaks in the balloon and syringe connector; retighten if necessary.
 - Discard the device if the balloon does not maintain pressure.
 - Check for air bubbles in the balloon. If air bubbles are visible, deflate the balloon, draw vacuum to remove bubbles and re-inflate.
 - Deflate the balloon and leave syringe at neutral. Do not lock.


THIS PAGE INTENTIONALLY LEFT BLANK

PROCEDURAL STEPS

STEP 1: POSITION BALLOON


- ◆ Insert the MYNX CONTROL VCD into the procedural sheath through the sheath valve. Advance the catheter until the sheath catch nears the hub of the sheath. Rotate the sheath catch as needed to hook onto the side port of the procedural sheath. **(Figure 2a)**
- ◆ Orient the device such that the tension indicator window on the handle assembly is facing upwards.
- ◆ Inflate the balloon until the black marker is fully visible on the inflation indicator and close stopcock. **(Figure 2b)**

Figure 2(a): Connect Device to Sheath; (b): Inflate Balloon

- ◆ Grasp the device handle and align the device with the tissue tract. Pull gently to retract the device until the black line in the tension indicator window aligns with the markers on the side, indicating the balloon is abutting the arteriotomy with the correct amount of tension. **(Figure 3)**

Figure 3: Tension Indicator


THIS PAGE INTENTIONALLY LEFT BLANK

STEP 2: DEPLOY SEALANT

- ◆ While maintaining tension **press button #1** until it is fully aligned with the handle, and the clock symbol ⌚ is visible in the tension indicator window. This deploys and compresses the sealant against the arteriotomy. **(Figure 4)**
- ◆ Lay the device down for **2 minutes**.


Figure 4: Deploy Sealant

STEP 3: REMOVE DEVICE

- ◆ Retract the syringe plunger to lock position. **(Figure 5)**


Figure 5: Lock Syringe Plunger

THIS PAGE INTENTIONALLY LEFT BLANK

- ◆ Apply light fingertip compression proximal to the insertion site and then lightly grasp the device at skin with thumb and forefinger and realign with the tissue tract. (Figure 6a)
- ◆ Open the stopcock to **deflate the balloon**. (Figure 6b)
 - To ensure complete balloon deflation, wait until air bubbles and fluid have stopped moving through the inflation tubing.
- ◆ Pick up the device handle and realign with the tissue tract. **Depress button #2** to pull the deflated balloon into the device. (Figure 6c)

Figure 6(a): Compress and Stabilize; (b): Deflate Balloon; (c): Retract Balloon

THIS PAGE INTENTIONALLY LEFT BLANK

- ◆ While maintaining fingertip compression on the skin, remove the device from the patient.
- ◆ Continue to apply fingertip compression for up to **1 minute** or as needed.
- ◆ Apply a sterile dressing once hemostasis is achieved.

It is recommended that the patient follow physician orders regarding patient ambulation and discharge. Refer to Patient Brochure for post-care instructions.

PRODUCT INFORMATION DISCLOSURE

The MYNX CONTROL Vascular Closure Device (VCD) will be operable during normal and proper use and prior to its expiration date. Cardinal Health shall have no liability to the purchaser or any other third party for product inoperability related to negligence, alteration, modification, abuse, misuse or storage of the MYNX CONTROL VCD by the purchaser or user which is inconsistent with the product specifications or other written instructions supplied by Cardinal Health. Cardinal Health shall have no liability for any matter beyond its control, including patient care, diagnosis, treatment, and surgical procedures that may directly affect the MYNX CONTROL VCD and the results obtained from its use. CARDINAL HEALTH MAKES NO OTHER WARRANTIES TO PURCHASER OR ANY USER WITH RESPECT TO THE MYNX CONTROL VCD, EXPRESS OR IMPLIED, AND CARDINAL HEALTH SPECIFICALLY DISCLAIMS ANY IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. UNDER NO CIRCUMSTANCES SHALL CARDINAL HEALTH'S LIABILITY TO PURCHASER, ITS CUSTOMERS AND USERS EXCEED THE PRICE PAID BY PURCHASER FOR THE MYNX CONTROL VCD SOLD BY CARDINAL HEALTH GIVING RISE TO THE CLAIM. CARDINAL HEALTH WILL NOT BE LIABLE FOR ANY INDIRECT, SPECIAL, INCIDENTAL, CONSEQUENTIAL OR PUNITIVE DAMAGES. Cardinal Health neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with the MYNX CONTROL VCD.

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