MYNXGRIP

MYNXGRIP VASCULAR CLOSURE DEVICE

INSTRUCTIONS FOR USE

EN	MynxGrip Vascular Closure Device
ES	Dispositivo de cierre vascular MynxGrip
DE	MynxGrip Gefäßverschluss-System
FR	Dispositif de fermeture vasculaire MynxGrip
IT	Dispositivo di occlusione vascolare MynxGrip
NL	MynxGrip Vascular Closure Device

MX6721	EN MynxGrip Vascular Closure Device 6F/7F ES Dispositivo de cierre vascular MynxGrip 6 F/7 F DE MynxGrip Gefäßverschluss-System 6 Fr/7 Fr FR Dispositif de fermeture vasculaire MynxGrip 6 Fr/7F IT Dispositivo di occlusione vascolare MynxGrip da 6 Fr/7 Fr NL MynxGrip Vascular Closure Device 6 Fr/7 Fr
MX5021	EN MynxGrip Vascular Closure Device 5F ES Dispositivo de cierre vascular MynxGrip 5 F DE MynxGrip Gefäßverschluss-System 5 Fr FR Dispositif de fermeture vasculaire MynxGrip 5F IT Dispositivo di occlusione vascolare MynxGrip da 5 Fr NL MynxGrip Vascular Closure Device 5 Fr

Instructions for Use

MynxGrip Vascular Closure Device

5F

6F/7F

To ensure proper deployment and use of this device and to prevent injury to patients, read all information contained in these instructions for use.

DEVICE DESCRIPTION

The MynxGrip Vascular Closure Device (MynxGrip) is designed to achieve femoral artery and femoral vein hemostasis via delivery of an extravascular, water-soluble synthetic sealant using a balloon catheter in conjunction with a standard procedural sheath. The sealant is made of a polyethylene glycol (PEG) material which expands upon contact with subcutaneous fluids to seal the arteriotomy or venotomy. The sealant is resorbed by the body within 30 days.

MynxGrip is supplied with a 10 ml locking syringe used for balloon inflation and deflation. There are no components manufactured from latex rubber.

The MynxGrip 5F device has a gray shuttle. The MynxGrip 6F/7F device has a green shuttle. Refer to Figure 1 in the Procedure and Device Preparation Section.

INDICATIONS FOR USE

MynxGrip is indicated for use to seal femoral arterial and femoral venous access sites while reducing times to hemostasis and ambulation in patients who have undergone diagnostic or interventional endovascular procedures utilizing a 5F, 6F or 7F procedural sheath.

CONTRAINDICATIONS

There are no known contraindications for MynxGrip.

WARNINGS

Do not use if components or packaging appear to be damaged or defective or if any portion of the packaging has been previously opened.

DO NOT REUSE OR RESTERILIZE. MynxGrip is for single use only. The balloon catheter is loaded with a single Hydrogel sealant. Reuse of the device would result in no delivery of Hydrogel sealant.

Do not use MynxGrip if the puncture site is located above the most inferior border of the inferior epigastric artery (IEA) (for arterial application) and/or above the inguinal ligament based upon osseous landmarks, since such a puncture site may result in a retroperitoneal hematoma/bleed. Perform a femoral angiogram or venogram to verify the location of the puncture site.

Do not use MynxGrip if the puncture is through the posterior wall or if there are multiple punctures, as such punctures may result in a retroperitoneal hematoma/bleed.

PRECAUTIONS

MynxGrip should only be used by a trained licensed physician or healthcare professional.

MynxGrip should not be used in patients with a known allergy to PEG.

SPECIAL PATIENT POPULATIONS

The safety and effectiveness of MynxGrip have not been established in the following patient populations:

- ◆ Pediatric patients or others with small common femoral arteries or veins (<5 mm in diameter)
- ◆ Patients with clinically significant peripheral vascular disease in the vicinity of the puncture
- ◆ Patients with prior surgical procedure, PTA, stent placement, or vascular graft in the common femoral artery or vein
- ◆ Patients with bleeding disorders such as thrombocytopenia (platelet count <100,000/mm³), hemophilia, von Willebrand's disease or anemia (Hgb <10 g/dL, Hct <30%)
- ◆ Patients with uncontrolled hypertension (systolic BP >180 mm Hg)
- ◆ Patients with morbid obesity (BMI >40 kg/m²)
- ◆ Patients who are pregnant or lactating
- ◆ Patients with documented INR >1.5 or patients currently receiving glycoprotein IIb/IIIa platelet inhibitors

ADVERSE EVENTS

The Mynx Vascular Closure Device (Mynx) was evaluated in a prospective multi-center, non-randomized clinical trial involving 190 patients to achieve femoral arterial access site hemostasis following diagnostic angiography (n=95) or interventional procedures (n=95). The control arm (standard compression) from the MATRIX VSG System Clinical Trial (MATRIX) was used as a historical control group. Table 1 is a report of the total number and the rate of major and minor complications in the Mynx Trial and the standard compression arm of the MATRIX Trial. Clinical data collected on the Mynx device is directly applicable to the MynxGrip.

**Table 1: Reported Major and Minor Complications (Intent-to-Treat)
All Patients**

Major Complications by Event All Patients	Mynx* (n=190)	Standard Compression (Historical Control) (n=164)	p-value†
Vascular Repair	0.0% (0/190)	0.0% (0/164)	N/A
Permanent Access Site-Related Nerve Injury	0.0% (0/190)	0.0% (0/164)	N/A
Surgery For Access Site-Related Nerve Injury	0.0% (0/190)	0.0% (0/164)	N/A
Access Site-Related Bleeding Requiring Transfusion	0.5% (1/190)	0.0% (0/164)	0.0021
New Ipsilateral Lower Extremity Ischemia Requiring Invasive/Non Invasive Intervention	0.0% (0/190)	0.0% (0/164)	N/A
Access Site-Related Infection – Major	0.0% (0/190)	0.0% (0/164)	N/A
Local Access Site Inflammatory Reaction – Major	0.0% (0/190)	0.0% (0/164)	N/A
Generalized Infection	0.0% (0/190)	0.0% (0/164)	N/A
Any Major Complication	0.5% (1/190)	0.0% (0/164)	0.0021
Minor Complications by Event All Patients	Mynx* (n=190)	Standard Compression (Historical Control) (n=164)	p-value†
Pseudoaneurysm – Treated With Thrombin Injection	0.5% (1/190)	0.0% (0/164)	0.0021
Pseudoaneurysm – Not Requiring Treatment‡	2.6% (5/190)	0.0% (0/164)	0.0831
AV Fistula	0.0% (0/190)*	0.0% (0/164)	N/A
Hematoma ≥6 cm	3.2% (6/190)*§	0.6% (1/164)	0.0853
Access Site-Related Bleeding Requiring >30 min to Achieve Hemostasis	0.0% (0/190)	0.6% (1/164)	0.0002
Late Access Site-Related Bleeding (Following Hospital Discharge)	0.0% (0/190)	0.0% (0/164)	N/A
Ipsilateral Lower Extremity Arterial Emboli	0.0% (0/190)	0.0% (0/164)	N/A
Transient Loss of Ipsilateral Lower Extremity Pulse	0.0% (0/190)	0.0% (0/164)	N/A
Ipsilateral Deep Vein Thrombosis	0.0% (0/190)*	0.0% (0/164)	N/A
Transient Access Site-Related Nerve Injury	0.0% (0/190)	0.0% (0/164)	N/A
Access Site-Related Vessel Laceration	0.0% (0/190)	0.0% (0/164)	N/A
Access Site Wound Dehiscence	0.0% (0/190)	0.0% (0/164)	N/A
Local Access Site Infection – Minor	0.0% (0/190)	0.0% (0/164)	N/A
Local Access Site Inflammatory Reaction – Minor	0.0% (0/190)	0.0% (0/164)	N/A
Any Minor Complication	3.7% (7/190)	1.2% (2/164)	0.0921

* Does not include four (4) non-device related events consisting of: AV Fistula (n=2), Ipsilateral Deep Vein Thrombosis (n=1), and Hematoma (n=1).

§ Does not include 1 patient with a pre-existing hematoma prior to deployment of Mynx.

‡ In a subset of Mynx patients, ultrasound assessment was performed prior to patient discharge which revealed some patients had developed small pseudoaneurysms which resolved spontaneously. The events are reported in the table above, but are not included in the overall rate of minor complications.

† p-value based on test for non-inferiority.

Table 2: Reported Major and Minor Complications (Intent-to-Treat) Diagnostic Patients

Major Complications by Event Diagnostic Patients	Mynx* Diagnostic (n=95)	Standard Compression (Historical Control) Diagnostic (n=83)	p-value[†]
Vascular Repair	0.0% (0/95)	0.0% (0/83)	N/A
Permanent Access Site-Related Nerve Injury	0.0% (0/95)	0.0% (0/83)	N/A
Surgery For Access Site-Related Nerve Injury	0.0% (0/95)	0.0% (0/83)	N/A
Access Site-Related Bleeding Requiring Transfusion	0.0% (0/95)	0.0% (0/83)	N/A
New Ipsilateral Lower Extremity Ischemia Requiring Invasive/Non Invasive Intervention	0.0% (0/95)	0.0% (0/83)	N/A
Access Site-Related Infection – Major	0.0% (0/95)	0.0% (0/83)	N/A
Local Access Site Inflammatory Reaction – Major	0.0% (0/95)	0.0% (0/83)	N/A
Generalized Infection	0.0% (0/95)	0.0% (0/83)	N/A
Any Major Complication	0.0% (0/95)	0.0% (0/83)	N/A
Minor Complications by Event Diagnostic Patients	Mynx* Diagnostic (n=95)	Standard Compression (Historical Control) Diagnostic (n=83)	p-value[†]
Pseudoaneurysm – Treated With Thrombin Injection	1.1% (1/95)	0.0% (0/83)	0.0459
Pseudoaneurysm – Not Requiring Treatment [‡]	1.1% (1/95)	0.0% (0/83)	0.0459
AV Fistula	0.0% (0/95)*	0.0% (0/83)	N/A
Hematoma ≥6 cm	2.1% (2/95)*§	1.2% (1/83)	0.0599
Access Site-Related Bleeding Requiring >30 min to Achieve Hemostasis	0.0% (0/95)	0.0% (0/83)	N/A
Late Access Site-Related Bleeding (Following Hospital Discharge)	0.0% (0/95)	0.0% (0/83)	N/A
Ipsilateral Lower Extremity Arterial Emboli	0.0% (0/95)	0.0% (0/83)	N/A
Transient Loss of Ipsilateral Lower Extremity Pulse	0.0% (0/95)	0.0% (0/83)	N/A
Ipsilateral Deep Vein Thrombosis	0.0% (0/95)	0.0% (0/83)	N/A
Transient Access Site-Related Nerve Injury	0.0% (0/95)	0.0% (0/83)	N/A
Access Site-Related Vessel Laceration	0.0% (0/95)	0.0% (0/83)	N/A
Access Site Wound Dehiscence	0.0% (0/95)	0.0% (0/83)	N/A
Local Access Site Infection – Minor	0.0% (0/95)	0.0% (0/83)	N/A
Local Access Site Inflammatory Reaction – Minor	0.0% (0/95)	0.0% (0/83)	N/A
Any Minor Complication	3.2% (3/95)	1.2% (1/83)	0.1458

* Excludes non-device related events consisting of: AV Fistula (n=2) and Hematoma (n=1).

§ Excludes 1 patient with pre-existing hematoma prior to deployment of Mynx.

‡ In a subset of Mynx patients, ultrasound assessment was performed prior to patient discharge which revealed some patients had developed small pseudoaneurysms which resolved spontaneously. The events are reported in the table above, but are not included in the overall rate of minor complications.

† p-value based on test for non-inferiority.

Table 3: Reported Major and Minor Complications (Intent-to-Treat) Interventional Patients

Major Complications by Event Interventional Patients	Mynx* Interventional (n=95)	Standard Compression (Historical Control) Interventional (n=81)	p-value†
Vascular Repair	0.0% (0/95)	0.0% (0/81)	N/A
Permanent Access Site-Related Nerve Injury	0.0% (0/95)	0.0% (0/81)	N/A
Surgery For Access Site-Related Nerve Injury	0.0% (0/95)	0.0% (0/81)	N/A
Access Site-Related Bleeding Requiring Transfusion	1.1% (1/95)	0.0% (0/81)	0.0459
New Ipsilateral Lower Extremity Ischemia Requiring Invasive/Non Invasive Intervention	0.0% (0/95)	0.0% (0/81)	N/A
Access Site-Related Infection – Major	0.0% (0/95)	0.0% (0/81)	N/A
Local Access Site Inflammatory Reaction – Major	0.0% (0/95)	0.0% (0/81)	N/A
Generalized Infection	0.0% (0/95)	0.0% (0/81)	N/A
Any Major Complication	1.1% (1/95)	0.0% (0/81)	0.0459
Minor Complications by Event Interventional Patients	Mynx* Interventional (n=95)	Standard Compression (Historical Control) Interventional (n=81)	p-value†
Pseudoaneurysm – Treated With Thrombin Injection	0.0% (0/95)	0.0% (0/81)	N/A
Pseudoaneurysm – Not Requiring Treatment‡	4.2% (4/95)	0.0% (0/81)	0.4819
AV Fistula	0.0% (0/95)	0.0% (0/81)	N/A
Hematoma ≥6 cm	4.2% (4/95)	0.0% (0/81)	0.4819
Access Site-Related Bleeding Requiring >30 min to Achieve Hemostasis	0.0% (0/95)	1.2% (1/81)	0.0060
Late Access Site-Related Bleeding (Following Hospital Discharge)	0.0% (0/95)	0.0% (0/81)	N/A
Ipsilateral Lower Extremity Arterial Emboli	0.0% (0/95)	0.0% (0/81)	N/A
Transient Loss of Ipsilateral Lower Extremity Pulse	0.0% (0/95)	0.0% (0/81)	N/A
Ipsilateral Deep Vein Thrombosis	0.0% (0/95)*	0.0% (0/81)	N/A
Transient Access Site-Related Nerve Injury	0.0% (0/95)	0.0% (0/81)	N/A
Access Site-Related Vessel Laceration	0.0% (0/95)	0.0% (0/81)	N/A
Access Site Wound Dehiscence	0.0% (0/95)	0.0% (0/81)	N/A
Local Access Site Infection – Minor	0.0% (0/95)	0.0% (0/81)	N/A
Local Access Site Inflammatory Reaction – Minor	0.0% (0/95)	0.0% (0/81)	N/A
Any Minor Complication	4.2% (4/95)	1.2% (1/81)	0.2950

* Excludes non-device related event consisting of: Ipsilateral Deep Vein Thrombosis (n=1).

‡ In a subset of Mynx patients, ultrasound assessment was performed prior to patient discharge which revealed some patients had developed small pseudoaneurysms which resolved spontaneously. The events are reported in the table above, but are not included in the overall rate of minor complications.

† p-value based on test for non-inferiority.

POTENTIAL ADVERSE EVENTS

In addition to the complications noted in the Mynx clinical trial, the following potential complications, which may be related to the endovascular procedure or the vascular closure, may occur: allergic reaction, ecchymosis, superficial vein thrombosis, foreign body/local reaction, retroperitoneal bleed, vessel occlusion, pulmonary embolism, or death.

CLINICAL TRIAL

Mynx was evaluated in a prospective, multi-center, non-randomized clinical trial designed to evaluate its safety and effectiveness in sealing femoral arterial access sites in patients following diagnostic or interventional catheterization procedures. Clinical data collected on the Mynx device is directly applicable to the MynxGrip device. The study was conducted in Germany at five institutions involving 190 patients. Patients eligible for participation included candidates for early ambulation and patients who were clinically indicated for a diagnostic or interventional procedure involving access through the femoral artery using a 5F, 6F or 7F sheath with an overall sheath length ≤ 15.7 cm.

The objective for the safety endpoints of the study was to demonstrate non-inferiority to a historical control group, and the objective for the effectiveness endpoints of the study was to demonstrate superiority to the historical group. The historical control group consisted of patients treated with standard compression who were enrolled in the control arm of the MATRIX Trial conducted under IDE# G030182.

Patients were required to be at least 18 years of age, to have signed an Informed Consent Form, and to have undergone a diagnostic or interventional procedure through the common femoral artery. Patients were excluded if they presented with clinically significant peripheral vascular disease, prior procedure in the ipsilateral common femoral artery ≤ 30 days before the Mynx study catheterization procedure, ipsilateral closure device, known allergy to contrast medium or device materials, a myocardial infarction with elevated ST segment ≤ 24 hours prior to procedure, uncontrolled hypertension, existing bleeding disorder, evidence of infection or local inflammation, chronic corticosteroid therapy ≥ 1 month duration, common femoral artery diameter < 5 mm, pre-existing bleeding around the procedural sheath, preexisting hematoma, intraluminal thrombus, pseudoaneurysm, AV fistula, any type of dissection, fibrotic, calcified, or $> 50\%$ stenotic femoral artery, arterial puncture outside the common femoral artery, ipsilateral venous sheath, multiple arterial sticks, suspected posterior femoral arterial wall puncture, antegrade puncture, intra-aortic balloon pump, or planned extended hospitalization.

The majority of the patients were male ($n=133$) vs. female ($n=57$) with all patients' ages ranging from 40 to 85 years. Of the 190 patients enrolled, 50% were diagnostic patients and 50% were interventional patients. The average activated clotting time (ACT) for interventional patients was 223 ± 84 seconds with a range of 108-634 seconds, and for diagnostic patients, 187 ± 81 seconds with a range of 133-280 seconds.

EFFECTIVENESS RESULTS

The study objective for the primary effectiveness endpoints was to demonstrate that the Mynx patient results were superior to the historical control arm (standard compression patients from the MATRIX Trial). The results of the effectiveness measures are summarized in Table 4 for all patients and in Table 5 for diagnostic and interventional patients.

**Table 4: Effectiveness and Secondary Endpoints
All Patients**

Effectiveness All Patients	Mynx	Standard Compression (Historical Control)	p-value[†]
Time to Hemostasis (minutes)*			
mean ± standard deviation (n)	1.3 ± 2.3 (183)	25.4 ± 16.2 (161)	<0.0001
median (Q1, Q3) (n)	0.5 (0.0, 2.0) (183)	20.0 (15.0, 30.0) (161)	<0.0001
range (min, max)	(0.0, 22.5)	(6.0, 120.0)	N/A
Time to Ambulation (hours)*			
mean ± standard deviation (n)	2.6 ± 2.6 (181)	7.4 ± 4.8 (160)	<0.0001
median (Q1, Q3) (n)	2.0 (1.8, 2.2) (181)	6.0 (4.5, 7.4) (160)	<0.0001
range (min, max)	(1.3, 20.0)	(1.6, 26.9)	N/A
Secondary Endpoints	Mynx	Standard Compression (Historical Control)	p-value[†]
Device Success	93.2% (177/190)	N/A	N/A
Procedure Success	99.5% (189/190)	100.0% (164/164)	1.0000

* The number of patients used to calculate effectiveness measures differ from overall study sample size due to missing values.

† p-value based on test for superiority.

Note: Times to ambulation varied across investigational sites due to individual hospital practices and protocol.

Table 5: Effectiveness and Secondary Endpoints Diagnostic and Interventional Patients

Effectiveness	Mynx Diagnostic	Standard Compression (Historical Control) Diagnostic	p-value†	Mynx Interventional	Standard Compression (Historical Control) Interventional	p-value†
	Time to Hemostasis (minutes)*					
mean ± standard deviation (n)	1.0 ± 1.3 (92)	23.6 ± 17.1 (83)	<0.0001	1.5 ± 2.9 (91)	27.3 ± 15.2 (78)	<0.0001
median (Q1, Q3) (n)	0.5 (0.0, 1.6) (92)	19.0 (14.0, 25.0) (83)	<0.0001	0.6 (0.0, 2.0) (91)	25.0 (19.0, 30.0) (78)	<0.0001
range (min, max)	(0.0, 6.0)	(6.0, 120.0)	N/A	(0.0, 22.5)	(10.0, 120.0)	N/A
Time to Ambulation (hours)*						
mean ± standard deviation (n)	2.5 ± 2.1 (91)	5.4 ± 2.7 (82)	<0.0001	2.8 ± 3.0 (90)	9.4 ± 5.6 (78)	<0.0001
median (Q1, Q3) (n)	2.0 (1.9, 2.3) (91)	5.2 (4.3, 6.1) (82)	<0.0001	1.9 (1.8, 2.2) (90)	7.1 (5.5, 11.8) (78)	<0.0001
range (min, max)	(1.4, 19.6)	(1.6, 26.9)	N/A	(1.3, 20.0)	(2.5, 22.3)	N/A
Secondary Endpoints	Mynx Diagnostic	Standard Compression (Historical Control) Diagnostic	p-value†	Mynx Interventional	Standard Compression (Historical Control) Interventional	p-value†
Device Success	93.7% (89/95)	N/A	N/A	92.6% (88/95)	N/A	N/A
Procedure Success	100% (95/95)	100% (83/83)	N/A	98.9% (94/95)	100% (81/81)	1.0000

* The number of patients used to calculate effectiveness measures differ from overall study sample size due to missing values.

† p-value based on test for superiority.

Note: Times to ambulation varied across investigational sites due to individual hospital practices and protocol.

Procedure success was defined as successfully achieving hemostasis using any method with freedom from major complications. Device success was defined as the ability to deploy the delivery system, deliver the sealant, and achieve hemostasis with Mynx at the femoral artery puncture site. The procedure success rate was 99.5% demonstrating successful hemostasis in all patients with only one major complication reported. The device success rate for Mynx was 93.2%.

Time to discharge was a secondary endpoint for the diagnostic patients only. The time to discharge was defined as the time from advancer tube removal to the time a patient was discharged from the hospital. The mean time to discharge was 32.3 ± 55.6 hours for the Mynx group compared to 20.1 ± 36.1 hours for the historical control standard compression group ($p=0.02$). In the subset of diagnostic patients, the mean time to discharge for Mynx patients was 35.3 ± 71.4 hours compared to 17.3 ± 47.6 hours for the historical control standard compression group ($p=0.06$). The mean time to discharge for Mynx interventional patients was 29.0 ± 31.1 hours compared to 23.0 ± 16.9 hours for the interventional historical control standard compression group ($p=0.15$). The differences in time to discharge may be due to the differences in standard of care for patients treated in Europe (Mynx patients) compared to patients treated in the United States (standard compression historical control group).

**Table 6:
Cumulative Time to Hemostasis
All Patients***

Time to Hemostasis (minutes)	Mynx (n=183)	Standard Compression (Historical Control) (n=161)
1	59.0% (108)	0.0% (0)
2	83.6% (153)	0.0% (0)
3	91.8% (168)	0.0% (0)
4	94.5% (173)	0.0% (0)
5	96.2% (176)	0.0% (0)
10	98.9% (181)	5.0% (8)
15	99.4% (182)	26.7% (43)
20	99.4% (182)	53.4% (86)
25	100% (183)	64.6% (104)
>30	100% (183)	100% (161)

**Table 7:
Cumulative Time to Ambulation
All Patients***

Time to Ambulation (hours)	Mynx (n=181)	Standard Compression (Historical Control) (n=160)
2	50.8% (92)	1.3% (2)
3	90.0% (163)	3.8% (6)
4	92.3% (167)	12.5% (20)
5	93.9% (170)	35.6% (57)
10	97.2% (176)	83.8% (134)
15	98.3% (178)	89.4% (143)
20	100% (181)	95.6% (153)
25	100% (181)	99.4% (159)
>30	100% (181)	100% (160)

*The number of patients differs from overall study sample size due to missing values

**Table 8: Cumulative Time to Discharge
Diagnostic Patients***

Time to Discharge (hours)	Mynx Diagnostic (n=79)	Standard Compression Diagnostic (n=82)
2	0.0% (0)	0.0% (0)
3	2.5% (2)	1.2% (1)
4	16.5% (13)	4.9% (4)
5	19.0% (15)	18.3% (15)
10	25.3% (20)	81.7% (67)
15	29.1% (23)	82.9% (68)
20	38.0% (30)	84.1% (69)
25	67.1% (53)	89.0% (73)
30	82.3% (65)	91.5% (75)
>30	100% (79)	100% (82)

*The number of patients differs from overall study sample size due to missing values

Overall, the primary effectiveness endpoints were successfully met in the Mynx study. Patients treated with Mynx had reduced times to hemostasis and ambulation compared to patients treated with standard compression. In addition, the Mynx procedures demonstrated a high level of both procedural and device success.

CONCLUSIONS

The results from this clinical trial demonstrate that patients who have undergone diagnostic or interventional procedures utilizing a 5F, 6F or 7F procedural sheath and treated with Mynx have superior times to hemostasis and ambulation compared to patients in the standard compression control arm of the Matrix Trial.

PROCEDURE AND DEVICE PREPARATION

The techniques and procedures described in these Instructions for Use do not represent all medically acceptable protocols, nor are they intended as a substitute for the physician's experience and judgment in treating any specific patient.

HOW SUPPLIED


The MynxGrip Vascular Closure Device is supplied sterile. Do not use if MynxGrip components or packaging appear to be damaged or defective or if any portion of the packaging has been previously opened.

DO NOT REUSE OR RESTERILIZE. MynxGrip is for single use only.

The MynxGrip Vascular Closure Device includes (**Figure 1**):

- (1) Balloon catheter with integrated sealant
- (1) 10 ml locking syringe

Figure 1: MynxGrip Device Balloon Catheter

The MynxGrip 5F device has a gray shuttle. The MynxGrip 6F/7F device has a green shuttle.

PROCEDURE PREPARATION

If using a MynxGrip 5F device, confirm that the procedural sheath is 5F with an overall length not exceeding 15.7 cm. When using a MynxGrip 6F/7F device, confirm that the procedural sheath is 6F or 7F with an overall length not exceeding 15.7 cm.

Do not attempt to use the MynxGrip 5F device to close access sites where a 6F or larger procedural sheath was used.

NOTE:

- ◆ If a patient has had a procedural sheath left in place for an extended period of time, consideration should be given to the use of prophylactic antibiotics before inserting MynxGrip.
- ◆ For arterial closure: If the puncture is at or below the femoral bifurcation, or is an antegrade puncture, the balloon may be prepped with a diluted contrast solution (50% contrast / 50% saline), in place of 100% saline, in order to visualize the balloon while pulling back to the arteriotomy, and to ensure that the balloon properly abuts the arteriotomy.

Confirm via femoral arteriogram or venogram:


- ◆ Common femoral artery or vein single wall puncture.
- ◆ Evidence of adequate flow.
- ◆ There is no evidence of significant PVD in the vicinity of the puncture.

DEVICE PREPARATION

- ◆ Flush the procedural sheath with sterile heparinized saline.
- ◆ Remove the syringe from the tray.
- ◆ Remove the balloon catheter from the tray by holding the shuttle and pulling the device out of the protective tubing.
- ◆ Fill the syringe with 2 to 3 ml of sterile saline, attach to the stopcock and draw vacuum.
 - Check the luer connector and tighten if necessary.
- ◆ Inflate the balloon until the black marker on the inflation indicator is fully visible (**Figure 2**).
 - Check for leaks in the balloon and the syringe connector; retighten if necessary.
 - Discard the device if the balloon does not maintain pressure.
 - Check for air bubbles in the balloon. If air bubbles are visible, deflate the balloon, draw vacuum to remove bubbles and re-inflate.
- ◆ Deflate the balloon and leave the syringe at neutral. Do not lock.
- ◆ Do not move the sealant sleeve.

Figure 2: MynxGrip Inflation Indicator

Black marker should be fully visible to indicate balloon is inflated


PROCEDURAL STEPS

STEP 1: POSITION BALLOON

- ◆ Insert MynxGrip into the procedural sheath up to the white shaft marker.
- ◆ Inflate the balloon until the black marker is fully visible on the inflation indicator (**Figure 2**).
- ◆ Close the stopcock.
- ◆ Grasp the handle and withdraw the balloon catheter until the balloon abuts the distal tip of the procedural sheath (**Figure 3a**).
- ◆ Continue to withdraw the balloon catheter until the balloon abuts the arteriotomy or venotomy (**Figure 3b**).

Figure 3a:
First point of tactile resistance

Figure 3b:
Second point of tactile resistance


The balloon is now providing temporary hemostasis and the tip of the introducer sheath is just above the vessel.

STEP 2: DEPLOY SEALANT


- ◆ Align the balloon catheter with the tissue tract.
- ◆ While pulling lightly on the device handle (to ensure the balloon is abutting the arteriotomy or venotomy), open the procedural sheath stopcock and confirm temporary hemostasis.
- ◆ Detach shuttle and advance in a continuous motion until a definitive stop is felt (**Figure 4a**).
- ◆ Immediately grasp the procedural sheath and withdraw it from the tissue tract. Continue retracting until the shuttle locks on the handle (**Figure 4b**).

Figures 4a and 4b:

Advance shuttle, grasp procedural sheath and retract until shuttle locks on handle

- ◆ While maintaining light tension to keep the balloon abutted against the arteriotomy or venotomy, immediately grasp the advancer tube at the skin and gently advance until the single marker is fully visible and then hold in place for **up to 30 seconds** (Figure 5).
- ◆ Lay the device down for **up to 90 seconds**.

STEP 3: REMOVE DEVICE

- ◆ Retract syringe plunger to lock position.
- ◆ Apply light fingertip compression proximal to the insertion site, then lightly grasp the advancer tube at skin with thumb and forefinger and realign with the tissue tract.
- ◆ Open the stopcock to ***deflate the balloon***.
- ◆ To ensure complete balloon deflation, wait until air bubbles and fluid have stopped moving through the inflation tubing.
- ◆ Slowly withdraw the balloon catheter through the advancer tube lumen (**Figure 6**).


NOTE: If unusual resistance is felt during balloon catheter withdrawal, pull the advancer tube and balloon catheter together through the tissue tract.

- ◆ While maintaining fingertip compression on the skin, remove the advancer tube from the tissue tract.
- ◆ Continue to apply fingertip compression for up to ***1 minute*** or as needed. If hemostasis is not achieved, apply additional compression as necessary. Apply a sterile dressing once hemostasis is achieved.

It is recommended that the patient follow physician orders regarding patient ambulation and discharge. Refer to Patient Brochure for post-care instructions.

Figure 6:

Grasp advancer tube at skin and remove balloon

PRODUCT INFORMATION DISCLOSURE

The MynxGrip will be operable during normal and proper use and prior to its expiration date. Cardinal Health shall have no liability to the purchaser or any other third party for product inoperability related to negligence, alteration, modification, abuse, misuse or storage of the MynxGrip by the purchaser or user which is inconsistent with the product specifications or other written instructions supplied by Cardinal Health. Cardinal Health shall have no liability for any matter beyond its control, including patient care, diagnosis, treatment, and surgical procedures that may directly affect the MynxGrip and the results obtained from its use. CARDINAL HEALTH MAKES NO OTHER WARRANTIES TO PURCHASER OR ANY USER WITH RESPECT TO THE MYNXGRIP, EXPRESS OR IMPLIED, AND CARDINAL HEALTH SPECIFICALLY DISCLAIMS ANY IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. UNDER NO CIRCUMSTANCES SHALL CARDINAL HEALTH'S LIABILITY TO PURCHASER, ITS CUSTOMERS AND USERS EXCEED THE PRICE PAID BY PURCHASER FOR THE MYNXGRIP SOLD BY CARDINAL HEALTH GIVING RISE TO THE CLAIM. CARDINAL HEALTH WILL NOT BE LIABLE FOR ANY INDIRECT, SPECIAL, INCIDENTAL, CONSEQUENTIAL OR PUNITIVE DAMAGES. Cardinal Health neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with MynxGrip.

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INSTRUCTIONS FOR USE - MynxGrip Vascular Closure Device

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Emergo Europe
Molenstraat 15
2513 BH, The Hague
The Netherlands

©2015 Cardinal Health. All Rights Reserved. CARDINAL HEALTH, the Cardinal Health LOGO, ESSENTIAL TO CARE, MYNX and MYNXGRIP are trademarks or registered trademarks of Cardinal Health.

©2015 Cardinal Health. Reservados todos los derechos. CARDINAL HEALTH, el logotipo de Cardinal Health, ESSENTIAL TO CARE, MYNX y MYNXGRIP son marcas comerciales o marcas comerciales registradas de Cardinal Health.

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This product and/or its method of manufacture are covered by one or more of the following US patents: 6,605,294; 6,887,974; 7,331,979; 7,335,220; 7,780,980; 7,790,192; 7,803,172; 7,806,856; 7,985,240; 8,105,622; 8,128,654; 8,382,797; and 8,382,798. Other US and foreign patents pending.

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