Transend™ 300

Guidewire with ICE™ Hydrophilic Coating

Directions for Use	2
Instrucciones de uso	4
Mode d'emploi	6
Gebrauchsanweisung	8
Istruzioni per l'uso	10
Gebruiksaanwijzing	12
Instruções de Utilização	14

Transend™ 300

Guidewire with ICE™ Hydrophilic Coating

Rx ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

WARNING

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

DEVICE DESCRIPTION

The Boston Scientific Transend 300 ES and Transend 300 Floppy guidewires are steerable guidewires available in a nominal diameter of 0.014 in (0.36 mm) and a length of 300 cm. The distal 2 cm of the guidewires are shapeable. The proximal portion of the guidewires is coated with Polytetrafluoroethylene (PTFE). All guidewires have a core wire and an ICE hydrophilic coated surface. The distal 3 cm tip of the wire is jacketed with a platinum coil to provide radiopacity.

Contents

- One (1) Insertion Tool
- One (1) Torque Device
- One (1) Guidewire

INTENDED USE/INDICATIONS FOR USE

The Transend 300 ES guidewire and Transend 300 Floppy guidewire are intended for general intravascular use, including the neuro and peripheral vasculature. The guidewires can be torqued to facilitate the selective placement of diagnostic or therapeutic catheters. These devices are not intended for use in coronary arteries. A torque device (pin vise) is included with the guidewire to facilitate directional manipulation of the guidewire.

CONTRAINDICATIONS

None known.

WARNINGS/ADVERSE EVENTS

Vessel trauma may result from the improper use of this device. Follow the instructions for use carefully. When the guidewire is in the body, it should be manipulated only under fluoroscopy. Do not attempt to move the wire without

observing the resultant tip response. Do not leave the wire in a prolapsed condition, as damage to the wire may occur. Never advance the guidewire against excessive resistance without first determining the reason for resistance under fluoroscopy. Excessive force against resistance may result in separation of the guidewire tip, damage to the catheter and/or vessel perforation.

Potential adverse events associated with endovascular procedures (including guidewire usage in the neuro and peripheral vasculature) include but are not limited to:

- Aneurysm rupture
- Access site complications including infection, hematoma and nerve injury
- Cerebral ischemia
- Death
- Embolism (catheter/device, air bubble, plaque, thrombus or char)
- Intracerebral/intracranial hemorrhage
- Pseudoaneurysm
- Seizure
- Stroke
- Transient ischemia attack
- Vasospasm
- Vessel perforation, dissection, trauma or damage

PRECAUTIONS

- This device should be used by physicians thoroughly trained in percutaneous, intravascular techniques and procedures.
- The Boston Scientific guidewire, the guidewire insertion tool, and the torque device are supplied STERILE and non-pyrogenic in unopened, undamaged packages. Verify that the sterility of the device has not been compromised by assuring package integrity has been maintained.
- If the integrity of the packaging has been compromised, **do not use or attempt to resterilize**. Contact your local Boston Scientific representative.
- Prior to a procedure, all equipment to be used for the procedure should be carefully examined to verify proper function and integrity.
- The outer diameter may reach up to 0.0145 in (0.37 mm) on the Transend 300 guidewire family
- Use a standard microcatheter and appropriate guiding catheter with either Transend 300 cm guidewire. Recommended guiding catheter specifications include a 90 cm length and a minimum ID of 0.050 in (1.27 mm).
- Check labeled diameter of diagnostic or therapeutic catheter and verify compatibility with the guidewire outer diameter prior to use.
- Due to the variations of certain catheter tip inner diameters, abrasion of the hydrophilic coating may occur during manipulation. If any resistance is felt during introduction of the catheter, use of a different catheter may be warranted.
- Inspect guidewire prior to use for any surface irregularities and bends or kinks. Damaged and/or irregular guidewires should not be used.
- Exercise care in handling a guidewire during a procedure to reduce the possibility of accidental breakage, bending or kinking.
- To avoid guidewire damage and possible shearing of plastic, do not withdraw or manipulate the guidewire through a metal needle cannula.

- Excessive tightening of the torque device onto the wire may result in abrasion of the coating of the wire.
- Neither the guidewire insertion tool, nor the torque device are intended to enter the body.

HOW SUPPLIED

Handling and Storage

- Store in a cool, dry, dark place.
- Do not use if labeling is incomplete or illegible.
- Do not use if package is opened and/or damaged.
- Use the device prior to the “Use By” date noted on the product label.

PREPARATION FOR USE

Before moving the guidewire from the carrier tube, flush the carrier tube with heparinized saline to hydrate the hydrophilic segment of the guidewire. Repeat injection if difficult removal of the wire occurs.

DIRECTIONS FOR USE

1. Prior to inserting the guidewire into a catheter, flush the catheter with heparinized saline. This will prime the catheter and also provide smooth movement of the guidewire within the catheter.
2. Gently remove the guidewire from the carrier tube and inspect the wire prior to use to verify that it is undamaged.
3. If desired, gently shape the guidewire flexible tip according to standard practices.
4. Carefully insert and advance the guidewire into the catheter. The guidewire insertion tool provided may be used to facilitate the introduction of guidewires into catheter hubs or hemostatic valves.
5. If desired, the guidewire can be used with the torque device provided. To use the torque device, slip the device over the proximal end of the guidewire. When the torque device is in the desired location on the guidewire, tighten the cap to secure the torque device in place. The torque device may be repositioned or removed by loosening the cap and retightening the cap.
6. To rotate the tip of the guidewire to facilitate vessel selection, slowly and carefully rotate the torque device attached to the guidewire.

WARRANTY

Boston Scientific Corporation (BSC) warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond BSC’s control directly affect the instrument and the results obtained from its use. BSC’s obligation under this warranty is limited to the repair or replacement of this instrument and BSC shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. BSC neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **BSC assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.**

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© 2011 Boston Scientific Corporation or its affiliates. All rights reserved.
Boston Scientific International S.A.
Boston Scientific (Australia) Pty Ltd