- EN**
8 **Tornado® Embolization Coils and Microcoils**
Instructions for Use
- CS**
10 **Embolizační spirálky a mikrospirálky**
Tornado®
Návod k použití
- DA**
12 **Tornado® emboliseringscoils og microcoils**
Brugsanvisning
- DE**
14 **Tornado® Embolisationsspiralen und**
Microcoils
Gebrauchsanweisung
- EL**
16 **Σπειράματα και Microcoils εμβολισμού**
Tornado®
Οδηγίες χρήσης
- ES**
19 **Espiraless y microespirales de embolización**
Tornado®
Instrucciones de uso
- FR**
21 **Coils et microcoils d'embolisation Tornado®**
Mode d'emploi
- HU**
24 **Tornado® embolizációs spirálok és**
mikrospirálok
Használati utasítás
- IT**
26 **Spirali e microspirali per embolizzazione**
Tornado®
Istruzioni per l'uso
- NL**
29 **Tornado® embolisatiecoils en Microcoils**
Gebruiksaanwijzing
- NO**
31 **Tornado® emboliseringsspiraler og**
mikrospiraler
Bruksanvisning
- PL**
33 **Spirale i mikrospirale embolizacyjne**
Tornado®
Instrukcja użycia
- PT**
36 **Espiraless e microcoils de embolização**
Tornado®
Instruções de utilização
- SV**
38 **Tornado® emboliseringscoils och Microcoils**
Bruksanvisning


.018 INCH MICROCOILS 0,018 INCH MICROCOILS MICROCOILS DE 0,018 INCH MICROCOILS VAN 0,018 INCH MICROCOILS DE 0,018 INCH
 MIKROSPIRÁLKY 0,018 INCH MICROCOILS 0,018 INCH 0,018 INCH MÉRETŰ MIKROSPIRÁLOK MIKROSPIRALER PÅ 0,018 INCH 0,018 INCH MICROCOILS
 0,018 INCH MICROCOILS MICROESPIRALES DE 0,018 INCH MICROSPIRALI DA 0,018 INCH MIKROSPIRALE O ŠREDNICY 0,018 INCH

1

3


2

4



.035 INCH / .038 INCH COILS
SPIRÁLKY 0,035 INCH / 0,038 INCH
0,035 INCH / 0,038 INCH COILS
0,035 INCH / 0,038 INCH SPIRALEN
ΣΠΕΙΡΑΜΑΤΑ 0,035 INCH / 0,038 INCH
ESPIRALES DE 0,035 INCH / 0,038 INCH
COILS DE 0,035 INCH/0,038 INCH
0,035 INCH / 0,038 INCH MÉRETŰ SPIRÁLOK
SPIRALI DA 0,035 INCH/0,038 INCH
COILS VAN 0,035 INCH/0,038 INCH
SPIRALER PÅ 0,035 INCH / 0,038 INCH
SPIRALE O ŠREDNICY 0,035 INCH / 0,038 INCH
ESPIRAIS DE 0,035 INCH / 0,038 INCH
0,035 INCH/0,038 INCH COILS

5

6

7

- 1 Weaving
- 2 Sheath
- 3 Catheter

- 1 Pohyb sem a tam
- 2 Sheath
- 3 Katetr

- 1 Snoning
- 2 Sheath
- 3 Kateter

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TORNADO® EMBOLIZATION COILS AND MICROCOILS

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

Tornado Embolization Coils (0.035" and 0.038") and 0.018" Microcoils are made of platinum with spaced synthetic fibers, and are supplied preloaded in a loading cartridge. They are designed to be delivered to the target vessel using a soft, straight wire guide through a standard angiographic catheter.

INTENDED USE

Tornado Embolization Coils and Microcoils are intended for arterial and venous embolization in the peripheral vasculature.

MRI SAFETY INFORMATION

Nonclinical testing has demonstrated that single and multiple Tornado Embolization Coils and Microcoils are **MR Conditional** according to ASTM F2503. A patient with this device may be safely scanned after placement under the following conditions.

- Static magnetic field of 3.0 tesla or 1.5 tesla only
- Maximum magnetic field spatial gradient of 1700 gauss/cm (17 T/m) or less
- Maximum MR system reported, whole-body-averaged specific absorption rate (SAR) of ≤ 2.0 W/kg (Normal Operating Mode)

Under the scan conditions provided above, the Tornado Embolization Coils and Microcoils are not expected to result in a temperature rise of more than 3.0°C after 15 minutes of continuous scanning.

The image artifact extends approximately 5.7 mm from the Tornado Embolization Coils and Microcoils as found during nonclinical testing when imaged with a gradient echo pulse sequence and a 3.0 Tesla MRI system.

For US Patients Only

Cook recommends that the patient register the MR conditions disclosed in this IFU with the MedicAlert Foundation. The MedicAlert Foundation can be contacted in the following manners:

Mail: MedicAlert Foundation International
2323 Colorado Avenue
Turlock, CA 95382

Phone: 888-633-4298 (toll free)
209-668-3333 from outside the US

Fax: 209-669-2450

Web: www.medicalert.org

CONTRAINDICATIONS

None known.

WARNINGS

- Positioning of Embolization Coils and Microcoils should be done with particular care. Coils should not be left too close to the inlets of arteries and should be intermeshed with previously placed coils if possible. A minimal but sufficient arterial blood flow should remain to hold the coils against the previously placed coils until a solid clot ensures permanent fixation. The purpose of these suggestions is to minimize the possibility of loose coils becoming dislodged and obstructing a normal and essential arterial channel.
- Tornado Embolization Coils and Microcoils are not recommended for use with polyurethane catheters or catheters with sideports. If a catheter with sideports is used, the embolus may lodge in the sideport or pass inadvertently through it. Use of a polyurethane catheter may also result in lodging of the embolus within the catheter.
- If difficulties occur when deploying the embolization coil, withdraw the wire guide, coil and angiographic catheter simultaneously as a unit.

PRECAUTIONS

- The product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of vascular access sheaths, angiographic catheters and wire guides should be employed.
- Perform an angiogram prior to embolization to determine correct catheter position.

- Prior to introduction of the embolization coil, flush the angiographic catheter with saline.
- **If using a .018 inch Tornado Embolization Microcoil, ensure that the delivery catheter has an internal diameter (ID) of .018 to .025 inch.**

PRODUCT RECOMMENDATIONS

The following table offers specific recommendations.

| Coil Size Diameter | Catheter Type & Size | Wire Guide Type & Size |
|--------------------|--|---|
| .018 inch | MCS-2.5 MCS-2.8 | STF-18 TSF-18 TSFB-18 |
| .035 inch | HNB4.0-35 HNB4.1-35 HNB[R]5.0-35 SCBR4.0-35 SCBR5.0-35 | TSF-35 TSFB-35 TSFBP-35 TSFNA-35 TSFNB-35 |
| .038 inch | HNB[R]5.0-38 HNB[R]6.0-38 SCBR4.0-38 SCBR5.0-38 | TSF-38 TSFB-38 TSFNA-38 TSFNB-38 |

Coil Delivery Technique and Coil Size Selection

Long-term occlusion depends on achieving cross-sectional occlusion of the blood vessel, and coaxial catheters provide the ability to control placement of coils and permanent occlusion. The combination of the coaxial technique and either the anchor or scaffold technique significantly enhances stability of coil deployment.

- **Coaxial technique:** The use of an outer guiding sheath/catheter is the most important step in preventing coil elongation and uncertain long-term occlusion. The outer guiding sheath/catheter provides support, and the inner catheter provides finer selective maneuvers. (**Fig. 7**)
- **Anchor technique:** The anchor technique provides safe and distal occlusion when there is a question about instability of coils. At least 2 cm of a coil is advanced into the side branch, which is normally sacrificed. The rest of the coil is then deployed just proximal to that side branch, and additional coils are packed. (**Fig. 8**) **NOTE:** Use a coil longer than 10 cm for anchor technique.
- **Scaffold technique:** The scaffold technique is used for high-flow vessels when there is concern about migration of a softer coil. A high radial force coil is placed initially. Then, several Inconel coils or platinum coils (soft coils) may be packed within the scaffold. (**Fig. 9**)

In general, the first coil selected should have a diameter that is 20% larger, or at least 2 mm oversized, than the vessel that is being occluded.

INSTRUCTIONS FOR USE

For .018 Inch Microcoils:

1. Perform an angiogram prior to embolization to determine optimal catheter position.
2. Push the loading cartridge completely into the delivery catheter. (**Fig. 1**)
3. Advance the Luer lock connector fitting toward the catheter hub and lock into place. (**Fig. 2**)
4. Use the pusher stylet to load the Tornado Embolization Microcoil into the delivery catheter. (**Fig. 3**) **NOTE:** The stylet must be pushed as far as possible into the loading cartridge to ensure proper loading. Remove the pusher stylet and loading cartridge.
5. To obtain secure placement of the Tornado Embolization Microcoil, we recommend the use of the pusher stylet to push the Microcoil further into the delivery catheter. (**Fig. 4**)
6. Once coil is positioned into delivery catheter, verify position of catheter tip prior to deployment. Coil design permits delivery into the target vessel by saline/contrast flush or by push technique using an appropriately sized wire guide or pusher.
7. Perform final angiogram to confirm coil position within target vessel.

For .035 and .038 Inch Coils:

1. Perform an angiogram prior to embolization to determine optimal catheter position.
2. Firmly grasp the loading cartridge between thumb and forefinger. Introduce the metal end of the loading cartridge into the base of the catheter hub. Lock loading cartridge onto catheter hub by turning Luer lock adapter clockwise. (**Fig. 5**)
3. Maintaining position of the cartridge, advance the stiff portion of the wire guide into the loading cannula. Push the coil into the first 20 to 30 centimeters of the angiographic catheter. (**Fig. 6**) Remove the wire guide and loading cartridge.
4. With the flexible tip of the wire guide, advance the embolization coil to the tip of the catheter. Verify position of the angiographic catheter prior to deployment.

5. Deploy the coil by advancing the wire guide past the tip of the catheter.
6. Perform final angiogram to confirm coil position within target vessel.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2017
Bloomington, IN 47404 U.S.A.