
# Retrievable & Permanent



**OPTEASE®** Retrievable Vena Cava Filter






### **OPTEASE®**

#### Retrievable Vena Cava Filter

The OPTEASE® Vena Cava Filter was designed to perform predictably and dependably. It allows you to focus on the effective treatment of your patients rather than the risk of caval perforation, migration or strut embolization.




Filter Features			
Material	Nitinol		
Max Caval Diameter	30 mm		
Low Profile	6F Sheath		
Ease of Use	Femoral, Jugular and Antecubital Vein Placement Options		
<b>Retrieval Direction</b>	Femoral		
Retrieval Period (per IFU)	Can be retrieved up to and including 12 days after placement		

WARNING: Implant of the OPTEASE® Vena Cava Filter with the hook oriented in the cranial direction can result in life threatening or serious injury including, but not limited to dissection, vessel perforation, migration of the filter with secondary damage to cardiac structures and ineffective pulmonary embolism prevention.

## **Engineered to Perform.**

Built to Stay.



## Two of the **Most Tested Filters** on the Market

#### **OPTEASE®**

Retrievable Vena Cava Filter

**PROOF Trial:** Protection from Pulmonary Embolism with the Cordis OPTEASE® Retrievable Vena Cava Filter<sup>1</sup>

Prospective, multi-center, single filter trial with 150 patients enrolled

## PE Protection Symptomatic PE 1 Month (n=111) 0.0% 6 Month (n=70) 0.0%

## Low Complications Caval Perforation

1 Month (n=111) 6 Month (n=70)	
Filter Migration 1 Month (n=111)	

#### Symptomatic Filter Thrombosis

	IVIOITUI	$(\Pi = \Pi \Pi \Pi)$ .	 	0.	<b>b</b> 70 (1 Case)
)	Month	(n=70)	 	No	New Cases

#### **TRAPEASE®**

Permanent Vena Cava Filter


#### Cordis TRAPESE® Permanent Vena Cava Filter: Experience in 751 Patients<sup>2</sup>

The largest single filter, single center VCF with CT follow up

#### PE Protection

#### **Low Complications**

0	Migration	0.0%
0	Caval Penetration	0.0%
0	Tilting	0.0%
0	Embolization of Filter Fragments	0.0%
	Asymptomatic Near Caval Oschusions	0.7%



### **OPTEASE®**

Retrievable Vena Cava Filter

## Retrievable & Permanent

PROVEN & RELIABLE

#### Cordis OPTEASE® Retrievable Vena Cava Filter Ordering Information

Description	Access Site	Catalog Numbers
Cordis OPTEASE® Vena Cava Filter and Introduction Kit (55 cm)	Femoral	466-F210AF
Cordis OPTEASE® Vena Cava Filter and Introduction Kit (55 cm)	Jugular	466-F210AJ
Cordis OPTEASE® Vena Cava Filter and Introduction Kit (90 cm)	Antecubital, Jugular	466-F210BJ

#### Cordis TRAPEASE® Permanent Vena Cava Filter Ordering Information

Description	Access Site	Catalog Numbers
Cordis TRAPEASE® Vena Cava Filter and Introduction Kit (55 cm)	Jugular, Femoral	466-P306A
Cordis TRAPEASE® Vena Cava Filter and Introduction Kit (90 cm)	Antecubital, Jugular, Femoral	466-P306B

<sup>&</sup>lt;sup>1</sup>Ziegler, JW et al. PROOF Trial: Protection from Pulmonary Embolism with the OPTEASE® Filter. JVIR 2008; 19:1165-1170. 
<sup>2</sup>Kalva, SP et al. TRAPEASE® Vena Cava Filter: Experience in 751 Patients. J Endovasc Ther 2006; 13:365-372.

**WARNING**: Implant of the OPTEASE® Vena Cava Filter with the hook oriented in the cranial direction can result in life threatening or serious injury including, but not limited to dissection, vessel perforation, migration of the filter with secondary damage to cardiac structures and ineffective pulmonary embolism prevention.

#### IMPORTANT INFORMATION:

Prior to use, refer to the full "Instructions For Use" supplied with these devices for more information on indications, contraindications, suggested procedures, warnings and precautions. Contact your Cordis sales representative for availability and ordering. As part of the Cordis policy of continuous product development, we reserve the right to change product specifications without prior notification. EMEA For Healthcare Professionals only. EU790 2/16.

© 2016 Cardinal Health. All rights reserved. CORDIS, the Cordis LOGO, and OPTEASE® are trademarks or registered trademarks of Cardinal Health.

• Tel: +32 2 746 30 00 • Fax: +32 2 746 34 05 - 2E-800-2123-1