INSTRUCTIONS FOR USE FOR:

GORE

MOLDING & OCCLUSION

B A L L O O N C A T H E T E R

en

English

hu

Magyar

cs

Čeština

it

Italiano

da

Dansk

no

Norsk

nl

Nederlands

pl

Polski

fi

Suomi

pt

Português

fr

Français

es

Español

de

Deutsch

sv

Svenska

el

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tr

Türkçe

INSTRUCTIONS FOR USE FOR GORE® MOLDING AND OCCLUSION BALLOON CATHETER

INTRODUCTION

The GORE® Molding and Occlusion Balloon Catheter is designed to assist in the expansion of self-expanding stent grafts as well as for the occlusion of large diameter vessels. Ballooning of a stent graft may improve modeling of the graft material and fixation of the stent graft to the vessel wall. There are risks involved with any medical procedure.

Both physician and patient should fully understand those risks associated with surgery, and additional new risks associated specifically with the use of this endoluminal device.

(See CAUTION/ PRECAUTION NOTICES, WARNINGS, COMPLICATIONS, and CONTRAINDICATIONS, herein.)

NOTE: Read the entire GORE Molding and Occlusion Balloon Catheter Instructions for Use prior to performing ballooning procedure. Failure to do so may result in complications

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

WARNING: Only physicians trained in vascular surgery, interventional radiology or cardiology, and who have completed training or have experience with stent grafts and balloon catheters should consider using this device.

DESCRIPTION

The GORE® Molding and Occlusion Balloon Catheter (MOB) is a sterile, single use, single-lobed polyurethane balloon catheter (see **Figure 1**). The compliant polyurethane balloon is mounted on the leading end of a 3 lumen catheter shaft (two inflation lumens and one guidewire lumen). Termination of the leading end of the catheter is an atraumatic catheter leading tip for smooth transition from the guidewire to catheter transition. Radiopaque markers (approximately 40 mm apart) indicate the proximal and distal end of the balloon which aid in proper balloon placement under fluoroscopy. Both of the (2) inflation lumens are in communication with each end of the balloon to facilitate balloon catheter preparation, inflation and deflation. At the trailing end of the balloon catheter is a dual port (balloon inflation and guidewire) y-arm. The inflation port of the y-arm is in communication with both of the balloon inflation lumens and is affixed with a luer lock and three way stopcock via an extension tube. The guidewire lumen of the y-arm allows introduction of a 0.035" (0.89 mm) diameter guidewire for over-the-wire access. The trailing end of the guidewire lumen is affixed with a flushing / guidewire port with luer lock used for flushing the guidewire lumen. The balloon catheter proximal y-arm is provided within a housing which contains markings of the balloon length and inflation range diameter.

The GORE® Molding and Occlusion Balloon Catheter is available in a single catheter length of 90 cm. The balloon (MOB37) can be inflated to diameters of 10 mm to 37 mm (see **Table 1** for recommended balloon inflation volume). The GORE® Molding and Occlusion Balloon Catheter has been tested to temporarily occlude vessels up to 37 mm in diameter. The balloon catheter profile is 10 Fr GORE® DrySeal Flex introducer sheath compatible.

The GORE® Molding and Occlusion Balloon Catheter consists of the following components:

- 3 Lumen Catheter Shaft
- Compliant polyurethane balloon
- Radiopaque marker bands
- Strain relief
- Balloon Inflation Extension tube with 3-way stopcock
- Proximal Y-arm hub and cover with balloon length and inflation range diameter markings

The GORE® Molding and Occlusion Balloon Catheter has a compliant polyurethane balloon with a maximum inflation diameter of 37 mm.

Catalogue Number	Catheter Effective Length
MOB37	90 cm

Figure 1. Illustration of Device


Table 1. Inflation Table

Balloon Inflation Chart	
Diameter (mm)	Volume (ml)
10	2
20	7
30	16
37*	35

* Maximum balloon inflation diameter

CAUTION: The above balloon inflation chart volumes and diameters are only an approximate guide. Balloon inflation should only be performed while closely monitoring under fluoroscopy. DO NOT inflate the balloon to a diameter larger than the vessel it is being inflated within.

Do not exceed 37 mm balloon diameter, balloon rupture may occur.

Do not exceed 48 cc inflation volume at 37 mm constrained balloon diameter.

Balloon over-inflation may result in vessel damage or rupture.

INTENDED USE

The GORE® Molding and Occlusion Balloon Catheter is intended for temporary occlusion of large diameter vessels or to assist the expansion of self-expanding endovascular prostheses (stent grafts).

CONTRAINDICATIONS

The GORE® Molding and Occlusion Balloon Catheter is contraindicated in patients who:

- are contraindicated to contrast media or anticoagulants
- have an arterial entry site that cannot accommodate a 10 Fr introducer sheath
- are minors
- are pregnant

WARNINGS

- Excessive inflation volume may result in balloon rupture, embolization, vessel damage, vessel rupture, or patient death.
- Do not continue if resistance is felt during advancement or retraction of the balloon catheter. Stop and assess the cause of the resistance.
- Continued movement against resistance may result in prosthesis migration, vessel damage, and / or catheter damage.
- Do not advance or retract the balloon catheter while the balloon is inflated. Ensure the balloon is deflated by locking the balloon inflation stop cock while maintaining a vacuum on the syringe and using fluoroscopy to confirm balloon complete deflation before moving the balloon catheter.
- Do not inflate the balloon in areas of significant calcified plaque. Balloon rupture and / or vessel damage may occur.
- If the deflated balloon catheter gets caught on the leading edge of the introducer sheath during retraction, carefully push the balloon catheter forward until it is outside the introducer sheath, slightly re-inflate then deflate the balloon until a vacuum is applied, while maintaining a vacuum on the balloon / syringe, re-lock the balloon inflation stopcock. While keeping the introducer sheath stationary, carefully try to retract the balloon catheter into the introducer sheath. Alternative methods are to pull the guidewire back until the leading guidewire tip is well inside the balloon catheter and / or to pull the introducer sheath back slightly until the leading tip of the sheath and balloon is in a straighter section of the vessel. At this point, try retracting the balloon catheter from the sheath. If resistance is still felt, remove the balloon catheter and the introducer sheath together.
- Incorrect usage of the balloon catheter may result in prosthesis movement and may require additional prosthesis placement or surgical intervention.
- Observe inflation and deflation of the balloon under fluoroscopy to ensure proper functioning of the balloon catheter. Withdrawal of the balloon catheter prior to deflation may result in prosthesis migration, vessel damage, and / or catheter damage.
- Do not use high pressure injector.
- Do not use the GORE® Molding and Occlusion Balloon Catheter for more than 20 inflation / deflation cycles.

PRECAUTIONS

- Do not resterilize. The GORE® Molding and Occlusion Balloon Catheter is intended for single use only.
- Do not use after the labeled "use by" (expiration) date.
- Do not use device if damaged or if sterile barrier has been compromised.

ADVERSE EVENTS

Adverse events which may require intervention include, but are not limited to:

fever,
infection,
hemorrhage or bleeding at the site of catheter introduction,
hematoma,
local neurologic damage,
arteriovenous fistula,
pseudoaneurysm,
thrombosis,
embolization;
trauma to the vessel wall, including spasm, dissection, perforation or rupture,
renal insufficiency or failure,
allergic reaction to contrast media,
respiratory complications,
stroke,
arrhythmia,
myocardial infarction,
prosthesis movement or occlusion,
aneurysm rupture,
death

CONTENTS/HOW SUPPLIED

The GORE® Molding and Occlusion Balloon Catheter is supplied sterile and non-pyrogenic.

STORAGE AND HANDLING

Store in cool dry place.

REQUIRED MATERIALS

- 35 ml syringe (recommended)
- Heparinized saline
- 0.035" (0.89 mm) stiff guidewire, 180 cm or longer
- 10 Fr or larger introducer sheath
- Radiopaque contrast media

DIRECTIONS FOR USE

Balloon Catheter Preparation

1. Remove the GORE® Molding and Occlusion Balloon Catheter from its packaging and examine for possible damage.
2. Remove balloon sleeve.
3. Flush heparinized saline through the guidewire lumen.
4. Fill a 35 cc syringe with 30:70 radiopaque contrast media to saline solution (diluted contrast solution) and attach to the stopcock of the balloon inflation port.
5. Orient the balloon catheter with the leading tip downward. Aspirate air from the balloon and catheter shaft using the syringe. Next, partially inflate the balloon with diluted contrast solution. Repeat aspiration and inflation until all of the air has been evacuated from the catheter and balloon.

6. Inflate the balloon catheter with 10 cc of diluted contrast solution. Inspect the balloon catheter for the presence of leaks.

NOTE: If the balloon does not inflate, catheter leakage is noted during inflation, or deflation vacuum cannot be maintained on the balloon catheter, do not use balloon catheter.

Balloon Catheter Insertion

- Utilizing the 35 ml syringe, draw a vacuum on the balloon catheter and close the stopcock to maintain vacuum on the balloon catheter prior to insertion into the introducer sheath.
- To aid in insertion of the GORE® Molding and Occlusion Balloon Catheter through introducer sheath hemostatic valve, it is recommended to wet the balloon with heparinized saline prior to insertion into the introducer sheath valve.

Balloon Catheter Inflation(s)

- Maintain a vacuum on the balloon catheter, advance the balloon catheter over a 0.035" (0.89 mm) guidewire, through a 10 Fr or larger introducer sheath. To minimize the potential for balloon catheter kink, advance the balloon catheter in small (1 cm) increments until fully advanced through the introducer sheath.
- Using fluoroscopic guidance, position the balloon catheter using the radiopaque markers.
- The volume of diluted contrast solution used to inflate the balloon catheter will determine the approximate balloon diameter achieved (Table 1). Observe inflation of the balloon under fluoroscopy to determine if less or additional inflation volume is required to optimize balloon contact. To avoid vessel trauma, do not over-inflate the balloon in relation to the diameter of the vessel or other devices.

Balloon Catheter Withdrawal

- When the procedure has been completed, deflate the balloon by drawing a vacuum on the balloon catheter, locking the balloon inflation stopcock while under vacuum, and remove the balloon from vasculature / introducer sheath.
- Observe deflation and withdrawal of the balloon catheter under fluoroscopy.
- If the deflated balloon catheter gets caught on the leading edge of the introducer sheath during retraction, carefully push the balloon catheter forward until it is outside the introducer sheath, slightly re-inflate then deflate. Reassure there is vacuum applied to the balloon via the syringe and re-close the balloon inflation port stopcock. Carefully try to retract the balloon catheter into the introducer sheath. Alternative methods are to pull the guidewire back until the leading guidewire tip is well inside the balloon catheter and / or to pull the introducer sheath back slightly until the leading tip of the sheath is in a straighter section of the vessel. At this point try retracting the balloon catheter into the sheath. If resistance is still felt, remove the balloon catheter and the introducer sheath together.

BALLOON CATHETER RELATED ADVERSE EVENT REPORTING

Any adverse event involving the GORE® Stent Graft Molding Balloon Catheter should be reported to W. L. Gore & Associates immediately. To report an event in the US, call 800.437.8181. Outside the US, contact your local technical representative.

DEFINITIONS

 Authorised Representative in the European Community

 Batch Code

 Catalogue Number

 Caution


 Only CAUTION: USA Federal Law restricts the sale, distribution, or use of this device to, by, or on the order of a physician.

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W. L. GORE & ASSOCIATES B.V.
CREAGH MEDICAL LTD,