INSTRUCTIONS FOR USE FOR:

TRI-LOBE

BALLOON CATHETER

en

English

hu

Magyar

cs

Čeština

it

Italiano

da

Dansk

no

Norsk

nl

Nederlands

pl

Polski

fi

Suomi

pt

Português

fr

Français

es

Español

de

Deutsch

sv

Svenska

el

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tr

Türkçe

INSTRUCTIONS FOR USE FOR GORE® Tri-Lobe Balloon Catheter

CAUTION: USA Federal Law restricts the sale, distribution, or use of this device to, by, or on the order of a physician.

Carefully read all instructions prior to use. Observe all warnings and precautions noted throughout these instructions. Failure to do so may result in complications.

DESCRIPTION

The GORE® Tri-Lobe Balloon Catheter (balloon catheter) (Figure 1) is a compliant, tri-lobed polyurethane balloon catheter. The lobed design of the balloon catheter allows for inflation without complete blockage of aortic blood flow (see Figure 2 for balloon inflation volume). The three polyurethane balloons are mounted on the leading end of a multi-lumen catheter shaft. Radiopaque markers indicate the balloon edges. Each of the three inflation lumens is in communication with one of the balloons. The inflation port is in communication with all of the inflation lumens and is affixed with a luer lock. The guidewire lumen allows introduction of a 0.035" (0.89 mm) diameter guidewire for over-the-wire access. The trailing end of the guidewire lumen is affixed with a flushing / guidewire port with luer lock used for flushing the guidewire lumen. A Tuohy-Borst valve is integrated into the trailing end of the guidewire lumen.

The GORE® Tri-Lobe Balloon Catheter is available in two sizes. The smaller balloon (BCM1634) can be inflated to diameters of 16 mm to 32 mm and the larger balloon (BCL2645) can be inflated to diameters of 26 mm to 42 mm (Figure 2). Both balloon sizes are 18 Fr sheath compatible.

FIGURE 1: GORE® TRI-LOBE BALLOON CATHETER


FIGURE 2: GORE® TRI-LOBE BALLOON CATHETER INFLATION VOLUME GUIDE

BCM1634 – 10 mm Balloon Catheter	
Intended Balloon Diameter (mm)	Recommended Inflation Volume (mL)
16	4
20	6
24	8
28	10
32	13
BCL2645 – 16 mm Balloon Catheter	
Intended Balloon Diameter (mm)	Recommended Inflation Volume (mL)
26	12
30	14
34	18
38	21
42	25

The above balloon inflation chart volumes and diameters are only an approximate guide. Balloon inflation should only be performed while closely monitoring under fluoroscopy. Do not inflate the balloon to a diameter larger than the vessel it is being inflated within.

INTENDED USE

The GORE® Tri-Lobe Balloon Catheter is intended to assist in the dilatation of self-expanding endoprostheses in large diameter vessels.

CONTRAINDICATIONS

There are no known contraindications.

WARNINGS

- Excessive inflation volume may result in balloon rupture, embolization, vessel damage, vessel rupture, or patient death.
- Do not continue if resistance is felt during advancement or retraction of the balloon catheter. Stop and assess the cause of the resistance. Continued movement against resistance may result in prosthesis migration, vessel damage, and / or catheter damage.
- Do not advance or retract the balloon catheter while the balloons are inflated. Ensure the balloons are deflated using fluoroscopy before moving the balloon catheter.
- Do not inflate the balloons in areas of significant calcified plaque. Balloon rupture and / or vessel damage may occur.
- Do not use in the ascending aorta as W. L. Gore & Associates, Inc. (Gore) has insufficient information to support this application.
- If the deflated balloon catheter gets caught on the leading edge of the introducer sheath during retraction, carefully push the balloon catheter slightly forward, while keeping the introducer sheath stationary, and rotate the balloon catheter. Carefully try to retract the balloon catheter into the introducer sheath. Alternative methods are to pull the guidewire back until the guidewire tip is well inside the balloon catheter and / or to pull the introducer sheath back slightly until the tip of the sheath is in a straighter section of the vessel. At this point, try retracting the balloon catheter from the sheath. If resistance is still felt, remove the balloon catheter and the introducer sheath together.
- Incorrect usage of the balloon catheter may result in prosthesis movement and may require additional prosthesis placement or surgical intervention.
- Observe inflation and deflation of the balloons under fluoroscopy to ensure proper functioning of the balloon catheter. Withdrawal of the balloon catheter prior to deflation may result in prosthesis migration, vessel damage, and / or catheter damage.
- Do not use high pressure injector.

PRECAUTIONS

- The GORE® Medical Device is designed for single use only; do not reuse device. Gore does not have data regarding reuse of this device. Reuse may cause device failure or procedural complications including device damage, compromised device biocompatibility, and device contamination. Reuse may result in infection, serious injury, or patient death.
- Do not resterilize.
- Do not use after the labeled "use by" (expiration) date.
- Do not use device if damaged or if sterile barrier has been compromised.

ADVERSE EVENTS

Adverse events which may require intervention include, but are not limited to: fever; infection; hemorrhage or bleeding at the site of catheter introduction; hematoma; lymph fistula; local neurologic damage; arteriovenous fistula; pseudoaneurysm; thrombosis; embolization; trauma to the vessel wall, including spasm, dissection, perforation or rupture; renal insufficiency or failure; allergic reaction to contrast media; respiratory complications; stroke; arrhythmia; myocardial infarction; prosthesis movement or occlusion; aneurysm rupture and death.

CONTENTS / HOW SUPPLIED

The GORE® Tri-Lobe Balloon Catheter is supplied sterile and non-pyrogenic.

STORAGE AND HANDLING

Store in a cool dry place.

REQUIRED MATERIALS

- 60 cc or smaller syringe
- Heparinized saline
- 0.035" (0.89 mm) guidewire, 250 cm or longer
- 18 Fr or larger introducer sheath
- Radiopaque contrast media
- 3-way stopcock

DIRECTIONS FOR USE

Balloon Catheter Preparation

1. Remove the GORE® Tri-Lobe Balloon Catheter from its packaging and examine for possible damage.
2. Open the Tuohy-Borst valve on the flushing / guidewire port. Flush heparinized saline through the guidewire lumen.
3. Attach a 3-way stopcock to the inflation port.
4. Fill a 60 cc or smaller syringe with 30:70 radiopaque contrast media to saline solution (diluted contrast solution) and attach to the stopcock of the inflation port.
5. Orient the balloon catheter with the tip downward. Aspirate air from the balloons and shaft using the syringe. Next, partially inflate the balloons with diluted contrast solution. Repeat aspiration and inflation until all of the air has been evacuated from the balloons.
6. Inflate the balloon catheter with 15 cc of diluted contrast solution. Inspect the balloon catheter for the formation of three lobes and for the presence of leaks.
NOTE: If three lobes are not formed or if balloon catheter leakage is noted during inflation, do not use balloon catheter.
7. Draw a vacuum on the balloon catheter and close the stopcock to maintain the vacuum on the balloon catheter.

Balloon Catheter Inflation(s)

1. Maintaining a vacuum on the balloon catheter, advance the balloon catheter over a 0.035" (0.89 mm) guidewire, through an 18 Fr or larger introducer sheath.
2. Using fluoroscopic guidance, position the balloon catheter using the radiopaque markers.
3. The volume of diluted contrast solution used to inflate the balloon catheter will determine the balloon diameter achieved (Figure 2). Observe inflation of the balloons under fluoroscopy to determine if additional inflation volume is required to optimize balloon contact. To avoid vessel trauma, do not over-inflate the balloons in relation to the diameter of the artery or other devices.
4. Deflate the balloon catheter, rotate the balloon catheter approximately 60° and repeat the inflation.

Balloon Catheter Withdrawal

1. When the procedure has been completed, deflate the balloons by drawing a vacuum on the balloon catheter and remove from vasculature.
2. Observe deflation and withdrawal of the balloon catheter under fluoroscopy.
3. If the deflated balloon catheter gets caught on the leading edge of the introducer sheath during retraction, carefully push the balloon catheter slightly forward, while keeping the introducer sheath stationary, and rotate the balloon catheter. Carefully try to retract the balloon catheter into the introducer sheath. Alternative methods are to pull the guidewire back until the guidewire tip is well inside the balloon catheter and / or to pull the introducer sheath back slightly until the tip of the the sheath is in a straighter section of the vessel. At this point try retracting the balloon catheter from the sheath. If resistance is still felt, remove the balloon catheter and the introducer sheath together.

BALLOON CATHETER RELATED ADVERSE EVENT REPORTING

Any adverse event involving the GORE® Tri-Lobe Balloon Catheter should be reported to W. L. Gore & Associates immediately. To report an event in the US, call 800.437.8181. Outside the US, contact your local technical representative.

DEFINITIONS

 Authorised Representative in the European Community

 Catalogue Number

 Caution

 CAUTION: USA Federal Law restricts the sale, distribution, or use of this device to, by, or on the order of a physician.

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© 2002, 2003, 2005, 2008, 2010, 2012, 2014-2015, 2017, 2020 W. L. Gore & Associates, Inc.
GORE® and designs are trademarks of W. L. Gore & Associates.
W. L. GORE & ASSOCIATES B.V.