# INSTRUCTIONS FOR USE

## Echelon™ Micro Catheter

This page is intentionally left blank.

## TABLE OF CONTENTS

### *Echelon™ Micro Catheter*

---

English.....	5
Français.....	6
Deutsch.....	8
Italiano.....	9
Español.....	11
Svenska.....	12
Nederlands.....	14
Português.....	15
Suomi.....	17
Dansk.....	18
Ελληνικά.....	20
Česky.....	21
Magyar.....	23
Русский.....	24
Polski.....	26
Türkçe.....	27
Norsk.....	29
Slovenčina.....	30
Română.....	32
Български.....	33
한국어.....	35
عربي.....	37
<b>Symbol Glossary</b> .....	<b>38</b>

---

This page is intentionally left blank.

## Instructions for Use

### Echelon™ Micro Catheter

#### CAUTION

Federal (USA) law restricts this device to sale, distribution, and use by or on the order of a physician.

This device should be used only by physicians with a thorough understanding of angiography and/or percutaneous neurointerventional procedures.

#### CONTENTS

One Echelon™ Micro Catheter with tip shaping mandrel.

#### DESCRIPTION

The Echelon™ Micro Catheter is an endhole, single-lumen catheter designed to be introduced over a steerable guidewire into the vasculature. The proximal end of the catheter incorporates a standard luer adapter to facilitate the attachment of accessories. The catheter has a semi-rigid proximal shaft which transitions into the flexible distal shaft to facilitate the advancement of the catheter in the anatomy. Dual radiopaque markers at the distal end facilitate fluoroscopic visualization. The outer surface of the catheter is coated to increase lubricity.

#### INDICATIONS FOR USE

The Echelon™ Micro Catheter is intended to access peripheral and neuro vasculature for the controlled selective infusion of physician-specified therapeutic agents such as embolization materials and of diagnostic materials such as contrast media.

#### CONTRAINDICATIONS

- Not intended for use in the coronary vasculature.
- Contraindicated when, in the medical judgment of the physician, such a procedure may compromise the patient's condition.
- Contraindicated for neonatal and pediatric use.

#### POTENTIAL COMPLICATIONS

Potential complications include, but are not limited to:

- Hematoma at the site of entry
- Vessel perforation
- Vessel spasm
- Hemorrhage
- Pain and tenderness
- Vascular thrombosis
- Thrombolytic episodes
- Neurological deficits including stroke and death
- Ischemia

#### WARNINGS

- Infusion pressure with this device should not exceed 600 PSI (4137 kPa). Pressures in excess of 600 PSI (4137 kPa) may result in catheter rupture, which may result in patient injury.
- If flow through the micro catheter becomes restricted, do not attempt to clear the device by high pressure infusion. Either remove the micro catheter to determine the cause of the obstruction or replace it with a new micro catheter. Excessive pressure may cause catheter rupture, which may result in patient injury.
- Never advance or withdraw an intraluminal device against resistance until the cause of the resistance is determined by fluoroscopy. Excessive force against resistance may result in damage to the device or vessel perforation or other patient injuries.
- This device is supplied STERILE for single use only. Do not reprocess or re-sterilize. Reprocessing and re-sterilization increase the risks of patient infection and compromised device performance.

#### PRECAUTIONS

- Ensure embolic material compatibility with catheter prior to use.
- Prior to use, carefully examine the micro catheter and packaging to verify that they have not been damaged during shipment.
- Prior to use, all accessory devices and agents should be fully prepared according to the manufacturer's instructions.
- Always monitor infusion rates when using the micro catheter.
- When injecting contrast for angiography, ensure that the micro catheter is not kinked or occluded.
- The micro catheter has a lubricious hydrophilic coating on the outside of the catheter. It must be kept hydrated in order to be lubricious. This can be accomplished by attaching the Y-connector to a continuous saline drip.
- When the infusion catheter is in the body, it should be manipulated only under fluoroscopy. Do not attempt to move the micro catheter without observing the resultant tip response.
- Because the micro catheter may be easily advanced into narrow, selective vasculature, repeatedly verify (for example by slightly withdrawing the catheter) that the micro catheter has not been advanced so far as to interfere with its removal.

#### STORAGE

The micro catheter should be stored in a dry place at between 50°F (10°C) and 90°F (32°C).

#### DIRECTIONS FOR USE

Flow Rates @ 100 PSI (690 kPa) and 300 PSI (2,070 kPa) – mL/sec

	Saline		76% Contrast/ Saline (50/50)		60% Ionic Contrast		76% Ionic Contrast	
Pressure	100	300	100	300	100	300	100	300
<b>Echelon 14 (105-5092-150, 145-5092-150, 190-5092-150)</b>								
Flowrates	0.25	0.85	0.33	0.46	0.16	0.27	0.06	0.09
PVA Particles	Maximum 500µm							
Guide Catheter Minimum ID	0.050" (1.28 mm)							
<b>Echelon 10 (105-5091-150, 145-5091-150, 190-5091-150)</b>								
Flowrates	0.25	0.85	0.33	0.46	0.16	0.27	0.06	0.09
PVA Particles	Maximum 500µm							
Guide Catheter Minimum ID	0.035" (0.89mm)							

1. Place the appropriate guide catheter following recommended procedures. Connect a hemostatic side arm adapter to the guide catheter to prevent backflow of blood during insertion of micro catheter. Connect a one-way stopcock to the hemostatic side arm adapter to continuously flush the guide catheter with saline.
2. Model numbers 145-5091-150, 145-5092-150, 190-5091-150, and 190-5092-150 are provided in a tray with lid, and the catheter tip is held within the groove of a curve retainer. To open the tray, lift the lid at the corner of the tray with the catheter hub. The catheter and curve retainer may be removed from the tray together.
3. To free the catheter tip from the curve retainer, pinch two opposite corners back. This opens the groove. Lift the catheter tip out of the groove.
4. Hydrate the catheter outer surface.
  - a. If the catheter is provided in a coil dispenser, flush the coil with heparinized saline through the female luer attached to the coil.
  - b. Otherwise, the catheter may be hydrated in the tray.
5. Prior to use, flush the catheter lumen with heparinized saline by attaching a saline filled syringe to the catheter hub.
6. Remove the appropriate steerable guidewire from its packaging and inspect for damage. Follow manufacturer's instructions for preparing and using the guidewire.
7. Inspect the catheter prior to use to verify that it is undamaged. Retain the packaging for storage of the catheter if it is not used immediately in the procedure.
8. Carefully insert guidewire into the hub of the micro catheter and advance guidewire into the catheter lumen.
9. Model numbers 145-5091-150, 145-5092-150, 190-5091-150, and 190-5092-150 include a split introducer to aid in insertion into the hemostatic side arm adapter. To use, slide the Split Introducer from the proximal end until it covers the distal end of catheter.
10. Close the one-way stopcock.
11. Loosen the hemostatic valve.
12. Introduce the guidewire and micro catheter as a unit through the hemostatic port of the hemostatic sidearm adapter into the lumen of the guiding catheter. Advance guidewire/catheter assembly to the distal tip of the guiding catheter.
  - a. If using the Split Introducer slide the introducer back to the proximal end of the catheter next to the hub.
  - b. Or peel introducer away from the microcatheter.
13. Tighten the valve around the micro catheter to prevent backflow, but allowing some movement through the valve by the micro catheter.
14. Open the one-way stopcock.
15. Alternatively advance the guidewire and micro catheter until the desired site has been reached.
16. When ready to infuse, withdraw the guidewire completely from the micro catheter. Connect a syringe containing infusate to the hub of the micro catheter, and infuse as required.
17. After completing the procedure, withdraw the micro catheter and discard.

**NOTE:** To facilitate catheter handling, the proximal portion of the micro catheter is uncoated to ensure a non-slip grip.

## STEAM SHAPING MANDREL

In order to maintain catheter integrity and dimensional stability of the inner diameter, it is recommended that the user follow these instructions.

### WARNINGS

Shaping Mandrel is not intended for use in the human body.

Use only a steam source to shape the catheter tip. Do not use other heat sources.

Prior to use, inspect the catheter tip for any damage that may have resulted from shaping. Do not use a catheter that has been damaged in any way. Damaged catheters may rupture causing vessel trauma or tip detachment during steering maneuvers.

Remove shaping mandrel from card and insert into distal tip of the catheter.

18. Carefully bend catheter tip and shaping mandrel to desired shape. A slight over exaggeration of the shape may be required to accommodate for catheter relaxation.
19. Shape the catheter by holding the shaped portion approximately 1 inch (2.5 cm not less than 1 cm) from the steam source for about 20 seconds (DO NOT EXCEED 30 SECONDS.)
20. Allow catheter tip to cool in air or saline prior to removing mandrel.
21. Remove mandrel from catheter and discard. Multiple shaping is not recommended.
22. Inspect the tip for any damage that may have resulted from steam shaping the catheter tip. If any damage is found, do not use the catheter.

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Micro Therapeutics, Inc.