INSTRUCTIONS FOR USE

Navien™ A+ Intracranial Support Catheter

This page is intentionally left blank.

TABLE OF CONTENTS

Navien™ A+ Intracranial Support Catheter

English	4
Français.	5
Deutsch.	6
Italiano	7
Español	8
Svenska.	9
Nederlands.	10
Português	11
Suomi.	12
Dansk	13
Ελληνικά	14
Česky	15
Magyar	16
Русский	17
Polski	18
Türkçe	19
Norsk	20
Slovenčina	21
Română.	22
Български	23
한국어.	24
عربي	25
Symbol Glossary	26

Navien™ A+ Intracranial Support Catheter

DEVICE DESCRIPTION

The Navien™ A+ Intracranial Support Catheter is a single lumen, flexible, variable stiffness composite catheter. The catheter shaft has a hydrophilic coating to reduce friction during use. The Navien A+ Intracranial Support Catheter shaft is visible under fluoroscopy. The Navien A+ Intracranial Support Catheter dimensions are included on the individual device label. The Navien A+ Intracranial Support Catheter inner lumen can accommodate guide wires up to 0.038 inches in diameter to aid in placement of the catheter system.

The proximal end of the Navien A+ Intracranial Support Catheter has a luer fitting to allow attachment of accessories and infusion of liquids through the system. The Navien A+ Intracranial Support Catheter is offered in various sizes to accommodate physician preferences and anatomical variations. The catheter is provided sterile, non-pyrogenic, and is intended for single use only.

INDICATIONS FOR USE

The Navien A+ Intracranial Support Catheter is indicated for the introduction of interventional / diagnostic devices into the peripheral, coronary, and neuro vasculature.

The Navien A+ Intracranial Support Catheter is also indicated for the removal/ aspiration of fresh, soft emboli and thrombi from selected blood vessels in the arterial system, including neurovasculature.

CONTRAINDICATIONS

There are no known contraindications.

COMPATIBILITY

Refer to product label for device dimensions. Refer to labeling provided with other medical technologies to determine compatibility.

WARNINGS

- The Navien A+ Intracranial Support Catheter should only be used by physicians who have received appropriate training in interventional techniques.
- Do not reuse. The device is intended for single use only. Discard the catheter after one procedure. Structural integrity and/or function may be impaired through reuse or cleaning. Catheters are extremely difficult to clean after exposure to biological materials and may cause adverse patient reactions if reused.

PRECAUTIONS

- Store in a cool, dry, dark place.
- Do not use kinked or damaged devices.
- Do not use open or damaged packages.
- Use prior to the "Use By" date.
- Exposure to temperatures above 54°C (130°F) may damage device and accessories. Do not autoclave.
- Use the Navien A+ Intracranial Support Catheter in conjunction with fluoroscopic visualization.
- Inspect the catheter system before use to verify that its size and condition are suitable for the specific procedure.

- Do not advance or withdraw Navien A+ Intracranial Support Catheter against resistance without careful assessment of the cause using fluoroscopy. If the cause cannot be determined, withdraw the device. Moving or torquing the device against resistance may result in damage to the vessel or device.
- Maintain a constant infusion of appropriate flush solution.
- If flow through the catheter becomes restricted, do not attempt to clear the lumen by infusion. Remove and replace the device.
- Extreme care must be taken to avoid damage to the vasculature through which the catheter passes. The catheter may occlude smaller vessels. Care must be taken to avoid complete blood flow blockage.
- Torquing the catheter may cause damage which could result in kinking and possible separation along the catheter shaft. Should the system become severely kinked, withdraw the entire system (Navien A+ Intracranial Support Catheter, guide wire and catheter sheath introducer).

POTENTIAL ADVERSE EVENTS

Possible complications include, but are not limited to, the following:

- acute occlusion
- air embolism
- death
- distal embolization
- emboli
- false aneurysm formation
- hematoma or hemorrhage at puncture site
- infection
- intracranial hemorrhage
- ischemia
- neurological deficits including stroke
- vessel spasm, thrombosis, dissection or perforation

RECOMMENDED PREPARATION AND USE

1. Select an appropriate sized Navien A+ Intracranial Support Catheter based on the anatomy and length.
2. Flush package hoop with saline before removing catheter.
3. Gently remove the Navien A+ Intracranial Support Catheter and packaging card from the pouch by grasping the packaging card and slowly pulling the catheter and card out of the pouch.
4. Remove the Navien A+ Intracranial Support Catheter from the packaging card by carefully removing the hub from the card tabs before gently removing the Navien A+ Intracranial Support Catheter components from the protective sheath held by the card tabs.
5. Inspect the Navien A+ Intracranial Support Catheter for kinks or other damage. If any damage is observed, replace with new device.
6. Connect a rotating hemostasis valve (RHV) to the hub of the Navien A+ Intracranial Support Catheter and flush the inner lumen with heparinized saline.
7. Introduce the catheter into the vasculature through the catheter sheath introducer, and over a guide wire up to 0.038 inches (0.97 mm) in diameter using a percutaneous entry technique of choice. A split sheath may be provided in the package to provide support and facilitate the introduction of the Navien A+ Intracranial Support Catheter distal tip into the catheter sheath introducer.
8. Under fluoroscopic guidance, advance the system over the guide wire until the desired position is attained.
9. Remove the guide wire prior to the introduction of other intravascular devices.

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Micro Therapeutics, Inc.