CONQUEST® 40

PTA Dilatation Catheter

ENGLISH

Instructions For Use

CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

Device Description

The CONQUEST® 40 PTA Dilatation Catheter is a high performance balloon catheter consisting of an over-the-wire catheter with a balloon fixed at the distal tip. The proprietary ultra non-compliant, low profile balloon is designed to provide consistent balloon diameters and lengths even at high pressures. Two radiopaque markers delineate the working length of the balloon and aid in balloon placement. The coaxial catheter includes a tapered atraumatic tip to facilitate advancement of the catheter to and through the stenosis. The proximal portion of the catheter includes a female luer lock hub connected to the inflation lumen, and a female luer-lock hub connected to the guidewire lumen. The over-the-wire catheter is compatible with 0.035" guidewires and is available in 50cm and 75cm working lengths.

Packaged with every product is a profile reducing sheath that is positioned over the balloon for protection before use. A re-wrapping tool is also provided on the catheter shaft. A stylet is placed into the tip of the catheter to aid in re-wrap/refolding of the balloon. This product is not manufactured with any latex.

Indications For Use

CONQUEST® 40 PTA Dilatation Catheter is recommended for use in Percutaneous Transluminal Angioplasty of the femoral, iliac, and renal arteries and for the treatment of obstructive lesions of native or synthetic arteriovenous dialysis fistulae. This device is also recommended for post-dilatation of stents and stent grafts in the peripheral vasculature. This catheter is not for use in coronary arteries.

Contraindications

None known.

Warnings

1. Contents supplied **STERILE** using ethylene oxide (EO). Non-Pyrogenic. Do not use if sterile barrier is opened or damaged. Single patient use only. Do not reuse, reprocess or re-sterilize.
2. This device has been designed for single use only. Reusing this medical device bears the risk of cross-patient contamination as medical devices – particularly those with long and small lumina, joints, and/or crevices between components – are difficult or impossible to clean once body fluids or tissues with potential pyrogenic or microbial contamination have had contact with the medical device for an indeterminate period of time. The residue of biological material can promote the contamination of the device with pyrogens or microorganisms which may lead to infectious complications.
3. Do not re-sterilize. After re-sterilization, the sterility of the product is not guaranteed because of an indeterminate degree of potential pyrogenic or microbial contamination which may lead to infectious complications. Cleaning, reprocessing and/or re-sterilization of the present medical device increases the probability that the device will malfunction due to potential adverse effects on components that are influenced by thermal and/or mechanical changes.
4. To reduce the potential for vessel damage, the inflated diameter and length of the balloon should approximate the diameter and length of the vessel just proximal and distal to the stenosis.
5. To reduce the potential for stent or stent graft damage and/or vessel damage from the stent or stent graft, the diameter of the balloon should be no greater than the diameter of the stent or

stent graft. Refer to the stent or stent graft IFU for safety information including the **WARNINGS, PRECAUTIONS** and potential **ADVERSE EFFECTS** regarding the use of balloon post-dilatation.

6. When the catheter is exposed to the vascular system, it should be manipulated while under high-quality fluoroscopic observation. Do not advance or retract the catheter unless the balloon is fully deflated. If resistance is met during manipulation, determine the cause of the resistance before proceeding. Applying excessive force to the catheter can result in tip breakage or balloon separation.
7. Do not exceed the RBP recommended for this device. Balloon rupture may occur if the RBP rating is exceeded. To prevent over-pressurization, use of a pressure monitoring device is recommended or the use of indicated syringes.
8. After use, this product may be a potential biohazard. Handle and dispose of in accordance with acceptable medical practices and applicable local, state and federal laws and regulations.

Precautions

1. Carefully inspect the catheter prior to use to verify that catheter has not been damaged during shipment and that its size, shape and condition are suitable for the procedure for which it is to be used. Do not use if product damage is evident.
2. The CONQUEST® 40 catheter shall only be used by physicians trained in the performance of Percutaneous Transluminal Angioplasty.
3. The minimal acceptable sheath French size is printed on the package label. Do not attempt to pass the PTA catheter through a smaller size sheath introducer than indicated on the label.
4. Do not remove the guidewire in situ to shoot contrast through the wire lumen or perform a wire exchange. If the wire is removed while the balloon catheter is situated in tortuous anatomy, the risk of kinking the catheter is increased.
5. Use the recommended balloon inflation medium (a range of 30-50% contrast medium/a range of 50-70% sterile saline solution). It has been shown that a 30/70% contrast / saline ratio has yielded faster balloon inflation / deflation times.
6. Never use air or other gaseous medium to inflate the balloon.
7. If resistance is felt during post procedure withdrawal of the catheter through the introducer sheath, determine if contrast is trapped in the balloon with fluoroscopy. If contrast is present, push the balloon out of the sheath and then completely evacuate the contrast before proceeding to withdraw the balloon.
8. If resistance is still felt during post procedure withdrawal of the catheter, it is recommended to remove the balloon catheter and guidewire/ introducer sheath as a single unit.
9. Do not continue to use the balloon catheter if the shaft has been bent or kinked.
10. Prior to re-insertion through the introducer sheath, the balloon should be wiped clean with gauze, rinsed with sterile normal saline, and refolded with the balloon re-wrap tool. Balloon re-wrapping should only occur while the balloon catheter is supported with a guidewire or stylet.

Potential Adverse Reactions

The complications which may result from a peripheral balloon dilatation procedure include:

- Additional intervention
- Allergic reaction to drugs or contrast medium
- Aneurysm or pseudoaneurysm
- Arrhythmias
- Embolization
- Hematoma
- Hemorrhage, including bleeding at the puncture site
- Hypotension/hypertension
- Inflammation
- Occlusion
- Pain or tenderness
- Pneumothorax or hemothorax
- Sepsis/infection

- Shock
- Short term hemodynamic deterioration
- Thrombosis
- Vessel dissection, perforation, rupture, or spasm

Directions For Use

Handling & Storage

Store in a cool, dry, dark place. Do not store near radiation or ultra-violet light sources. Rotate inventory so that the catheters and other dated products are used prior to the "Use By" date. Do not use if packaging is damaged or opened.

Equipment For Use

- Contrast medium
- Sterile saline solution
- Inflation device with manometer (10ml or larger) OR
- 1 cc, 3 cc and/or 10 cc luer lock syringe (based on labeling)
- 3-way stopcock
- Appropriate introducer sheath and dilator set
- 0.035" guidewire

*Not all syringes and stopcocks are designed to hold pressures that are generated with angioplasty. Acrylic devices may avoid leakage when pressurized.

Dilatation Catheter Preparation

1. Remove Catheter from package.
2. Verify the balloon size is suitable for the procedure and the selected accessories accommodate the catheter as labeled.
3. Remove the balloon guard and stylet by grasping the balloon catheter just proximal to the balloon and with the other hand, gently grasp the balloon protector and slide distally off of the balloon catheter.
4. Slide the re-wrap tool to the proximal end of the catheter shaft.
5. Prior to use, the air in the balloon catheter should be removed. To facilitate purging, select a syringe or inflation device with a 10 ml or larger capacity and fill approximately half of it with the appropriate balloon inflation medium (30-50% contrast medium/50-70% sterile saline solution). Do not use air or any gaseous medium to inflate the balloon.
6. Connect a stopcock to the balloon inflation female luer hub on the dilatation catheter.
7. Connect the syringe to the stopcock.
8. Hold the syringe with the nozzle pointing downward, open the stopcock and aspirate for approximately 15 seconds. Release the plunger.
9. Repeat step #8 two more times or until bubbles no longer appear during aspiration (negative pressure). Once completed, evacuate all air from the barrel of the syringe/inflation device.
10. Prepare the wire lumen of the catheter by attaching a syringe to the wire lumen hub and flushing with sterile saline solution.

Use of the CONQUEST® 40 PTA Dilatation Catheter with an Inflation Device

1. Backload the distal tip of the CONQUEST® 40 balloon catheter over the pre-positioned guidewire and advance the tip to the introduction site.
2. Advance the catheter through the introducer sheath and over the wire to the site of inflation. If the stenosis cannot be crossed with the desired dilatation catheter, use a smaller diameter catheter to pre-dilate the lesion to facilitate passage of a more appropriately sized dilatation catheter.
3. Position the balloon relative to the lesion to be dilated, ensure the guidewire is in place and inflate the balloon to the appropriate pressure.
4. If greater pressure is required than the inflation device offers, close the stopcock to the inflation device and open the stopcock to the indicated 1 cc or 3 cc syringe. Apply pressure. Close the stopcock to hold the pressure if needed.
5. Open the stopcock to the inflation device and apply negative pressure to fully evacuate fluid from the balloon. Confirm that the balloon is fully deflated under fluoroscopy.
6. While maintaining negative pressure and the position of the guidewire, withdraw the deflated dilatation catheter over the wire through the introducer sheath. Use of a gentle clockwise motion may be used to help facilitate catheter removal through the introducer sheath.

Use of the CONQUEST® 40 PTA Dilatation Catheter with Syringe Assembly

1. Backload the distal tip of the CONQUEST® 40 balloon catheter over the pre-positioned guidewire and advance the tip to the introduction site.
2. Advance the catheter through the introducer sheath and over the wire to the site of inflation. If the stenosis cannot be crossed with the desired dilatation catheter, use a smaller diameter catheter to pre-dilate the lesion to facilitate passage of a more appropriately sized dilatation catheter.
3. Position the balloon relative to the lesion to be dilated, ensure the guidewire is in place.
4. Prepare syringe assembly by half filling a 10 cc and the indicated 1 cc or 3 cc syringe with dilute contrast (30-50% contrast medium/50-70% sterile saline solution) .
5. Attach the 10 cc syringe to the in-line port of a three way stopcock. Attach the 1 cc or 3 cc syringe to the side port.
6. Attach the three-way stopcock to the balloon.
7. Apply pressure with the 10 cc syringe to inflate the balloon.
8. Turn the valve of the stopcock to access the 1 cc or 3 cc syringe and inflate to the appropriate pressure.
9. Apply negative pressure with the 10 cc syringe to fully evacuate fluid from the balloon. Confirm the balloon is fully deflated under fluoroscopy.
10. While maintaining negative pressure and the position of the guidewire, withdraw the deflated dilatation catheter over the wire through the introducer sheath. Use of a gentle counterclockwise motion may be used to help facilitate catheter removal through the introducer sheath.

Balloon Reinsertion

Precaution: Do not continue to use the balloon catheter if the shaft has been bent or kinked.

Precaution: Prior to re-insertion through the introducer sheath, the balloon should be wiped clean with gauze, rinsed with sterile normal saline, and refolded with the balloon re-wrap tool. Balloon re-wrapping should only occur while the balloon catheter is supported with a guidewire or stylet.

1. Load the balloon catheter onto a guidewire or stylet.
2. Advance the balloon re-wrap tool over the catheter to the proximal end of the balloon.
3. Grasp the catheter shaft just proximal to the balloon with one hand, and with the other hand gently slide the rewrap tool over the balloon to the catheter tip and then back over the balloon to the catheter.
4. Slide the re-wrap tool to the proximal end of the catheter shaft and if necessary, remove the stylet.
5. Advance the balloon catheter over the pre-positioned guidewire to the introduction site and through the introducer sheath. If resistance is encountered, replace the previously used balloon catheter with a new balloon.
6. Continue the procedure according to the "Use of the CONQUEST® 40 Catheter" section herein.

Warning: After use, this product may be a potential biohazard. Handle and dispose of in accordance with acceptable medical practices and applicable local, state and federal laws and regulations.

Warranty

Bard Peripheral Vascular warrants to the first purchaser of this product, that this product will be free from defects in materials and workmanship for a period of one year from the date of first purchase and liability under this limited product warranty will be limited, to repair or replacement of the defective product, in Bard Peripheral Vascular's sole discretion, or refunding your net price paid. Wear and tear from normal use or defects resulting from misuse of this product are not covered by this limited warranty.

TO THE EXTENT ALLOWABLE BY APPLICABLE LAW, THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL BARD PERIPHERAL VASCULAR BE LIABLE TO YOU FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM YOUR HANDLING OR USE OF THIS PRODUCT.

Some countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your country.

An issue or revision date and revision number for these instructions are included for the user's information on the last page of this booklet.

In the event 36 months have elapsed between this date and product use, the user should contact Bard Peripheral Vascular to see if additional product information is available.

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**Do Not Use if the Product Sterile Barrier System or its Packaging is Compromised**

Ne pas utiliser si le système de barrière stérile du produit ou son emballage sont ou paraissent endommagés.
 Bei beeinträchtigtem Sterilschutzsystem oder beschädigter Verpackung das Produkt nicht verwenden
 Non utilizzare se il sistema della barriera sterile del prodotto o la confezione sono compromessi
 No utilizar si el sistema de barrera estéril del producto o su envase están dañados
 Niet gebruiken wanneer het steriele-barrièresysteem van het product of de verpakking is aangetast
 Não utilizar se o sistema de barreira estéril do produto ou respectiva embalagem estiverem comprometidos
 Μη χρησιμοποιείτε το προϊόν εάν έχει παραβιαστεί το σύστημα φραγμού αποστείρωσης ή η συσκευασία του
 Má ikke anvendes, hvis produktets sterilbarrieresystem eller emballagen er kompromitteret
 Använd inte om produktens sterila barriär eller förpackning är skadad
 Ei saa käyttää, jos tuotteen steriilisuojaus tai pakkaus on vaurioitunut.
 Skal ikke brukes hvis produktets sterile barrieresystem eller pakningen er forringet
 Nie używać, jeśli naruszono sterylność produktu lub jego opakowanie
 Ne használja, ha a termék védő steril lezárási rendszer vagy a csomagolás sérült!
 Pokud je narušena sterilni bariéra nebo poškozen obal, výrobek nepoužívejte
 Ürünün Steril Bariyer Sistemi ya da Ambalajı Bozulmuşsa Ürünü Kullanmayın
 Запрещается применять изделие, если его стерильная упаковка или внешняя упаковка повреждены
 如果產品消毒屏障系統或包裝受損，請勿使用

**Bard and Conquest are trademarks and/or registered trademarks of C. R. Bard Inc. or an affiliate**

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**Copyright ©2014 C. R. Bard, Inc. All rights reserved. Printed in the USA.**

Copyright ©2014 C. R. Bard, Inc. Tous droits réservés. Imprimé aux États-Unis.
 Copyright ©2014 C. R. Bard, Inc. Alle Rechte vorbehalten. Gedruckt in USA.
 Copyright ©2014 C. R. Bard, Inc. Tutti i diritti riservati. Stampato negli Stati Uniti.
 Copyright ©2014, C. R. Bard, Inc. Reservados todos los derechos. Impreso en EE.UU.
 Copyright ©2014 C. R. Bard, Inc. Alle rechten voorbehouden. Gedrukt in de VS.
 Copyright ©2014, C. R. Bard, Inc. Todos os direitos reservados. Impresso nos E.U.A.
 Πνευματικά Δικαιώματα ©2014 C. R. Bard, Inc. Με την επιφύλαξη παντός δικαιώματος. Εκτυπώθηκε στις Η.Π.Α.
 Copyright ©2014 C. R. Bard, Inc. Alle rettigheder forbeholdes. Trykt i USA.
 Copyright ©2014 C. R. Bard, Inc. Alla rättigheter förbehållna. Tryckt i USA.
 Copyright ©2014 C. R. Bard, Inc. Kaikki oikeudet pidätetään. Painettu Yhdysvalloissa.
 Copyright ©2014 C. R. Bard, Inc. Alle rettigheter forbeholdt. Trykt i USA.
 Copyright ©2014 C. R. Bard, Inc. Wszelkie prawa zastrzeżone. Wydrukowano w USA.
 Copyright ©2014, C. R. Bard, Inc. Minden jog fenntartva. Nyomtatva az USA-ban.
 Copyright ©2014 C. R. Bard, Inc. Všechna práva vyhrazena. Vytisťeno v USA.
 Yapıt Hakkı ©2014 C. R. Bard, Inc. Tüm hakları saklıdır. ABD'de basılmıştır.
 Авторские права ©2014 Компания C. R. Bard, Inc. Все права защищены. Напечатано в США.
 Copyright ©2014 C. R. Bard, Inc. 保留所有權利。在美國印製。
 ©2014 C. R. Bard, Inc.

**Authorised Representative in the European Community**

Représentant autorisé au sein de la Communauté européenne
 Bevollmächtigter in der Europäischen Gemeinschaft
 Rappresentante autorizzato nella Comunità Europea
 Representante autorizado en la Comunidad Europea
 Gemachtigde vertegenwoordiger binnen de Europese Gemeenschap
 Representante autorizado na Comunidade Europeia
 Εξουσιοδοτημένος Αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα
 Autoriseret repræsentant i Det Europæiske Fællesskab
 Auktoriserad representant i EG
 Valtuutettu edustaja Euroopan yhteisössä
 Autorisert representant i EU
 Autoryzowany przedstawiciel w Unii Europejskiej
 Bejegyzett képviselő az Európai Közösségben
 Autorizovaný zástupce pro Evropské společenství
 Avrupa Topluluğu Yetkili Temsilcisi
 Уполномоченный представитель в Европейском сообществе
 歐洲共同體的授權代表

Manufacturer:
Bard Peripheral Vascular, Inc.
1625 West 3rd Street
Tempe, AZ 85281 USA
Tel: 1-480-894-9515
1-800-321-4254
Fax: 1-480-966-7062
1-800-440-5376
bardpv.com

**Authorised Representative
in the European Community**
Bard Limited
Forest House
Tilgate Forest Business Park
Brighton Road, Crawley
West Sussex
RH11 9BP, UK

**BARID | PERIPHERAL
VASCULAR**