# VASCU<sup>TRAK</sup><sup>®</sup>

## PTA Dilatation Catheter

**ENGLISH**

### INSTRUCTIONS FOR USE

**CAUTION:** Federal (U.S.A.) Law restricts this device to sale by or on the order of a physician.

#### Device Description

The VASCU<sup>TRAK</sup><sup>®</sup> PTA Dilatation Catheter consists of flexible shaft with a semi-compliant balloon fixed at the distal end and a 0.018" or 0.014" guidewire lumen through the distal tip. The catheter shaft contains an inflation lumen that begins at the proximal female luer lock hub and ends in the proximal portion of the balloon. The proximal portion of the shaft is comprised of a stainless steel tube while the distal portion of the shaft contains a stainless steel core wire that parallels the inflation lumen, exits the shaft proximal to the balloon and terminates distal to the balloon. Two radiopaque marker bands located on the outer core wire delineate the working length of the balloon to aid in balloon placement. The VASCU<sup>TRAK</sup><sup>®</sup> PTA Dilatation Catheter is available in two catheter lengths, 140 cm and 80 cm. The stated catheter length is the usable length of the catheter. The 140 cm catheter shaft includes two non-radiopaque markers located 90 cm and 100 cm from the balloon catheter tip to help confirm when the balloon catheter tip exits guiding catheters or sheaths. The 80 cm catheter shaft does not contain visual markers on the shaft.

Packaged with every product is a profile reducing sheath that is positioned over the balloon for protection before use.

These products are not made with natural rubber latex.

**Balloon Compliance Table**

Pressure ATM	Balloon Diameter (mm)							
	2.0	2.5	3.0	3.5	4.0	5.0	6.0	7.0
2.0	1.71	1.93	2.32	2.60	3.39	4.25	5.01	5.88
4.0	1.94	2.13	2.53	2.92	3.89	4.61	5.41	6.33
6.0	2.08	2.31	2.71	3.15	4.17	4.96	5.85	6.83
8.0	2.17	2.49	2.89	3.39	4.33	5.28	6.24	7.29
10.0	2.23	2.63	3.04	3.60	4.44	5.58	6.54	7.62
*12.0	2.28	2.74	3.16	3.77	4.49	5.81	6.81	7.96

**Nominal Balloon Pressure of 8.0 ATM for balloon diameter sizes 2.0 – 3.5 mm and 6.0 ATM for balloon diameters 4.0 – 7.0 mm.**

\*Rated Burst Pressure.

#### Indications for Use

The VASCU<sup>TRAK</sup><sup>®</sup> PTA Dilatation Catheter is intended to dilate stenoses in the iliac, femoral, ilio-femoral, popliteal, infra-popliteal and renal arteries and for the treatment of obstructive lesions of native or synthetic arteriovenous dialysis fistulae. This device is also recommended for post dilatation of balloon expandable stents, self-expanding stents, and stent grafts in the peripheral vasculature.

#### Contraindications

The VASCU<sup>TRAK</sup><sup>®</sup> PTA Catheter is contraindicated:

- where there is the inability to cross the target lesion with a guidewire
- for use in the coronary or neuro vasculature

#### Warnings

1. Contents supplied **STERILE** using ethylene oxide (EO). **Non-Pyrogenic. Do not use if sterile barrier is opened or damaged. Do not reuse, reprocess or re-sterilize. Use the catheter prior to the "Use By" date specified on the package label.**
2. This device has been designed for single use only. Reusing this medical device bears the risk of cross-patient contamination as medical devices – particularly those with long and small lumina, joints, and/or crevices between components – are difficult or impossible to clean once body fluids or tissues with potential pyrogenic or microbial contamination have had contact with the medical device for an indeterminable period of time. The residue of biological material can promote the contamination of the device with pyrogens or microorganisms which may lead to infectious complications.
3. Do not resterilize. After resterilization, the sterility of the product is not guaranteed because of an indeterminable degree of potential pyrogenic or microbial contamination which may lead to infectious complications. Cleaning, reprocessing and/or resterilization of the present medical device increases the probability that the device will malfunction due to potential adverse effects on components that are influenced by thermal and/or mechanical changes.
4. To reduce the potential for vessel damage or difficulty in deflating, the inflated diameter and length of the balloon should approximate the diameter and length of the vessel just proximal and distal to the stenosis.

5. When the catheter is exposed to the vascular system, it should be manipulated while under high-quality fluoroscopic observation. Do not advance or retract the catheter unless the balloon is fully deflated. If resistance is met during manipulation, determine the cause of the resistance before proceeding. Applying excessive force to the catheter can result in tip breakage or balloon separation.
6. Do not exceed the RBP recommended for this device. Balloon rupture or difficulty in deflation may occur if the RBP rating is exceeded. To prevent over pressurization, use of a pressure monitoring device is recommended.
7. After use, this product may be a potential biohazard. Handle and dispose of in accordance with acceptable medical practices and applicable local, state and federal laws and regulations.

#### Precautions

1. Carefully inspect the catheter prior to use to verify that catheter has not been damaged during shipment and that its size, shape and condition are suitable for the procedure for which it is to be used. Do not use if product damage is evident.
2. The VASCU<sup>TRAK</sup><sup>®</sup> PTA Balloon Dilatation Catheter should only be used by physicians experienced in the performance of percutaneous transluminal angioplasty.
3. In order to activate the hydrophilic coating, wet the VASCU<sup>TRAK</sup><sup>®</sup> balloon and catheter with sterile saline or wipe the balloon catheter with sterile saline saturated gauze immediately prior to its insertion in the body. Do not wipe the balloon catheter with dry gauze.
4. When back loading the catheter onto the guidewire, support the catheter and ensure that the guidewire tip does not snag or come into contact with the balloon.
5. The minimal acceptable sheath French size is printed on the package label. Do not attempt to pass the PTA catheter through a smaller size introducer sheath than indicated on the label.
6. Use the recommended balloon inflation medium (50% contrast medium/50% sterile saline solution). It has been shown that a 25/75% contrast / saline ratio has yielded faster balloon inflation / deflation times. Never use air or other gaseous medium to inflate the balloon.
7. VASCU<sup>TRAK</sup><sup>®</sup> PTA Dilatation Catheters should be used with caution for procedures involving calcified lesions, stents or synthetic vascular grafts due to the abrasive nature of these lesions.
8. Never attempt to move the guidewire when the balloon is inflated.
9. Fully evacuate the balloon prior to withdrawing the system. Larger sizes of VASCU<sup>TRAK</sup><sup>®</sup> Balloons may exhibit slower deflation times. If the balloon does not deflate, advance a sheath or catheter over the proximal portion of the balloon to straighten out the transition from connection of the balloon to the inflation lumen.
10. If resistance is felt during post procedure withdrawal of the catheter through the introducer sheath, determine if contrast medium is trapped in the balloon with fluoroscopy. If contrast is present, push the balloon out of the sheath and then completely evacuate the contrast before proceeding to withdraw the balloon.
11. If resistance is still felt during post procedure withdrawal of the catheter, it is recommended to remove the balloon catheter and guidewire/introducer sheath as a single unit.
12. Do not continue to use the balloon catheter if the shaft has been bent or kinked.
13. Prior to re-insertion through the introducer sheath, re-activate the hydrophilic coating and clean the balloon catheter by wiping the balloon catheter with sterile saline saturated gauze and rinsing with sterile saline. Do not wipe the balloon catheter with dry gauze.

#### Potential Adverse Reactions

The complications that may result from a peripheral balloon dilatation procedure include:

- Additional intervention
- Allergic reaction to drugs or contrast medium
- Aneurysm or pseudoaneurysm
- Arrhythmias
- Embolization
- Hematoma
- Hemorrhage, including bleeding at the puncture site
- Hypotension/hypertension
- Inflammation
- Occlusion
- Pain or tenderness
- Pneumothorax or hemothorax
- Sepsis/infection
- Shock
- Short term hemodynamic deterioration
- Stroke
- Thrombosis
- Vessel dissection, perforation, rupture, or spasm

#### Directions for Use

##### Handling & Storage

Store in a cool, dry, dark place. Do not store near radiation or ultra-violet light sources.

Rotate inventory so that the catheters and other dated products are used prior to the "Use By" date. Do not use if packaging is damaged or opened.

#### Equipment Required

- Contrast medium
- Sterile saline solution
- Luer lock syringe/inflation device with manometer (10 ml or larger)
- Appropriate introducer sheath and dilator set
- 0.014" or 0.018" guidewire

#### Dilatation Catheter Preparation

1. Remove catheter from package. Verify the balloon size is suitable for the procedure and the selected accessories accommodate the catheter as labeled.
2. Remove the balloon guard by grasping the balloon catheter just proximal to the balloon and with the other hand, gently grasp the balloon guard and slide distally off of the balloon catheter.
3. Prior to use, the air in the balloon catheter should be removed. To facilitate purging, select a syringe or inflation device with a 10 ml or larger capacity and fill approximately half of it with the appropriate balloon inflation medium (50% contrast medium/50% sterile saline solution). Do not use air or any gaseous medium to inflate the balloon. Note: it has been shown that a 25/75% contrast / saline ratio has yielded faster balloon inflation / deflation times.
4. Connect a stopcock to the balloon inflation female luer hub on the dilatation catheter.
5. Connect the syringe to the stopcock.
6. Hold the syringe with the nozzle pointing downward, open the stopcock and aspirate for approximately 15 seconds. Release the plunger.
7. Repeat step #6 two more times or until bubbles no longer appear during aspiration (negative pressure). Once completed, evacuate all air from the barrel of the syringe/inflation device.
8. In order to activate the coating, wet the balloon catheter with a sterile saline or wipe the balloon catheter with saline saturated gauze immediately prior to its insertion in the introducer sheath.

**Precaution:** Do not wipe the balloon catheter with dry gauze.

#### Use of the VASCU TRAK® PTA Balloon Dilatation Catheters

1. Backload the distal tip of the VASCU TRAK® PTA Balloon Dilatation Catheters over the pre-positioned guidewire and advance the tip to the introduction site by grasping the distal tip.

**Precaution:** When back loading the catheter onto the guidewire, support the catheter and ensure that the guidewire tip does not snag or come into contact with the balloon.

2. Using short advances, advance the catheter through the introducer sheath and over the wire to the site of inflation. (Note: If using a hydrophilic guidewire, ensure that it is kept hydrated with sterile normal saline, at all times.) If the stenosis cannot be crossed with the desired dilatation catheter, use a smaller diameter catheter to pre-dilate the lesion to facilitate passage of a more appropriately sized dilatation catheter.
3. Position the balloon relative to the lesion to be dilated, ensure the guidewire is in place and **slowly** inflate the balloon with an inflation device. It is recommended that the balloon pressure be increased 1 atmosphere every 30 seconds until the lesion is effaced or the rated burst pressure is reached.

**Precaution:** Do not inflate the balloon past its labeled rated burst pressure.

4. **Slowly** apply negative pressure to fully evacuate fluid from the balloon. Confirm that the balloon is fully deflated under fluoroscopy and that no contrast is left in the balloon.

**Precaution:** Fully evacuate the balloon prior to withdrawing the system. Larger sizes of VASCU TRAK® Balloons may exhibit slower deflation times. If the balloon does not deflate, advance a sheath or catheter over the proximal portion of the balloon to straighten out the connection of the balloon to the inflation lumen.

5. While maintaining negative pressure and the position of the guidewire, grasp the balloon catheter just outside the sheath and withdraw the deflated dilatation catheter over the wire through the introducer sheath.

#### Balloon Reinsertion

**Precaution:** Do not continue to use the balloon catheter if the shaft has been bent or kinked.

**Precaution:** Prior to re-insertion through the introducer sheath, re-activate the hydrophilic coating and clean the balloon catheter by wiping the balloon catheter with sterile saline-soaked gauze and rinsing with sterile saline. Do not wipe the balloon catheter with dry gauze.

1. Load the balloon catheter onto a guidewire ensuring that the guidewire tip does not snag or come into contact with the balloon.
2. Advance the balloon catheter over the pre-positioned guidewire to the introduction site and through the introducer sheath. If resistance is encountered, replace the previously used balloon catheter with a new balloon.
3. Continue the procedure according to the "Use of the VASCU TRAK® PTA Balloon Dilatation Catheter" section herein.

**Warning: After use, this product may be a potential biohazard. Handle and dispose of in accordance with acceptable medical practices and applicable local, state and federal laws and regulations.**

#### Warranty

Bard Peripheral Vascular warrants to the first purchaser of this product, that this product will be free from defects in materials and workmanship for a period of one year from the date of first purchase and liability under this limited product warranty will be limited, to repair or replacement of the defective product, at Bard Peripheral Vascular's sole discretion, or refunding your net price paid. Wear and tear from normal use or defects resulting from misuse of this product is not covered by this limited warranty.

**TO THE EXTENT ALLOWABLE BY APPLICABLE LAW, THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL BARD PERIPHERAL VASCULAR BE LIABLE TO YOU FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM YOUR HANDLING OR USE OF THIS PRODUCT.**

Some countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your country.

An issue or revision date and revision number for these instructions are included for the user's information on the last page of this booklet. In the event 36 months have elapsed between this date and product use, the user should contact Bard Peripheral Vascular to see if additional product information is available.

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**Do Not Use if the Product Sterilization Barrier or its Packaging is Compromised**

Ne pas utiliser si la barrière de stérilisation ou l'emballage du produit est endommagé  
Nicht verwenden, wenn die Sterilisationsabdeckung oder die Verpackung des Produkts beschädigt ist  
Non utilizzare se la barriera sterile del prodotto o la confezione è compromessa  
No utilice el producto si la barrera de esterilización o su envase no están en perfecto estado  
Niet gebruiken wanneer de sterilisatiebarrière van het product of de verpakking is aangetast  
Não utilizar se a barreira de esterilização do produto ou respectiva embalagem estiverem comprometidas  
Μη χρησιμοποιείτε το προϊόν εάν έχει παραβιαστεί ο φραγμός αποστείρωσης ή η συσκευασία του  
Må ikke anvendes, hvis produktets steriliseringsbarriere eller emballagen er kompromitteret  
Använd inte produkten om sterilsyddet är skadat eller förpackningen är bruten  
Nåä käytät tuotetta, jos sen steriloitinsuojaus tai pakkaus on vaurioitunut  
Bruk ikke produktet hvis den sterile barrieren eller emballasjen er brutt  
Nie używać, jeżeli naruszono sterylność produktu lub jego opakowanie  
Ne használja, ha a termékét védő steril zár vagy a csomagolás sérült  
Nepoužívejte, pokud je porušena sterilní bariéra nebo obal produktu  
Úrün Steril Bariyeri veya Ambalajı zarar görmüşse kullanmayın  
如果产品灭菌防护膜或其包装已损坏, 请勿使用

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**Bard and Vascutrak are trademarks and/or registered trademarks of C. R. Bard, Inc. or an affiliate.**

Bard et Vascutrak sont des marques et (ou) des marques déposées de C. R. Bard, Inc. ou d'une de ses filiales.  
Bard und Vascutrak sind Marken und/oder eingetragene Marken von C. R. Bard, Inc., oder einer Tochtergesellschaft.  
Bard e Vascutrak sono marchi commerciali e/o registrati di C. R. Bard, Inc. o di una sua affiliata.  
Bard y Vascutrak son marcas comerciales y/o comerciales registradas de C. R. Bard, Inc. o de una filial.  
Bard en Vascutrak zijn handelsmerken en/of gedeponeerde handelsmerken van C. R. Bard, Inc. of een daarmee geaffilieerd bedrijf.  
Bard e Vascutrak são marcas comerciais e/ou registradas da C. R. Bard, Inc. ou de uma empresa afiliada.  
Οι ονομασίες Bard και Vascutrak είναι εμπορικά σήματα ή και σήματα κατατεθέντα της C. R. Bard, Inc. ή κάποιας θυγατρικής της.  
A Bard és a Vascutrak er varemerker og/eller registrerede varemerker tilhørende C. R. Bard, Inc. eller et tilknyttet selskab.  
Bard och Vascutrak är varumärken/inregistrerade varumärken som tillhör C. R. Bard, Inc. eller affilierat företag.  
Bard ja Vascutrak ovat C. R. Bard, Inc:n tai sen tytäryhtiön tavaramerkkejä ja/tai rekisteröityjä tavaramerkkejä.  
Bard og Vascutrak er varemerker og/eller registrerte varemerker som tilhører C. R. Bard, Inc. eller et tilknyttet selskap.  
Bard i Vascutrak są znakami towarowymi i/lub zarejestrowanymi znakami towarowymi firmy C. R. Bard, Inc. lub jej spółek stowarzyszonych.  
A Bard és a Vascutrak a C. R. Bard, Inc. illetve társállalatának védjegye és/vagy bejegyzett védjegye.  
Bard a Vascutrak jsou ochranné známky a/nebo registrované ochranné známky společnosti C. R. Bard, Inc. nebo její pobočky.  
Bard ve Vascutrak, C. R. Bard, Inc.'in veya bağlı kuruluşunun ticari markaları ve/veya tescilli ticari markalarıdır.  
Bard和Vascutrak是C. R. Bard, Inc.及其附属机构的商标和/或注册商标。  
Bard Vascutrak C. R. Bard, Inc.  
Bard a Vascutrak sú obchodné značky a/alebo zaregistrované obchodné značky spoločnosti C. R. Bard, Inc. alebo jej partnera.

**Copyright © 2017 C. R. Bard, Inc. All rights reserved. Printed in the USA.**

Copyright © 2017 C. R. Bard, Inc. Tous droits réservés. Imprimé aux États-Unis.  
Copyright © 2017 C. R. Bard, Inc. Alle Rechte vorbehalten. Gedruckt in USA.  
Copyright © 2017 C. R. Bard, Inc. Tutti i diritti riservati. Stampato negli Stati Uniti.  
Copyright © 2017 C. R. Bard, Inc. Reservados todos los derechos. Impreso en EE.UU.  
Copyright © 2017 C. R. Bard, Inc. Alle rechten voorbehouden. Gedrukt in de VS.  
Copyright © 2017 C. R. Bard, Inc. Todos os direitos reservados. Impresso nos EUA.  
Πνευματικά Δικαιώματα © 2017 C. R. Bard, Inc. Με την επιφύλαξη παντός δικαιώματος. Εκτυπώθηκε στις Η.Π.Α.  
Copyright © 2017 C. R. Bard, Inc. Alle rettigheder forbeholdes. Trykt i USA.  
Copyright © 2017 C. R. Bard, Inc. Alla rättigheter förbehållna. Tryckt i USA.  
Copyright © 2017 C. R. Bard, Inc. Kaikki oikeudet pidätetään. Painettu Yhdysvalloissa.  
Copyright © 2017 C. R. Bard, Inc. Alle rettigheter forbeholdt. Trykt i USA.  
Copyright © 2017 C. R. Bard, Inc. Wszelkie prawa zastrzeżone. Wydrukowano w USA.  
Copyright © 2017 C. R. Bard, Inc. Minden jog fenntartva. Nyomatva az USA-ban.  
Copyright © 2017 C. R. Bard, Inc. Všechna práva vyhrazena. Vytlačeno v USA.  
Yapıt Hakkı © 2017 C. R. Bard, Inc. Tüm hakları saklıdır. ABD'de basılmıştır.  
版权所有 © 2017 C. R. Bard, Inc. 保留一切权利。于美国印制。  
© 2017 C. R. Bard, Inc.  
Copyright © 2017 C. R. Bard, Inc. Všetky práva vyhrazené. Vytlačené v USA.

**Authorised Representative in the European Community**

Représentant autorisé au sein de la Communauté européenne  
Bevollmächtigter in der Europäischen Gemeinschaft  
Rappresentante autorizzato nella Comunità Europea  
Representante autorizado en la Comunidad Europea  
Gemachtigde binnen de Europese Gemeenschap  
Representante autorizado na Comunidade Europeia  
Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα  
Autoriseret repræsentant i Det Europæiske Fællesskab  
Auktoriserad representant inom EG  
Valtuutettu edustaja Euroopan yhteisössä  
Autorisert representant i EU  
Autoryzowany przedstawiciel na terenie Unii Europejskiej  
Meghatalmazott képviselő az Európai Közösségben  
Autorizovaný zástupce v Evropské unii  
Avrupa Topluluğu Yetkilii Temsilcisi  
欧共同体授权代表

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**Manufacturer:**  
**Bard Peripheral Vascular, Inc.**  
1625 West 3rd Street  
Tempe, AZ 85281  
USA

TEL: 1-480-894-9515  
1-800-321-4254  
FAX: 1-480-966-7062  
1-800-440-5376  
[www.bardpv.com](http://www.bardpv.com)

**Authorised Representative  
in the European Community**  
Bard Limited  
Forest House  
Tilgate Forest Business Park  
Brighton Road, Crawley  
West Sussex, RH11 9BP, UK