ULTRAVERSE® 014/018

PTA Balloon Dilatation Catheter

ENGLISH

INSTRUCTIONS FOR USE

CAUTION: Federal (U.S.A.) Law restricts this device to sale by or on the order of a physician.

Device Description

The ULTRAVERSE® 014 and 018 PTA Balloon Dilatation Catheters are small vessel balloon catheters consisting of an over the wire co-axial lumen catheter with a balloon fixed at the distal tip. For all balloon lengths, radiopaque markers delineate the working length of the balloon and aid in balloon placement. For balloon lengths of 100mm and greater, two radiopaque markers are positioned on the distal portion of the balloon and one radiopaque marker is positioned on the proximal portion of the balloon to differentiate between the distal and proximal ends of the balloon. The catheters include an atraumatic tip to facilitate advancement of the catheter to and through the stenosis. To facilitate catheter advancement through the vasculature and the vessel stenosis, the ULTRA-CROSS® Dual Layer Hydrophilic Coating is present on the distal segment of the shaft and the balloon. ULTRAVERSE® 014 and 018 are available in multiple shaft lengths. ULTRAVERSE® 014 is compatible with 0.014" guidewires, while ULTRAVERSE® 018 is compatible with 0.014" and 0.018" guidewires. The proximal portion of the catheters include a female luer lock hub connected to the inflation lumen, and a female luer lock hub connected to the guidewire lumen.

Packaged with every product is a profile reducing sheath that is positioned over the balloon for protection before use. A stylet is placed into the tip of the catheter to aid in rewrap/refolding of the balloon. These products are not made with natural rubber latex.


Figure 1

The GEOALIGN® Marking System is a non-radiopaque ruler on the catheter shaft measured from the distal tip. The GEOALIGN® markings are designated on the catheter shaft by 1cm increment bands with an accuracy within ± 1 mm. The distance from the distal catheter tip is labeled in 10cm increments. Thicker bands denote the midway point (5cm) between the labeled distances. The GEOALIGN® Marking System is designed to be used as a tool to externally measure the intravascular advancement and/or retraction of the catheter. This can provide an intravascular reference regarding the location of the distal tip of the catheter or an approximate intravascular length measurement between two points. The GEOALIGN® Marking System may also facilitate geographic alignment of an adjunctive therapy that includes the same GEOALIGN® Marking System.

NOTE: The GEOALIGN® Marking System provides an approximation that may not be an exact representation of the actual distance traveled intravascularly and should be confirmed under fluoroscopy.


Figure 2

NOTE: The GEOALIGN® Marking System includes non-radiopaque white markings and are designed to be utilized outside the sheath.

Indications for Use

ULTRAVERSE® 014 and ULTRAVERSE® 018 PTA Dilatation Catheters are recommended for use in percutaneous transluminal angioplasty (PTA) of the renal, popliteal, tibial, femoral, and peroneal arteries. These catheters are not for use in coronary arteries.

Contraindications

None known.

Warnings

1. Contents supplied **STERILE** using ethylene oxide (EO). Non-pyrogenic. Do not use if sterile barrier is opened or damaged. Do not reuse, reprocess or sterilize.
2. This device has been designed for single use only. Reusing this medical device bears the risk of cross-patient contamination as medical devices – particularly those with long and small lumina, joints, and/or crevices between components – are difficult or impossible to clean once body fluids or tissues with potential pyrogenic or microbial contamination have had contact with the medical device for an indeterminate period of time. The residue of biological material can promote the contamination of the device with pyrogens or microorganisms which may lead to infectious complications.

3. Do not resterilize. After resterilization, the sterility of the product is not guaranteed because of an indeterminable degree of potential pyrogenic or microbial contamination which may lead to infectious complications. Cleaning, reprocessing and/or resterilization of the present medical device increases the probability that the device will malfunction due to potential adverse effects on components that are influenced by thermal and/or mechanical changes.
4. To reduce the potential for vessel damage, the inflated diameter and length of the balloon should approximate the diameter and length of the vessel just proximal and distal to the stenosis.
5. When the catheter is exposed to the vascular system, it should be manipulated while under high-quality fluoroscopic observation. Do not advance or retract the catheter unless the balloon is fully deflated. If resistance is met during manipulation, determine the cause of the resistance before proceeding. Applying excessive force to the catheter can result in tip breakage or balloon separation.
6. Do not exceed the RBP recommended for this device. Balloon rupture may occur if the RBP rating is exceeded. To prevent over pressurization, use of a pressure monitoring device is recommended.
7. After use, this product may be a potential biohazard. Handle and dispose of in accordance with acceptable medical practices and applicable local, state and federal laws and regulations.
8. The safety and effectiveness of the device has not been established, or is unknown, in vascular regions other than those specifically indicated.

Precautions

1. Refer to accessory IFU for potential access site warnings, precautions, and adverse events.
2. Carefully inspect the catheter prior to use to verify that catheter has not been damaged during shipment and that its size, shape and condition are suitable for the procedure for which it is to be used. Do not use if product damage is evident.
3. The ULTRAVERSE® 014 and ULTRAVERSE® 018 PTA Balloon Dilatation Catheters shall only be used by physicians trained in the performance of percutaneous transluminal angioplasty.
4. The minimal acceptable sheath French size is printed on the package label. Do not attempt to pass the PTA catheter through a smaller size introducer sheath than indicated on the label.
5. Do not remove the guidewire in situ to shoot contrast through the wire lumen or perform a wire exchange. If the wire is removed while the balloon catheter is situated in tortuous anatomy, the risk of kinking the catheter is increased.
6. Use the recommended balloon inflation medium (25% contrast medium/75% sterile saline solution). It has been shown that a 25/75% contrast / saline ratio has yielded faster balloon inflation / deflation times. Never use air or other gaseous medium to inflate the balloon.
7. If resistance is felt during post procedure withdrawal of the catheter through the introducer sheath, determine if contrast is trapped in the balloon with fluoroscopy. If contrast is present, push the balloon out of the sheath and then completely evacuate the contrast before proceeding to withdraw the balloon.
8. If resistance is still felt during post procedure withdrawal of the catheter, it is recommended to remove the balloon catheter and guidewire / introducer sheath as a single unit, and replace the previously used balloon catheter with a new balloon. Exercise caution when removing the device.
9. Do not continue to use the balloon catheter if the shaft has been bent or kinked. Do not excessively bend, twist, or alter the shape of the device as it may compromise the integrity of the hydrophilic coating.
10. Prior to re-insertion through the introducer sheath, the balloon should be wiped clean with wet gauze and rinsed with sterile normal saline. Avoid excessively wiping the coated portions of the device, or wiping with dry gauze, as this may damage the hydrophilic coating.
11. Balloon re-wrapping should only occur while the balloon catheter is supported with a guidewire or stylet.
12. This device is coated with a hydrophilic coating at the distal segment of the shaft and the balloon. Please refer to the Directions for Use section for further information on how to prepare and use this device to ensure it performs as intended. Failure to abide by the warnings in this labeling might result in damage to the hydrophilic coating, which may require intervention or result in serious adverse events.
13. In order to activate the hydrophilic coating, it is recommended to wet the ULTRAVERSE® catheter with sterile saline solution immediately prior to its insertion in the body. Using different media other than the recommended solution could affect the hydrophilic coating and its performance.
14. The GEOALIGN® Marking System is designed to be used as an additional reference tool to accompany the interventionalist standard operation procedure. The use of fluoroscopic imaging is recommended following positioning of the catheter to the target lesion and prior to balloon deployment.
15. Avoid using alcohol, antiseptic solutions, or other solvents to pre-treat the device because this may cause unpredictable changes in the hydrophilic coating which could affect the device safety and performance.
16. Avoid pre-soaking devices for extended periods, as this may impact the hydrophilic coating performance.
17. It is recommended to consider the use of anti-coagulants, anti-platelet agents, and/or vasodilators in conformance with the accepted standard of practice or institutional guidelines surrounding peripheral endovascular procedures.

Potential Adverse Reactions

The complications which may result from a peripheral balloon dilatation procedure include:

- Additional intervention
- Allergic reaction to drugs or contrast medium
- Aneurysm or pseudoaneurysm
- Arrhythmias
- Compartment Syndrome
- Embolization
- Hematoma
- Hemorrhage, including bleeding at the puncture site
- Hypotension/hypertension
- Inflammation
- Occlusion
- Pain or tenderness
- Pneumothorax or hemothorax
- Sepsis/infection
- Shock
- Short term hemodynamic deterioration
- Stroke
- Thrombosis
- Vessel dissection, perforation, rupture, or spasm

Directions for Use

Handling & Storage

Store in a cool, dry, dark place. Do not store near radiation or ultra-violet light sources.

Rotate inventory so that the catheters and other dated products are used prior to the 'Use By' date. Do not use if packaging is damaged or opened.

Equipment Required

- Contrast medium
- Sterile saline solution
- Luer lock syringe/inflation device with manometer (10 ml or larger)
- Appropriate introducer sheath and dilator set
- 0.014" or 0.018" guidewire

Dilatation Catheter Preparation

NOTE: Remove the stylet prior to removing the protective sheath as part of the preparation for use.

1. Remove Catheter from package. Verify the balloon size is suitable for the procedure and the selected accessories accommodate the catheter as labeled.
2. Remove the balloon guard and stylet by grasping the balloon catheter just proximal to the balloon and with the other hand, gently grasp the balloon protector and slide distally off of the balloon catheter.
3. Prior to use, the air in the balloon catheter should be removed. To facilitate purging, select a syringe or inflation device with a 10 ml or larger capacity and fill approximately half of it with the appropriate balloon inflation medium (25% contrast medium/75% sterile saline solution). Do not use air or any gaseous medium to inflate the balloon.
4. Connect a stopcock to the balloon inflation female luer hub on the dilatation catheter.
5. Connect the syringe to the stopcock.
6. Hold the syringe with the nozzle pointing downward, open the stopcock and aspirate for approximately 15 seconds. Release the plunger.
7. Repeat step #6 two more times or until bubbles no longer appear during aspiration (negative pressure). Once completed, evacuate all air from the barrel of the syringe/inflation device.
8. Prepare the wire lumen of the catheter by attaching a syringe to the wire lumen hub and flushing with sterile saline solution.
9. In order to activate the hydrophilic coating, wet the balloon catheter with sterile saline immediately prior to its insertion into the introducer sheath.

Use of the Ultraverse® 014 and Ultraverse® 018 PTA Balloon Dilatation Catheters

1. Backload the distal tip of the ULTRVERSE® 014 or ULTRVERSE® 018 PTA Balloon Dilatation Catheters over the pre-positioned guidewire and advance the tip to the introduction site.
2. Advance the catheter through the introducer sheath and over the wire to the site of inflation. (Note: If using a hydrophilic guidewire, ensure that it is kept hydrated with sterile normal saline, at all times.) If the stenosis cannot be crossed with the desired dilatation catheter, use a smaller diameter catheter to pre-dilate the lesion to facilitate passage of a more appropriately sized dilatation catheter.
3. Position the balloon relative to the lesion to be dilated, ensure the guidewire is in place and inflate the balloon to the appropriate pressure.

NOTE: Balloon lengths 100mm and longer have two marker bands located at the distal portion of the working length of the balloon and one marker band at the proximal portion. The working length of the balloon is measured at the proximal edge of the proximal marker band to the distal edge of the most distal marker band.

4. Apply negative pressure to fully evacuate fluid from the balloon. Confirm that the balloon is fully deflated under fluoroscopy and that no contrast is left in the balloon.

5. While maintaining negative pressure and the position of the guidewire, grasp the balloon catheter just outside the sheath and withdraw the deflated dilatation catheter over the wire through the introducer sheath. Use of a gentle counterclockwise motion may be used to help facilitate catheter removal through the introducer sheath.

Optional: When used as a location reference tool, utilize the GEOALIGN® Marking System to assist with measuring how far the catheter has advanced. The GEOALIGN® Marking System provides an approximation that may not be an exact representation of the actual distance traveled intravascularly and should be confirmed under fluoroscopy (Reference Figure 3).


Figure 3

Optional: When used for intravascular length measurement, advance or retract the catheter between Points A and B utilizing the GEOALIGN® Marking System to measure the difference. The GEOALIGN® Marking System provides an approximation that may not be an exact representation of the actual distance traveled intravascularly and should be confirmed under fluoroscopy (Reference Figure 4).


Figure 4

Optional: To facilitate repeat catheter alignment with an additional device with the GEOALIGN® Marking System, ensure the GEOALIGN® markings are at the same location outside the sheath as the initial device.

Balloon Reinsertion

Precaution: Do not continue to use the balloon catheter if the shaft has been bent or kinked.

Precaution: Prior to re-insertion through the introducer sheath, the balloon should be wiped clean with gauze and rinsed with sterile normal saline. Balloon rewinding should only occur while the balloon catheter is supported with a guidewire or stylet.

1. Load the balloon catheter onto a guidewire.
2. Advance the balloon catheter over the prepositioned guidewire to the introduction site and through the introducer sheath. If resistance is encountered, replace the previously used balloon catheter with a new balloon.
3. Continue the procedure according to the "Use of the ULTRVERSE® 014 and ULTRVERSE® 018 PTA Balloon Dilatation Catheter" section herein.

Warning: After use, this product may be a potential biohazard. Handle and dispose of in accordance with acceptable medical practices and applicable local, state and federal laws and regulations.

Warranty

Bard Peripheral Vascular warrants to the first purchaser of this product, that this product will be free from defects in materials and workmanship for a period of one year from the date of first purchase and liability under this limited product warranty will be limited, to repair or replacement of the defective product, at Bard Peripheral Vascular's sole discretion, or refunding your net price paid. Wear and tear from normal use or defects resulting from misuse of this product is not covered by this limited warranty.

TO THE EXTENT ALLOWABLE BY APPLICABLE LAW, THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL BARD PERIPHERAL VASCULAR BE LIABLE TO YOU FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM YOUR HANDLING OR USE OF THIS PRODUCT.

Some states/countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your state/country.

An issue or revision date and revision number for these instructions are included for the user's information on the last page of this booklet. In the event 36 months have elapsed between this date and product use, the user should contact Bard Peripheral Vascular to see if additional product information is available.

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**BD, the BD logo, Bard, GeoAlign, Ultra-Cross and Ultraverse are trademarks and/or registered trademarks of Becton, Dickinson and Company.**
**Copyright © 2019 BD and its subsidiaries. All rights reserved.**
Copyright © 2019 BD et ses filiales. Tous droits réservés.
Copyright © 2019 BD und dessen Tochtergesellschaften. Alle Rechte vorbehalten.
Copyright © 2019 BD e le sue sussidiarie. Tutti i diritti riservati.
Copyright © 2019 BD y subsidiarias. Todos los derechos reservados.