

# SOFIA<sup>®</sup> EX

Intracranial Support Catheter

NA

## EXcellent Trackability. EXtra Stability.

### EXcellent Trackability for distal navigation and positioning.<sup>1</sup>

- Soft distal tip designed for easy navigation to the desired target location




### EXtra Stability for reliable device delivery.<sup>2</sup>

- Exhibits minimal tip movement when deploying and re-sheathing treatment devices<sup>2</sup>
- Kink resistant construction provides consistent support for device delivery in tortuous anatomy


Distal Soft Length: 9 cm

Proximal Length: 96 cm



Distal Soft Length: 9 cm

Proximal Length: 106 cm



### SOFIA® EX Intracranial Support Catheter

Catalog Number	Distal/Proximal ID		Distal OD			Proximal OD			Usable Length (cm)	Distal Soft Length (cm)	Tip Shape
	(in)	(mm)	(F)	(in)	(mm)	(F)	(in)	(mm)			
ISC5105ST	0.058	1.47	5.2	0.068	1.73	5.4	0.071	1.80	105	9	Straight
ISC5115ST									115		

1. TR19-145A; n=3
2. TR19-073A; n=3

#### INDICATIONS FOR USE

The SOFIA® EX Catheter is indicated for general intravascular use, including the neuro and peripheral vasculature. The SOFIA® EX Catheter can be used to facilitate introduction of diagnostic agents or therapeutic devices. The SOFIA® EX Catheter is not intended for use in coronary arteries.

RX Only: Federal (USA) law restricts this device to sale by or on the order of a physician.

This device should be used only by physicians trained in percutaneous, intravascular techniques and procedures at medical facilities with the appropriate fluoroscopy equipment.

MICROVENTION and SOFIA are registered trademarks of MicroVention, Inc.  
© 2019 MicroVention, Inc. MM761 NA 06/19

**MicroVention Worldwide  
Innovation Center** 1.714.247.8000  
35 Enterprise  
Aliso Viejo, CA 92656, USA  
**Customer Service** 1.800.990.8368  
**Web** microvention.com

