

Use the MReye[®] coil as an anchor in large and high-flow vessels.



MReye[®] EMBOLIZATION COIL

High radial force

The high radial force of the MReye coil provides strong purchase and stability in large or high-flow vessels.

Extra-long fibers

Fibers on the highly visible Inconel[®] MReye coil allow you to achieve rapid occlusion in high-flow vessels.

Inconel is a registered trademark of the Huntington Alloys Corporation.

Multiple configurations

Designed for .035 and .038 inch diameter systems, the MReye coil comes in a variety of lengths (2-16 cm) and diameters (2-45 mm) so that you can occlude large or high-flow vessels.



www.cookmedical.com

MREYE EMBOLIZATION COIL

Order Number	Reference Part Number	Recommended Catheter ID/ End Hole Diameter inch	Extended Embolus Length cm	Coiled Embolus Diameter mm	Approximate Number of Loops
.035 inch					
G42406	IMWCE-35-2-3	.035	2	3	2.1
G47900	IMWCE-35-3-2	.035	3	2	4.8
G42408	IMWCE-35-3-3	.035	3	3	3.1
G42409	IMWCE-35-3-4	.035	3	4	2.4
G42410	IMWCE-35-3-5	.035	3	5	1.9
G42414	IMWCE-35-4-3	.035	4	3	4.2
G20098	IMWCE-35-4-4	.035	4	4	3.2
G42419	IMWCE-35-5-3	.035	5	3	5.3
G36405	IMWCE-35-5-5	.035	5	5	3.1
G20291	IMWCE-35-5-6	.035	5	6	2.7
G19907	IMWCE-35-5-7	.035	5	7	2.3
G42423	IMWCE-35-5-8	.035	5	8	2.0
G19908	IMWCE-35-5-10	.035	5	10	1.6
G19957	IMWCE-35-8-5	.035	8	5	5.1
G20948	IMWCE-35-8-8	.035	8	8	3.2
G20207	IMWCE-35-8-10	.035	8	10	2.5
G20452	IMWCE-35-10-8	.035	10	8	4.0
G20035	IMWCE-35-10-10	.035	10	10	3.2
G47887	IMWCE-35-12-10	.035	12	10	3.8
G47894	IMWCE-35-14-8	.035	14	8	5.6
G47890	IMWCE-35-14-12	.035	14	12	3.7
G19929	IMWCE-35-15-10	.035	15	10	4.8
G47895	IMWCE-35-15-12	.035	15	12	4.0
G20055	IMWCE-35-15-15	.035	15	15	3.2
G19949	IMWCE-35-15-20	.035	15	20	2.4

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

MREYE EMBOLIZATION COIL

Order Number	Reference Part Number	Recommended Catheter ID/ End Hole Diameter inch	Extended Embolus Length cm	Coiled Embolus Diameter mm	Approximate Number of Loops
.038 inch					
G42448	IMWCE-38-2-3	.038	2	3	2.1
G20338	IMWCE-38-3-3	.038	3	3	3.2
G42451	IMWCE-38-3-5	.038	3	5	1.9
G42454	IMWCE-38-4-3	.038	4	3	4.2
G20013	IMWCE-38-4-4	.038	4	4	3.2
G42458	IMWCE-38-5-3	.038	5	3	5.3
G42460	IMWCE-38-5-5	.038	5	5	3.2
G42463	IMWCE-38-5-8	.038	5	8	2.0
G42465	IMWCE-38-5-10	.038	5	10	1.6
G42466	IMWCE-38-5-12	.038	5	12	1.3
G42467	IMWCE-38-5-15	.038	5	15	1.1
G42472	IMWCE-38-8-5	.038	8	5	5.1
G42475	IMWCE-38-8-10	.038	8	10	2.5
G42477	IMWCE-38-8-15	.038	8	15	1.7
G19920	IMWCE-38-10-10	.038	10	10	3.2
G47914	IMWCE-38-12-10	.038	12	10	3.8
G47917	IMWCE-38-12-14	.038	12	14	2.7
G19886	IMWCE-38-15-10	.038	15	10	4.8
G42486	IMWCE-38-15-12	.038	15	12	4.0
G42487	IMWCE-38-15-15	.038	15	15	3.2
G19887	IMWCE-38-15-20	.038	15	20	2.4
G42489	IMWCE-38-15-25	.038	15	25	1.9
G52983	IMWCE-38-15-30	.038	15	30	1.6
G52984	IMWCE-38-16-45	.038	16	45	1.1

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.



MReye®

EMBOLIZATION COIL

- MReye coils are used for arterial and venous embolization in the peripheral vasculature.
- MReye coils maintain higher radial force than platinum varieties.
- The radiodensity of the coils is higher than stainless steel. This makes the coils more visible under fluoroscopy.
- MReye coils are MR conditional up to 3T, SAR of 2.0 W/kg.
- The coils come in .035 and .038 inch diameters.
- MReye coils are available in extended embolus lengths of 2-16 cm and coiled embolus diameters of 3-45 mm.
- The approximate number of loops on the coils ranges from 1.1 to 5.6.

Customer Service

EMEA: EDI - www.cookmedical.com/edi.do
Distributors: +353 61239240, ssc.distributors@cookmedical.com
Austria: +43 179567121, oe.orders@cookmedical.com
Belgium: +32 27001633, be.orders@cookmedical.com
Denmark: +45 38487607, da.orders@cookmedical.com
Finland: +358 972519996, fi.orders@cookmedical.com
France: +33 171230269, fr.orders@cookmedical.com
Germany: +49 6950072804, de.orders@cookmedical.com
Hungary: +36 17779199, hu.orders@cookmedical.com
Ireland: +353 61239252, ie.orders@cookmedical.com
Italy: +39 0269682853, it.orders@cookmedical.com
Netherlands: +31 202013367, nl.orders@cookmedical.com
Norway: +47 23162968, no.orders@cookmedical.com
Spain: +34 912702691, es.orders@cookmedical.com
Sweden: +46 858769468, se.orders@cookmedical.com
Switzerland - French: +41 448009609, fr.orders@cookmedical.com
Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
Switzerland - German: +41 448009609, de.orders@cookmedical.com
United Kingdom: +44 2073654183, uk.orders@cookmedical.com

www.cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do
Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
E-mail: orders@cookmedical.com

Australia:
Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283
E-mail: cau.custserv@cookmedical.com

