

# VSI Guidewire

0.18" and 0.014" Guidewires

**Coil Tips** .....  
available in soft, flexible  
stainless steel or highly  
radiopaque tungsten

**50cm 0.014"**  
for specialty use

**Stainless Steel or Nitinol Mandrels** .....  
designed for flexibility and pushability

**Separately Packaged**  
0.018" guidewires available  
in a variety of lengths

# Ordering Information

## 0.018" Guidewires

### Stainless Steel Mandrel - Stainless Steel Soft Tip

| MODEL | WIRE DIAMETER | WIRE LENGTH | SPECIAL FEATURES | BOX QUANTITY |
|-------|---------------|-------------|------------------|--------------|
| 7300  | 0.018"        | 40 cm       | -                | 10           |
| 7301  | 0.018"        | 60 cm       | -                | 10           |

### Nitinol Mandrel - Stainless Steel Soft Tip

| MODEL | WIRE DIAMETER | WIRE LENGTH | SPECIAL FEATURES | BOX QUANTITY |
|-------|---------------|-------------|------------------|--------------|
| 7302  | 0.018"        | 40 cm       | -                | 10           |
| 7303  | 0.018"        | 60 cm       | -                | 10           |

### Stainless Steel Mandrel - Tungsten Soft Tip

| MODEL | WIRE DIAMETER | WIRE LENGTH | SPECIAL FEATURES | BOX QUANTITY |
|-------|---------------|-------------|------------------|--------------|
| 7304  | 0.018"        | 40 cm       | -                | 5            |
| 7305  | 0.018"        | 60 cm       | -                | 5            |
| 7326  | 0.018"        | 80 cm       | PTFE-coated      | 5            |
| 7327  | 0.018"        | 100 cm      | PTFE-coated      | 5            |

### Nitinol Mandrel - Tungsten Soft Tip

| MODEL | WIRE DIAMETER | WIRE LENGTH | SPECIAL FEATURES | BOX QUANTITY |
|-------|---------------|-------------|------------------|--------------|
| 7306  | 0.018"        | 40 cm       | Angled Tip       | 5            |
| 7342  | 0.018"        | 40 cm       | -                | 5            |
| 7307  | 0.018"        | 60 cm       | -                | 5            |
| 7324  | 0.018"        | 80 cm       | -                | 5            |
| 7325  | 0.018"        | 130 cm      | Marked*          | 5            |

Devices packaged in the quantities listed above. Individual units are not sold separately.

\*Reference markers at 5 cm intervals between distal 20 cm and 65 cm.

The VSI Guidewire is indicated for percutaneous entry of peripheral vessels using the Seldinger Technique. This device is not intended for use in the coronary or cerebral vasculature.

## 0.014" Guidewires

### Nitinol Mandrel - Stainless Steel Soft Tip

| MODEL | WIRE DIAMETER | WIRE LENGTH | SPECIAL FEATURES | BOX QUANTITY |
|-------|---------------|-------------|------------------|--------------|
| 7343  | 0.014"        | 50 cm       | -                | 10           |

Packaged in quantities of 10 units per box. Individual units are not sold separately.

The VSI Guidewire is indicated for percutaneous entry of peripheral vessels using the Seldinger Technique. This device is not intended for use in the coronary or cerebral vasculature.

Please see Instructions for Use for a complete listing of the indications, contraindications, warnings and precautions.

CAUTION: Federal Law (U.S.A.) restricts this device to sale by or on the order of a physician.

Teleflex, and the Teleflex logo are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries. All other trademarks and registered trademarks are property of their respective owners.

Information in this material is not a substitute for the product Instructions for Use. Not all products may be available in all countries. Please contact your local representative. Revised: 02/2020.

© 2020 Teleflex Incorporated. All rights reserved. MC-006091