


# Navigate smoothly through the vessel with hydrophilic microwire technology.



## Approach<sup>®</sup> Hydro ST

MICROWIRE GUIDE

### Short Taper


The short taper is designed specifically for peripheral intervention. The wire guide's distal 13.5 cm is enhanced with a hydrophilic coating that provides smooth advancement through vessels.

### Visualization

Platinum coils in the distal tip give you the visibility you need to advance the wire guide with control and accuracy.

### Supportive Shaft

A supportive stainless-steel shaft helps you navigate support catheters smoothly through the vessel.



# Approach® Hydro ST

MICROWIRE GUIDE

Used to facilitate delivery of percutaneous catheters into the peripheral vasculature.

Order Number	Reference Part Number	Wire Guide Diameter inch	Wire Guide Length cm
G52937	HMW-14-135-ST	.014	135
G52938	HMW-14-190-ST	.014	190
G52939	HMW-14-300-ST	.014	300

## Pair with your Cook essentials.

### ACCESS

**Micropuncture®**  
Introducer Set

**Micropuncture**  
Pedal Access Set

### TARGET


**CXI®**  
Support Catheter

**Flexor® Ansel Guiding Sheath**  
with High-Flex Dilator

### TREAT

**Advance® 14LP**  
Low-Profile PTA Balloon Dilatation Catheter

**Advance Serenity™ 14**  
Hydrophilic PTA Balloon Dilatation Catheter



### Customer Service

EU Website: [cookmedical.eu](http://cookmedical.eu)  
 EDI: [cookmedical.eu/edi](http://cookmedical.eu/edi)  
 Distributors: +353 61239240, [ssc.distributors@cookmedical.com](mailto:ssc.distributors@cookmedical.com)  
 Austria: +43 179567121, [oe.orders@cookmedical.com](mailto:oe.orders@cookmedical.com)  
 Belgium: +32 27001702, [be.orders@cookmedical.com](mailto:be.orders@cookmedical.com)  
 Denmark: +45 38487607, [da.orders@cookmedical.com](mailto:da.orders@cookmedical.com)  
 Finland: +358 972519996, [fi.orders@cookmedical.com](mailto:fi.orders@cookmedical.com)  
 France: +33 171230269, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
 Germany: +49 6950072804, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
 Hungary: +36 17779199, [hu.orders@cookmedical.com](mailto:hu.orders@cookmedical.com)  
 Iceland: +354 8007615, [IS.orders@cookmedical.com](mailto:IS.orders@cookmedical.com)  
 Ireland: +353 61239252, [ie.orders@cookmedical.com](mailto:ie.orders@cookmedical.com)  
 Italy: +39 0269682853, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
 Netherlands: +31 202013367, [nl.orders@cookmedical.com](mailto:nl.orders@cookmedical.com)  
 Norway: +47 23162968, [no.orders@cookmedical.com](mailto:no.orders@cookmedical.com)  
 Spain: +34 912702691, [es.orders@cookmedical.com](mailto:es.orders@cookmedical.com)  
 Sweden: +46 858769468, [se.orders@cookmedical.com](mailto:se.orders@cookmedical.com)  
 Switzerland - French: +41 448009609, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
 Switzerland - Italian: +41 448009609, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
 Switzerland - German: +41 448009609, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
 United Kingdom: +44 2073654183, [uk.orders@cookmedical.com](mailto:uk.orders@cookmedical.com)

USA Website: [cookmedical.com](http://cookmedical.com)  
 EDI: [cookmedical.com/edi.do](http://cookmedical.com/edi.do)  
**Americas:**  
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335  
 E-mail: [customersupport@cookmedical.com](mailto:customersupport@cookmedical.com)  
**Australia:**  
 Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283  
 E-mail: [cau.custserv@cookmedical.com](mailto:cau.custserv@cookmedical.com)



AI, ESC, IR, OHNS, PI, RH, SUR-8.5X11