


CONFIDENTLY
CROSS



Intuition™
Advanced Workhorse Guidewire



Intuition™

Advanced Workhorse Guidewire

The Intuition™ guidewire's torque performance and tip weight give you precision and confidence as you control and direct the wire.

- Improved torque control
- Designed to provide steerability and control in tortuous anatomy

ACCU-CORE TECHNOLOGY OFFERS YOU:


TORQUE PERFORMANCE

Our core wire process is designed to provide torque performance to control and direct the wire.

TACTILE FEEL AND CROSSABILITY

The rounded tip and body provide control.

ACCU-CORE TECHNOLOGY




Centered core wire

PREVIOUS DESIGN



Uncentered core wire

Intuition™ Advanced Workhorse Guidewire



ADDITIONAL FEATURES

- Femoral (proximal) markers on wire, 90 cm and 100 cm from distal tip
- DOC™* compatible proximal ends on 180-cm wires

A **smooth**, responsive feel you can **depend** on

Intuition™ Advanced Workhorse Guidewire

A choice of two coatings to meet your needs — Available in Pro/Pel™ silicone coating or Hydro-Track™ hydrophilic coating, for a more lubricious wire without increasing tip stiffness.

Ordering Information — Contact your Medtronic sales representative.

Product	Product Code	Length (cm)	Coating	Tip
Intuition™ Guidewire	INTU180HS	180	Hydro-Track™	Straight tip
	INTU180SS	180	Pro/Pel™	Straight tip
	INTU300HS	300	Hydro-Track™	Straight tip
	INTU300SS	300	Pro/Pel™	Straight tip

Note: All guidewires are 0.014-in wires and 180-cm wires are DOC™* compatible.

Medtronic

Medtronic
3576 Unocal Place
Santa Rosa, CA 95403
USA
Tel: 707.525.0111

LifeLine Customer Support
Tel: 877.526.7890
Tel: 763.526.7890

Product Services
Tel: 888.283.7868
Fax: 800.838.3103

Europe
Medtronic Intl. Trading SARL
Route du Molliat 31
Case Postale
CH-1131 Tolochenaz
Switzerland
Tel: 41.21.802.7000
Fax: 41.21.802.7900

Canada
Medtronic of Canada Ltd.
99 Hereford Street
Brampton, Ontario L6Y 0R3
Canada
Tel: 905.460.3800
Fax: 905.460.3998
Toll-free: 800.268.5346

Asia Pacific
Medtronic Intl. Ltd.
49 Changi South Avenue 2
Nasaco Tech Centre
Singapore 486056
Singapore
Tel: 65.6436.5000
Fax: 65.6776.6335

Latin America
Medtronic
9850 Doral Blvd., 4th Floor
Miami, FL 33178
USA
Tel: 786.709.4200
Fax: 786.709.4244

© 2018 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. ™* Third-party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company. For distribution only in markets where the Intuition™ crossing guidewire has been approved. Not for distribution in Japan or France. UC201707399IE 12/18