



TZ Medical
Sparked by your ideas

ACCUFLATE INFLATION DEVICE

The TZ Medical AccuFlate is a 20ml disposable device designed to accurately control and monitor pressure during inflation and deflation of dilation balloon catheters. Our patented technology allows the Inflation Device to generate pressure up to 30 atm with accuracy of ± 1 atm.



Patented self-locking mechanism for easy one-handed locking control, as well as gradual and rapid pressure release ergonomic handle.



Precision manometer with a maximum pressure of 30 atm.



Optional, male rotating luer lock, 3-way stopcock with "ON" or "OFF" tap handle.

Item Number	Description	UOM
-------------	-------------	-----

IDD-0001	Balloon Inflation Device	Box of 5
----------	--------------------------	----------

FIND OUT MORE



tzmedical.com



800.944.0187



info@tzmedical.com

19337 R1.0