


Address transjugular procedures confidently with an array of specialized products.



PRODUCTS FOR TRANSJUGULAR
LIVER ACCESS AND BIOPSY

Cook Medical Transjugular Liver Access and Biopsy

For transjugular procedures, you need the flexibility to reach your objective, and the precision to make your access count.

Cook Medical's specialized sets are purposefully designed for the specific patient needs and physician preferences associated with transjugular liver procedures. The controlled maneuverability of Flexor® introducer sheaths, the precision of Quick-Core® biopsy needles, and a full range of accessory products work together to help you meet your objective.

Transjugular Liver Access

Rösch-Uchida Technique

Rösch-Uchida Transjugular Liver Access Set 5

Colapinto Technique

Ring Transjugular Intrahepatic Access Set 6

Transjugular Liver Biopsy

Core Biopsy

Liver Access and Biopsy Set 9

Accessory Products

Access Needles and Sets 11

Introducers 11


Wire Guides 12-13


Dilators 13


Adapters and Stopcocks 14

Connecting Tube 14

Transjugular Liver Access









Rösch-Uchida

TRANSJUGULAR LIVER ACCESS SET

- Stiffening cannula provides added support during puncture.
- Trocar stylet facilitates vascular access within the liver for subsequent wire guide, catheter, and sheath placement.
- The Flexor® sheath design provides maximum flexibility without kinking or compression.
- Radiopaque band on introducer increases visibility for positioning accuracy.

Order Number	Reference Part Number	Introducer Fr/Length cm	Trocar Stylet Diameter inch	Trocar Stylet Length cm	Stiffening Cannula gage	Stiffening Cannula Length cm
Rösch-Uchida Transjugular Liver Access Set						
G06929	RUPS-100	10/40	.038	62.5	14	51.5

Determine the maximum diameter of the instrument or catheter to be introduced to ensure its passage through the sheath.







Ring

TRANSJUGULAR INTRAHEPATIC ACCESS SET


- Ross Modified Colapinto Needle facilitates vascular access within the liver for subsequent wire guide, catheter, and sheath placement.
- The Flexor® sheath design provides maximum flexibility without kinking or compression.
- Radiopaque band on introducer increases visibility for precise positioning.
- Availability in both 9 and 10 Fr sizes allows physician to address a variety of patient needs.


Order Number	Reference Part Number	Needle gage	Needle Length cm	Introducer Fr	Introducer Length cm
Ring Transjugular Intrahepatic Access Set					
G06541	RTPS-100	16	50.5	9	38.5
G29769	RTPS-100-10.0	16	50.5	10	38.5

Determine the maximum diameter of the instrument or catheter to be introduced to ensure its passage through the sheath.



Transjugular Liver Biopsy






LABS

LIVER ACCESS AND BIOPSY SET


- Includes a full complement of components for transjugular liver access and tissue biopsy.
- Indicated for adult and pediatric use.
- Precurved catheter and cannula ideal for hepatic vein access.
- Single-throw mechanical needle passes easily into tissue.
- Check-Flo® Performer adapter with sidearm tubing prevents blood reflux and allows the injection of contrast media or the aspiration of blood to confirm positioning.

Order Number	Reference Part Number	Needle gage	Needle Length cm	Needle Throw Length mm	Sheath Fr
Liver Access and Biopsy Set					
G57485	LABS-100-J-01	18	60	20	7
G57487	LABS-200-J-01	19	60	20	7
Liver Access and Biopsy Needle					
G57482	RMT-14XT-50.5-LABS-100-01	18	60	20	7
G57483	RMT-14XT-50.5-LABS-200-01	19	60	20	7

Determine the maximum diameter of the instrument or catheter to be introduced to ensure its passage through the sheath.



Accessory Products



ACCESSORY PRODUCTS

Access Needles and Sets



Order Number	Reference Part Number	gage	Length cm
Percutaneous Entry Needle			
G00166	SDN-18-7.0	18	7
EchoTip® Percutaneous Entry Needle			
G04717	SDN-18-7.0-U	18	7



Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Wire Guide Diameter inch	Wire Guide Length cm	Needle gage/Length cm
Micropuncture Transitionless Access Set						
G43872	MPIS-505-SST	5	10	.018	40	21/7

Introducers



Order Number	Reference Part Number	Fr	Accepts Wire Guide Diameter inch	Length cm
Flexor® Guiding Sheath				
G32233	KCFW-10.0-38-40-RB	10	.038	40



Order Number	Reference Part Number	Fr	Accepts Wire Guide Diameter inch	Length cm
Flexor® RTPS Guiding Sheath				
G13160	KCFW-9.0-35-38.5-RB-RTPS-100	9	.035	38.5
G13081	KCFW-10.0-35-RB-RTPS-100	10	.035	38.5

ACCESSORY PRODUCTS

Wire Guides



Order Number	Reference Part Number	Diameter inch	Length cm	Tip Configuration
Roadrunner® The Firm™ Hydrophilic-Coated Wire Guide				
G09534	RPC-35-80	.035	80	angled
G06979	RPC-35-145	.035	145	angled
G07518	RPC-35-180	.035	180	angled
Roadrunner® Uniglide® Hydrophilic-Coated Wire Guide – Standard Shaft				
G56172	HPWA-35-80	.035	80	angled
G56173	HPWA-35-150	.035	150	angled
G56174	HPWA-35-180	.035	180	angled



Order Number	Reference Part Number	Diameter inch	Length cm	Flexible Tip Length cm	Tip Curve Radius mm
Amplatz Extra-Stiff Support Wire Guide					
G04003	THSF-35-80-AES	.035	80	7	straight
G03095	THSF-35-145-AES	.035	145	7	straight
G03330	THSF-35-180-AES	.035	180	7	straight
G04237	THSCF-35-80-3-AES	.035	80	7	3
G03203	THSCF-35-145-3-AES	.035	145	7	3
G03565	THSCF-35-180-3-AES	.035	180	7	3
Amplatz Ultra-Stiff Support Wire Guide					
G28787	THSF-35-90-AUS2	.035	90	7	straight
G19852	THSF-35-145-AUS2	.035	145	7	straight
G20131	THSF-35-180-AUS2	.035	180	7	straight
G28761	THSCF-35-90-3-AUS2	.035	90	7	3
G19853	THSCF-35-180-3-AUS2	.035	180	7	3


ACCESSORY PRODUCTS

Wire Guides *(continued)*



Order Number	Reference Part Number	Diameter inch	Length cm	Tip Curve Radius mm
Tefcor Moveable Core Wire Guide				
G01730	TCMTNA-35-145-1.5	.035	145	1.5
G01678	TCMTNA-35-145-3	.035	145	3


Dilators



Order Number	Reference Part Number	Fr	Accepts Wire Guide Diameter inch	Length cm
Dilator				
G00892	JCD7.0-35-20	7	.035	20
G00973	JCD7.0-38-20	7	.038	20
G00984	JCD8.0-35-20	8	.035	20
G00980	JCD8.0-38-20	8	.038	20
G00992	JCD9.0-38-20	9	.038	20
G00993	JCD10.0-38-20	10	.038	20
G00994	JCD11.0-38-20	11	.038	20
G00995	JCD12.0-38-20	12	.038	20
AQ® Hydrophilic-Coated Dilator				
G10208	JCD7.0-38-20-HC	7	.038	20
G10289	JCD8.0-38-20-HC	8	.038	20
G10290	JCD9.0-38-20-HC	9	.038	20
G10207	JCD10.0-38-20-HC	10	.038	20
G10292	JCD12.0-38-20-HC	12	.038	20


COMPLEMENTARY PRODUCTS

Adapters and Stopcocks



Order Number	Reference Part Number
Check-Flo® Hemostasis Assembly	
G06002	CFM-100
One-Way Plastic Stopcock	
G00213	POWS-FLL-MLL
Three-Way Plastic Stopcock	
G00219	PTWS-2FLL-MLL-R
High-Pressure Three-Way Stopcock	
G06540	PTWSC-2FLL-MLL-R-HP

Connecting Tube



Order Number	Reference Part Number	Fr	Length cm
Medium-Pressure Connecting Tube			
G00043	DPT9.5-30-P-FM	9.5	30

Determine the maximum diameter of the instrument or catheter to be introduced to ensure its passage through the sheath.

Customer Service

EU Website: cookmedical.eu
EDI: cookmedical.eu/edi
Distributors: +353 61239240, ssc.distributors@cookmedical.com
Austria: +43 179567121, oe.orders@cookmedical.com
Belgium: +32 27001702, be.orders@cookmedical.com
Denmark: +45 38487607, da.orders@cookmedical.com
Finland: +358 972519996, fi.orders@cookmedical.com
France: +33 171230269, fr.orders@cookmedical.com
Germany: +49 6950072804, de.orders@cookmedical.com
Hungary: +36 17779199, hu.orders@cookmedical.com
Iceland: +354 8007615, IS.orders@cookmedical.com
Ireland: +353 61239252, ie.orders@cookmedical.com
Italy: +39 0269682853, it.orders@cookmedical.com
Netherlands: +31 202013367, nl.orders@cookmedical.com
Norway: +47 23162968, no.orders@cookmedical.com
Spain: +34 912702691, es.orders@cookmedical.com
Sweden: +46 858769468, se.orders@cookmedical.com
Switzerland - French: +41 448009609, fr.orders@cookmedical.com
Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
Switzerland - German: +41 448009609, de.orders@cookmedical.com
United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com

EDI: cookmedical.com/edi.do

Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

E-mail: customersupport@cookmedical.com

Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283

E-mail: cau.custserv@cookmedical.com



AI, ESC, IR, OHNS, PI, RH, SUR-8.5X11