



Compression. Comfort. Convenience.



why safeguard radial?

The Safeguard Radial Compression Device is designed to help you achieve PATENT HEMOSTASIS at the conclusion of your radial procedure. The device is COMFORTABLE for your patients, CONVENIENT for hospital staff, and EASY to use.



patent hemostasis

Focuses compression only on the radial artery where you want it

Puncture site clearly visible for visual assessment

Stable positioning of bulb for optimal compression, not occlusion



comfort

Non-woven textile band for patient comfort

No rigid plastic edges

Low profile design



convenient and easy

One size fits all, easily adjusted to any wrist size

Uses standard Luer lock syringe, no need for special syringe

Easy to monitor distal and proximal pulses to ensure patent hemostasis is achieved.

Easy inflation and deflation



Service Shear And Shear An



ORDER INFORMATION

Catalog Number: Description: Length

SGR26 (Former catalog number 0684-00-0572)

SAFEGUARD Radial Compression Device

26 cm

Quantity: 10 units per box, packaged individually Contact Merit Medical for more information.

MERITATION CA.



Merit Medical Systems, Inc. 1600 West Merit Parkway South Jordan, Utah 84095 1-801-253-1600 1-800-35-MERIT www.merit.com

Merit Medical Europe, Middle East, and Africa (EMEA)

Amerikalaan 42, 6199 AE Maastricht-Airport The Netherlands

Tel: +31 43 358 82 22

Merit Medical Ireland Ltd.

Parkmore Business Park West Galway, Ireland +353 (0) 91 703 733

Free phone for specific country

Austria 0800 295 374 | Belgium 0800 72 906 (Dutch) 0800 73 172 (French) | Denmark 80 88 00 24 France 0800 91 60 30 | Finland 0800 770 586 | Germany 0800 182 0871 | Ireland (Republic) 1800 553 163 Italy 800 897 005 | Luxembourg 8002 25 22 Netherlands 0800 022 81 84 \(\text{\text{\$\ I}} \) Norway 800 11629 Sweden 020 792 445 | UK 0800 973 115

ID 112113 400508001/A