




PERIPHERAL INTERVENTION

DRAINAGE SOLUTIONS

Everything you
Need for Drainage
Procedures



LOCKING DRAINAGE CATHETER

Designed for abscess, nephrostomy and miscellaneous drainage procedures, the ReSolve locking drainage catheter comes packaged with a metal stiffening cannula, a plastic stiffening cannula, and a trocar stylet.



- Catheter tip positioned within pigtail curve to reduce tissue irritation.
- Depth markers located on the catheter shaft aid accurate placement and act as a reference for correct positioning.
- Large drainage holes to maximise drainage.
- Tapered tip and hydrophilic coating for a smooth introduction to reduce patient discomfort.
- Colour-coded hubs that match standard French sizing charts for easy identification.
- Flexible polyurethane material to prevent kinking and ensure excellent catheter rebound.
- Pigtail straightener allows for easier insertion of stiffening cannulas.
- Repositioning Tool (included) for unlocking and repositioning the catheter at any time.
- Highly flexible radiopaque marker band, integral to the catheter shaft, forged from platinum.

ORDERING INFORMATION

Catalogue Number	French Size	Working Length	Max Guide Wire	Hub Colour	Drainage Holes
NON-MARKER BAND					
RLC-6-038	6.5F	25 cm	0.038" (0.97 mm)	Green	8
RLC-7-038	7.5F	25 cm	0.038" (0.97 mm)	Red	8
RLC-8-038	8.5F	25 cm	0.038" (0.97 mm)	Blue	8
RLC-10-038	10F	25 cm	0.038" (0.97 mm)	Purple	8
RLC-12-038	12F	25 cm	0.038" (0.97 mm)	Yellow	9
RLC-14-038	14F	25 cm	0.038" (0.97 mm)	Orange	9
MARKER BAND					
RLC-6-038MB	6.5F	25 cm	0.038" (0.97 mm)	Green	8
RLC-7-038MB	7.5F	25 cm	0.038" (0.97 mm)	Red	8
RLC-8-038MB	8.5F	25 cm	0.038" (0.97 mm)	Blue	8
RLC-10-038MB	10F	25 cm	0.038" (0.97 mm)	Purple	8
RLC-12-038MB	12F	25 cm	0.038" (0.97 mm)	Yellow	9
RLC-14-038MB	14F	25 cm	0.038" (0.97 mm)	Orange	9

Packaged 5 per box

RC-RT00L Repositioning Tool Packaged 15 per box

DRAINAGE CATHETER TRAYS

Conveniently containing the ReSolve Locking Drainage Catheter, MAK-NV™ Introducer System, Merit Drainage Depot™ Drainage Bag, and dilators, the ReSolve trays offer all you need from stick to fixation.



ORDERING INFORMATION

Catalogue Number	French Size	20 cm Dilator French Size	Marker Band	Securement Device
RLC-NV-6	6.5F	-	No	-
RLC-NV-7	7.5F	-	No	-
RLC-NV-8	8.5F	8F	No	-
RLC-NV-10	10F	8, 10F	No	-
RLC-NV-12	12F	8, 10, 12F	No	-
RLC-NV-14	14F	8, 10, 12, 14F	No	-
RLCMB-NV-6	6.5F	-	Yes	-
RLCMB-NV-7	7.5F	-	Yes	-
RLCMB-NV-8	8.5F	8F	Yes	-
RLCMB-NV-10	10F	8, 10F	Yes	-
RLCMB-NV-12	12F	8, 10, 12F	Yes	-
RLCMB-NV-14	14F	8, 10, 12, 14F	Yes	-
RLC-NV-6-SFX	6.5F	-	No	StayFIX (small)
RLC-NV-7-SFX	7.5F	-	No	StayFIX (small)
RLC-NV-8-SFX	8.5F	8F	No	StayFIX (small)
RLC-NV-10-SFX	10F	8, 10F	No	StayFIX (small)
RLC-NV-12-SFX	12F	8, 10, 12F	No	StayFIX (small)
RLC-NV-14-SFX	14F	8, 10, 12, 14F	No	StayFIX (large)
RLCMB-NV-6-SFX	6.5F	-	Yes	StayFIX (small)
RLCMB-NV-7-SFX	7.5F	-	Yes	StayFIX (small)
RLCMB-NV-8-SFX	8.5F	8F	Yes	StayFIX (small)
RLCMB-NV-10-SFX	10F	8, 10F	Yes	StayFIX (small)
RLCMB-NV-12-SFX	12F	8, 10, 12F	Yes	StayFIX (small)
RLCMB-NV-14-SFX	14F	8, 10, 12, 14F	Yes	StayFIX (large)
RLC-NV-6-REV	6.5F	-	No	Revolution
RLC-NV-7-REV	7.5F	-	No	Revolution
RLC-NV-8-REV	8.5F	8F	No	Revolution
RLC-NV-10-REV	10F	8, 10F	No	Revolution
RLC-NV-12-REV	12F	8, 10, 12F	No	Revolution
RLC-NV-14-REV	14F	8, 10, 12, 14F	No	Revolution
RLCMB-NV-6-REV	6.5F	-	Yes	Revolution
RLCMB-NV-7-REV	7.5F	-	Yes	Revolution
RLCMB-NV-8-REV	8.5F	8F	Yes	Revolution
RLCMB-NV-10-REV	10F	8, 10F	Yes	Revolution
RLCMB-NV-12-REV	12F	8, 10, 12F	Yes	Revolution
RLCMB-NV-14-REV	14F	8, 10, 12, 14F	Yes	Revolution

Packaged as individual units



ReSolve+™ LOCKING DRAINAGE CATHETER WITH SECUREMENT DEVICE

Having your drainage catheter and securement device conveniently packaged together saves both time and inventory space.

ORDERING INFORMATION

Catalogue Number	French Size	Working Length	Marker Band	Securement Device
NON-MARKER BAND				
RLC-6-REV	6.5F	25 cm	No	Revolution™
RLC-6-SFX	6.5F	25 cm	No	StayFIX® (small)
RLC-7-REV	7.5F	25 cm	No	Revolution™
RLC-7-SFX	7.5F	25 cm	No	StayFIX® (small)
RLC-8-REV	8.5F	25 cm	No	Revolution™
RLC-8-SFX	8.5F	25 cm	No	StayFIX® (small)
RLC-10-REV	10F	25 cm	No	Revolution™
RLC-10-SFX	10F	25 cm	No	StayFIX® (small)
RLC-12-REV	12F	25 cm	No	Revolution™
RLC-12-SFX	12F	25 cm	No	StayFIX® (small)
RLC-14-REV	14F	25 cm	No	Revolution™
RLC-14-SFX	14F	25 cm	No	StayFIX® (large)
MARKER BAND				
RLCMB-6-REV	6.5F	25 cm	Yes	Revolution™
RLCMB-6-SFX	6.5F	25 cm	Yes	StayFIX® (small)
RLCMB-7-REV	7.5F	25 cm	Yes	Revolution™
RLCMB-7-SFX	7.5F	25 cm	Yes	StayFIX® (small)
RLCMB-8-REV	8.5F	25 cm	Yes	Revolution™
RLCMB-8-SFX	8.5F	25 cm	Yes	StayFIX® (small)
RLCMB-10-REV	10F	25 cm	Yes	Revolution™
RLCMB-10-SFX	10F	25 cm	Yes	StayFIX® (small)
RLCMB-12-REV	12F	25 cm	Yes	Revolution™
RLCMB-12-SFX	12F	25 cm	Yes	StayFIX® (small)
RLCMB-14-REV	14F	25 cm	Yes	Revolution™
RLCMB-14-SFX	14F	25 cm	Yes	StayFIX® (large)

Packaged 5 per box



NON-LOCKING DRAINAGE CATHETER

Designed for drainage of abscess, nephrostomy, and miscellaneous drainage procedures, the ReSolve comes in a J or straight tip and is packaged with a metal stiffening cannula and a trocar stylet.

- Large drainage holes to maximise drainage.
- Tapered tip and hydrophilic coating for a smooth introduction to reduce patient discomfort.
- Colour-coded hubs that match standard French sizing charts for easy identification.
- Flexible polyurethane material to prevent kinking and to ensure excellent catheter rebound.

ORDERING INFORMATION



Catalogue Number	French Size	Working Length	Max Guide Wire	Tip Config.	Drainage Holes
RNL-6-038S	6.5F	20 cm	0.038" (0.97 mm)	Straight	4
RNL-7-038S	7.5F	25 cm	0.038" (0.97 mm)	Straight	4
RNL-8-038S	8.5F	25 cm	0.038" (0.97 mm)	Straight	4
RNL-10-038S	10F	25 cm	0.038" (0.97 mm)	Straight	4
RNL-12-038S	12F	25 cm	0.038" (0.97 mm)	Straight	5
RNL-6-038J	6.5F	20 cm	0.038" (0.97 mm)	J Curve	4
RNL-7-038J	7.5F	25 cm	0.038" (0.97 mm)	J Curve	4
RNL-8-038J	8.5F	25 cm	0.038" (0.97 mm)	J Curve	4
RNL-10-038J	10F	25 cm	0.038" (0.97 mm)	J Curve	4
RNL-12-038J	12F	25 cm	0.038" (0.97 mm)	J Curve	5

Packaged 5 per box



VALVED ONE-STEP CENTESIS CATHETER

The self-closing valve prevents air infiltration and fluid leakage. When activated with a luer, the valve fully opens to allow maximum flow.



- Self-closing valve prevents air infiltration and fluid leakage.
- Four large holes spiral around catheter to maximise drainage.
- Tapered tip for a smooth introduction to reduce patient discomfort.
- Echo-enhanced needle for increased needle visibility.
- Bevel indicator gives a tactile and visual reference for bevel orientation.
- Clear hub increases flashback visibility with ergonomic design for easy gripping and handling.



ORDERING INFORMATION

Catalogue Number	French Size	Length	Introducer Needle	Drainage Holes
SLIP FIT				
OSC-4F-7T	4F	7 cm	20G	4
OSC-4F-10T	4F	10 cm	20G	4
OSC-5F-7T	5F	7 cm	19G	4
OSC-5F-10T	5F	10 cm	19G	4
OSC-5F-15T	5F	15 cm	19G	4
FIXED LUER				
OSC-4F-7LT	4F	7 cm	20G	4
OSC-4F-10LT	4F	10 cm	20G	4
OSC-5F-7LT	5F	7 cm	19G	4
OSC-5F-10LT	5F	10 cm	19G	4
OSC-5F-15LT	5F	15 cm	19G	4

Packaged 5 per box

ONE-STEP CENTESIS CATHETER

For ease of insertion, the One-Step centesis catheter has a smooth transition for decreased resistance and an echo-enhanced needle for more accurate placement.



- Available with luer locking or slip fit needle.
- Four large holes spiral around catheter to maximise drainage.
- Tapered tip for a smooth introduction to reduce patient discomfort.
- FEP material resists buckling for a smooth insertion.
- Unique hub design for easy gripping and handling during procedure.
- Echo-enhanced needle for increased needle visibility.

ORDERING INFORMATION

Catalogue Number	French Size	Length	Introducer Needle	Drainage Holes
SLIP FIT				
OSC-4F-7S	4F	7 cm	20G	4
OSC-4F-10S	4F	10 cm	20G	4
OSC-5F-7S	5F	7 cm	19G	4
OSC-5F-10S	5F	10 cm	19G	4
OSC-5F-15S	5F	15 cm	19G	4
FIXED LUER				
OSC-4F-7L	4F	7 cm	20G	4
OSC-4F-10L	4F	10 cm	20G	4
OSC-5F-7L	5F	7 cm	19G	4
OSC-5F-10L	5F	10 cm	19G	4
OSC-5F-15L	5F	15 cm	19G	4
FIXED LUER NO SIDE HOLES				
OSC-4F-7LNH	4F	7 cm	20G	-
OSC-4F-10LNH	4F	10 cm	20G	-
OSC-5F-7LNH	5F	7 cm	19G	-
OSC-5F-10LNH	5F	10 cm	19G	-
OSC-5F-15LNH	5F	15 cm	19G	-

Packaged 5 per box

MERIT MAK-NV MINI ACCESS KIT NON-VASCULAR

Excellent transitions for a smooth introduction.

- Radiopaque marker tip clearly shows tip of introducer.
- Superior tip design decreases buckling and introducer peel back during insertion.
- Echo-enhanced needle increases visualisation for ultrasound placement.
- 4F introducer ID allows 0.018" guide wire to be left in place for safety during larger wire placement.



ORDERING INFORMATION

MAK-NV-001	0.038" (0.97 mm) X 150 cm double ended heavy duty 3 mm J guide wire 0.018" (0.46 mm) X 60 cm stainless steel guide wire with platinum tip 6F (2.0 mm) X 20 cm coaxial introducer with dilator and stiffening cannula 21G X 15 cm introducer needle with trocar stylet
MAK-NV-002	0.038" (0.97 mm) X 150 cm double ended heavy duty 3 mm J guide wire 0.018" (0.46 mm) X 60 cm nitinol guide wire with platinum tip 6F (2.0 mm) X 20 cm coaxial introducer with dilator and stiffening cannula 21G X 15 cm introducer needle with trocar stylet
MAK-NV-003	0.018" (0.46 mm) X 60 cm stainless steel guide wire with platinum tip 6F (2.0 mm) X 20 cm coaxial introducer with dilator and stiffening cannula 21G X 15 cm introducer needle with trocar stylet
MAK-NV-004	0.018" (0.46 mm) X 60 cm nitinol guide wire with platinum tip 6F (2.0 mm) X 20 cm coaxial introducer with dilator and stiffening cannula 21G X 15 cm introducer needle with trocar stylet
MAK-NV-005*	21G X 15 cm introducer needle with trocar stylet
MAK-NV-006*	6F (2.0mm) X 20 cm coaxial introducer with dilator and stiffening cannula 21G X 15 cm introducer needle with trocar stylet
MAK-NV-007*	6F (2.0 mm) X 20 cm coaxial introducer with dilator and stiffening cannula
MAK-NV-008	0.038" (0.97 mm) X 80 cm double ended heavy duty 3 mm J guide wire 0.018" (0.46 mm) X 60 cm stainless steel guide wire with platinum tip 6F (2.0 mm) X 20 cm coaxial introducer with dilator and stiffening cannula 21G X 15 cm introducer needle with trocar stylet
MAK-NV-009	0.038" (0.97 mm) X 80 cm double ended heavy duty 3mm J guide wire 0.018" (0.46 mm) X 60 cm nitinol guide wire with platinum tip 6F (2.0 mm) X 20 cm coaxial introducer with dilator and stiffening cannula 21G X 15 cm introducer needle with trocar stylet

Packaged 10 per box *Packaged 5 per box



ORDERING INFORMATION

MAK-NV-010	0.038" (0.97 mm) X 150 cm double ended heavy duty 3 mm J guide wire 0.018" (0.46 mm) X 60 cm stainless steel guide wire with platinum tip 6F (2.0 mm) X 20 cm coaxial introducer with dilator and stiffening cannula 21G X 15 cm introducer needle with Chiba stylet
MAK-NV-011	0.038" (0.97 mm) X 150 cm double ended heavy duty 3 mm J guide wire 0.018" (0.46 mm) X 60 cm nitinol guide wire with platinum tip 6F (2.0 mm) X 20 cm coaxial introducer with dilator and stiffening cannula 21G X 15 cm introducer needle with Chiba stylet
MAK-NV-012	0.018" (0.46 mm) X 60 cm stainless steel guide wire with platinum tip 6F (2.0 mm) X 20 cm coaxial introducer with dilator and stiffening cannula 21G X 15 cm introducer needle with Chiba stylet
MAK-NV-013	0.018" (0.46 mm) X 60 cm nitinol guide wire with platinum tip 6F (2.0 mm) X 20 cm coaxial introducer with dilator and stiffening cannula 21G X 15 cm introducer needle with Chiba stylet
MAK-NV-014*	21G X 15 cm introducer needle with Chiba stylet
MAK-NV-015*	6F (2.0 mm) X 20 cm coaxial introducer with dilator and stiffening cannula 21G X 15 cm introducer needle with Chiba stylet
MAK-NV-016	0.038" (0.97 mm) X 80 cm double ended heavy duty 3 mm J guide wire 0.018" (0.46 mm) X 60 cm stainless steel guide wire with platinum tip 6F (2.0 mm) X 20 cm coaxial introducer with dilator and stiffening cannula 21G X 15 cm introducer needle with Chiba stylet
MAK-NV-017	0.038" (0.97 mm) X 80 cm double ended heavy duty 3 mm J guide wire 0.018" (0.46 mm) X 60 cm nitinol guide wire with platinum tip 6F (2.0 mm) X 20 cm coaxial introducer with dilator and stiffening cannula 21G X 15 cm introducer needle with Chiba stylet
MAK-NV-018	0.038" (0.97 mm) X 150 cm double ended heavy duty 3 mm J guide wire 0.018" (0.46 mm) X 60 cm stainless steel guide wire with platinum tip 6F (2.0 mm) X 20 cm coaxial introducer with dilator and stiffening cannula 21G X 20 cm introducer needle with trocar stylet
MAK-NV-019	0.038" (0.97 mm) X 150 cm double ended heavy duty 3 mm J guide wire 0.018" (0.46 mm) X 60 cm nitinol guide wire with platinum tip 6F (2.0 mm) X 20 cm coaxial introducer with dilator and stiffening cannula 21G X 20 cm introducer needle with trocar stylet
MAK-NV-020*	21G X 20 cm introducer needle with trocar stylet
MAK-NV-021*	21G X 20 cm introducer needle with Chiba stylet
MAK005-60	0.018" (0.46 mm) X 60 cm stainless steel guide wire with platinum tip
MAK001N-60	0.018" (0.46 mm) X 60 cm nitinol guide wire with platinum tip

Packaged 10 per box *Packaged 5 per box



REVOLUTION™

CATHETER SECUREMENT DEVICE

Saves time, provides convenience and enhances compliance with safety regulations for securing catheters with a needleless procedure.

- Secures non vascular catheters from 6F-24F.
- Hydrocolloid is absorbent around puncture site, and conforms well to body anatomy.
- Provides safe and secure fixation for up to 7 days.
- Included shower cap covers insertion site for added protection while showering.

CATALOGUE NUMBER

CFG624

Packaged 25 per box

StayFIX®

FIXATION DEVICE

StayFIX fixation device design secures catheters in place without the use of sutures, reducing the exposure to sharps and sharps-related accidents.

- Secures non-vascular catheters from 5F - 22F.
- Provides safe and secure fixation for up to 7 days.

CATALOGUE NUMBER

680ME 5F-14F - Packaged 25 per box

685ME 12F-22F - Packaged 20 per box

* *Distribution agreement limits sales to selected countries. Contact Customer Service.*

DRAINAGE BAGS

These bags come with a VELCRO® brand strap and are designed with patient mobility and comfort in mind.

- Adjustable length tubing to fit a variety of patients
- Cloth backing for patient comfort
- Twist valve for easy disposal of fluids
- Easy-to-read markings to identify fluid volume
- Internal anti-reflux valve prevents backflow

MDD600 600 mL bag

MDD600PML 600 mL bag

FZ624 600 mL cloth-backed bag

Packaged 20 per box

K08-01907 Jackson-Pratt with 100 mL bulb

K10-03450 MDD600 drainage bag with Jackson-Pratt® 100 mL bulb designed for vacuum drainage

Packaged 10 per box

ABSCCESS DRAINAGE KIT

K10-00230B OSHA-compliant closed waste system with 1400 mL MDD bag, tubing and CRV Valve
Designed for quick and safe drainage elimination.

Packaged 15 per box

DRAINAGE CONNECTOR TUBING

K10-04129 30 cm connecting tubing

K10-04130 42 cm connecting tubing

K10-04131 30 cm connecting tubing with one-way stopcock

K10-04132 42 cm connecting tubing with one-way stopcock

Packaged 25 per box

DRAINAGE TUBING

FAS5248 48" large bore tubing with pinch clamp and non-vented spike. Packaged 25 units per box.

K08-02027 Tubing kit with barbed male luer. Packaged 25 units per box.

K10-04750 48" Large bore tubing with roller clamp and 16G x 1.5" (38.1 mm) needle.
Packaged 15 per box

VALVE ADAPTER

VA-40 Designed to be back loaded over a guide wire and attached to the MAK-NV to facilitate aspiration or infusion. Accepts up to a 4F device

Packaged 10 per box

LARGE BORE STOPCOCKS

Clear body with large lumen designed to increase flow rate.

L3SNC 3-way large bore stopcock

L4SNC 4-way large bore stopcock

Packaged 50 per box

VacLok® SYRINGES

- Patented clear body design with multiple locking positions.
- Maintains negative pressure while reducing physician strain.

VAC120 20 mL VacLok Syringe

VAC130 30 mL VacLok Syringe


Packaged 20 per box

VAC110 10 mL VacLok Syringe

VAC160 60 mL VacLok Syringe

Packaged 25 per box

VELCRO® brand strap is a registered trademark of Velcro Industries B.V.
Jackson-Pratt® is a registered trademark of Cardinal Health.



BackStop®

BACKSTOP MBS100 AND BackStop+™ MBSP100

The OSHA-compliant BackStop and BackStop Plus waste basins are designed with particulate-free materials, a seal tight-lid and a 500 mL fluid capacity. Designed for quick and safe drainage disposal. *Packaged 10 per box*



MiniStop®

MINISTOP MMS100 AND MiniStop+™ MMSP100


The OSHA-compliant MiniStop and MiniStop Plus waste basins are designed with particulate-free materials, a seal tight lid and a 100 mL fluid capacity. Designed for smaller drainage procedures. *Packaged 20 per box*



DugOut®

DUGOUT DUG100 AND DUG1SS

The OSHA-compliant DugOut waste basin is designed with particulate-free materials, a secondary compartment for storage of ancillary products or saline and a seal tight lid with a 1000 mL fluid capacity. Designed for quick and safe disposal of the largest volume drainage procedures. *Packaged 6 per box*



InQwire®

GUIDE WIRES

InQwire® guide wires are precoated, creating a smooth, consistent surface for drainage catheters to track over. 80 cm length reduces excess wire and risk for contamination.

CATALOGUE NUMBER


IQ35F80B	Bentson 0.035"(0.89 mm)
6680-61	Super Stiff 0.035"(0.89 mm)
IQ35F80J1O5RS	Rosen (1.5 mm J) 0.035"(0.89 mm)
IQ35F80J3	Rosen (3 mm J-tip) 0.035"(0.89 mm)



SAFETY PARACENTESIS PROCEDURE TRAY*

Includes necessary products to complete a paracentesis procedure, including three safety needles to ensure a safe working environment.

* Not available in countries where CE mark is required.



TRAY COMPONENTS

1 each	One-Step™ centesis catheter (either 5F 7 cm or 5F 10 cm)
1 each	Fenestrated Drape w/ Adhesive Backing
1 each	Absorbent Paper Towel
1 each	ChloraPrep®
1 each	20 mL Lidocaine HCl Injection USP 1%
1 each	18G x 1" (25.4 mm) Shielded Needle
1 each	25G x 1" (25.4 mm) Shielded Needle
1 each	22G x 1.5" (38.1 mm) Shielded Needle
1 each	10 mL Luer Lock Syringe
10 each	Gauze 4"x4"
1 each	#11 Futura® Safety Scalpel
1 each	60 mL BD Luer Lock Syringe
2 each	Specimen Vials (10 mL)
1 each	Non-Vented Female Cap
1 each	Large Bore Stopcock (3-way)
1 each	Tubing Set, 48" (122cm) with Pinch Clamp and Non-Vented Spike
1 each	Bandage Dressing (1" x 3")
1 each	CSR Wrap
1 each	10 mL Medallion® Luer Lock Syringe
1 each	Lidocaine Label



CATALOGUE NUMBER

SPPT-100	(Without centesis catheter)
SPPT-5F-7S	(Includes 5F, 7 cm One-Step centesis catheter with slip fit)
SPPT-5F-10S	(Includes 5F, 10 cm One-Step centesis catheter with slip fit)
SPPT-5F-7L	(Includes 5F, 7 cm One-Step centesis catheter with fixed luer)
SPPT-5F-10L	(Includes 5F, 10 cm One-Step centesis catheter with fixed luer)


Packaged 5 per box

ChloraPrep® is a registered trademark of CareFusion 2200, Inc.

THORACENTESIS AND PARACENTESIS SET*

Includes necessary products to complete a thoracentesis and paracentesis procedure, including three safety needles to ensure a safe working environment.

* Not available in countries where CE mark is required.



TRAY COMPONENTS

- 1 each One-Step™ Valved Centesis Catheter (either 5F 7 cm or 5F 10 cm)
- 1 each Fenestrated Drape w/ Adhesive Backing
- 1 each Absorbent Paper Towel
- 1 each ChloraPrep®
- 1 each 20 mL Lidocaine HCl Injection USP 1%
- 1 each 18G x 1" (25.4 mm) Shielded Needle
- 1 each 25G x 1" (25.4 mm) Shielded Needle
- 1 each 22G x 1.5" (38.1 mm) Shielded Needle
- 1 each 10 mL Luer Lock Syringe
- 10 each Gauze 4"x4"
- 1 each #11 Futura® Safety Scalpel
- 1 each 60 mL BD Luer Lock Syringe
- 2 each Specimen Vials (10 mL)
- 1 each Non-Vented Female Cap
- 1 each Large Bore Stopcock (3-way)
- 1 each Tubing Set, 48" (122 cm) with Pinch Clamp and Non-Vented Spike
- 1 each Tubing Set, 48" (122 cm) with Roller Clamp and 16G x 1.5" needle
- 1 each Bandage Dressing (1" x 3")
- 1 each CSR Wrap
- 1 each 10 mL Medallion® Luer Lock Syringe
- 1 each Lidocaine Label

CATALOGUE NUMBER

TAPS-100	(Without centesis catheter)
TAPS-5F-7T	(Includes 5F, 7 cm One-Step valved centesis catheter with slip fit)
TAPS-5F-10T	(Includes 5F, 10 cm One-Step valved centesis catheter with slip fit)
TAPS-5F-7LT	(Includes 5F, 7 cm One-Step valved centesis catheter with fixed luer)
TAPS-5F-10LT	(Includes 5F, 10 cm One-Step valved centesis catheter with fixed luer)

Packaged 5 per box

ChloraPrep® is a registered trademark of CareFusion 2200, Inc.







Merit Medical Europe, Middle East & Africa

Amerikalaan 42, 6199 AE Maastricht-Airport, Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.

Parkmore Business Park West, Galway, Ireland
+35 39 170 37 33

Merit Medical France

136 Ave. Joseph Kessel, 78960 Voisins Le Bretonneux, France
+33 (0)1 30 57 24 24

Merit Medical Systems, Inc.

1600 West Merit Parkway, South Jordan, UT, 84095 USA
+1 801 253 1600

Asia, Pacific Rim, Latin America, Caribbean & Canada

+1 801 208 4438

Free phone for specific country

Austria

0800 295 374

Belgium

0800 72 906 (Dutch)
0800 73 172 (French)

Denmark

80 88 00 24

Finland

0800 770 586

France

0800 91 60 30

Germany

0800 182 0871

Ireland (Republic)

1800 553 163

Italy

800 897 005

Luxembourg

8002 25 22

Netherlands

0800 022 81 84

Norway

800 11629

Sweden

020 792 445

Sweden

020 792 445

Switzerland

00800 9000 5060

UK

0800 973 115

400797002/C ID 140613