




Solitaire™ **X**
Revascularization Device
Portfolio



3MM
NOW AVAILABLE

Medtronic



SOLITAIRE™ X REVASCULARIZATION DEVICE PORTFOLIO INFORMATION¹

Model	Recommended Vessel Diameter ^A (mm)		Microcatheter ID Range	Push Wire Length	Stent Diameter	Usable Length ^B	Stent Length	Length from Distal Tip to Fluoroscope Marker	Radiopaque Markers		Radiopaque Stent Markers Spacing
	(min)	(max)							(min - max)	Distal	
SFR4-3-20-10	1.5	3.0	0.017" – 0.027" 0.43mm – 0.69mm	200	3.0	20.0	30.6	< 150	3	1	10
SFR4-4-20-05	1.5	4.0	0.021" – 0.027" 0.53mm – 0.69mm	200	4.0	20.0	31.0	<130	3	1	5
SFR4-4-20-10	1.5	4.0	0.021" – 0.027" 0.53mm – 0.69mm	200	4.0	20.0	31.0	<130	3	1	10
SFR4-4-40-10	1.5	4.0	0.021" – 0.027" 0.53mm – 0.69mm	200	4.0	40.0	50.0	<130	3	1	10
SFR4-6-20-10	2.0	5.5	0.021" – 0.027" 0.53mm – 0.69mm	200	6.0	20.0	31.0	<130	4	1	10
SFR4-6-24-06	2.0	5.5	0.021" – 0.027" 0.53mm – 0.69mm	200	6.0	24.0	37.0	<130	4	1	6
SFR4-6-40-10	2.0	5.5	0.021" – 0.027" 0.53mm – 0.69mm	200	6.0	40.0	47.0	<130	4	1	10

A. Based on the smallest vessel diameter at thrombus site.

B. Usable length that is at least as long as the length of the thrombus.

Up to 3 retrieval passes¹

Medtronic

9775 Toledo Way
Irvine, CA 92618
USA
Tel 877.526.7890
Fax 763.526.7888

medtronic.com/solitaire

1. M003592CDOC2

CAUTION: Federal (USA) law restricts these devices to sale distribution and use by or on order of a physician. Indications, contraindications, warnings and instructions for use for Solitaire™ X Revascularization Device can be viewed at www.medtronic.com/manuals. Indications, contraindications, warnings and instructions for all other products can be found in the product labeling supplied with each device.

© 2021 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. UC201911494aEN