

SYNTERVENTION

SWASHER[®]

Ultra Low-Lint Surgical Towel



The SWASHER Ultra Low-Lint Surgical Towel

The SWASHER is a sterile, biocompatible surgical towel designed to reduce lint and particulate contamination in interventional and surgical procedures.

89% Less Lint


- 89% less lint than standard cotton surgical towels.
- Place medical instruments on the SWASHER without worrying about contaminant transfer.




Simple lint roller test performed on cotton towel.



Simple lint roller test performed on SWASHER.



Cotton towel fibers found on guide wire.



Cotton towel fibers trapped on stent strut.

SWASHER has 89% Less Lint than Standard Cotton Surgical Towels*

Product	10 µm	25 µm	50 µm	100 µm
Cotton Towel	333,740	23,200	11,400	5,000
SWASHER	38,600	730	130	0

*Independent LPC testing. Data on file at Syntervention, Inc.


Economical

- Weighs less to save on shipping and medical waste disposal costs.
- Durable and more absorbent for extended use.
- Use our calculators to estimate how much your facility could save. Visit syntervention.com/calculators.

Better Absorbency

- Excellent fluid capacity with SWASHER.
- See the SWASHER in action at syntervention.com/video-library.

SWASHER Weighs Less and Absorbs More





SWASHER Order Information

Number	Description	Size	Case Quantity
5010B	SWASHER, Blue, 2 per pack	17 x 26	100 units/case
5011B	SWASHER, Blue, 4 per pack	17 x 26	50 units/case

Remit To	Syntervention, Inc., 1444 Jeffreys Road, Suite 133 Rocky Mount, NC 27804
Contact Us	Phone: 888.505.8802 Fax: 888.853.6559 Email: info@syntervention.com
Tax ID	45-3609564
Payment Terms	Net 30 days.
Shipment	Items are shipped ground unless otherwise specified.
Product Storage	Store in a cool, dry place. Avoid exposure to moisture and excessive heat.
Product Sterilization	Sterilized with EO or gamma radiation. Non-sterile products available upon request.
Shelf Life	See package label for expiration date. 1-3 years.
Latex and DEHP Status	Latex, BPA and DEHP free.
Returned Goods Policy	A Returned Goods Authorization (RGA) number is required for all returns and replacements. Authorization can be obtained by calling Syntervention Customer Service at 888.505.8802.

To view calculators visit www.syntervention.com/calculators.

To view the Learning Center Video Library visit www.syntervention.com/video-library.

Syntervention and SWASHER, are trademarks or registered trademarks of Syntervention, Inc.



SYNTERVENTION

Syntervention, Inc.
1444 Jeffreys Road, Suite 133
Rocky Mount, NC 27804
Customer Service: 888.505.8802
Fax: 888.853.6559
www.syntervention.com