

HEMO-CATH[®]



HEMO-CATH[®] LT LONG-TERM HEMODIALYSIS CATHETER



Circle "C" Internal Lumen Design



Step Tip Design



Straight & Pre-Curved
Configurations



8F & 12.5F Available

HEMO-CATH[®] ST SHORT-TERM HEMODIALYSIS CATHETER



Double "D" Internal Lumen Design



Tapered Tip Design



Straight Configuration



8F & 11.5F Available

HEMO-CATH® LT

LONG-TERM HEMODIALYSIS CATHETER

HEMO-CATH® ST

SHORT-TERM HEMODIALYSIS CATHETER

8F BASIC SET

DSL18P	8F X 18CM HEMO-CATH® LT SET	5/BOX
DSL24P	8F X 24CM HEMO-CATH® LT SET	5/BOX

Set Contents: (1) Catheter (1) Introducer Needle (1) Guidewire (1) Tunneler (1) Peelable Introducer (1) Adhesive Wound Dressing (1) 10cc Syringe (1) Hemo-Cath® Hub Clip (2) End Caps (1) Patient ID Card

12.5F BASIC SET

DMC101241	12.5F X 15CM HEMO-CATH® LT STRAIGHT	5/BOX
DMC101242	12.5F X 18CM HEMO-CATH® LT STRAIGHT	5/BOX
DMC101243	12.5F X 24CM HEMO-CATH® LT STRAIGHT	5/BOX
DSL28	12.5F X 28CM HEMO-CATH® LT STRAIGHT	5/BOX
DSL32	12.5F X 32CM HEMO-CATH® LT STRAIGHT	5/BOX
DSL28PC	12.5F X 28CM PC HEMO-CATH® LT PRE-CURVED	5/BOX
DSL32PC	12.5F X 32CM PC HEMO-CATH® LT PRE-CURVED	5/BOX

Set Contents: (1) Catheter (1) Introducer Needle (1) Guidewire (1) Tunneler (1) Scalpel (1) Peelable Introducer (1) Adhesive Wound Dressing (1) 10cc Syringe (1) Hemo-Cath® Hub Clip (2) End Caps

12.5F FULL TRAY

MCDLTSL28	12.5F X 28CM HEMO-CATH® LT STRAIGHT	5/BOX
MCDLTSL32	12.5F X 32CM HEMO-CATH® LT STRAIGHT	5/BOX
MCDLTSL28PC	12.5F X 28CM PC HEMO-CATH® LT PRE-CURVED	5/BOX
MCDLTSL32PC	12.5F X 32CM PC HEMO-CATH® LT PRE-CURVED	5/BOX

Set Contents: (1) Catheter (1) Introducer Needle (1) Guidewire (1) Tunneler (1) Peelable Introducer (1) Adhesive Wound Dressing (1) Chloraprep (2) 5cc Syringe (2) 10cc Syringes (3) 4x4 Gauze (3) Safety Hypodermic Needles (1) 5mL Lidocaine (1) Filter Straw (1) Scalpel (1) Surgical Gloves (1) 2-0 Silk Sutures (1) Fenestrated Drape (1) Hemostat (1) Lidocaine/Saline Sticker Sheet (1) Hemo-Cath® Hub Clip (2) End Caps

IMPORTANT RISK INFORMATION

LONG TERM

Indications for Use: The Medcomp® Hemo-Cath® LT Silicone Double Lumen Catheter can be utilized for long term implantation as well as temporary access for hemodialysis, hemoperfusion, or apheresis therapy. The cannula may be inserted via the Seldinger Technique due to the inner Teflon stylet, increasing linear strength. The stylet is removed after insertion, leaving the soft silicone cannula in the body. The flexible silicone make-up conforms well to the vessel anatomy, resulting in higher patient tolerance during extended use.

Contraindications: The Subclavian Approach is NOT recommended for use with the Medcomp® Double Lumen Subclavian-Femoral Catheter in Hemodialysis or Hemoperfusion Procedures used for the management of acute poisoning or other situations in which a ventilator might be used due to risk of traumatic pneumothorax posing a dangerous complication for the patient.

SHORT TERM

Indications for Use: The Medcomp® Hemo-Cath® ST Silicone Double Lumen Catheter can be utilized for temporary access for hemodialysis, hemoperfusion, or apheresis therapy. The cannula may be inserted via the Seldinger technique due to the inner Teflon stylet, increasing linear strength. The stylet is removed after insertion, leaving the soft silicone cannula in the body. The flexible silicone make-up conforms well to the vessel anatomy, resulting in higher patient tolerance during extended use.

Contraindications: The Subclavian Approach is NOT recommended for use with the Medcomp® Hemo-Cath® ST Double Lumen Subclavian-Femoral Catheter in Hemodialysis or Hemoperfusion Procedures used for the management of acute poisoning or other situations in which a ventilator might be used due to risk of traumatic pneumothorax posing a dangerous complication for the patient.

Refer to Instructions for Use provided with the product for complete instructions, warnings, precautions, and contraindications. Observe all instructions for use prior to using products. Failure to do so may result in patient complications.

8F BASIC SET

DSL12P	8F X 12CM HEMO-CATH® ST	5/BOX
--------	-------------------------	-------

Set Contents: (1) Catheter (1) Introducer Needle (1) Guidewire (1) Dilator (1) Peelable Introducer (1) Adhesive Wound Dressing (1) 10cc Syringe (1) Hemo-Cath® Hub Clip (2) End Caps (1) Patient ID Card

11.5F BASIC SET

DSL12	11.5F X 12CM HEMO-CATH® ST SET	5/BOX
DSL15	11.5F X 15CM HEMO-CATH® ST SET	5/BOX
DSL20	11.5F X 20CM HEMO-CATH® ST SET	5/BOX
DSL24	11.5F X 24CM HEMO-CATH® ST SET	5/BOX
DSL27	11.5F X 27CM HEMO-CATH® ST SET	5/BOX

Set Contents: (1) Catheter (1) Introducer Needle (1) Guidewire (1) Dilator (1) Peelable Introducer (1) Adhesive Wound Dressing (1) Hemo-Cath® Hub Clip (2) End Caps (1) Patient ID Card

11.5F FULL TRAY

MCDLTSL15	11.5F X 15CM HEMO-CATH® ST SET	5/BOX
MCDLTSL20	11.5F X 20CM HEMO-CATH® ST SET	5/BOX
MCDLTSL24	11.5F X 24CM HEMO-CATH® ST SET	5/BOX
MCDLTSL27	11.5F X 27CM HEMO-CATH® ST SET	5/BOX

Set Contents: (1) Catheter (1) Introducer Needle (1) Guidewire (1) Dilator (1) Peelable Introducer (1) Adhesive Wound Dressing (2) Chloraprep (2) 5cc Syringe (2) 10cc Pre-Filled Saline Syringes (3) 4x4 Gauze (3) Safety Hypodermic Needles (1) 5mL Lidocaine (1) Filter Straw (1) Scalpel (1) Surgical Gloves (1) 2-0 Silk Sutures (1) Fenestrated Drape (1) Hemostat (1) Needle Stick Pad (1) Lidocaine/Saline Sticker Sheet (1) Hemo-Cath® Hub Clip (2) End Caps