

# Arrow Edge® Hemodialysis Catheterization Kit

## Contents:

- 1 Catheter: 15 Fr.
- 1 Arrow® SmartSeal® Hemostatic Peelable Dialysis Sheath: 16 Fr.
- 1 Spring Wire Guide: .038" (.97 mm) x 39-1/2" (100 cm) straight stiff tip on one end — "J" tip on other with Arrow Advancer™ Guidewire Advancer
- 1 Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm)
- 1 Pre-Loaded Tunneler (metal): with tunneling sheath
- 2 Luer-Lock Caps
- 1 Dressing: Tegaderm® 10 cm x 12 cm
- 1 Tissue Dilator: 12 Fr.
- 1 Tissue Dilator: 14 Fr.
- 1 Safety Scalpel: #11
- 1 SharpsAway® II Locking Disposal Cup



Arrow Edge® Hemodialysis  
Catheterization Kit (CS-15282-I)

To learn more about Arrow Edge® Catheters:  
contact your Teleflex Representative

Teleflex is a global provider of medical technologies designed to improve the health and quality of people's lives. We apply purpose-driven innovation – a relentless pursuit of identifying unmet clinical needs – to benefit patients and healthcare providers. Our portfolio is diverse, with solutions in the fields of vascular and interventional access, surgical, anesthesia, cardiac care, urology, emergency medicine and respiratory care. Teleflex employees worldwide are united in the understanding that what we do every day makes a difference. For more information, please visit [teleflex.com](http://teleflex.com).

Teleflex is the home of Arrow®, Deknatel®, Hudson RCI®, LMA®, Pilling®, Rüscher®, UroLift®, and Weck® – trusted brands united by a common sense of purpose.

## Corporate Office

Phone +1 610 225 6800, 550 E. Swedesford Road, Suite 400, Wayne, PA 19087, USA

## Regional Offices

**United States:** Phone +1 919 544 8000, Toll Free 866 246 6990, [cs@teleflex.com](mailto:cs@teleflex.com), 3015 Carrington Mill Boulevard, Morrisville, NC 27560, USA

**Latin America:** Phone +1 919 433 4999, [la.cs@teleflex.com](mailto:la.cs@teleflex.com), 3015 Carrington Mill Boulevard, Morrisville, NC 27560, USA

**International:** Phone +353 (0)9 06 46 08 00, [orders.intl@teleflex.com](mailto:orders.intl@teleflex.com), Teleflex Medical Europe Ltd., IDA Business and Technology Park, Dublin Road, Athlone, Co Westmeath, Ireland

Does not contain natural rubber latex.  
All Chronic Catheters sold individually.

Caution: U.S. federal law limits this device to sale by or on order of a physician. Contents of unopened, undamaged package are sterile. Disposable. Refer to Instructions for Use (IFU) for current warnings, indications, contraindications, precautions and Instructions For Use.

Contact Teleflex for reference information on the above.


## REFERENCES:

- 1 Alcohol, alcohol-based solutions (e.g., Hibiclens®, Chloraprep®), iodine-based solutions (Povidone-Iodine), PEG-based ointments (e.g., Bactroban®), hydrogen peroxide or ExSept Plus® are accepted for use with this catheter.
- 2 In vitro test performed by an independent laboratory (Citech, Plymouth Meeting, PA).


Tegaderm is a registered trademark of 3M Company. Bactroban is a registered trademark of GlaxoSmithKline LLC. Hibiclens is a registered trademark of Regent Medical Limited. Chloraprep is a registered trademark of Carefusion 2200, Inc. ExSept Plus is a registered trademark of Angelini Pharma Inc.

Teleflex, the Teleflex logo, Arrow, Arrow Advancer, Arrow Edge, Arrow Simplicity, SharpsAway and SmartSeal are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries. Information in this material is not a substitute for the product Instructions for Use. Not all products may be available in all countries. Please contact your local representative. Revised: 05/2021.

© 2021 Teleflex Incorporated. All rights reserved. MC-006568 Rev 1




## Arrow Edge® Antegrade-Tunneled Chronic Hemodialysis Catheter




# Arrow Edge® Antegrade-Tunneled Chronic Hemodialysis Catheter

## High Flow, Optimal Performance




**Arrow® SmartSeal®**  
Hemostatic peelable dialysis sheath-featured in kits



Arrow® SmartSeal® Hemostatic Peelable Dialysis Sheath

Conventional valve

**Arrow Edge® Catheter: Pressure vs. Flow Rate<sup>2</sup>**



**Delivers high flow rates up to 400 cc/min. at low pressures<sup>2</sup>**

Crescent-shaped distal arterial tip allows blood flow even if catheter tip rests against vessel wall.

## Ordering Information

### Arrow Edge® Hemodialysis Catheterization Kit with SmartSeal® Hemostatic Sheath and Arrow Simplicity® Color-Coded Micro-Introducer

COLOR CODE	PRODUCT NO.	DESCRIPTION	TIP-TO-CUFF INSERTION LENGTH (CM)	CE
19	CS-15242-IM	Arrow 15 Fr. x 9.5" (24 cm) Radiopaque Polyurethane Catheter 21 Ga. Micro-Introducer Set	19	N
23	CS-15282-IM	Arrow 15 Fr. x 11" (28 cm) Radiopaque Polyurethane Catheter 21 Ga. Micro-Introducer Set	23	N
27	CS-15322-IM	Arrow 15 Fr. x 12.5" (32 cm) Radiopaque Polyurethane Catheter 21 Ga. Micro-Introducer Set	27	N
31	CS-15362-IM	Arrow 15 Fr. x 14" (36 cm) Radiopaque Polyurethane Catheter 21 Ga. Micro-Introducer Set	31	N

### Arrow Edge® Hemodialysis Catheterization Kit with SmartSeal® Hemostatic Sheath

COLOR CODE	PRODUCT NO.	DESCRIPTION	TIP-TO-CUFF INSERTION LENGTH (CM)	CE
19	CS-15242-I	Arrow 15 Fr. x 9.5" (24 cm) Radiopaque Polyurethane Catheter	19	Y
23	CS-15282-I	Arrow 15 Fr. x 11" (28 cm) Radiopaque Polyurethane Catheter	23	Y
27	CS-15322-I	Arrow 15 Fr. x 12.5" (32 cm) Radiopaque Polyurethane Catheter	27	Y
31	CS-15362-I	Arrow 15 Fr. x 14" (36 cm) Radiopaque Polyurethane Catheter	31	Y
	CS-15552-I	Arrow 15 Fr. x 22" (55 cm) Radiopaque Polyurethane Catheter	50	Y

### Arrow Simplicity® Micro-Introducer Set

COLOR CODE	PRODUCT #	DESCRIPTION	QTY/BOX
	MP-04450-AS	Color-Coded Micro-Introducer Set	5 ea.

### Arrow Edge® Hemodialysis Catheterization Set with SmartSeal® Hemostatic Sheath

COLOR CODE	PRODUCT NO.	DESCRIPTION	TIP-TO-CUFF INSERTION LENGTH (CM)
19	CSD-15242-I	Edge, Antegrade, HEMODIALYSIS SET: 2-L15 FRX 19 cm	19
23	CSD-15282-I	Edge, Antegrade, HEMODIALYSIS SET: 2-L15 FRX 23 cm	23
27	CSD-15322-I	Edge, Antegrade, HEMODIALYSIS SET: 2-L15 FRX 27 cm	27
31	CSD-15362-I	Edge, Antegrade, HEMODIALYSIS SET: 2-L15 FRX 31 cm	31
	CSD-15552-I	Edge, Antegrade, HEMODIALYSIS SET: 2-L15 FRX 50 cm	50



Colored bands on wire guide correspond to color-coded catheter packaging. Arrow Simplicity Micro-Introducer Set can only be used with Arrow Chronic Hemodialysis Catheters.

### Arrow Simplicity® Micro-Introducer Set Contents:

- 1 Introducer Needle: 21 Ga. x 2.75" (7 cm) thin wall echogenic tip needle with bevel indicator and green hub
- 1 Coaxial Sheath/Dilator: 4 Fr. x 4" (10 cm)
- 1 Micro-Introducer Wire Guide: 0.018" (0.46 mm) x 45 cm stainless steel guidewire with platinum coil and color band markings