

Vascular Stent Systems

PROVEN Long Term Primary Patency through **10 years** with the **S.M.A.R.T.**[®] Vascular Stent Systems

Long Term Outcomes of S.M.A.R.T. Stent Implantation in Fem-Pop¹

	Years After SMART Implantation				
	1	3	5	7	10
Primary Patency					
Overall	87.9	79.4	68.1	53.7	50.3
TASC IIA/B	92.3	84.9	76.5	64.7	62.0
TASC IIC/D	80.0	70.4	58.9	34.4	NA
Run-Off ≥ 2	91.2	83.7	69.8	59.6	59.6
Run-Off ≤ 2	85.0	76.2	67.0	48.9	42.5
Secondary Patency	96.5	91.7	85.0	73.8	67.7

- 319 limbs in 269 patients
- Critical Limb Ischemia (CLI) observed in 14.1% of limbs
- Mean target lesion length was 118.5±81.0 mm
- 57.7% of lesions had poor distal runoff
- 37.3% classified as TASC IIC/D
- 41.1% CTO

SMART Vascular Stent System and SMART Control Vascular Stent System



Ordering Information

Product Codes		Stent Specifications			
80 cm Usable Length	120 cm Usable Length	Indication	Stent Diameter (mm)	Stent Length (mm)	Recommended Vessel Size (mm)
C06020SL	C06020ML	Iliac SFA	6	20	4-5
C06030SL	C06030ML	Iliac SFA	6	30	4-5
C06040SL	C06040ML	Iliac SFA	6	40	4-5
C06060SL	C06060ML	Iliac SFA	6	60	4-5
C06080SL	C06080ML	Iliac SFA	6	80	4-5
C06100SL	C06100ML	Iliac SFA	6	100	4-5
	C06120ML	SFA	6	120	4-5
	C06150ML	SFA	6	150	4-5
C07020SL	C07020ML	Iliac SFA	7	20	5-6
C07030SL	C07030ML	Iliac SFA	7	30	5-6
C07040SL	C07040ML	Iliac SFA	7	40	5-6
C07060SL	C07060ML	Iliac SFA	7	60	5-6
C07080SL	C07080ML	Iliac SFA	7	80	5-6
C07100SL	C07100ML	Iliac SFA	7	100	5-6
	C07120ML	SFA	7	120	5-6
	C07150ML	SFA	7	150	5-6
C08020SL	C08020ML	Iliac SFA	8	20	6-7
C08030SL	C08030ML	Iliac SFA	8	30	6-7
C08040SL	C08040ML	Iliac SFA	8	40	6-7
C08060SL	C08060ML	Iliac SFA	8	60	6-7
C08080SL	C08080ML	Iliac SFA	8	80	6-7
C08100SL	C08100ML	Iliac SFA	8	100	6-7
	C08120ML	SFA	8	120	6-7
	C08150ML	SFA	8	150	6-7
C09020SL	C09020ML	Iliac	9	20	7-8
C09030SL	C09030ML	Iliac	9	30	7-8
C09040SL	C09040ML	Iliac	9	40	7-8
C09060SL	C09060ML	Iliac	9	60	7-8
C10020SL	C10020ML	Iliac	10	20	8-9
C10030SL	C10030ML	Iliac	10	30	8-9
C10040SL	C10040ML	Iliac	10	40	8-9
C10060SL	C10060ML	Iliac	10	60	8-9

For customer service, call **1.800.327.7714**.

For more information, visit **cordis.com**

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

IMPORTANT INFORMATION: Prior to use, refer to the "Instructions for Use" supplied with these devices for indications, contraindications, side effects, suggested procedures, warnings and precautions.

CORDIS, Cordis LOGO, and S.M.A.R.T. are trademarks of Cardinal Health and may be registered in the US and/or in other countries.

© 2018 Cardinal Health. All Rights Reserved.

100517766 6/2018