



Merit Laureate[®]

Hydrophilic Guide Wire




Merit Laureate[®]


Hydrophilic Guide Wire

EXPERIENCE THE LATEST IN HYDROPHILIC TECHNOLOGY

The Merit Laureate hydrophilic guide wire provides an optimum level of lubricity to cross the most difficult lesions and facilitates successful clinical outcomes in the most challenging procedures.



The resilient hydrated gel coating provides a constant fluid layer on the wire surface for optimum lubricity and long lasting durability.





Nitinol core wire and radiopaque jacket are designed to provide enhanced visibility.

Designed to provide true 1:1 torque response to facilitate navigation and rapid vessel selection.




MARKET LEADER

MERIT LAUREATE




The market leader shows a slow torque response followed by rapid exaggerated movement.

The Merit Laureate shows definitive turn for turn responsiveness (1:1 torque).



Comparative torque testing showed the Merit Laureate guide wire consistently met the 1:1 torque target.

*Bench test results may not necessarily be indicative of clinical performance.



Catalogue Number	Configuration	Tip Shape	Diameter	Length	Description
0.035" STANDARD ANGLED					
LWSTDA3580	Standard	Angled	0.035" (0.89 mm)	80 cm	
LWSTDA35150	Standard	Angled	0.035" (0.89 mm)	150 cm	
LWSTDA35180	Standard	Angled	0.035" (0.89 mm)	180 cm	
LWSTDA35220	Standard	Angled	0.035" (0.89 mm)	220 cm	
LWSTDA35260EX	Standard	Angled	0.035" (0.89 mm)	260 cm	Exchange
0.035" STANDARD STRAIGHT					
LWSTDS3580	Standard	Straight	0.035" (0.89 mm)	80 cm	
LWSTDS35150	Standard	Straight	0.035" (0.89 mm)	150 cm	
LWSTDS35180	Standard	Straight	0.035" (0.89 mm)	180 cm	
LWSTDS35220	Standard	Straight	0.035" (0.89 mm)	220 cm	
LWSTDS35260EX	Standard	Straight	0.035" (0.89 mm)	260 cm	Exchange
0.035" STIFF ANGLED					
LWSTFA3580	Stiff	Angled	0.035" (0.89 mm)	80 cm	
LWSTFA35150	Stiff	Angled	0.035" (0.89 mm)	150 cm	
LWSTFA35180	Stiff	Angled	0.035" (0.89 mm)	180 cm	
LWSTFA35220	Stiff	Angled	0.035" (0.89 mm)	220 cm	
LWSTFA35260EX	Stiff	Angled	0.035" (0.89 mm)	260 cm	Exchange
0.035" STIFF STRAIGHT					
LWSTFS3580	Stiff	Straight	0.035" (0.89 mm)	80 cm	
LWSTFS35150	Stiff	Straight	0.035" (0.89 mm)	150 cm	
LWSTFS35180	Stiff	Straight	0.035" (0.89 mm)	180 cm	
LWSTFS35220	Stiff	Straight	0.035" (0.89 mm)	220 cm	
LWSTFS35260EX	Stiff	Straight	0.035" (0.89 mm)	260 cm	Exchange
0.038" STANDARD ANGLED					
LWSTDA3880	Standard	Angled	0.038" (0.97 mm)	80 cm	
LWSTDA38120	Standard	Angled	0.038" (0.97 mm)	120 cm	
LWSTDA38150	Standard	Angled	0.038" (0.97 mm)	150 cm	
LWSTDA38180	Standard	Angled	0.038" (0.97 mm)	180 cm	
LWSTDA38260EX	Standard	Angled	0.038" (0.97 mm)	260 cm	Exchange
0.038" STANDARD STRAIGHT					
LWSTDS38150	Standard	Straight	0.038" (0.97 mm)	150 cm	
LWSTDS38180	Standard	Straight	0.038" (0.97 mm)	180 cm	
LWSTDS38260EX	Standard	Straight	0.038" (0.97 mm)	260 cm	Exchange
0.038" STIFF ANGLED					
LWSTFA3880	Stiff	Angled	0.038" (0.97 mm)	80 cm	
LWSTFA38150	Stiff	Angled	0.038" (0.97 mm)	150 cm	
LWSTFA38180	Stiff	Angled	0.038" (0.97 mm)	180 cm	
0.038" STIFF STRAIGHT					
LWSTFS38150	Stiff	Straight	0.038" (0.97 mm)	150 cm	
LWSTFS38180	Stiff	Straight	0.038" (0.97 mm)	180 cm	
0.025" STANDARD ANGLED					
LWSTDA25150	Standard	Angled	0.025" (0.64 mm)	150 cm	
LWSTDA25180	Standard	Angled	0.025" (0.64 mm)	180 cm	
LWSTDA25260EX	Standard	Angled	0.025" (0.64 mm)	260 cm	Exchange
0.025" STANDARD STRAIGHT					
LWSTDS25150	Standard	Straight	0.025" (0.64 mm)	150 cm	
LWSTDS25180	Standard	Straight	0.025" (0.64 mm)	180 cm	
LWSTDS25260EX	Standard	Straight	0.025" (0.64 mm)	260 cm	Exchange
0.018" STANDARD ANGLED					
LWSTDA1880	Standard	Angled	0.018" (0.46 mm)	80 cm	
LWSTDA18150	Standard	Angled	0.018" (0.46 mm)	150 cm	
LWSTDA18180	Standard	Angled	0.018" (0.46 mm)	180 cm	
LWSTDA18260EX	Standard	Angled	0.018" (0.46 mm)	260 cm	Exchange
0.018" STANDARD STRAIGHT					
LWSTDS18150	Standard	Straight	0.018" (0.46 mm)	150 cm	
LWSTDS18180	Standard	Straight	0.018" (0.46 mm)	180 cm	
LWSTDS18260EX	Standard	Straight	0.018" (0.46 mm)	260 cm	Exchange



Understand. Innovate. Deliver.™

Merit Medical France
136 Ave. Joseph Kessel,
78960 Voisins Le Bretonneux, France
+33 (0)1 30 57 24 24

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
+1 801-208-4300
+1 800-35-MERIT

**Merit Medical Europe, Middle East,
& Africa (EMEA)**
Amerikalaan 42, 6199 AE Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733

Austria
0800 295 374

Belgium
0800 72 906 (Dutch)
0800 73 172 (Français)

Denmark
80 88 00 24

Finland
0800 770 586

France
0800 91 60 30

Germany
0800 182 0871

Ireland (Republic)
1800 553 163

Italy
800 897 005

Luxembourg
8002 25 22

Netherlands
0800 022 81 84

Norway
800 11629

Portugal
308 801 034

Russia
+7 495 221 89 02

Spain
+34 911238406

Sweden
020 792 445

Switzerland
(Deutsch)

+41 225180252
(Français)

+41 225948000
(Italiano)

+41 225180035
UK
0800 973 115