



# Merit Advance<sup>®</sup>

Angiographic Safety Needle

Sticking with Safety.



Customers rely on Merit Medical when  
needlestick prevention matters most.

# Merit Advance®

Angiographic Safety Needle

The unique design of the Merit Advance Angiographic Safety Needle provides an easy to activate, low profile safety mechanism, which shields the bevel of the needle after use.

## FEATURES & BENEFITS

### LOW PROFILE SAFETY MECHANISM

Less likely to change or affect clinical technique.

### COLOR CODED SAFETY MECHANISM

For fast identification of needle gauge size.

### ERGONOMIC HUB DESIGN

Facilitates grip, stability, and needle control during use.

### TRANSPARENT POLYCARBONATE HUB

For instant bleed-back confirmation.

### SHARP CANNULA BEVEL

Provides quick access to the vasculature and minimizes trauma to surrounding tissue.



## ORDERING INFORMATION

| Catalog Number | Gauge | Length        | Maximum Guide Wire | Hub           |
|----------------|-------|---------------|--------------------|---------------|
| SL18T71W       | 18    | 7.0cm / 2.75" | 0.038 (0.97mm)     | without wings |
| SL18T71WS      | 18    | 7.0cm / 2.75" | 0.038 (0.97mm)     | with wings    |
| SL21T71W       | 21    | 7.0cm / 2.75" | 0.021 (0.53mm)     | without wings |
| SL21T71WS      | 21    | 7.0cm / 2.75" | 0.021 (0.53mm)     | with wings    |

Packaged 25 per box – 4 boxes per case.  
Latex free.



Understand. Innovate. Deliver.™

Merit Medical Systems, Inc.  
1600 West Merit Parkway  
South Jordan, Utah 84095  
1.801.208.4300  
1.800.35.MERIT

Merit Medical Europe, Middle East, & Africa (EMEA)  
Amerikalaan 42, 6199 AE Maastricht-Airport  
The Netherlands  
+31 43 358 82 22

Merit Medical Ireland Ltd.  
Parkmore Business Park West  
Galway, Ireland  
+353 (0) 91 703 733

Merit.com

Austria  
0800 295 374

Belgium  
0800 72 906 (Dutch)  
0800 73 172 (Français)

Denmark  
80 88 00 24

Finland  
0800 770 586

France  
0800 91 60 30

Germany  
0800 182 0871

Ireland (Republic)  
1800 553 163

Italy  
800 897 005

Luxembourg  
8002 25 22

Netherlands  
0800 022 81 84

Norway  
800 11629

Portugal  
308 801 034

Russia  
+7 495 221 89 02

Spain  
+34 911238406

Sweden  
020 792 445

Switzerland  
(Deutsch)  
+41 225180252  
(Français)  
+41 225948000  
(Italiano)  
+41 225180035  
UK  
0800 973 115