



M·CATH

Excellent control in CTO


M-CATH™

Microcatheter

- Small entry profile: 1.6 Fr
- Unique shaft design
- High push proximal shaft
- Enhanced torquability
- High trackability




Tapered and robust tip for enhanced crossability



Hybrid shaft design for maximum pushability and torquability



Ergonomic design for high comfort and control



Ordering Information / Technical Information

Reference	Useable Length	Entry Profile	Dist. Outer Diameter	Max. Outer Diameter	Dist. Inner Diameter	Prox. Inner Diameter	Max. Guidewire Size
750135730	135 cm	1.6 Fr 0.53 mm	2.25 Fr 0.75 mm	3.3 Fr 1.10 mm	0.016" 0.40 mm	0.018" 0.45 mm	0.014" 0.36 mm

Excellent control in CTO