



### Catheter Specifications

Part Number	Inner Diameter	Outer Diameter Proximal Portion	Outer Diameter Distal Portion	Usable Length	Hydrophilic Coating Length	Max Guide Wire	Maximum Injection Pressure	Distance From Tip To 1 <sup>st</sup> Markerband	Distance Between 1 <sup>st</sup> And 2 <sup>nd</sup> Markerbands	Max. Microsphere Size	Max. Coil Size	Cath. Vol.
MIV-20500	0.021" (0.54 mm)	2.9F (0.97 mm)	2.4F (0.80 mm)	125 cm	80 cm	≤ 0.018" (0.46 mm)	6,900 kPa (1,000 psi)	0.5 mm	13.5 mm	700 μm	0.018"	0.49 mL

Distal Tip Cross Section



Close-up of Steerable Tip



### Microsphere Compatibility

#### QuadraspHERE® Microspheres

Dry Particle Size Range	Hydrated Size Range	Color Code	Product Number
30-60 μm	120-240 μm	Orange	V225QS
50-100 μm	200-400 μm	Green	V325QS
100-150 μm	400-600 μm	Blue	V525QS



#### Embosphere® Microspheres

Size Range	Color Code	Code pre-filled syringe 2 mL	Microspheres Color
40-120 μm	Orange	S120GH	-
100-300 μm	Yellow	S220GH	-
300-500 μm	Blue	S420GH	-
500-700 μm	Red	S620GH	-



#### EmboGold® Microspheres

40-120 μm	Orange	S120EG	red
100-300 μm	Yellow	S220EG	red
300-500 μm	Blue	S420EG	red
500-700 μm	Red	S620EG	red



#### Bearing nsPVA® Embolization Particles

Vial Size Range	Color Code	Vial Product Number
45-150 μm	Yellow	V100EP
150-250 μm	Purple	V200EP
250-355 μm	Blue	V300EP
355-500 μm	Green	V400EP
500-710 μm	Orange	V600EP



Understand. Innovate. Deliver.™

Merit Medical Systems, Inc.  
1600 West Merit Parkway  
South Jordan, Utah 84095  
1.801.208.4300  
1.800.35.MERIT

Merit Medical Europe, Middle East, & Africa (EMEA)  
Amerikalaan 42, 6199 AE Maastricht-Airport  
The Netherlands  
+31 43 358 82 22

Merit Medical Ireland Ltd.  
Parkmore Business Park West  
Galway, Ireland  
+353 (0) 91 703 733

merit.com

Austria  
0800 295 374

Belgium  
0800 72 906 (Dutch)  
0800 73 172 (Français)

Denmark  
80 88 00 24

Finland  
0800 770 586

France  
0800 91 60 30

Germany  
0800 182 0871

Ireland (Republic)  
1800 553 163

Italy  
800 897 005

Luxembourg  
8002 25 22

Netherlands  
0800 022 81 84

Norway  
800 11629

Portugal  
308 801 034

Russia  
+7 495 221 89 02

Spain  
+34 911238406

Sweden  
020 792 445

Switzerland  
(Deutsch)  
+41 22 518 02 30  
(Français)  
+41 22 518 02 52  
(Italiano)  
+41 22 518 00 35  
UK  
0800 973 115