



EN Snare®

Catalog Number	Loop Working Diameter	Snare Length	French Size	Catheter Length
----------------	-----------------------	--------------	-------------	-----------------

STANDARD SNARE KIT

EN2006010	6 - 10 mm	120 cm	6F	100 cm
EN2006015	9 - 15 mm	120 cm	6F	100 cm
EN2006020	12 - 20 mm	120 cm	6F	100 cm
EN2007030	18 - 30 mm	120 cm	7F	100 cm
EN2007045	27 - 45 mm	120 cm	7F	100 cm

MINI SNARE KIT

EN1003004	2 - 4 mm	175 cm	3.2F	150 cm
EN1003008	4 - 8 mm	175 cm	3.2F	150 cm



ONE Snare®

Catalog Number	Loop Diameter	Snare Length	French Size	Catheter Length	Catalog Number	Loop Diameter	Snare Length	French Size	Catheter Length
----------------	---------------	--------------	-------------	-----------------	----------------	---------------	--------------	-------------	-----------------

STANDARD SNARE KIT

ONE500	5 mm	120 cm	4F	100 cm
ONE1000	10 mm	120 cm	4F	100 cm
ONE1500	15 mm	120 cm	6F	100 cm
ONE2000	20 mm	120 cm	6F	100 cm
ONE2500	25 mm	120 cm	6F	100 cm
ONE3000	30 mm	120 cm	6F	100 cm
ONE3500	35 mm	120 cm	6F	100 cm

PETITE SNARE KIT

ONE1001	10 mm	65 cm	4F	48 cm
ONE2501	25 mm	65 cm	6F	48 cm

MICROSNARE KIT

ONE200	2 mm	175 cm	2.3-3F	150 cm
ONE201	2 mm	200 cm	2.3-3F	175 cm
ONE400	4 mm	175 cm	2.3-3F	150 cm
ONE401	4 mm	200 cm	2.3-3F	175 cm
ONE700	7 mm	175 cm	2.3-3F	150 cm
ONE701	7 mm	200 cm	2.3-3F	175 cm

SNARE REPLACEMENT CATHETERS

ONE4000	-	-	4F	100 cm
ONE6000	-	-	6F	100 cm