



SLX®

**Indications for Use:** The Medcomp® SLX Double Lumen Catheter can be utilized for temporary access for hemodialysis, hemoperfusion, or apheresis therapy. The cannula may be inserted via the Seldinger technique due to the inner Teflon stylet, increasing linear strength. The stylet is removed after insertion, leaving the soft silicone cannula in the body. The flexible silicone make-up conforms well to the vessel anatomy, resulting in higher patient tolerance during extended use.

**Contraindications:** The Subclavian Approach is NOT recommended for use with the Medcomp® SLX Double Lumen Subclavian-Femoral Catheter in Hemodialysis or Hemoperfusion Procedures used for the management of acute poisoning or other situations in which a ventilator might be used due to risk of traumatic pneumothorax posing a dangerous complication for the patient.

Refer to Instructions for Use provided with the product for complete instructions, warnings, precautions, and contraindications. Observe all instructions for use prior to using products. Failure to do so may result in patient complications.





DOUBLE "D" INTERNAL  
LUMEN DESIGN

**Silicone**  
Material

**Tapered Tip**  
Tip Design

**14F**  
French Size

**Straight**  
Configuration



14F LUMEN



TAPERED TIP



SILICONE MATERIAL





BASIC SET		
MC061402	14F X 15CM SLX CATHETER	5/BOX
MC061403	14F X 20CM SLX CATHETER	5/BOX
MC061404	14F X 24CM SLX CATHETER	5/BOX
MC061405	14F X 30CM SLX CATHETER	5/BOX

**BASIC SET CONTENTS:**

- (1) Catheter
- (1) Introducer Needle
- (1) Guidewire
- (2) Dilators
- (1) Scalpel
- (1) Peelable Introducer
- (1) Adhesive Wound Dressing
- (2) End Caps
- (1) Hemo-Cath® Hub Clip
- (1) 2-0 Suture

MAX BARRIER TRAY		
MC061422MB	14F X 15CM SLX CATHETER	4/BOX
MC061423MB	14F X 20CM SLX CATHETER	4/BOX
MC061424MB	14F X 24CM SLX CATHETER	4/BOX
MC061425MB	14F X 30CM SLX CATHETER	4/BOX

**MAX BARRIER TRAY CONTENTS:**

- (1) Catheter with Stylets
- (1) Introducer Needle
- (1) Guidewire
- (2) Dilators
- (1) Scalpel
- (1) Peelable Introducer
- (1) Adhesive Wound Dressing
- (2) End Caps
- (1) Hemo-Cath Hub Clip
- (3) 4" x 4" Gauze
- (3) 2" x 2" Gauze
- (2) Chloraprep
- (1) 5cc Syringe
- (2) Pre-Filled Saline Syringes
- (1) 2-0 Silk Suture
- (3) Safety Hypodermic Needles
- (1) Hemostat
- (1) Fenestrated Drape
- (1) 5cc Lidocaine
- (1) Filter Straw
- (1) Surgical Gloves
- (1) Needle Stick Pad
- (1) Steri-Strips
- (2) Face Mask
- (1) Lidocaine/Saline Stickers
- (1) Surgical Tape
- (1) XL Gown
- (1) Bouffant Cap
- (1) Surgical Mask



[Refer to the Table of Contents for Important Risk Information regarding this device.](#)