




Metacross[®]

PTA Balloon Dilatation Catheter

Crosstella[®]

PTA Balloon Dilatation Catheter

Crosterio[®]

PTA Balloon Dilatation Catheter

MAXIMIZE
YOUR OPTIONS
FOR **RX AND OTW**
PERFORMANCE

 **TERUMO**
INTERVENTIONAL
SYSTEMS

CHOOSE FROM A FULL LINE OF RX AND OTW OPTIONS

Designed for optimal performance

METACROSS® PTA Balloon Dilatation Catheter (0.035")

CROSSTELLA® PTA Balloon Dilatation Catheter (0.018")

CROSPERIO® PTA Balloon Dilatation Catheter (0.014")

Providing end-to-end performance with less effort and more options, our rapid exchange (RX) and over-the-wire (OTW) PTA balloon dilatation catheters let you choose the optimal solution for you and your patients.

- 3 guidewire sizes available (0.014", 0.018", and 0.035")¹⁻⁶
- Smaller balloon profiles help ease crossability^{7,a}
 - Supports re-crossability within the same procedure¹⁻⁶
- Well-designed shaft and stiff tapered core wire engineered for maximizing pushability^b
 - Effective transition of force from proximal to distal end



Balloon Catheter	Guidewire Compatibility (in)	Usable Shaft Length (cm)	
Metacross® RX	0.035	135	
Metacross® OTW	0.035	65/90/135	
Crosstella® RX	0.018	150	Prox
Crosstella® OTW	0.018	150	
Crosporio® RX	0.014	150	
Crosporio® OTW	0.014	65/150	

* Diameter range:
 METACROSS® RX/OTW Balloon Catheter (3-12 mm);
 CROSSTELLA® RX/OTW Balloon Catheter (2-6 mm);
 CROSPERIO® RX/OTW Balloon Catheter (1.5-4 mm).
^b RX design feature only.
^c METACROSS® RX/OTW Balloon Catheter illustrations shown.

PUSHING BOUNDARIES

Terumo Interventional Systems is **committed to your success** with innovative procedural solutions and ongoing support for your most challenging cases.

We are relentlessly seeking new ways to help you apply effective solutions and achieve **better outcomes for more patients.**



CHOOSE FROM A FULL LINE OF RX AND OTW OPTIONS

Designed for optimal performance

- METACROSS®** PTA Balloon Dilatation Catheter (0.035")
- CROSSTELLA®** PTA Balloon Dilatation Catheter (0.018")
- CROSPERIO®** PTA Balloon Dilatation Catheter (0.014")

Providing end-to-end performance with less effort and more options, our rapid exchange (RX) and over-the-wire (OTW) PTA balloon dilatation catheters let you choose the optimal solution for you and your patients.

- 3 guidewire sizes available (0.014", 0.018", and 0.035")¹⁻⁶
- Smaller balloon profiles help ease crossability^{7,a}
 - Supports re-crossability within the same procedure¹⁻⁶
- Well-designed shaft and stiff tapered core wire engineered for maximizing pushability^b
 - Effective transition of force from proximal to distal end



RAPID EXCHANGE: BETTER HANDLING BY DESIGN

Achieve optimal results with Rapid Exchange

- Enables faster, more convenient exchange of devices
- Supports the use of shorter guidewires
- Encourages procedural efficiency, reducing fluoroscopy time and radiation dose
- Optimizes control, allowing a single operator to mitigate lost wire position

Metacross® RX


PTA Balloon Dilatation Catheter

ENGINEERED FOR PERFORMANCE

Excellent crossability with innovative catheter design⁷

- Designed to provide best-in-class pushability⁷
- Innovative shaft provides robust strength constructed to cross the lesion⁷
- Smallest balloon profile among tested semi-compliant balloons⁷

METACROSS RX Balloon Catheter has excellent crossability compared to select competitors



Sample size: METACROSS RX Balloon Catheter n=10, competitors n=5; size used for comparison: 6 mm x 40 mm.

Balloon Catheter	Guidewire Compatibility (in)	Usable Shaft Length (cm)	Shaft Diameter		Proximal Tapered	Proximal Marker (mm)	Center (mm)	Distal Marker (mm)	Distal Tapered (mm)	Distal (mm)
Metacross® RX	0.035	135	Proximal 3.9 Fr (1.30 mm)	Distal 5.7 Fr (1.90 mm)	1.7	1.7	1.6	1.7	1.4	1.02
Metacross® OTW	0.035	65/90/135	5 Fr compatible: 5.4 Fr (1.80 mm) 6/7 Fr compatible: 5.7 Fr (1.90 mm)		1.7	1.7	1.6	1.7	1.4	1.07
Crosstella® RX	0.018	150	Proximal 3.4 Fr (1.12 mm) or 4.1 Fr (1.35 mm)	Distal 4 Fr (1.32 mm)	1.28	1.32	1.28	1.36	1.29	0.62
Crosstella® OTW	0.018	150	4 Fr compatible: 3.9 Fr (1.30 mm) 5 Fr compatible: 4.3 Fr (1.43 mm)		1.19	1.30	1.26	1.30	1.01	0.58
Crosperio® RX	0.014	150	Proximal 3.5 Fr (1.16 mm)	Distal 3.5 Fr (1.16 mm)	0.85	0.85	0.82	0.86	0.83	0.47
Crosperio® OTW	0.014	65/150	3.8 Fr (1.28 mm)		0.90	0.87	0.81	0.83	0.76	0.46

* Diameter range:
METACROSS® RX/OTW Balloon Catheter (3-12 mm);
CROSSTELLA® RX/OTW Balloon Catheter (2-6 mm);
CROSPERIO® RX/OTW Balloon Catheter (1.5-4 mm).
^a RX design feature only.
^c METACROSS® RX/OTW Balloon Catheter illustrations shown.

METACROSS® RX

CROSSTELLA® RX

CROSPERIO® RX

See inserts for ordering information.

Metacross®

PTA Balloon Dilatation Catheter

Crosstella®

PTA Balloon Dilatation Catheter

Crosperio®

PTA Balloon Dilatation Catheter



FIND OUT MORE  Phone: 800.862.4143  terumo.com  Fax: 800.411.5870

Indications

CROSPERIO® RX PTA Balloon Dilatation Catheter, CROSPERIO® OTW PTA Balloon Dilatation Catheter, CROSSTELLA® RX PTA Balloon Dilatation Catheter, CROSSTELLA® OTW PTA Balloon Dilatation Catheter, METACROSS® RX PTA Balloon Dilatation Catheter, and METACROSS® OTW PTA Balloon Dilatation Catheter are intended to dilate stenosis in the iliac, femoral, ilio-femoral, popliteal, infra-popliteal, and renal arteries, and for the treatment of obstructive lesions of native and synthetic arteriovenous dialysis fistulae. This device is also indicated for stent post-dilatation in the peripheral vasculature.

RX ONLY. Refer to the product labels and package insert for complete warnings, precautions, potential complications, and instructions for use.

References:

1. METACROSS® RX PTA Balloon Dilatation Catheter IFU, January 2016. **2.** METACROSS® OTW PTA Balloon Dilatation Catheter IFU, March 2016. **3.** CROSSTELLA® RX PTA Balloon Dilatation Catheter IFU, February 2016. **4.** CROSSTELLA® OTW PTA Balloon Dilatation Catheter IFU, May 2016. **5.** CROSPERIO® RX PTA Balloon Dilatation Catheter IFU, March 2016. **6.** CROSPERIO® OTW PTA Balloon Dilatation Catheter IFU, March 2016. **7.** Data on file.

©2021 Terumo Medical Corporation. All rights reserved. The third-party trademarks used herein are trademarks of their respective owners. All other brand names are trademarks or registered trademarks of Terumo. PM-04467

 **TERUMO**
INTERVENTIONAL
SYSTEMS

Metacross® RX

PTA Balloon Dilatation Catheter

ORDERING INFORMATION

METACROSS® RX 0.035" PTA Balloon Dilatation Catheter – 135 cm								
BALLOON DIAMETER (mm)	BALLOON LENGTH (mm)							
	20	40	60	80	100	120	150	200
3	BD-L30020CR	BD-L30040CR	BD-L30060CR	BD-L30080CR	BD-L30100CR	BD-L30120CR	BD-L30150CR	BD-L30200CR
4	BD-L40020CR	BD-L40040CR	BD-L40060CR	BD-L40080CR	BD-L40100CR	BD-L40120CR	BD-L40150CR	BD-L40200CR
5	BD-L50020CR	BD-L50040CR	BD-L50060CR	BD-L50080CR	BD-L50100CR	BD-L50120CR	BD-L50150CR	BD-L50200CR
6	BD-L60020CR	BD-L60040CR	BD-L60060CR	BD-L60080CR	BD-L60100CR	BD-L60120CR	BD-L60150CR	BD-L60200CR
7	BD-L70020CR	BD-L70040CR	BD-L70060CR	BD-L70080CR	BD-L70100CR	BD-L70120CR	BD-L70150CR	BD-L70200CR
8	BD-L80020CR	BD-L80040CR	BD-L80060CR	BD-L80080CR	—	—	—	—
9	BD-L90020CR	BD-L90040CR	BD-L90060CR	BD-L90080CR	—	—	—	—
10	BD-LX0020CR	BD-LX0040CR	BD-LX0060CR	BD-LX0080CR	—	—	—	—
12	BD-LX2020CR	BD-LX2040CR	—	—	—	—	—	—

6Fr

7Fr



Metacross® OTW

PTA Balloon Dilatation Catheter

ORDERING INFORMATION

METACROSS® OTW 0.035" PTA Balloon Dilatation Catheter – 65 cm								
BALLOON DIAMETER (mm)	BALLOON LENGTH (mm)							
	20	40	60	80	100	120	150	200
3	—	BD-L30040SW	BD-L30060SW	—	—	—	—	—
4	BD-L40020SW	BD-L40040SW	—	—	—	—	—	—
5	—	BD-L50040SW	BD-L50060SW	—	—	—	BD-L50150SW	—
6	BD-L60020SW	BD-L60040SW	BD-L60060SW	—	—	—	BD-L60150SW	—
7	—	BD-L70040SW	BD-L70060SW	BD-L70080SW	—	—	—	—
8	BD-L80020SW	BD-L80040SW	BD-L80060SW	BD-L80080SW	—	—	—	—
9	BD-L90020SW	BD-L90040SW	BD-L90060SW	—	—	—	—	—
10	—	BD-LX0040SW	BD-LX0060SW	BD-LX0080SW	—	—	—	—
12	—	BD-LX2040SW	—	—	—	—	—	—

5Fr

6Fr

7Fr

METACROSS® OTW 0.035" PTA Balloon Dilatation Catheter – 90 cm								
BALLOON DIAMETER (mm)	BALLOON LENGTH (mm)							
	20	40	60	80	100	120	150	200
3	BD-L30020MW	BD-L30040MW	BD-L30060MW	BD-L30080MW	BD-L30100MW	BD-L30120MW	BD-L30150MW	BD-L30200MW
4	BD-L40020MW	BD-L40040MW	BD-L40060MW	BD-L40080MW	BD-L40100MW	BD-L40120MW	BD-L40150MW	BD-L40200MW
5	BD-L50020MW	BD-L50040MW	BD-L50060MW	BD-L50080MW	BD-L50100MW	BD-L50120MW	BD-L50150MW	BD-L50200MW
6	BD-L60020MW	BD-L60040MW	BD-L60060MW	BD-L60080MW	BD-L60100MW	BD-L60120MW	BD-L60150MW	BD-L60200MW
7	BD-L70020MW	BD-L70040MW	BD-L70060MW	BD-L70080MW	BD-L70100MW	BD-L70120MW	BD-L70150MW	BD-L70200MW
8	BD-L80020MW	BD-L80040MW	BD-L80060MW	BD-L80080MW	—	—	—	—
9	BD-L90020MW	BD-L90040MW	BD-L90060MW	BD-L90080MW	—	—	—	—
10	BD-LX0020MW	BD-LX0040MW	BD-LX0060MW	BD-LX0080MW	—	—	—	—
12	BD-LX2020MW	BD-LX2040MW	—	—	—	—	—	—

5Fr

6Fr

7Fr

METACROSS® OTW 0.035" PTA Balloon Dilatation Catheter – 135 cm								
BALLOON DIAMETER (mm)	BALLOON LENGTH (mm)							
	20	40	60	80	100	120	150	200
3	BD-L30020CW	BD-L30040CW	BD-L30060CW	BD-L30080CW	BD-L30100CW	BD-L30120CW	BD-L30150CW	BD-L30200CW
4	BD-L40020CW	BD-L40040CW	BD-L40060CW	BD-L40080CW	BD-L40100CW	BD-L40120CW	BD-L40150CW	BD-L40200CW
5	BD-L50020CW	BD-L50040CW	BD-L50060CW	BD-L50080CW	BD-L50100CW	BD-L50120CW	BD-L50150CW	BD-L50200CW
6	BD-L60020CW	BD-L60040CW	BD-L60060CW	BD-L60080CW	BD-L60100CW	BD-L60120CW	BD-L60150CW	BD-L60200CW
7	BD-L70020CW	BD-L70040CW	BD-L70060CW	BD-L70080CW	BD-L70100CW	BD-L70120CW	BD-L70150CW	BD-L70200CW
8	BD-L80020CW	BD-L80040CW	BD-L80060CW	BD-L80080CW	—	—	—	—
9	BD-L90020CW	BD-L90040CW	BD-L90060CW	BD-L90080CW	—	—	—	—
10	BD-LX0020CW	BD-LX0040CW	BD-LX0060CW	BD-LX0080CW	—	—	—	—
12	BD-LX2020CW	BD-LX2040CW	—	—	—	—	—	—

5Fr

6Fr

7Fr

FIND OUT MORE ☎ Phone: 800.862.4143 🌐 terumo.com 📠 Fax: 800.411.5870

RX ONLY. Refer to the product labels and package insert for complete warnings, precautions, potential complications, and instructions for use.

©2021 Terumo Medical Corporation. All rights reserved. The third-party trademarks used herein are trademarks of their respective owners. All other brand names are trademarks or registered trademarks of Terumo. PM-04467

TERUMO
INTERVENTIONAL
SYSTEMS

Crosstella® RX


PTA Balloon Dilatation Catheter

ORDERING INFORMATION

CROSSTELLA® RX 0.018" PTA Balloon Dilatation Catheter – 150 cm								
BALLOON DIAMETER (mm)	BALLOON LENGTH (mm)							
	20	40	60	80	100	120	150	200
2	BD-F20020LR	BD-F20040LR	BD-F20060LR	BD-F20080LR	BD-F20100LR	BD-F20120LR	BD-F20150LR	BD-F20200LR
2.5	BD-F25020LR	BD-F25040LR	BD-F25060LR	BD-F25080LR	BD-F25100LR	BD-F25120LR	BD-F25150LR	BD-F25200LR
3	BD-F30020LR	BD-F30040LR	BD-F30060LR	BD-F30080LR	BD-F30100LR	BD-F30120LR	BD-F30150LR	BD-F30200LR
4	BD-F40020LR	BD-F40040LR	BD-F40060LR	BD-F40080LR	BD-F40100LR	BD-F40120LR	BD-F40150LR	BD-F40200LR
5	BD-F50020LR	BD-F50040LR	BD-F50060LR	BD-F50080LR	BD-F50100LR	BD-F50120LR	BD-F50150LR	BD-F50200LR
6	BD-F60020LR	BD-F60040LR	BD-F60060LR	BD-F60080LR	BD-F60100LR	BD-F60120LR	BD-F60150LR	BD-F60200LR

4Fr

5Fr



Crosstella® OTW

PTA Balloon Dilatation Catheter

ORDERING INFORMATION

CROSSTELLA® OTW 0.018" PTA Balloon Dilatation Catheter – 150 cm

BALLOON DIAMETER (mm)	BALLOON LENGTH (mm)							
	20	40	60	80	100	120	150	200
2	BD-F20020LW	BD-F20040LW	BD-F20060LW	BD-F20080LW	BD-F20100LW	BD-F20120LW	BD-F20150LW	BD-F20200LW
2.5	BD-F25020LW	BD-F25040LW	BD-F25060LW	BD-F25080LW	BD-F25100LW	BD-F25120LW	BD-F25150LW	BD-F25200LW
3	BD-F30020LW	BD-F30040LW	BD-F30060LW	BD-F30080LW	BD-F30100LW	BD-F30120LW	BD-F30150LW	BD-F30200LW
4	BD-F40020LW	BD-F40040LW	BD-F40060LW	BD-F40080LW	BD-F40100LW	BD-F40120LW	BD-F40150LW	BD-F40200LW
5	BD-F50020LW	BD-F50040LW	BD-F50060LW	BD-F50080LW	BD-F50100LW	BD-F50120LW	BD-F50150LW	BD-F50200LW
6	BD-F60020LW	BD-F60040LW	BD-F60060LW	BD-F60080LW	BD-F60100LW	BD-F60120LW	BD-F60150LW	BD-F60200LW

4Fr

5Fr



FIND OUT MORE  Phone: 800.862.4143  terumo.com  Fax: 800.411.5870

RX ONLY. Refer to the product labels and package insert for complete warnings, precautions, potential complications, and instructions for use.

©2021 Terumo Medical Corporation. All rights reserved. The third-party trademarks used herein are trademarks of their respective owners. All other brand names are trademarks or registered trademarks of Terumo. PM-04467

TERUMO
INTERVENTIONAL
SYSTEMS


Crosperio® RX

PTA Balloon Dilatation Catheter

ORDERING INFORMATION

CROSPERIO® RX 0.014" PTA Balloon Dilatation Catheter – 150 cm						
BALLOON DIAMETER (mm)	BALLOON LENGTH (mm)					
	20	40	80	120	150	200
1.5	BD-B15020LR	BD-B15040LR	BD-B15080LR	BD-B15120LR	—	—
2	—	BD-B20040LR	BD-B20080LR	BD-B20120LR	BD-B20150LR	BD-B20200LR
2.5	—	BD-B25040LR	BD-B25080LR	BD-B25120LR	BD-B25150LR	BD-B25200LR
3	—	BD-B30040LR	BD-B30080LR	BD-B30120LR	BD-B30150LR	BD-B30200LR
3.5	—	BD-B35040LR	BD-B35080LR	BD-B35120LR	BD-B35150LR	BD-B35200LR
4	—	BD-B40040LR	BD-B40080LR	BD-B40120LR	BD-B40150LR	BD-B40200LR

4Fr



Crosperio® OTW

PTA Balloon Dilatation Catheter

ORDERING INFORMATION

CROSPERIO® OTW 0.014" PTA Balloon Dilatation Catheter – 65 cm						
BALLOON DIAMETER (mm)	BALLOON LENGTH (mm)					
	20	40	80	120	150	200
1.5	—	—	—	—	—	—
2	—	—	BD-B20080SW	—	—	BD-B20200SW
2.5	—	—	BD-B25080SW	—	—	BD-B25200SW
3	—	—	BD-B30080SW	—	—	BD-B30200SW
3.5	—	—	BD-B35080SW	—	—	BD-B35200SW
4	—	—	BD-B40080SW	—	—	BD-B40200SW

4Fr

CROSPERIO® OTW 0.014" PTA Balloon Dilatation Catheter – 150 cm						
BALLOON DIAMETER (mm)	BALLOON LENGTH (mm)					
	20	40	80	120	150	200
1.5	BD-B15020LW	BD-B15040LW	BD-B15080LW	BD-B15120LW	—	—
2	—	BD-B20040LW	BD-B20080LW	BD-B20120LW	BD-B20150LW	BD-B20200LW
2.5	—	BD-B25040LW	BD-B25080LW	BD-B25120LW	BD-B25150LW	BD-B25200LW
3	—	BD-B30040LW	BD-B30080LW	BD-B30120LW	BD-B30150LW	BD-B30200LW
3.5	—	BD-B35040LW	BD-B35080LW	BD-B35120LW	BD-B35150LW	BD-B35200LW
4	—	BD-B40040LW	BD-B40080LW	BD-B40120LW	BD-B40150LW	BD-B40200LW

4Fr

FIND OUT MORE  Phone: 800.862.4143  terumo.com  Fax: 800.411.5870

RX ONLY. Refer to the product labels and package insert for complete warnings, precautions, potential complications, and instructions for use.

©2021 Terumo Medical Corporation. All rights reserved. The third-party trademarks used herein are trademarks of their respective owners. All other brand names are trademarks or registered trademarks of Terumo. PM-04467

